

Sterile closed system bags

The 4l sterile closed system bag intended for the collection of urine. It is primarily used post operatively following transurethral resection (TUR) in order to rinse the bladder.

The product is intended for single use only.



4L Irrigation Bag with Non-return Valve (NRV)

| REF* | Tube length cm | Capacity ml | Standard sample port (for use with needle) | KombiKon™ needle-free sample port | Tube clamp | Drip chamber | Integrated hanger | Bottom outlet/ cross valve | With pump |
|--------------|-------------------|----------------|---|---|---------------|-----------------|----------------------|----------------------------------|--------------|
| 48 010 1S-xx | 100 | 4,000 | Yes | – | Yes | Yes | Yes | Yes | – |
| 48 011 1S-xx | 127 | 4,000 | Yes | – | Yes | Yes | Yes | Yes | Yes |
| 48 110 1S-xx | 100 | 4,000 | – | Yes | Yes | Yes | Yes | Yes | – |
| 48 111 1S-xx | 127 | 4,000 | – | Yes | Yes | Yes | Yes | Yes | Yes |

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

| | |
|--------------------------------------|--|
| Contains latex? | No |
| Contains DEHP? | Yes |
| Front, back and NRV film/foil | Polyvinyl chloride |
| Tubes | Polyvinyl chloride |
| Front filter | Oleophobic ePTFE membrane and Polyester |
| Bottom outlet | Polypropylene and High density polyethylene |
| Drip chamber | Polyvinyl chloride, Acrylonitrile-butadiene-styrene and hydrophobic membrane |
| Pumpball | Soft - Polyvinyl chloride |
| Sample port/KombiKon™ | Polyvinyl chloride, Polypropylene, Thermoplastic elastomer |
| Cap | Polyethylene |
| Tube Clamp | Polypropylene |
| Hanger | Polypropylene and Acetal (POM) copolymer |

| | |
|---|--------------|
| CE-marking? | Yes (0086) |
| Product classification | Class I (EU) |
| Quality system certified to | ISO 13485 |
| Global Medical Device Nomenclature (GMDN) | 58921 |

| | |
|----------------------|----------------|
| Sterile? | Yes |
| Sterilized by | Ethylene oxide |
| Shelf life | 5 years |

4L Irrigation bag

Sterile closed system bags

Packaging:

4L Irrigation bag

| REF* | Shelf box | | Shipping carton | | | | MOQ, pcs. |
|--------------|-----------|----------------|-----------------|-----------------|---------------|------------------------|-----------|
| | Qty., pcs | Dimensions, mm | Qty., pcs | Dimensions, mm | Gross Wt., kg | Volume, m ³ | |
| 48 010 1S-xx | * | * | 20 | 582 x 250 x 234 | 4.4 | 0.034 | 840 |
| 48 011 1S-xx | * | * | 20 | 582 x 378 x 203 | 5.1 | 0.045 | 640 |
| 48 110 1S-xx | * | * | 20 | 582 x 250 x 234 | 4.4 | 0.034 | 840 |
| 48 111 1S-xx | * | * | 20 | 582 x 378 x 203 | 5.1 | 0.045 | 840 |

* = No shelf box

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

