



Bacterial & Viral Filter and Heat & Moisture Exchanger

Intended use:

Bacterial & Viral Filter and Heat & Moisture Exchanger provide filtration and moisture return required during anesthesia and critical care. The product is intended for single use only.



Product range:

	REF*	Item	Description	Filter Efficiency	Weight	Pressure Drop	Connectors
	93-764MM	Heat & Moisture Exchanger	Heat & Moisture Exchanger with Bacterial & Viral Filter, Luer lock sample port	>99.999%	27g	<1.5 cm H ₂ O @ 30 LPM	22F/15M – 22M/15F
	93-864MM	Bacterial & Viral Filter	Bacterial & Viral Filter	>99.999%	23g	<1.5 cm H ₂ O @ 30 LPM	22F/15M – 22M/15F

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains phthalates? No
 Heat & moisture exchanger Polyurethane
 Bacterial & viral filter. Propene polymer
 Connectors Propene polymer
 Luer lock sample port. Propene polymer

Certifications:

CE-marking? Yes (CE 0120)
 Product classification Class IIA (EU)
 Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? No
 Shelf life 5 years

Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Qty, pcs
93-764MM	1	10	50	350 x 300 x 200	1.94	1,600
93-864MM	1	10	50	350 x 300 x 200	1.69	1,600

EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-764MM	(01)06945099316011	(01)26945099316015	(01)46945099316019
93-864MM	(01)06945099316028	(01)26945099316022	(01)46945099316026

Bacterial & Viral Filter and Heat & Moisture Exchanger







Unomedical

Packaging material:

Primary package	Poly bag - Low density polyethylene / Insert - 80g paper
Secondary package	Poly bag - Low density polyethylene / Sticker - Coated paper
Tertiary package	Carton - Corrugated board / Sticker - Coated paper

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec

ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-018735-MM, ver. FEB 2018