

Leg Straps

Elasticated Velcro, Deluxe

Intended use:

Leg straps for securing of leg bags on patient's leg.


List of product range:

Leg straps	Size, mm	REF*
Deluxe leg strap	570	63 078 0-xx
Elasticated Velcro leg strap	570	63 077 4-xx

Material specifications:

Contains latex? No
Contains DEHP? No
Velcro leg strap Silicone, Polyamide, Polyester, Lycra
Deluxe strap Silicone, Polyamide, Polyester, Lycra

Certification:

CE-marking? Yes
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 37685

Sterilisation:

Sterile? No
Shelf life 5 years

Packaging:

REF*	Inner carton		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
63 077 4-xx	40	270 x 180 x 56	320	388 x 286 x 242	2.25	0.027	320
63 078 0-xx	na	na	150	375 x 280 x 278	4.05	0.029	900

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:


With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.