


Naso-Fix™

Securement device

Intended use:

Naso-Fix™ securement device is intended for the securement of oxygen catheters, duodenal and nasal tubes and to prevent accidental displacement of the catheter or tube from the nose. The product is intended for single use only.



List of product range:

	REF*	Size, cm
Naso-Fix™, Large	625M-L	8.0 x 8.7
Naso-Fix™, Small	625M-S	7.0 x 7.6
Naso-Fix™, Infant	625M-I	2.8 x 3.7

Material specifications:

Contains latex? No
 Contains DEHP? No
 Carrier Non-woven polypropylene
 Adhesive Medical grade acrylic
 Backing paper Release paper 67g/m² (one side)

Certification:

CE-marking? Yes
 Product classification Class I (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 36053

Sterilisation:

Sterile? No
 Shelf life 3 years

Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
Naso-Fix™, Large	625M-L	100	102 x 72 x 41	1,600	420 x 140 x 100	1.25	0.00588	1,600
Naso-Fix™, Small	625M-S	100	102 x 72 x 41	1,600	420 x 140 x 100	1.11	0.00588	1,600
Naso-Fix™, Infant	625M-I	100	102 x 72 x 41	1,600	420 x 140 x 100	0.75	0.00588	1,600

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.