

UnoMeter™ Safeti™ Plus System

Hourly Diuresis Monitoring System (with Change Bag and KombiKon™ Sampling Port)



Intended use:

Urine meter with chamber for accurate measuring of hourly diuresis.
Fitted with bottom outlet. Designed to prevent infection.

Standardisation:

CE-marking **CE 0086**
Quality system certified to . . **ISO 13485**
GMDN **32072**

List of product range:

Sterilisation:

Sterile **Yes**
Sterilised by **Ethylene Oxide**
Shelf life **5 years from the date of sterilisation**

REF	Description
158100710190	1.5L Change Bag, straps and round hooks, 150 cm, EasiMT™
158101110190	1.5L Change Bag, straps and round hooks, 110 cm, EasiMT™
158100610190	1.5L Change Bag, straps, 150 cm, EasiMT™
158101010190	1.5L Change Bag, straps, 110 cm, EasiMT™
30 810 1S - 90	1.5L Sterile single packed
47 810 1S - 90	3.5L Sterile single packed

Material:

Contains latex **No**
Contains paper **Yes**
Contains cotton **No**
UnoMeter™ Safeti™ Plus . . . **Polyvinyl Chloride, ABS, Polypropylene, Polyethylene, Styrene Acrylonitrile**
Airfilter **Teflon, non-woven Polyester backing**
KombiKon™ **Polyvinyl Chloride, Polypropylene, Polyethylene, TPE**

UnoMeter™ Safeti™ Plus System

Hourly Diuresis Monitoring System (with Change Bag and KombiKon™ Sampling Port)

UnoMeter™ Safeti™ Plus

Packaging:

REF	Shipping carton				MOQ pcs.
	Qty. Pcs.	Dimensions mm	Gross weight kg	Volume M ³	
158100710190	10	400x360x290	4.3	0.042	10
158101110190	10	400x360x290	4.3	0.042	10
158100610190	10	400x360x290	4.3	0.042	10
158101010190	10	400x360x290	4.3	0.042	10
30 810 1S - 90	100	375x280x180	3.08	0.319	100
47 810 1S - 90	100	375x280x180	5.62	0.029	100

Instruction for use:

Yes.....

Storage and handling:

