


A2

Basic drainage bags

Intended use:

Urine drainage bag for the collection and monitoring of urine output.
The product is intended for single use only.



List of product range:

Drainage bags with non-return valve (NRV)

REF*	Capacity, ml	Tube length, cm
34 001 0-xx	1,500	10
34 003 0-xx	1,500	35
34 005 0-xx	1,500	90
36 030 3-xx	2,000	90
36 105 0-xx	2,000	90

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Film Polyvinyl chloride
 Tube Polyvinyl chloride
 Connector Polyvinyl chloride
 Cap Low density polyethylene, High density polyethylene

Certification:

CE-marking? Yes
 Product classification Class I (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 58918

Sterilisation:

Sterile? No
 Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume m ³	
34 001 0-xx	500	375 x 280 x 278	9.0	0.030	24,000
34 003 0-xx	250	375 x 280 x 180	5.8	0.019	9,000
34 005 0-xx	250	388 x 286 x 242	8.0	0.027	14,000
36 030 3-xx	250	375 x 280 x 278	9.8	0.030	12,000
36 105 0-xx	250	388 x 286 x 242	8.7	0.027	14,000

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.