




3-in-1 Emergency Mask

Intended use:

Emergency masks can be used for delivery of oxygen in emergency medicine. They contain components to be used as high concentration non-rebreather, high concentration rebreather or medium concentration mask. The products are intended for single use only.



Product range:

	REF*	Item	Size	Description	Side valves	Tube length
 93-182MM	93-182MM	3-in-1 Emergency Mask Kit	Adult	Adult 3-in-1 Emergency Mask, Elongated oxygen mask with nose clip and elastic strap, 2 removable side valves, 1 removable check valve, Preattached reservoir bag, Oxygen adapter, 2.1m (7') No-crush oxygen tube, Fits-All connector	2	2.1m (7')
	93-182MMV (Vanilla scented mask)					
 93-282MM	93-282MM	3-in-1 Emergency Mask Kit	Pediatric	Pediatric 3-in-1 Emergency Mask, Elongated oxygen mask with nose clip and elastic strap, 2 removable side valves, 1 removable check valve, Preattached reservoir bag, Oxygen adapter, 2.1m (7') No-crush oxygen tube, Fits-All connector	2	2.1m (7')
	93-282MMV (Vanilla scented mask)					
 93-181MM	93-181MM	3-in-1 Emergency Mask Kit	Adult	Adult 3-in-1 Emergency Mask, Elongated oxygen mask with nose clip and elastic strap, 1 removable side valve, 1 removable check valve, Preattached reservoir bag, Oxygen adapter, 2.1m (7') No-crush oxygen tube, Fits-All connector	1	2.1m (7')
	93-281MM	3-in-1 Emergency Mask Kit	Pediatric	Pediatric 3-in-1 Emergency Mask, Elongated oxygen mask with nose clip and elastic strap, 1 removable side valve, 1 removable check valve, Preattached reservoir bag, Oxygen adapter, 2.1m (7') No-crush oxygen tube, Fits-All connector	1	2.1m (7')

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex? No
- Contains phthalates? No
- Contains bisphenols? No
- Oxygen nipple connector High density polyethylene
- Mask Polyvinyl chloride
- Elastic strap Polyester
- Nose clip Aluminum
- Side valves Silicone
- Check valve Silicone
- Oxygen connector on reservoir bag Polypropylene
- Reservoir bag Polyvinyl chloride
- Additional oxygen connector Polypropylene
- Oxygen tube Polyvinyl chloride
- Tube connector Polyvinyl chloride

Certifications:

- CE-marking? Yes (CE 0120)
- Product classification Class IIA (EU)
- Quality system certified to ISO 13485

Sterilization:

- Sterile? No
- Shelf life 5 years

Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-182MM	1	10	50	400 x 380 x 260	4.49	900
93-182MMV	1	10	50	400 x 380 x 260	4.49	900
93-282MM	1	10	50	400 x 380 x 260	3.69	900
93-282MMV	1	10	50	400 x 380 x 260	3.69	900
93-181MM	1	10	50	400 x 380 x 260	4.49	900
93-281MM	1	10	50	400 x 380 x 260	3.69	900





EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-182MM	(01)06945099301222	(01)26945099301226	(01)46945099301220
93-182MMV	(01)06945099301390	(01)26945099301394	(01)46945099301398
93-282MM	(01)06945099301239	(01)26945099301233	(01)46945099301237
93-282MMV	(01)06945099301406	(01)26945099301400	(01)46945099301404
93-181MM	(01)06945099306036	(01)26945099306030	(01)46945099306034
93-281MM	(01)06945099306043	(01)26945099306047	(01)46945099306041

Packaging material:

Primary package	Poly bag - Low density polyethylene / Insert - 80g paper
Secondary package	Poly bag - Low density polyethylene / Sticker - Coated paper
Tertiary package	Carton - Corrugated board / Sticker - Coated paper

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-019548-MM, ver. SEP 2018