


# PVC intermittent catheters

**Intended use:**

PVC intermittent catheters are intended to be used to provide an intermittent pathway for short term drainage of fluids from the bladder for either male (Nelaton and Tiemann) and female (Female) or for intravesical instillation. Products are designed for transient use only, i.e. max. 10-20 minutes. The product is intended for single use only.



**List of product range:**

CH size	Outer diameter, mm	Nelaton catheter	Tiemann catheter	Female catheter
		Length 360 mm	Length 360 mm	Length 147 mm
		REF*	REF*	REF*
06	2.0	0100118x	0300518x	0201218x
08	2.7	0100418x	0300618x	0201418x
10	3.3	0100718x	0300718x	0201518x
12	4.0	0100818x	0300818x	0201618x
14	4.7	0100918x	0300918x	0201718x
16	5.3	0101318x	0301118x	0201918x
18	6.0	0101518x	0301318x	0202018x
20	6.7	0101718x	0301418x	
22	7.4	0101918x		
24	8.0	0102018x		

\* Here noted as x. reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

Contains latex? . . . . . No  
 Contains DEHP? . . . . . Yes  
 Catheter . . . . . Polyvinyl chloride  
 Funnel . . . . . Polyvinyl chloride

**Certification:**

CE-marking? . . . . . Yes (CE 0086)  
 Product classification . . . . . Class I sterile (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . . 45603

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Ethylene oxide  
 Shelf life . . . . . 5 years

# PVC intermittent catheters

## Packaging:

REF*	Shelf box		Shipping carton				MOQ. pcs.
	Qty., pcs	Dimensions. mm	Qty., pcs	Dimensions. mm	Gross Wt., kg	Volume m <sup>3</sup>	
100118x	100	145 x 74 x 576	1000	291 x 304 x 585	6.7	0.072	2,000
100418x	100	145 x 74 x 576	1000	291 x 304 x 585	7.5	0.072	2,000
100718x	100	145 x 74 x 576	1000	391 x 304 x 585	7.5	0.072	2,000
100818x	100	145 x 93 x 446	1000	395 x 308 x 462	8.3	0.056	1,600
100918x	100	145 x 93 x 446	800	395 x 308 x 462	7.6	0.056	1,600
101318x	100	145 x 93 x 446	800	395 x 308 x 462	8.7	0.056	1,600
101518x	100	145 x 93 x 446	800	391 x 304 x 455	9.5	0.056	1,600
101718x	100	145 x 124 x 446	600	391 x 304 x 455	8.9	0.056	1,600
101918x	100	145 x 124 x 446	600	395 x 308 x 462	9.7	0.056	1,800
102018x	100	145 x 124 x 446	600	395 x 308 x 462	11.0	0.056	1,800
201218x	100	145 x 74 x 234	1000	395 x 308 x 249	4.2	0.031	3,000
201418x	100	145 x 74 x 234	1000	395 x 308 x 249	4.4	0.031	2,000
201518x	100	145 x 74 x 234	1000	395 x 308 x 249	4.5	0.031	2,000
201618x	100	145 x 93 x 234	800	395 x 308 x 249	4.6	0.031	1,600
201718x	100	145 x 74 x 234	1000	395 x 308 x 249	5.1	0.031	2,000
201918x	100	145 x 93 x 234	800	395 x 308 x 249	4.8	0.031	1,600
202018x	100	145 x 93 x 234	800	395 x 308 x 249	5.0	0.031	1,600
300518x	100	145 x 74 x 576	1000	391 x 304 x 585	6.3	0.072	2,000
300618x	100	145 x 74 x 576	1000	391 x 304 x 585	6.7	0.072	2,000
300718x	100	145 x 74 x 576	1000	391 x 304 x 585	7.5	0.072	2,000
300818x	100	145 x 93 x 446	800	395 x 308 x 462	7.6	0.056	1,600
300918x	100	145 x 93 x 446	800	395 x 308 x 462	8.3	0.056	1,600
301118x	100	145 x 93 x 446	800	395 x 308 x 462	8.7	0.056	1,600
301318x	100	145 x 93 x 446	800	395 x 308 x 462	9.5	0.056	1,600
301418x	100	145 x 124 x 446	600	395 x 308 x 462	8.9	0.056	1,200

\* Here noted as x. reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

