


OP-Flex™ : Sump Yankauer

Suction handles and sets

Intended use:

OP-Flex™ Sump Yankauer are suction handles and sets intended for intra operative removal of blood and fluids during open surgery. The Sump Yankauer is designed with a diffuser tip for delicate tissue. Must be attached to a canister via a tube before connecting to a vacuum source of max. 400 mmHg. The product is intended for single use only.



List of product range:

	REF*	Connecting tube Shore A 74		Inner/outer diameter of suction tip mm	Suction handle, length mm	Vacuum control	CH/FG (handle)	
		Length mm	CH/FG					
OP-Flex™ Sump Yankauer, diffuser tip, single bend	3408218x	2,100	24	4.6 / 6.5	155	No	20	Medium
	3408318x	3,500	24	4.6 / 6.5	155	No	20	Medium
	3304318x	-	-	4.6 / 6.5	155	No	20	Medium
	3407918x	2,100	24	4.6 / 6.5	155	Yes	20	Medium
	3408118x	3,500	24	4.6 / 6.5	155	Yes	20	Medium
	3304218x	-	-	4.6 / 6.5	155	Yes	20	Medium

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Handle Acrylonitrile-butadiene-styrene
 Sump tip Polyvinyl chloride
 Tubing Polyvinyl chloride
 Connectors Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device
 Nomenclature (GMDN) Handles: 35917/Sets: 36778

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

OP-Flex™ : Sump Yankauer

Suction handles and sets

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
3304218x	50	294 x 125 x 360	300	391 x 304 x 743	9.10	0.092	600
3304318x	50	294 x 125 x 360	300	391 x 304 x 743	9.10	0.092	600
3407918x	*	*	50	392 x 306 x 375	7.40	0.045	600
3408118x	*	*	50	392 x 306 x 375	8.50	0.045	150
3408218x	*	*	50	392 x 306 x 375	7.40	0.045	600
3408318x	*	*	50	392 x 306 x 375	8.50	0.045	150

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

