

Siliconized Foley catheter

Standard, 3-way, soft valve/hard valve

Intended use:

The siliconized Foley catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus.

The siliconized Foley catheter is intended for short-, medium- and long-term drainage of the bladder.

The product is intended for single use only.



List of product range:

Siliconized Foley catheter, Standard, 3-way, soft valve/hard valve

CH size	3-way SV 5-10 ml/cc balloon	3-way SV 30 ml/cc balloon	3-way HV 5-10 ml/cc balloon	3-way HV 30 ml/cc balloon
16	MM55111610	MM55111630	MM55121610	MM55121630
18	MM55111810	MM55111830	MM55121810	MM55121830
20	MM55112010	MM55112030	MM55122010	MM55122030
22	MM55112210	MM55112230	MM55122210	MM55122230
24	MM55112410	MM55112430	MM55122410	MM55122430
26	MM55112610	MM55112630	MM55122610	MM55122630

Material specifications:

Contains latex? Yes
 Contains DEHP (Hard valve)? Yes
 Contains DEHP (Soft valve)? No
 Catheter Latex
 Hard valve Polypropylene
 Soft valve Rubber
 Sleeve (Hard valve) Polyvinyl chloride
 Sleeve (Soft valve) Rubber
 Catheter surface Silicone emulsion

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Siliconized Foley catheter

Standard, 3-way, soft valve/hard valve

Packaging:

Siliconized Foley catheter, Standard, 3-way, soft valve/hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM55111610	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55111630	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55111810	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55111830	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112010	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112030	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112210	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112230	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112410	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112430	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112610	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55112630	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55121610	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55121630	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55121810	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55121830	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122010	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122030	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122210	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122230	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122410	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122430	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122610	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400
MM55122630	5	505 x 97 x 35	400	575 x 495 x 522	18.6	0.17	400

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

