

eakin **Cohesive**[®] seals

Absorbency
matters **most.**




EAKIN COHESIVE® SEALS FOR ABSORPTION

The damaging effect to the skin of prolonged exposure to moisture is a **common problem**^{1,2,3} experienced by up to **70%** of ostomates.^{4,5}

As it is not possible to control the source of excessive moisture, **the use of an absorptive ostomy seal could help prevent MASD.**⁶

In *in vitro* testing, eakin **Cohesive®** seals demonstrated **greater absorbency** of water than other hydrocolloid seals tested.⁷

Absorption with combined fecal enzyme (elastase and lipase) inhibition of ostomy seals.



eakin Cohesive® seals

Choose eakin
Cohesive® seals
to help create a
watertight
appliance

Reasons to use:

1. Smoothing out irregular skin and filling dips and creases^{7,8,9,10}
2. Encircling the stoma to prevent corrosive output from damaging the appliance flange^{9,10}
3. Adding protection for people with a high output ileostomy^{7,9,10}
4. Protecting against contact dermatitis^{7,8,9,10}
5. Helping to prevent repeated stoma problems^{9,10}

Case study submitted by:

Louise Tremblay, ET Nurse, Cité de la Santé, Canada

Patient History:

- 45 year old male
- History of Crohn's Disease
- Permanent Ileostomy since 1991



Source: TG Eakin LTD: Permission granted for use



Source: TG Eakin LTD: Permission granted for use



Source: TG Eakin LTD: Permission granted for use



Source: TG Eakin LTD: Permission granted for use

Nurse review:

- Stoma size varies between 25 - 28 mm, profile varies from flush to protruding up to 6 mm
- Abdomen is flat
- Hyper granulation
- Pouch weartime is 2 days

Problems:

- Output is very liquid and the pouch needs to be emptied 8-10 times/day
- Also experiencing burning skin

Nurse recommends:

- Recommended use of eakin **Cohesive SLIMS®**

Outcome:

- Skin has completely healed
- Patient no longer complains about burning skin
- Weartime has increased significantly to 5 days
- Patient has commented on the level of comfort that he has now
- The appliance is watertight with no more leaks, thanks to using eakin **Cohesive SLIMS®**

Clinical support for choosing eakin Cohesive® seals:

- Superior absorption⁷
- Helps to protect against enzymes^{7,10}
- Contributes to combatting irritants^{7,10}

A SEAL FOR EVERY SITUATION



eakin **Cohesive** SLIMS® *Your Everyday Seal*

- For everyday skin protection and security
- For frequent pouch changes

Product Code: **839005**



Large eakin **Cohesive**® seals *Greatest Coverage*

- For covering a large area of peristomal skin
- For complex skin problems or uneven skin surfaces

Product Code: **839001**



Small eakin **Cohesive**® seals *The Extra Protector*

- For maximum skin protection
- For those with high fluid content

Product Code: **839002**



eakin **Cohesive** StomaWrap™ *Your Time Saver*

- For large, oval or irregular shaped stomas
- For those with limited dexterity

Product Code: **839006**

REFERENCES

1. Voegeli, D., Moisture-associated skin damage: aetiology, prevention and treatment. Br J. Nurs. 2012;21(9):517-520
2. Voegeli, D., Moisture-associated skin damage: an overview for community nurses. Br J of Community Nursing 2013;18(1):6-10
3. Colwell C., et al, MASD Part 3: Peristomal Moisture Associated Dermatitis and Periwound Moisture-Associated Dermatitis. A Consensus J WOCN 2011;38(5):541-553
4. Colwell C., Goldberg M, Carmel J. The state of the standard diversion. J Wound Ostomy Continence Nurs. 2001;28:6-17. [PubMed]
5. Salvadoleña GD. Incidence of complications of the stoma and peristomal skin among individuals with colostomy, ileostomy and urostomy. J WOCN 2008;35(6):506-507
6. Gray M, Black JM, Baharestani MM et al. Moisture-associated skin damage: overview and pathophysiology. J Wound Ostomy Continence Nurs. 2011;38(3):233-41. <https://doi.org/10.1097/won.0b013e318215f798>
7. Mc Groggan G, Haughey S and McDowell K (2018) An absorbent, enzyme-inhibiting seal reduces peristomal skin complications. Gastrointestinal Nursing 16:1 42-49
8. Boyles A (2010) Keeping up to date with stoma accessories: enabling informed choice. Gastrointestinal Nursing 8:6 28-42
9. Evans S, Burch J (2017) An overview of stoma care accessory products for protecting peristomal skin. Gastrointestinal Nursing 15:7 25-34
10. O'Flynn K (2016) Protecting peristomal skin: a guide to conditions and treatment. Gastrointestinal Nursing 14:7 14-19



ConvaTec

ConvaTec Canada

100-1425, Route Transcanadienne,
Dorval, QC H9P 2V3

1-800-465-6302

convatec.canada@convatec.com

COH.LB.19 (CA.EN)

© 2019 ConvaTec Inc.

All trademarks are property of their respective owners.

Permission has been granted for use of photographs throughout.

AP-020580-CA

eakin® 