



Effective September 1, 2023

# Product Catalog 2023 2024

Advanced Wound Care  
Continence Care  
Critical Care  
Ostomy Care



1-800-465-6302 | [www.convatec.ca](http://www.convatec.ca)

**Connect your patients with a life they'll love living.**



Ron Nicholson  
Thriving with an ostomy since 2015

**The me+™ support program allows ostomy patients to:**

- ✓ Connect with products that are right for them
- ✓ Connect with ostomy nurses and product specialists
- ✓ Connect with peer stories for community support

**A good connection makes a healthy bond.**

**Scan to sign your patients up for the me+™ program**



# Convatec Customer Relations Center: What is it?

The Customer Relations Center (CRC) is Convatec's living communication hub of ostomy, wound, and continence and critical care expertise.

## For healthcare professionals:

- Clinical information and technical data on our products
- New product information
- Samples for your patients
- Patient literature on ostomy, wound, continence and critical care
- Advanced wound care expertise

## For retailers:

- Technical assistance and product questions
- Samples for your customers
- Contact with Territory Managers

## For people living with an ostomy:

- Free product samples
- Helpful literature
- Expert advice on how to care for their stoma

**The CRC team: Registered Nurses,  
NSWOC Nurse, and Product Specialists**

**1-800-465-6302**

Monday to Friday from 8:00 AM to 6:00 PM (EST)

Toll-free fax: 1-877-437-1777

You can also reach us by email at

[convatec.canada@convatec.com](mailto:convatec.canada@convatec.com)

or visit our website at [www.convatec.ca](http://www.convatec.ca)

## Table of Contents

<b>GENERAL INFORMATION (Terms and conditions of sale)</b> .....	6
<b>ORDERING INFORMATION</b> .....	7
<b>OSTOMY CARE ACCESSORIES</b>	
Diamonds® Gelling and Odor Control Sachets .....	8
AllKare® wipes .....	8
Niltac™ Sting-Free Adhesive Remover .....	8
Sillesse™ Sting-Free Skin Barrier .....	9
Esenta™ Sting-Free Adhesive Remover/Sting-Free Skin Barrier .....	8-9
Stomahesive® Skin Barrier, Powder and Paste .....	10
Eakin Cohesive® Paste .....	10
Stomahesive® Strips .....	10
Eakin Cohesive® Seals, StomaWrap™ and Skin Barriers .....	11
DuoLock® and ConvaTec® Tail Closures .....	11
Orabase® Oral Protective Paste .....	11
Ostomy Appliance Belt .....	12
ConvaTec® Night Drainage Container Set and Tubing .....	12
ConvaTec® Loop Ostomy Rod Package .....	12
Visi-Flow® .....	13
<b>OSTOMY CARE PRODUCTS</b>	
Natura® / Natura®+ Two-Piece Ostomy Systems .....	14
Esteem <i>synergy</i> ®/ Esteem <i>synergy</i> ®+ Ostomy Systems .....	41
Esteem®+ One-Piece Ostomy Systems .....	49
ActiveLife® One-Piece Ostomy Systems .....	61
<b>PEDIATRIC CARE PRODUCTS</b>	
Little Ones® Pediatric Ostomy Systems .....	67
<b>EAKIN WOUND POUCHES™</b>	
Eakin® Fistula and Wound Pouches .....	71
<b>CRITICAL CARE PRODUCTS</b>	
Flexi-Seal® Fecal Management Systems .....	78
<b>CONTINENCE CARE PRODUCTS</b>	
GentleCath™ Uncoated Intermittent Urinary Catheters .....	81
GentleCath™ Glide Hydrophilic Intermittent Urinary Catheters .....	82

# Table of Contents

## WOUND FILLERS AND WOUND HYDRATION

DuoDerm® Hydroactive™ Gel .....	86
---------------------------------	----

## SPECIALTY PRIMARY DRESSINGS

CarboFlex® Odor Control Dressing .....	86
--	----

## MOISTURE-RETENTIVE DRESSINGS

AQUACEL® Foam .....	87
---------------------	----

AQUACEL® Ag Foam .....	88
------------------------	----

AQUACEL® Foam PRO .....	88
-------------------------	----

FoamLite™ Convatec .....	89
--------------------------	----

DuoDerm® Extra Thin .....	90
---------------------------	----

DuoDerm Signal® .....	90
-----------------------	----

DuoDerm® CGF™ .....	91
---------------------	----

## EXUDATE MANAGEMENT DRESSINGS

AQUACEL® Extra™ .....	93
-----------------------	----

AQUACEL® Surgical .....	94
-------------------------	----

AQUACEL® Ag Surgical .....	94
----------------------------	----

KALTOSTAT® .....	95
------------------	----

## EXUDATE AND BIOFILM MANAGEMENT DRESSINGS

AQUACEL® Ag+ Extra™ .....	95
---------------------------	----

## NEGATIVE PRESSURE WOUND THERAPY

Avelle™ Negative Pressure Wound Therapy System .....	98
--	----

REFERENCES .....	100
------------------	-----

USEFUL RESOURCES .....	102
------------------------	-----

OSTOMY GLOSSARY OF TERMS .....	103
--------------------------------	-----

WOUND CARE GLOSSARY OF TERMS .....	104
------------------------------------	-----

PRODUCT CODE INDEX .....	105
--------------------------	-----

AUTHORIZATION TO RETURN MERCHANDISE FORM .....	111
--	-----

EXPLANATION OF PACKAGE CODE NUMBERING SYSTEM .....	114
--	-----

**Note:** In order to avoid misunderstanding and to protect the interest of our customers, we respectfully request careful reading of the "Terms and Conditions of Sale", as outlined. Non-compliance with these terms may result in suspension of your account.



## General Information

### TERMS AND CONDITIONS OF SALE

#### PRICES AND TERMS OF SALE

Prices and terms listed herein supersede those set forth on all previous price lists and all contracts and agreements that are expired. Due to market fluctuation, prices shown are subject to change without notice.

**RESTRICTION OF PURCHASE & SALES:** Purchase of Convatec products must be done through the Canadian order desk at 1-866-987-1471 or indirectly through a distributor who purchases directly from Convatec Canada. These products are not for distribution or resale outside of Canada or to anyone who - to the customers' knowledge - is susceptible to this practice.

#### PAYMENT TERMS

TERMS ARE NET 30 DAYS FROM DATE OF INVOICE.

Electronic Fund Transfers (EFT) are possible, please contact the Order Desk.

#### CHANGES IN PRICES

In all cases, merchandise will be invoiced at prices in effect on the date the order is received. Any pricing discrepancy MUST be reported within 60 days.

#### MINIMUM SHIPMENTS AND TRANSPORTATION COSTS

All orders amounting to a minimum of \$250 before tax for retail accounts and for health care facilities will be shipped prepaid via seller's choice of carrier. Any orders of value less than the required minimum will be subject to a \$50 service charge. Extra shipping costs incurred because of special services requested by a customer will be added to the invoice.

#### SHIPMENTS

Title to merchandise passes to the customer upon delivery of the goods to the customer.

#### SHORTAGE, DAMAGE OR LOSS IN SHIPMENT

Merchandise released for shipment leaves the depot in good order. Claims for damage, defacement or loss occurring in transit should be filed promptly with the transportation company. Customers should contact the Order Desk for processing of such claims and in supplying proof of delivery, providing these requests are received within 21 DAYS following the date of shipment. ANY COSTS INCURRED WILL BE PASSED ALONG TO THE CUSTOMER IF PROOF OF DELIVERY IS ESTABLISHED.

#### RETURNED GOODS POLICY

All returns MUST have written authorization from the Order Desk (telephone: 1-866-987-1471). This authorization will be confirmed by the representative and the credit will only be issued after inspection of the returned goods by Convatec. RETURNED GOODS WILL BE CREDITED AT ZERO IF NOT ACCOMPANIED BY A RETURN AUTHORIZATION NUMBER OR DOCUMENT. Order desk may ask for copy of purchase order for compliance purposes.

- All merchandise returned for credit must be prepaid and received in the original unopened containers. Partially full containers, those from which the seal, carton or label have been defaced or removed, will not be credited under any circumstances. Goods for return should be packaged carefully before sending to Convatec; products and containers received with any signs of damage will not be credited. Convatec will not ship back products that have been sent to us, as such, please ensure you respect the terms and conditions before sending products to Convatec.
- Products shipped in error by Convatec Canada will be fully credited if returned within 21 days of invoice, freight collect.
- Products ordered in error will be fully credited if returned within 21 days of invoice, freight prepaid.
- There will be a restocking service charge equal to 15% of the invoice price or the lowest price indicated for that product in the current price list.
- All returned goods must have been purchased in the last 12 months.
- All returned goods must have at least 6 months of shelf life remaining.
- All returned goods must be received by Convatec within 21 days after the written authorization was issued.
- Products that have been deleted from the line for a period of 6 months or more after date of deletion notification from Convatec, including special offers, will not be accepted for credit or exchange.
- Merchandise that has been involved in a fire, flood, sacrifice or bankruptcy sale or returned by other than original purchaser will not be accepted for credit or exchange.

#### CONSUMER SATISFACTION

We are quality driven and we want the consumers to be completely satisfied with our products. In the event of a product complaint from a consumer, the retailer should refer this consumer to the Convatec Customer Relations Center at **1-800-465-6302**.

#### TAXES

Provincial sales tax and/or GST, HST, QST will be charged, if applicable.

THE ABOVE TERMS CANNOT BE MODIFIED BY OUR REPRESENTATIVES.

# Ordering Information

PRODUCT ORDER DESK: TEL: 1-866-987-1471 - FAX: 1-866-987-1472  
TO REACH CONVATEC CUSTOMER RELATIONS CENTER:  
TEL: 1-800-465-6302  
FAX: 1-877-437-1777  
E-MAIL: [convatec.canada@convatec.com](mailto:convatec.canada@convatec.com)  
WEBSITE: [www.convatec.ca](http://www.convatec.ca)

## Merchandise returns

Returns of merchandise must be sent to:  
Convatec Canada  
RETURNS DEPARTMENT  
1375 Newton  
Boucherville, QC  
J4B 5H2

## Account payment

In order to expedite and facilitate the payment of your account, we maintain a Lock Box System. Please forward your cheque to the Post Office Box below. This is for cheques only.

### TORONTO

**Convatec Canada Ltd.**  
**PO Box 15000**  
**Station A**  
**Toronto, ON**  
**M5W 1C1**



1425, TransCanada, Suite 100 - Dorval, QC H9P 2V3 Tel: (514)-822-5900 Fax: (514)-822-5930  
[www.convatec.ca](http://www.convatec.ca)

## Ostomy Care Accessories



### Diamonds® Gelling and Odor Control Sachets

Improved Gelling performance. Provides a flatter, more discreet pouch.

Product Code	Box of
420791	100 Sachets per jar



### AllKare® Adhesive Remover Wipes

For the gentle removal of all adhesives, including skin barriers, tapes, and hydrocolloid dressings. The AllKare® Adhesive Remover Wipe formula effectively removes tape, skin barrier/wafer, and adhesive residues from the skin.

Product Code	Box of
37436	50
37443	100



### Niltac™ Sting-Free Adhesive Remover Wipes

Non-alcohol, silicone based formula that offers fast removal of adhesive residue while designed to minimize painful skin stripping on removal.<sup>6</sup>

Product Code	Box of
Non-Sterile 420788	30



### Niltac™ Sting-Free Adhesive Releaser Spray

Non-alcohol, silicone based formula that offers fast release of adhesives from the skin while designed to minimize painful skin stripping on removal.<sup>6</sup>

Product Code	Size	Box of
Non-Sterile 420787	50 ml	1

Coming soon



### Esenta™ Sting-Free Medical Adhesive Remover Wipes

Non-alcohol, silicone based formula that offers fast release of adhesives from the skin while designed to minimize painful skin stripping on removal.<sup>6</sup>

Product Code	Box of
Non-Sterile 423281	30

Coming soon



### Esenta™ Sting-Free Medical Adhesive Remover Spray

Non-alcohol, silicone based formula that offers fast removal of adhesive residue while designed to minimize painful skin stripping on removal.<sup>6</sup>

Product Code	Size	Box of
Non-Sterile 423289	50 ml	1



### AllKare® Protective Barrier Wipes

Provides a barrier film layer on the skin under tapes, skin barriers/wafers, adhesives, and hydrocolloid dressings to help protect against irritation and excoriation.

Product Code	Box of
37439	50
37444	100





**Sillesse™ Sting-Free Barrier Wipes**

Non-alcohol, silicone based formula for skin protection. Provides a barrier film layer on the skin that protects it from the damaging effects of adhesives, body waste and skin irritation caused by enzymatic stomal output.

Product Code	Box of
Non-Sterile 420789	30



**Sillesse™ Sting-Free Barrier Spray (non-aerosol)**

Non-alcohol, silicone based formula for skin protection. Provides a barrier film layer on the skin that protects it from the damaging effects of adhesives, body waste and skin irritation caused by enzymatic stomal output.

Product Code	Size	Box of
Non-Sterile 420790	50 ml	1

Coming soon



**Esenta™ Sting-Free Skin Barrier Wipes**

Non-alcohol, silicone based formula for skin protection. Provides a barrier film layer on the skin that protects it from the damaging effects of adhesives, body waste and skin irritation caused by enzymatic stomal output.

Product Code	Box of
Non-Sterile 423282	30

Coming soon



**Esenta™ Sting-Free Skin Barrier Spray**

Non-alcohol, silicone based formula for skin protection. Provides a barrier film layer on the skin that protects it from the damaging effects of adhesives, body waste and skin irritation caused by enzymatic stomal output.

Product Code	Size	Box of
Non-Sterile 423288	50 ml	1



**Stomahesive® Paste**

(Pectin-based Protective Skin Barrier)

Protective skin barrier and filler. Can lead to a longer wear time.

Product Code	Size	Box of
183910	Tube of 56.7 g (2 oz)	1 tube



**Eakin Cohesive® Paste**

(Alcohol-free)

Use to fill crevices, skin folds and scars to help create an even pouching surface.

This improves the likelihood of good adhesion and helps reduce the risk of leakage.

Product Code	Size	Box of
839010	Tube of 60 g (2.1 oz)	1 tube



**Stomahesive® Protective Powder**

This preparation absorbs excess moisture.

Product Code	Size	Box of
25510	Bottle of 28.3 g (1 oz)	1 bottle



**Stomahesive® Skin Barrier**

A multipurpose skin barrier designed to be skin friendly.

Product Code	Size	Box of
Non-Sterile 21712	approx. 10 cm X 10 cm (4" x 4")	5
Non-Sterile 21715	approx. 20 cm X 20 cm (8" x 8")	3



**Stomahesive® Strips**

(Moldable Adhesive)

Double sided Stomahesive® Strips fills uneven skin surfaces to create a secure platform for improved adhesion of ostomy appliances.

Product Code	Size	Box of
25542	approx. 15 mm X 120 mm (1" X 5")	15



### Eakin Cohesive® Seals

A donut-shaped ostomy accessory that can be molded to a variety of shapes. It acts as a gasket or barrier, and provides enhanced protection for hard-to-fit stomas.

Product Code	Size	Box of
839002	Small 5 cm (2")	20
839001	Large 10 cm (4")	10
839005	SLIMS™ 5 cm (2")	10



### Eakin Cohesive® StomaWrap™

Designed for large or oval shaped stomas and for people with limited dexterity.

Product Code	Size	Box of
839006	8.5 cm (3")	10



### Eakin Cohesive® Skin Barrier

A moldable and moisture-absorbing skin barrier. It is designed to help protect the skin from contact with body fluids.

Product Code	Size	Box of
839004	Small 10 cm X 10 cm (4" x 4")	5
839003	Large 10 cm X 20 cm (4" x 8")	5



### DuoLock® Curved Tail Closures

For use with Convatec Drainable Pouches

Product Code	Box of
175652	10



### ConvaTec® Tail Closures

For use with Convatec Drainable Pouches

Product Code	Box of
175635	10

## Other Accessories



### Orabase® Oral Protective Paste

For temporary relief of minor discomfort and protection of irritated areas of the mouth and gums.

Product Code	Size	Box of
62436	7.5 g tube	12 tubes

Note: Sold only in minimum quantities of one dozen.



**Ostomy Appliance Belt**

Adjustable up to 107 cm (42")

Product Code	Box of
175507	1



**ConvaTec® Night Drainage Container Set**

Set contains:

- 1 ConvaTec® Night Drainage Container with tubing
- 1 ConvaTec® Night Drainage Container Cover
- 1 universal adapter

Product Code	Box of
27060	1 set



**ConvaTec® Night Drainage Container Tubing**

Product Code	Box of
27062	1



**ConvaTec® Loop Ostomy Rod Package - Sterile**

For use with ConvaTec post-operative flange size 70 mm (2 3/4")

Product Code	Size	Box of
22355	65 mm	10

For use with ConvaTec post-operative flange size 100 mm (4")

Product Code	Size	Box of
22356	90 mm	10

## Visi-Flow® Irrigation Products



### Visi-Flow® Irrigator

Designed for individuals with a colostomy for whom irrigation is indicated. The Visi-Flow® Irrigator features a front-fill reusable bag with an exclusive “paddle wheel” flow indicator, thumb-control flow regulator, bold-print volume graduations and a flexible stoma cone attachment.





Product Code	Size	Box of
401989	2000 mL capacity (67 oz)	1



### Visi-Flow® Irrigation Sleeve

Transparent

Compatible with Natura® Skin Barriers. One tail closure per box.

Product Code	Color Match	Flange Size	Box of
401911		38 mm (1 1/2")	5
401912		45 mm (1 3/4")	5
401913		57 mm (2 1/4")	5
401914		70 mm (2 3/4")	5



### Visi-Flow® Stoma Cone

Product Code	Box of
22736	1



### Visi-Flow® Irrigation Sleeve Tail Closure

Product Code	Box of
175650	10



### Visi-Flow® Irrigation Adapter Faceplate

For clients using a one-piece system, this adaptor is required to perform irrigations.

Place the adaptor around the stoma, hold in place with the ostomy belt, so that the irrigation sleeve can be attached to perform the irrigation.

Product Code	Flange Size	Box of
401918	45 mm (1 3/4")	1
401919	70 mm (2 3/4")	1



Two-Piece Ostomy Systems - Natura®



# Natura® Two-Piece Ostomy Systems

An innovative line of pouches and skin barriers for every lifestyle

## Natura® Color Key and Sizing Guides

Match the color of the dot on the skin barrier package with the color of the dot on the pouch package.

Natura® Flange Sizes	
● 32 mm (1 1/4")	● 57 mm (2 1/4")
● 38 mm (1 1/2")	● 70 mm (2 3/4")
● 45 mm (1 3/4")	● 100 mm (4")

Natura® Moldable Skin Barrier Sizing (Non-Accordion Flange)			
STOMA SIZE	SKIN BARRIER SIZE	COLOR MATCH	FLANGE SIZE
13 mm - 22 mm (1/2" - 7/8")	Small	●	45 mm (1 3/4")
22 mm - 33 mm (7/8" - 1 5/16")	Medium	●	45 mm (1 3/4")
33 mm - 45 mm (1 5/16" - 1 3/4")	Large	●	57 mm (2 1/4")
Available in Flat only: 45 mm - 56 mm (1 3/4" - 2 1/8")	Extra-Large (with OVAL starter hole)	●	70 mm (2 3/4")

Natura® Moldable Skin Barrier with Accordion Flange Sizing		
STOMA SIZE	COLOR MATCH	FLANGE SIZE
13 mm - 22 mm (1/2" - 7/8")	●	57 mm (2 1/4")
22 mm - 33 mm (7/8" - 1 5/16")	●	57 mm (2 1/4")
33 mm - 45 mm (1 5/16" - 1 3/4")	●	70 mm (2 3/4")

Natura® Cut-To-Fit Skin Barrier Sizing (Non-Accordion Flange)		
CUTTING RANGE	COLOR MATCH	FLANGE SIZE
13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")
13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")
13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")
13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")
13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")
13 mm - 90 mm (1/2" - 3 1/2")	●	100 mm (4")

Natura® Cut-to-Fit Skin Barrier with Accordion Flange Sizing		
CUTTING RANGE	COLOR MATCH	FLANGE SIZE
13 mm - 21 mm (1/2" - 13/16")	●	45 mm (1 3/4")
13 mm - 33 mm (1/2" - 1 5/16")	●	57 mm (2 1/4")
13 mm - 45 mm (1/2" - 1 3/4")	●	70 mm (2 3/4")

Natura® Cut-to-Fit Convex Skin Barrier with Accordion Flange Sizing		
CUTTING RANGE	COLOR MATCH	FLANGE SIZE
13 mm - 22 mm (1/2" - 7/8")	●	45 mm (1 3/4")
13 mm - 35 mm (1/2" - 1 3/8")	●	57 mm (2 1/4")
13 mm - 48 mm (1/2" - 1 7/8")	●	70 mm (2 3/4")

# Natura® Two-Piece Products



The Natura® Two-Piece Ostomy System is Convatec's most extensive two-piece product line. The Natura® Flange System signals a secure, accurate closure with a series of audible "clicks." A simple snap-off, snap-on design allows users to change the pouch or pouch style quickly to match activities without removing the skin barrier from the skin.

When it comes to Natura® ostomy products, we recommend you take a look at the innovative pouch and skin barrier features that have been combined with our Stomahesive® and Durahesive® adhesives. The Natura® ostomy system offers special moldable skin barriers for a personalized, custom fit around the stoma to help protect the skin. Natura® pouches are available in drainable, closed-end, urostomy, transparent and opaque options. In addition, they are available with a filter and with the InvisiClose® clipless tail closure.



# Natura® Stomahesive® and Durahesive® Skin Barriers

(compatible with Natura® and Natura®+ pouches)

**Stomahesive® skin barriers** are ideal for people whose stoma output is semi-formed to formed stool. This reliable skin barrier provides a special formula that enables the skin barrier to adhere to both dry and moist skin. It holds firmly to skin while helping to protect against irritating output. Another feature of the Stomahesive® skin barrier is an outer protective coating which helps keep it in place. Plus, Stomahesive® skin barriers are easy to fit, apply, and remove.

**Durahesive® skin barriers** are designed for people whose stoma output is mostly liquid. When exposed to liquid output, Durahesive® skin barriers/wafers swell up to protect the stoma. This special effect, called turtlenecking, helps create a secure seal without harming the skin. Durahesive® skin barriers are also easy to fit, apply, and remove.



## Natura® Stomahesive® Skin Barrier with ConvaTec Moldable Technology™ and Accordion Flange

Hydrocolloid (HC) tape collar, tan

Product Code	Stoma Size	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")				
421033	small	●	57 mm (2 1/4")	10
	13 mm - 22 mm (1/2" - 7/8")			
421034	medium	●	57 mm (2 1/4")	10
	22 mm - 33 mm (7/8" - 1 5/16")			
13 cm x 13 cm (approx. 5.25" x 5.25")				
421035	large	●	70 mm (2 3/4")	10
	33 mm - 45 mm (1 5/16" - 1 3/4")			



## Natura® Durahesive® Skin Barrier with ConvaTec Moldable Technology™ and Accordion Flange

Hydrocolloid (HC) tape collar, tan

Product Code	Stoma Size	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")				
421039	small	●	57 mm (2 1/4")	10
	13 mm - 22 mm (1/2" - 7/8")			
421040	medium	●	57 mm (2 1/4")	10
	22 mm - 33 mm (7/8" - 1 5/16")			
13 cm x 13 cm (approx. 5.25" x 5.25")				
421041	large	●	70 mm (2 3/4")	10
	33 mm - 45 mm (1 5/16" - 1 3/4")			



**Natura® Stomahesive® Skin Barrier with ConvaTec Moldable Technology™**  
Hydrocolloid (HC) or Acrylic (AC) tape collar, white

Collar	Product Code	Stoma Size	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")					
HC	411801	small	●	45 mm (1 3/4")	10
AC	413421	13 mm - 22 mm (1/2" - 7/8")			
HC	411803	medium	●	45 mm (1 3/4")	10
AC	413422	22 mm - 33 mm (7/8" - 1 5/16")			
12.5 cm x 12.5 cm (approx. 5" x 5")					
HC	411805	large	●	57 mm (2 1/4")	10
AC	413423	33 mm - 45 mm (1 5/16" - 1 3/4")			
13 cm x 13 cm (approx. 5.25" x 5.25")					
HC	411807	extra-large	●	70 mm (2 3/4")	10
AC	413424	33 mm x 45 mm by 45 mm x 56 mm (1 5/16" x 1 3/4" by 1 3/4" x 2 1/8") Oval shape			




**Natura® Durahesive® Skin Barrier with ConvaTec Moldable Technology™**  
Hydrocolloid (HC) or Acrylic (AC) tape collar, white

Collar	Product Code	Stoma Size	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")					
HC	411800	small	●	45 mm (1 3/4")	10
AC	413417	13 mm - 22 mm (1/2" - 7/8")			
HC	411802	medium	●	45 mm (1 3/4")	10
AC	413418	22 mm - 33 mm (7/8" - 1 5/16")			
12.5 cm x 12.5 cm (approx. 5" x 5")					
HC	411804	large	●	57 mm (2 1/4")	10
AC	413419	33 mm - 45 mm (1 5/16" - 1 3/4")			
13 cm x 13 cm (approx. 5.25" x 5.25")					
HC	411806	extra-large	●	70 mm (2 3/4")	10
AC	413420	33 mm x 45 mm by 45 mm x 56 mm (1 5/16" x 1 3/4" by 1 3/4" x 2 1/8") Oval shape			



**Natura® Durahesive® Convex Skin Barrier with ConvaTec Moldable Technology™**  
Acrylic (AC) tape collar, white


Product Code	Stoma Size	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")				
404592	small	●	45 mm (1 3/4")	10
	13 mm - 22 mm (1/2" - 7/8")			
404593	medium	●	45 mm (1 3/4")	10
	22 mm - 33 mm (7/8" - 1 5/16")			
12.5 cm x 12.5 cm (approx. 5" x 5")				
404594	large	●	57 mm (2 1/4")	10
	33 mm - 45 mm (1 5/16" - 1 3/4")			





### Natura® Durahesive® Skin Barrier with Pre-Cut Opening, CONVEX-IT®

Acrylic (AC) tape collar, white




CONVEX

Product Code	Stoma Size	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")				
413177	13 mm (1/2")	●	45 mm (1 3/4")	10
413178	16 mm (5/8")	●	45 mm (1 3/4")	10
413179	19 mm (3/4")	●	45 mm (1 3/4")	10
413180	22 mm (7/8")	●	45 mm (1 3/4")	10
413181	25 mm (1")	●	45 mm (1 3/4")	10
413182	28 mm (1 1/8")	●	45 mm (1 3/4")	10
413183	32 mm (1 1/4")	●	45 mm (1 3/4")	10
413184	35 mm (1 3/8")	●	45 mm (1 3/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
413185	38 mm (1 1/2")	●	57 mm (2 1/4")	10
413186	41 mm (1 5/8")	●	57 mm (2 1/4")	10
413187	45 mm (1 3/4")	●	57 mm (2 1/4")	10
413188	50 mm (2")	●	57 mm (2 1/4")	10



### Natura® Durahesive® Convex Skin Barrier with Cut-to-Fit Opening and Accordion Flange

Hydrocolloid (HC) or Acrylic (AC) tape collar, tan



CONVEX

Collar	Product Code	Cutting Range	Color Match	Flange Size	Box of
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")					
AC	421638	13 mm - 22 mm	●	45 mm	10
HC	421639	(1/2" - 7/8")		(1 3/4")	
AC	421640	13 mm - 35 mm	●	57 mm	10
HC	421641	(1/2" - 1 3/8")		(2 1/4")	
13 cm x 13 cm (approx. 5.25" x 5.25")					
AC	421642	13 mm - 48 mm	●	70 mm	10
HC	421643	(1/2" - 1 7/8")		(2 3/4")	



### Natura® Stomahesive® Skin Barrier with Cut-to-Fit Opening and Accordion Flange

Acrylic (AC) tape collar, white

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
421456	small 13 mm - 21 mm (1/2" - 13/16")	●	45 mm (1 3/4")	10
421460	medium 13 mm - 33 mm (1/2" - 1 5/16")	●	57 mm (2 1/4")	10
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")				
421464	large 13 mm - 45 mm (1/2" - 1 3/4")	●	70 mm (2 3/4")	10



### Natura® Durahesive® Skin Barrier with Cut-to-Fit Opening and Accordion Flange

Acrylic (AC) tape collar, white

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
421454	small 13 mm - 21 mm (1/2" - 13/16")	●	45 mm (1 3/4")	10
421458	medium 13 mm - 33 mm (1/2" - 1 5/16")	●	57 mm (2 1/4")	10
11.5 cm x 11.5 cm (approx. 4.5" x 4.5")				
421462	large 13 mm - 45 mm (1/2" - 1 3/4")	●	70 mm (2 3/4")	10



### Natura® Stomahesive® Skin Barrier with Cut-to-Fit Opening

Acrylic (AC) tape collar, tan

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
125262	13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")	10
125263	13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")	10
125264	13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
125265	13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")	10
125266	13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")	10



### Natura® Durahesive® Flexible Skin Barrier with Cut-to-Fit Opening

Acrylic (AC) tape collar, tan

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
413164	13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")	10
413165	13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")	10
413166	13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
413167	13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")	10
413168	13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")	10



### Natura® Stomahesive® Flexible Skin Barrier with Cut-to-Fit Opening

Acrylic (AC) tape collar, white

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
125257	13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")	10
125258	13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")	10
125259	13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
125260	13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")	10
125261	13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")	10



### Natura® Durahesive® Flexible Skin Barrier with Cut-to-Fit Opening

Acrylic (AC) tape collar, white

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
413159	13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")	10
413160	13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")	10
413161	13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
413162	13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")	10
413163	13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")	10



### Natura® Stomahesive® Solid Skin Barrier with Cut-to-Fit Opening

No tape collar

Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
401573	13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")	10
401574	13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")	10
401575	13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")	10
401576	13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
401577	13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")	10
15 cm x 15 cm (approx. 6" x 6")				
401906	13 mm - 90 mm (1/2" - 3 1/2")	●	100 mm (4")	5



### Natura® Durahesive® Solid Skin Barrier with Cut-to-Fit Opening

No tape collar


Product Code	Cutting Range	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
413153	13 mm - 19 mm (1/2" - 3/4")	●	32 mm (1 1/4")	10
413154	13 mm - 25 mm (1/2" - 1")	●	38 mm (1 1/2")	10
413155	13 mm - 32 mm (1/2" - 1 1/4")	●	45 mm (1 3/4")	10
413156	13 mm - 45 mm (1/2" - 1 3/4")	●	57 mm (2 1/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
413157	13 mm - 57 mm (1/2" - 2 1/4")	●	70 mm (2 3/4")	10
15 cm x 15 cm (approx. 6" x 6")				
401905	13 mm - 90 mm (1/2" - 3 1/2")	●	100 mm (4")	5



## Natura® Stomahesive® Flexible Skin Barrier with Pre-Cut Opening

Acrylic (AC) tape collar, tan

Product Code	Stoma Size	Color Match	Flange Size	Box of
10 cm x 10 cm (approx. 4" x 4")				
125267	13 mm (1/2")	●	45 mm (1 3/4")	10
125268	16 mm (5/8")	●	45 mm (1 3/4")	10
125269	19 mm (3/4")	●	45 mm (1 3/4")	10
125270	22 mm (7/8")	●	45 mm (1 3/4")	10
125271	25 mm (1")	●	45 mm (1 3/4")	10
125272	28 mm (1 1/8")	●	45 mm (1 3/4")	10
125273	32 mm (1 1/4")	●	45 mm (1 3/4")	10
125274	35 mm (1 3/8")	●	45 mm (1 3/4")	10
12.5 cm x 12.5 cm (approx. 5" x 5")				
125275	38 mm (1 1/2")	●	57 mm (2 1/4")	10
125276	41 mm (1 5/8")	●	57 mm (2 1/4")	10
125277	45 mm (1 3/4")	●	57 mm (2 1/4")	10



Convex Insert Sizing			
STOMA SIZE NO WIDER THAN	COLOR MATCH	FLANGE SIZE (Outer Diameter)	CONVEX INSERT (Inner Diameter)
13 mm (1/2")	●	38 mm (1 1/2")	19 mm (3/4")
16 mm (5/8")	●	38 mm (1 1/2")	22 mm (7/8")
19 mm (3/4")	●	38 mm (1 1/2")	25 mm (1")
22 mm (7/8")	●	45 mm (1 3/4")	28 mm (1 1/8")
25 mm (1")	●	45 mm (1 3/4")	32 mm (1 1/4")
28 mm (1 1/8")	●	45 mm (1 3/4")	35 mm (1 3/8")
32 mm (1 1/4")	●	57 mm (2 1/4")	38 mm (1 1/2")
35 mm (1 3/8")	●	57 mm (2 1/4")	41 mm (1 5/8")

**Natura® Disposable Convex Inserts\***

Select an insert that is at least 7 mm (1/4") larger than the maximum diameter of the stoma.

Product Code	Stoma Size	Color Match	Box of
404006	For stoma size no wider than 13 mm (1/2") for use with 38 mm (1 1/2") Natura® wafers	●	5
404007	For stoma size no wider than 16 mm (5/8") for use with 38 mm (1 1/2") Natura® wafers	●	5
404008	For stoma size no wider than 19 mm (3/4") for use with 38 mm (1 1/2") Natura® wafers	●	5
404009	For stoma size no wider than 22 mm (7/8") for use with 45 mm (1 3/4") Natura® wafers	●	5
404010	For stoma size no wider than 25 mm (1") for use with 45 mm (1 3/4") Natura® wafers	●	5
404011	For stoma size no wider than 28 mm (1 1/8") for use with 45 mm (1 3/4") Natura® wafers	●	5
404012	For stoma size no wider than 32 mm (1 1/4") for use with 57 mm (2 1/4") Natura® wafers	●	5
404013	For stoma size no wider than 35 mm (1 3/8") for use with 57 mm (2 1/4") Natura® wafers	●	5

\* For Natura® Skin Barriers only. Not for use with Skin Barriers with Accordion Flange. Natura® Disposable Convex Inserts are not available for stoma sizes wider than 35 mm (1 3/8").



# Comfort. Confidence. Discretion.



\* Only on drainable pouches

+ pouch options can be found in these product lines:

NATURA®  ESTEEM SYNERGY®  ESTEEM® 



convatec

# THE ONLY THING THAT GETS CUT OUT IS THE WORRY.

No cutting or guessing means your patients say yes to everything else.



## ConvaTec Moldable Technology™ Skin Barrier

### ✔ LEAK-LOCK™

First-of-its-kind system provides a snug “turtleneck” effect for maximum leakage control and added security

### ✔ PATENTED REBOUNDED MEMORY TECHNOLOGY™

Customizable so that it molds to the stoma for a perfect fit every time

### ✔ STAYS ON

Barriers are made with our acclaimed adhesives, to keep them secure and give your patients peace of mind

Our moldable skin barriers are clinically proven to prevent peristomal skin irritation<sup>1</sup>

**95%** OF NEW USERS KEPT THEIR SKIN HEALTHY<sup>1</sup>



For more information, visit [convatec.ca](http://convatec.ca) or contact us at 1-800-465-6302 or [convatec.canada@convatec.com](mailto:convatec.canada@convatec.com)

Model portrayal. ©2021 ConvaTec Inc.  
All trademarks are property of their respective owners.  
AP-034394-MRL-CA

<sup>1</sup>Maria Teresa Szewczyk, MD, PhD; Grazyna Majewska, RN, ETN; Mary V. Cabral, MS, FNP-BC, CWOCN-AP; and Karin Holzel-Plontek, RN; The Effects of Using a Moldable Skin Barrier on Peristomal Skin Condition in Persons with an Ostomy: Results of a Prospective, Observational, Multinational Study, *Ostomy Wound Management* 2014;60(12):16-26.

Nothing protects like ConvaTec



**convatec**



## Natura® + Pouches

(Compatible with Natura® skin barriers)

Convatec developed these innovative pouches to be noticeably more comfortable and more discreet. The Natura®+ pouching system offers the latest technologies from Convatec including:



- Filter with anti-clogging film layer
- Soft and quiet materials
- InvisiClose® tail closure
- Lock-it Pocket™ that securely hides away the tail giving the comfort of a closed-end pouch

These pouches combined with Convatec advanced skin barrier adhesives offer security, comfort and discretion.

### Natura® + Drainable Pouch with InvisiClose® Tail Closure System, no Filter

Transparent with 1-sided comfort panel  
Pouch length: 30.5 cm (12")



Product Code	Flange Size	Color Match	Box of
416418	45 mm (1 3/4")	●	10
416421	57 mm (2 1/4")	●	10
416424	70 mm (2 3/4")	●	10
416472*	100 mm (4")	●	10

\*416472 pouch length is 35.6 cm (14").

### Natura® + Drainable Pouch with InvisiClose® Tail Closure System and Filter

Transparent with 1-sided comfort panel  
Pouch length: 30.5 cm (12"). 10 filter covers per box.



Product Code	Flange Size	Color Match	Box of
416416	45 mm (1 3/4")	●	10
416419	57 mm (2 1/4")	●	10
416422	70 mm (2 3/4")	●	10

### Natura® + Drainable Pouch with InvisiClose® Tail Closure System and Filter

Opaque with 2-sided comfort panel  
Pouch length: 30.5 cm (12"). 10 filter covers per box.



Product Code	Flange Size	Color Match	Box of
416415	38 mm (1 1/2")	●	10
416417	45 mm (1 3/4")	●	10
416420	57 mm (2 1/4")	●	10
416423	70 mm (2 3/4")	●	10



**Natura® + Drainable Pouch with InvisiClose® Tail Closure System, no Filter**

Opaque with 2-sided comfort panel

Pouch length: 30.5 cm (12")

Product Code	Flange Size	Color Match	Box of
421738	45 mm (1 3/4")	●	10
421739	57 mm (2 1/4")	●	10
421740	70 mm (2 3/4")	●	10



**Natura® + Drainable Pouch with InvisiClose® Tail Closure System, Window and Filter**

Opaque with 2-sided comfort panel

Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Flange Size	ColorMatch	Box of
421746	38 mm (1 1/2")	●	10
421747	45 mm (1 3/4")	●	10
421748	57 mm (2 1/4")	●	10
421749	70 mm (2 3/4")	●	10



### Natura®+ Closed-End Pouch with Filter

Opaque with 2-sided comfort panel

Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Flange Size	Color Match	Box of
416400	32 mm (1 1/4")	●	30
416403	38 mm (1 1/2")	●	30
416406	45 mm (1 3/4")	●	30
416409	57 mm (2 1/4")	●	30
416412	70 mm (2 3/4")	●	30



### Natura®+ Standard Closed-End Pouch with Window, no Filter

Opaque with 2-sided comfort panel

Pouch length: 20.8 cm (8.2")

Product Code	Flange Size	Color Match	Box of
421679	38 mm (1 1/2")	●	30
421681	45 mm (1 3/4")	●	30
421683	57 mm (2 1/4")	●	30
421893	70 mm (2 3/4")	●	30



### Natura®+ Closed-End Pouch with Window and Filter

Opaque with 2-sided comfort panel

Pouch length: 20.8 cm (8.2"). 30 filter covers per box.

Product Code	Flange Size	Color Match	Box of
421676	32 mm (1 1/4")	●	30
421678	38 mm (1 1/2")	●	30
421680	45 mm (1 3/4")	●	30
421682	57 mm (2 1/4")	●	30
421798	70 mm (2 3/4")	●	30



### Natura® + Mini Closed-End Pouch with Window and Filter

Opaque with 2-sided comfort panel

Pouch length: 16.2 cm (6.4"). 30 filter covers per box.

Product Code	Flange Size	Color Match	Box of
421796	45 mm (1 3/4")	●	30
421797	57 mm (2 1/4")	●	30



### Natura® + Urostomy Pouch with Soft-Tap Closure

Transparent with 1-sided comfort panel

Pouch length: 25.4 cm (10"). 2 Soft-Tap adapters per box.

Product Code	Flange Size	Color Match	Box of
413437	45 mm (1 3/4")	●	10
413438	57 mm (2 1/4")	●	10

# NATURA<sup>®</sup>



## Accordion Cut-to-fit Convex



**Lifts easily to provide generous finger room  
so coupling is easy and comfortable.**

- Comfortable and discreet
- Simple to measure and cut
- Compatible with a range of Natura<sup>®</sup> pouch options
- Available in Durahesive<sup>®</sup> with acrylic or hydrocolloid tape collars

## Natura® Drainable Pouches

(Compatible with Natura® skin barriers)

Natura® drainable pouches are fast drying. With streamlined body contour panels, low-profile design, and minimal headroom, they distribute pouch content evenly and away from the stoma. All drainable pouches have a wide tail opening for easy emptying.



### Natura® Drainable Pouch with Filter, Clipped Closure

Transparent with 1-sided comfort panel

Pouch length: 30.5 cm (12"). 2 tail closures and 10 filter covers per box.

Product Code	Flange Size	Color Match	Box of
411265	45 mm (1 3/4")	●	10
411266	57 mm (2 1/4")	●	10
411267	70 mm (2 3/4")	●	10



### Natura® Drainable Pouch with Filter, Clipped Closure

Opaque with 2-sided comfort panel

Pouch length: 30.5 cm (12"). 2 tail closures and 20 filter covers per box.

Product Code	Flange Size	Color Match	Box of
411490	38 mm (1 1/2")	●	20
411491	45 mm (1 3/4")	●	20
411492	57 mm (2 1/4")	●	20
411493	70 mm (2 3/4")	●	20



### Natura® Drainable Pouch with Filter, Clipped Closure

Opaque with 2-sided comfort panel

Pouch length: 25.4 cm (10"). 2 tail closures and 20 filter covers per box.

Product Code	Flange Size	Color Match	Box of
411289	38 mm (1 1/2")	●	20
411290	45 mm (1 3/4")	●	20
411291	57 mm (2 1/4")	●	20





### Natura® Drainable Pouch, Clipped Closure, no Filter

Transparent with 1-sided comfort panel  
Pouch length: 30.5 cm (12"). 2 tail closures per box.

Product Code	Flange Size	Color Match	Box of
401510	32 mm (1 1/4")	●	10
401511	38 mm (1 1/2")	●	10
401512	45 mm (1 3/4")	●	10
401513	57 mm (2 1/4")	●	10
401514	70 mm (2 3/4")	●	10



### Natura® Drainable Pouch, Clipped Closure, no Filter

Opaque with 2-sided comfort panel  
Pouch length: 30.5 cm (12"). 2 tail closures per box.

Product Code	Flange Size	Color Match	Box of
404025	32 mm (1 1/4")	●	20
404026	38 mm (1 1/2")	●	20
404027	45 mm (1 3/4")	●	20
404028	57 mm (2 1/4")	●	20
404029	70 mm (2 3/4")	●	20
401500	32 mm (1 1/4")	●	10
401501	38 mm (1 1/2")	●	10
401502	45 mm (1 3/4")	●	10
401503	57 mm (2 1/4")	●	10
401504	70 mm (2 3/4")	●	10



### Natura® Drainable Pouch, Clipped Closure, no Filter

Opaque with 1-sided comfort panel  
Pouch length: 30.5 cm (12"). 2 tail closures per box.

Product Code	Flange Size	Color Match	Box of
404014	32 mm (1 1/4")	●	20
404015	38 mm (1 1/2")	●	20
404016	45 mm (1 3/4")	●	20
404017	57 mm (2 1/4")	●	20
404018	70 mm (2 3/4")	●	20



**Natura® Drainable Pouch, Clipped Closure, no Filter**

Transparent with 1-sided comfort panel  
Pouch length: 35.5 cm (14"). 2 tail closures per box.

Product Code	Flange Size	Color Match	Box of
401515	100 mm (4")	●	10



**Natura® Drainable Pouch, Clipped Closure, no Filter**

Opaque with 2-sided comfort panel  
Pouch length: 25.4 cm (10"). 2 tail closures per box.

Product Code	Flange Size	Color Match	Box of
404030	32 mm (1 1/4")	●	20
404031	38 mm (1 1/2")	●	20
404032	45 mm (1 3/4")	●	20
404033	57 mm (2 1/4")	●	20
404034	70 mm (2 3/4")	●	20
401505	32 mm (1 1/4")	●	10
401506	38 mm (1 1/2")	●	10
401507	45 mm (1 3/4")	●	10
401508	57 mm (2 1/4")	●	10
401509	70 mm (2 3/4")	●	10



**Natura® High Output Drainable Pouch with Filter**

Opaque with 1-sided comfort panel, anti-reflux valve, replaceable filter, and outlet with spout and cap for high liquid output  
Pouch length: 35.5 cm (14"). 1 tail closure, 1 night drainage adapter and 5 replacement filters per box.

Product Code	Flange Size	Color Match	Box of
401557	45 mm (1 3/4")	●	5
401558	57 mm (2 1/4")	●	5
401559	70 mm (2 3/4")	●	5



**Natura® High Output Drainable Pouch with Filter**

Transparent with 1-sided comfort panel, anti-reflux valve, replaceable filter, and outlet with spout and cap for high liquid output  
Pouch length: 35.5 cm (14"). 1 tail closure, 1 night drainage adapter, and 5 replacement filters per box.

Product Code	Flange Size	Color Match	Box of
420695	45 mm (1 3/4")	●	5
420696	57 mm (2 1/4")	●	5
420697	70 mm (2 3/4")	●	5

## Natura® Closed-End Pouches

(Compatible with Natura® Skin Barriers)



Natura® Closed-End pouches are designed for one-time usage, and may be used everyday, for special occasions or intimate moments. There is no draining or uncomfortable clip to worry about. When you are ready, you simply remove the pouch, place it into a disposable bag, then discard it.



### Natura® Mini Pouch, no Filter

Opaque with 1-sided comfort panel

Pouch length: 12.5 cm (5")

Product Code	Flange Size	Color Match	Box of
401529	32 mm (1 1/4")	●	20
401530	38 mm (1 1/2")	●	20
401531	45 mm (1 3/4")	●	20
401532	57 mm (2 1/4")	●	20



### Natura® Flange Cap with Filter

Opaque with 1-sided comfort panel. 25 filter covers per box.

Product Code	Flange Size	Color Match	Box of
401909	45 mm (1 3/4")	●	25

# Natura® Urostomy Pouches

(Compatible with Natura® Skin Barriers)



Natura® urostomy pouches have an extended anti-reflux valve, comfort panels and a choice of two types of taps: the Accuseal® tap with valve and the bendable tap that is designed for easy manipulation and emptying.



## Natura® Urostomy Pouch with Bendable Tap

Transparent with 1-sided comfort panel

Pouch length: 23 cm (9"). 2 night drainage adapters per box.

Product Code	Flange Size	Color Match	Box of
401538	32 mm (1 1/4")	●	10
401539	38 mm (1 1/2")	●	10
401540	45 mm (1 3/4")	●	10
401541	57 mm (2 1/4")	●	10



## Natura® Urostomy Pouch with Bendable Tap

Transparent with 1-sided comfort panel

Pouch length: 25.4 cm (10"). 2 night drainage adapters per box.

Product Code	Flange Size	Color Match	Box of
401533	32 mm (1 1/4")	●	10
401534	38 mm (1 1/2")	●	10
401535	45 mm (1 3/4")	●	10
401536	57 mm (2 1/4")	●	10
401537	70 mm (2 3/4")	●	10



## Natura® Urostomy Pouch with Accuseal® Tap with Valve

Transparent with 1-sided comfort panel

Pouch length: 25.4 cm (10"). 2 Accuseal® night drainage adapters per box.

Product Code	Flange Size	Color Match	Box of
401542	32 mm (1 1/4")	●	10
401543	38 mm (1 1/2")	●	10
401544	45 mm (1 3/4")	●	10
401545	57 mm (2 1/4")	●	10
401546	70 mm (2 3/4")	●	10



### Natura® Urostomy Pouch with Accuseal® Tap with Valve

Opaque with 1-sided comfort panel

Pouch length: 23 cm (9"). 2 Accuseal® night drainage adapters per box.

Product Code	Flange Size	Color Match	Box of
401547	32 mm (1 1/4")	●	10
401548	38 mm (1 1/2")	●	10
401549	45 mm (1 3/4")	●	10
401550	57 mm (2 1/4")	●	10



### Natura® Urostomy Pouch with Accuseal® Tap with Valve

Opaque with 1-sided comfort panel

Pouch length: 25.4 cm (10"). 2 Accuseal® night drainage adapters per box.

Product Code	Flange Size	Color Match	Box of
401551	32 mm (1 1/4")	●	10
401552	38 mm (1 1/2")	●	10
401553	45 mm (1 3/4")	●	10
401554	57 mm (2 1/4")	●	10
401555	70 mm (2 3/4")	●	10

# Natura® Non-Sterile Post-Operative/Surgical Kits



**NEW** premier line of two-piece non-sterile post-operative/surgical kits - **Coming Soon in 2023**  
Simple, post-operative solutions for skin-friendly and secure product connections.

Featuring our innovative Convatec technologies:

- Unique, three-layer ConvaTec Moldable Technology™ skin barrier with Patented Rebounding Memory Technology™ molds to the stoma for a custom fit
- Durahesive® barrier suitable for people with liquid to semi-liquid outputs
- Accordion flange for low pressure support
- Natura®+ pouches designed with comfort in mind, including odour control and security features

Plus:

- 100 mm (4") skin barrier size available
- 5 individually labeled and packaged post-operative/surgical kits in a box
- Sterile Loop Rod component available separately



## Natura® Skin Barrier and Natura®+ Drainable Pouch Post-Operative/Surgical Kit - Non-Sterile

Durahesive® (Extended Wear) Solid Skin Barrier with Cut-to-Fit Opening, No Tape Collar, Drainable Pouch with InvisiClose® Tail Closure, Transparent with 1-sided comfort panel, no Filter  
Pouch Length: 30.5 cm (12")

Product Code	Flange Size	Cutting Range	Color Match	Box of
423556	57 mm (2 1/4")	13 - 45 mm (1/2" - 1 3/4")	●	5 Kits
423555	70 mm (2 3/4")	13 - 57 mm (1/2" - 2 1/4")	●	5 Kits
423554	100 mm (4")	13 - 90 mm (1/2" - 3 1/2")	●	5 Kits
Kit Components				
Kit	Skin Barrier in Kit		Pouch in Kit	
423556	413156		416421	
423555	413157		416424	
423554	401905		416472*	

\*416472 pouch length is 35.6 cm (14")



## Natura® Skin Barrier and Natura®+ Drainable Pouch Post-Operative/Surgical Kit - Non-Sterile

Durahesive® (Extended Wear) Skin Barrier with Cut-to-Fit Opening, Accordion Flange, Acrylic (AC) Tape Collar (White), Drainable Pouch with InvisiClose® Tail Closure, Transparent with 1-sided comfort panel, no Filter  
Pouch Length: 30.5 cm (12")

Product Code	Flange Size	Cutting Range	Color Match	Box of
423548	70 mm (2 3/4")	13 - 45 mm (1/2" - 1 3/4")	●	5 Kits
Kit Components				
Kit	Skin Barrier in Kit		Pouch in Kit	
423548	421462		416424	



### Natura® Skin Barrier and Natura®+ Drainable Pouch Post-Operative/Surgical Kit - Non-Sterile

Durahesive® (Extended Wear) Skin Barrier with ConvaTec Moldable Technology™, Hydrocolloid (HC) Tape Collar (White), Drainable Pouch with InvisiClose® Tail Closure, Transparent with 1-sided comfort panel, no Filter  
Pouch Length: 30.5 cm (12")

Product Code	Flange Size	Stoma Size	Color Match	Box of
423559	45 mm (1 3/4")	22 - 33 mm (7/8" - 1 5/16")	●	5 Kits
423558	57 mm (2 1/4")	33 - 45 mm (1 5/16" - 1 3/4")	●	5 Kits
423557	70 mm (2 3/4")	33 x 45 mm by 45 x 56 mm (1 5/16" x 1 3/4" by 1 3/4" x 2 1/8") Oval shape	●	5 Kits
Kit Components				
Kit	Skin Barrier in Kit		Pouch in Kit	
423559	411802		416418	
423558	411804		416421	
423557	411806		416424	



### Natura® Skin Barrier and Natura®+ Drainable Pouch Post-Operative/Surgical Kit - Non-Sterile

Durahesive® (Extended Wear) Skin Barrier with ConvaTec Moldable Technology™, Accordion Flange, Hydrocolloid (HC) Tape Collar (Tan), Drainable Pouch with InvisiClose® Tail Closure, Transparent with 1-sided comfort panel, no Filter  
Pouch Length: 30.5 cm (12")

Product Code	Flange Size	Stoma Size	Color Match	Box of
423549	70 mm (2 3/4")	33 - 45 mm (1 5/16" - 1 3/4")	●	5 Kits
Kit Components				
Kit	Skin Barrier in Kit		Pouch in Kit	
423549	421041		416424	



### Natura® Skin Barrier and Natura® Urostomy Pouch Post-Operative/Surgical Kit - Non-Sterile

Durahesive® (Extended Wear) Skin Barrier with Cut-to-Fit Opening, Accordion Flange, Acrylic (AC) Tape Collar (White), Urostomy Pouch with Accuseal® Tap with Valve, Transparent with 1-sided comfort panel  
Pouch Length: 25.4 cm (10")

Product Code	Flange Size	Cutting Range	Color Match	Box of
423551	57 mm (2 1/4")	13 - 33 mm (1/2" - 1 5/16")	●	5 Kits
423550	70 mm (2 3/4")	13 - 45 mm (1/2" - 1 3/4")	●	5 Kits
Kit Components				
Kit	Skin Barrier in Kit		Pouch in Kit	
423551	421458		401545	
423550	421462		401546	



### Natura® Skin Barrier and Natura® Urostomy Pouch Post-Operative/Surgical Kit - Non-Sterile

Durahesive® (Extended Wear) Skin Barrier with ConvaTec Moldable Technology™, Accordion Flange, Hydrocolloid (HC) Tape Collar (Tan), Urostomy Pouch with Accuseal® Tap with Valve, Transparent with 1-sided comfort panel  
 Pouch Length: 25.4 cm (10")

Product Code	Flange Size	Stoma Size	Color Match	Box of
423553	57 mm (2 1/4")	22 - 33 mm (7/8" - 1 5/16")	●	5 Kits
423552	70 mm (2 3/4")	33 - 45 mm (1 5/16" - 1 3/4")	●	5 Kits
Kit Components				
Kit	Skin Barrier in Kit		Pouch in Kit	
423553	421040		401545	
423552	421041		401546	





# Two-Piece Ostomy Systems




Esteem synergy®+, Convatec's lowest profile  
two-piece solution with Adhesive Coupling  
Technology™






## Color Key and Sizing Guides

Match the color of the dot on the skin barrier package with the color of the dot on the pouch package.





### Esteem synergy® Skin Barrier Sizing

 Small
 Medium
 Large
 Extra-Large


### Moldable Skin Barrier Sizing

Stoma Size	Color Match	Skin Barrier
13 mm - 22 mm (1/2" to 7/8")		Medium
22 mm - 32 mm (7/8" to 1 1/4")		Medium
32 mm - 45 mm (1 1/4" to 1 3/4")		Large

### Cut-To-Fit Skin Barrier Sizing

Cutting Range	Color Match	Skin Barrier
Up to 35 mm (1 3/8")		Small
Up to 48 mm (1 7/8")		Medium
Up to 61 mm (2 3/8")		Large
Up to 89 mm (3 1/2")		Extra-Large

# Esteem synergy® Two-Piece Ostomy Systems




The Esteem synergy® ostomy system utilizes an innovative Adhesive Coupling Technology™ that connects the pouch to the skin barrier and eliminates the need for plastic rings that must be snapped or locked together. Esteem synergy® ostomy system is Convatec's innovative product that provides the skin friendliness and versatility of a two-piece system, as well as the low profile, comfort, and flexibility of a conventional one-piece system.

The Esteem synergy® ostomy system has a flat, flexible adhesive foam ring on the pouch that securely adheres to a clear, flat, flexible plastic "Landing Zone" flange on the skin barrier.

# Esteem synergy® Skin Barriers

(Compatible with Esteem synergy®+ and Esteem synergy® pouches)




The landing zone flange on the skin barrier is square and is slightly larger in size compared to the round adhesive ring on the pouch.



## Esteem synergy® Durahesive® Convex Skin Barrier with ConvaTec Moldable Technology™

Acrylic (AC) tape collar, white



Product Code	Stoma Size	Color Match	Barrier Overall Dimensions	Box of
409268	small: 13 - 22 mm (1/2" - 7/8")	●	10 cm x 10 cm (4" x 4")	10
409269	medium: 22 - 33 mm (7/8" - 1 5/16")	●	11.5 cm x 11.5 cm (4.5" x 4.5")	10
409270	large: 33 - 45 mm (1 5/16" - 1 3/4")	●	12 cm x 12 cm (5" x 5")	10



## Esteem synergy® Stomahesive® Skin Barrier with ConvaTec Moldable Technology™

Hydrocolloid (HC) tape collar, white

Product Code	Stoma Size	Color Match	Barrier Overall Dimensions	Box of
411657	small: 13 - 22 mm (1/2" - 7/8")	●	10 cm x 10 cm (4" x 4")	10
411658	medium: 22 - 33 mm (7/8" - 1 5/16")	●	11.5 cm x 11.5 cm (4.5" x 4.5")	10
411659	large: 33 - 45 mm (1 5/16" - 1 3/4")	●	12 cm x 12 cm (5" x 5")	10



## Esteem synergy® Stomahesive® Skin Barrier with Cut-to-Fit Opening

Acrylic (AC) tape collar, white

Product Code	Cutting Range	Color Match	Barrier Overall Dimensions	Box of
405456	small: up to 35 mm (1 3/8")	●	10 cm x 10 cm (4" x 4")	10
405457	medium: up to 48 mm (1 7/8")	●	11.5 cm x 11.5 cm (4.5" x 4.5")	10
405458	large: up to 61 mm (2 3/8")	●	12 cm x 12 cm (5" x 5")	10

Note: Starter hole 13 mm (1/2")




### Esteem synergy® Stomahesive® Skin Barrier with Cut-to-Fit Opening

No tape collar

Product Code	Cutting Range	Color Match	Barrier Overall Dimensions	Box of
405468	medium: up to 48 mm (1 7/8")	●	11.5 cm x 11.5 cm (4.5" x 4.5")	10
403948	large: up to 61 mm (2 3/8")	●	12 cm x 12 cm (5" x 5")	10
405483	extra-large: up to 89 mm (3 1/2")	●	15.2 cm x 15.2 cm (6" x 6")	10

Note: Starter hole 13 mm (1/2")



### Esteem synergy® Stomahesive® Skin Barrier with Pre-Cut Opening

Acrylic (AC) tape collar, white

Product Code	Stoma Size	Color Match	Barrier Overall Dimensions	Box of
405474	small: 19 mm (3/4")	●	10 cm x 10 cm (4" x 4")	10
405475	small: 22 mm (7/8")	●	10 cm x 10 cm (4" x 4")	10
405476	small: 25 mm (1")	●	10 cm x 10 cm (4" x 4")	10
405477	small: 28 mm (1 1/8")	●	10 cm x 10 cm (4" x 4")	10
405478	small: 32 mm (1 1/4")	●	10 cm x 10 cm (4" x 4")	10
405479	small: 35 mm (1 3/8")	●	10 cm x 10 cm (4" x 4")	10
405480	medium: 38 mm (1 1/2")	●	11.5 cm x 11.5 cm (4.5" x 4.5")	10
405482	medium: 45 mm (1 3/4")	●	11.5 cm x 11.5 cm (4.5" x 4.5")	10

## Esteem synergy® + Pouches

(Compatible with Esteem synergy® Skin Barriers)



Convatec developed these innovative pouches to be noticeably more comfortable and more discreet than our previous pouches.

The Esteem synergy®+ pouching system includes:

- Filter with anti-clogging film layer
- Soft and quiet materials
- Streamlined rounded shape skin barriers
- InvisiClose® tail closure
- Lock-it Pocket™ that securely hides away the tail giving the comfort of a closed-end pouch

Convatec's pouches combined with Esteem synergy® advanced skin barrier adhesives offer security, comfort and discretion.

### Esteem synergy®+ Drainable Pouch with InvisiClose® Tail Closure System, no Filter

Transparent with 1-sided comfort panel

Pouch length: 30.5 cm (12")

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416900	Medium pouch 48 mm (1 7/8")	●	10
416903	Large pouch 61 mm (2 3/8")	●	10

### Esteem synergy®+ Drainable Pouch with InvisiClose® Tail Closure and Filter

Transparent with 1-sided comfort panel

Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416795	Small pouch 35 mm (1 3/8")	●	10
416798	Medium pouch 48 mm (1 7/8")	●	10
416901	Large pouch 61 mm (2 3/8")	●	10

Opaque with 2-sided comfort panel

Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416796	Small pouch 35 mm (1 3/8")	●	10
416799	Medium pouch 48 mm (1 7/8")	●	10
416902	Large pouch 61 mm (2 3/8")	●	10





### Esteem synergy® + Closed-End Pouch with Filter

Transparent with 1-sided comfort panel  
Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416774	Small pouch 35 mm (1 3/8")	●	30
416778	Medium pouch 48 mm (1 7/8")	●	30

Opaque with 2-sided comfort panel  
Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416775	Small pouch 35 mm (1 3/8")	●	30
416779	Medium pouch 48 mm (1 7/8")	●	30
416783	Large pouch 61 mm (2 3/8")	●	30

### Esteem synergy® + Closed-End Mini Pouch with Filter

Transparent with 1-sided comfort panel  
Pouch length: 15.2 cm (6"). 30 filter covers per box.

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416780	Medium pouch 48 mm (1 7/8")	●	30



Opaque with 2-sided comfort panel  
Pouch length: 15.2 cm (6"). 30 filter covers per box.

Product Code	Matches to Skin Barrier with:	Color Match	Box of
416781	Medium pouch 48 mm (1 7/8")	●	30

## Esteem synergy® Urostomy Pouches

(Compatible with Esteem synergy® Skin Barriers)

The Esteem synergy® Urostomy Pouch with Accuseal® tap with valve is designed for easy manipulation and rapid emptying. The built-in anti-reflux valve resists the back flow of urine. The pouch has a body-side comfort panel.



### Esteem synergy® Urostomy Pouch with Accuseal® Tap with Valve

Transparent with 1-sided comfort panel  
Pouch length: 26.2 cm (10.3"). 2 Accuseal® night drainage adapters per box.



Product Code	Matches to Skin Barrier with:	Color Match	Box of
405451	Small pouch 35 mm (1 3/8")	●	10

# For People On the Go



ESTEEM SYNERGY® 

## Esteem synergy®+, Convatec's lowest profile two-piece solution with Adhesive Coupling Technology™

The skin barrier and pouch adhere securely together via a soft, foam adhesive ring on the pouch.

### Product Highlights

- **Innovative Adhesive Coupling Technology™**—foam adhesive ring on the pouch adheres to a landing zone on the skin barrier.
- **Comfort and discretion**—with soft, quiet, fabric.
- **Flexible**—moves with you and conforms to body contours.



Try the Esteem synergy+ family of products today!  
With cut-to-fit, pre-cut and ConvaTec Moldable Technology™.





# One-Piece Ostomy Systems

Esteem<sup>®</sup>+, Convatec's low-profile one-piece solution that combines the skin barrier and pouch in a single unit





**convatec**

Ostomy Care

**Shaping Lives with Complete  
Ostomy Products, Care and Support**

**Esteem®+ | Natura®+**

One-piece and two-piece solutions with convexity



**Are your patients  
concerned about  
leakage?**

**Consider a Convatec Convex solution.**

Our convexity solutions are curved to fit a variety of body types and designed to help maintain a good seal around the stoma.

**To find out more, visit [convatec.ca](https://convatec.ca)**

# Esteem®+ One-Piece Ostomy System



Flexible and discreet, the Esteem®+ one-piece system is the all-in-one solution that combines the skin barrier and pouch in a single unit. Available in multiple adhesive technologies that gently adhere to skin and protect skin from breakdown.

The Esteem®+ one piece pouching system includes:

- Filter with an anti-clogging film layer
- InvisiClose® clipless tail closure with a strong hook-on-hook fastener
- The Lock-It Pocket™ which securely hides away the tail of the drainable pouch to provide the comfort of a closed end pouch
- Comfort panel pouch material that is soft and smooth against the skin

## Esteem®+ One-Piece Drainable Pouch with Stomahesive® ConvaTec Moldable Technology™ Skin Barrier and Filter

Transparent with 1-sided comfort panel and no tape collar  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.



Product Code	Stoma Size	Box of
413515	Medium 20 - 30 mm (13/16" - 1 3/16")	10
413516	Large 30 - 40 mm (1 3/16" - 1 9/16")	10

Opaque with 2-sided comfort panel and window, no tape collar.  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Stoma Size	Box of
413521	Medium 20 - 30 mm (13/16" - 1 3/16")	10
413522	Large 30 - 40 mm (1 3/16" - 1 9/16")	10
413523	Extra-Large 40 - 50 mm (1 9/16" - 1 15/16")	10

## Esteem®+ One-Piece Drainable Pouch with Durahesive Plus® Cut-to-Fit Skin Barrier, InvisiClose® Tail Closure and Filter

Opaque with 2-sided comfort panel, Acrylic (AC) tape collar, and belt tabs  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.



Product Code	Cutting Range	Box of
416975	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 64 mm (2 1/2")	10



**Esteem®+ One-Piece Drainable Pouch with Durahesive Plus® Cut-to-Fit Skin Barrier, InvisiClose® Tail Closure and Filter**

Transparent with 1-sided comfort panel, Acrylic (AC) tape collar, and belt tabs  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Cutting Range	Box of
416976	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 64 mm (2 1/2")	10



**Esteem®+ One-Piece Drainable Pouch with Clipped Closure, Modified Stomahesive®, Cut-to-Fit Skin Barrier and Filter**

Transparent with 1-sided comfort panel, no belt tabs  
 Pouch length: 30.5 cm (12"). 10 filter covers and 2 tail closures per box.

Product Code	Cutting Range	Box of
421848	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	10



**Esteem®+ One-Piece Drainable Pouch with Modified Stomahesive®, Cut-to-Fit Skin Barrier, InvisiClose® Tail Closure System, Window, no Filter**

Opaque with 2-sided comfort panel, no belt tabs  
 Pouch length: 30.5 cm (12")

Product Code	Cutting Range	Box of
421853	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	10



**Esteem®+ One-Piece Drainable Pouch with Modified Stomahesive®, Cut-to-Fit Skin Barrier, InvisiClose® Tail Closure System, Window and Filter**

Opaque with 2-sided comfort panel, no belt tabs  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Cutting Range	Box of
421858	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	10



**Esteem®+ One-Piece Drainable Pouch with Modified Stomahesive® Cut-to-Fit Skin Barrier and InvisiClose® Tail Closure, no Filter**

Transparent with 1-sided comfort panel and no tape collar  
Pouch length: 30.5 cm (12")

Product Code	Cutting Range	Box of
416721	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	10



**Esteem®+ One-Piece Drainable Pouch with Modified Stomahesive® Cut-to-Fit Skin Barrier, InvisiClose® Tail Closure and Filter**

Transparent with 1-sided comfort panel and no tape collar  
Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Cutting Range	Box of
416718	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	10

Opaque with 2-sided comfort panel and no tape collar  
Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Cutting Range	Box of
416719	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	10



**Esteem®+ One-Piece Drainable Pouch with Stomahesive® Cut-to-Fit Skin Barrier\* and InvisiClose® Tail Closure, no Filter**

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar  
Pouch length: 35.5 cm (14")

Product Code	Cutting Range	Box of
416908	Can be enlarged to accommodate stoma sizes up to 100 mm (4")	10

\*Double-layer adhesive (Stomahesive® body-side, Durahesive® top-side) which combines the skin friendly aspects of Stomahesive® with the turtlenecking protection of Durahesive®.



**Esteem® + One-Piece Drainable Pouch with Modified Stomahesive® Pre-Cut Skin Barrier, InvisiClose® Tail Closure and Filter**

Transparent with 1-sided comfort panel and no tape collar  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Stoma Size	Box of
416724	25 mm (1")	10
416728	30 mm (1 3/16")	10
416732	35 mm (1 3/8")	10

Opaque with 2-sided comfort panel and no tape collar  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.


Product Code	Stoma Size	Box of
416725	25 mm (1")	10
416729	30 mm (1 3/16")	10
416733	35 mm (1 3/8")	10
416737	40 mm (1 9/16")	10



**Esteem® + One-Piece Drainable Pouch with Durahesive® Pre-Cut Convex Skin Barrier, InvisiClose® Tail Closure and Filter**


Transparent with 1-sided comfort panel, Acrylic (AC) tape collar, and belt tabs  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Stoma Size	Box of
416738	19 mm (3/4")	10
416740	22 mm (7/8")	10
416742	25 mm (1")	10
416744	28 mm (1 1/8")	10
416746	32 mm (1 1/4")	10
416748	35 mm (1 3/8")	10
416750	38 mm (1 1/2")	10



Opaque with 2-sided comfort panel. Acrylic (AC) tape collar, and belt tabs  
 Pouch length: 30.5 cm (12"). 10 filter covers per box.

Product Code	Stoma Size	Box of
416739	19 mm (3/4")	10
416741	22 mm (7/8")	10
416743	25 mm (1")	10
416745	28 mm (1 1/8")	10
416747	32 mm (1 1/4")	10
416749	35 mm (1 3/8")	10
416751	38 mm (1 1/2")	10



**Esteem®+ One-Piece Closed-End Pouch with Stomahesive® ConvaTec Moldable Technology™ Skin Barrier and Filter**

Opaque with 2-sided comfort panel and no tape collar  
 Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Stoma Size	Box of
413509	Medium 20 - 30 mm (13/16" - 1 3/16")	30
413510	Large 30 - 40 mm (1 3/16" - 1 9/16")	30
413511	Extra Large 40 - 50 mm (1 9/16" - 1 15/16")	30



**Esteem®+ One-Piece Closed-End Pouch with Modified Stomahesive® Cut-to-Fit Skin Barrier and Filter**

Transparent with 1-sided comfort panel and no tape collar  
 Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Cutting Range	Box of
416700	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	30

Opaque with 2-sided comfort panel and no tape collar  
 Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Cutting Range	Box of
416701	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	30



**Esteem®+ One-Piece Closed-End Pouch with Modified Stomahesive® Cut-to-Fit Skin Barrier, Window and Filter**

Opaque with 2-sided comfort panel and no tape collar  
 Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Cutting Range	Box of
421821	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	30




**Esteem®+ One-Piece Closed-End Pouch with Modified Stomahesive®  
Pre-Cut Skin Barrier and Filter**

Opaque with 2-sided comfort panel and no tape collar  
Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Stoma Size	Box of
416704	25 mm (1")	30
416707	30 mm (1 3/16")	30
416710	35 mm (1 3/8")	30
416713	40 mm (1 9/16")	30
416715	50 mm (2")	30

Transparent with 1-sided comfort panel and no tape collar  
Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Stoma Size	Box of
416703	25 mm (1")	30
416706	30 mm (1 3/16")	30
416709	35 mm (1 3/8")	30
416712	40 mm (1 9/16")	30



**Esteem®+ One-Piece Closed-End Pouch with Modified Stomahesive®  
Pre-Cut Skin Barrier, Window and Filter**

Opaque with 2-sided comfort panel and no tape collar  
Pouch length: 20.3 cm (8"). 30 filter covers per box.

Product Code	Stoma Size	Box of
421687	25 mm (1")	30
421688	30 mm (1 3/16")	30
421689	35 mm (1 3/8")	30
421822	40 mm (1 9/16")	30
421823	50 mm (2")	30





**Esteem<sup>®</sup>+ One-Piece Mini Closed-End Pouch with Modified Stomahesive<sup>®</sup> Cut-to-Fit Skin Barrier, Window and Filter**

Opaque with 2-sided comfort panel and no tape collar  
 Pouch length: 15.2 cm (6"). 30 filter covers per box.

Product Code	Cutting Range	Box of
421816	Can be enlarged to accommodate stoma sizes 20 mm (3/4") to 70 mm (2 3/4")	30



**Esteem<sup>®</sup>+ One-Piece Closed-End Mini Pouch with Modified Stomahesive<sup>®</sup> Pre-Cut Skin Barrier and Filter**

Opaque with 2-sided comfort panel and no tape collar  
 Pouch length: 15.2 cm (6"). 30 filter covers per box.

Product Code	Stoma Size	Box of
416705	25 mm (1")	30



**Esteem<sup>®</sup>+ One-Piece Urostomy Pouch with Durahesive<sup>®</sup> Cut-to-Fit Skin Barrier and Accuseal<sup>®</sup> Tap with Valve**

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar. Includes belt tabs  
 Pouch length: 25.4 cm (10"). 2 Accuseal<sup>®</sup> night drainage adapters per box.

Product Code	Cutting Range	Box of
421631	13 mm - 45 mm (1/2" - 1 3/4")	10

# Post-Operative/Surgical Non-Sterile Kit



**Esteem®+ One-Piece Drainable Pouch with Stomahesive® Cut-to-Fit Skin Barrier\* and InvisiClose® Tail Closure Post-Operative/ Surgical Kit - Non-Sterile**

Transparent with 1-sided comfort panel and Acrylic tape collar (AC)  
 Pouch length: 35.5 cm (14")

Product Code	Cutting Range	Box of
416915	Can be enlarged to accommodate stoma sizes up to 100 mm (4")	5 individually wrapped pouches

\*Double-layer adhesive (Stomahesive® body-side, Durahesive® top-side) which combines the skin friendly aspects of Stomahesive® with the turtlenecking protection of Durahesive®.

# ESTEEM® Flex Convex


The Esteem® + Flex Convex is a one-piece system that fits the contours of the body, moving with you to give you the confidence to live life the way you want.

## Key benefits include:

- Moves with you - flexible convexity fits to body contours for comfort and easy positioning.
- Skin health and security - hydrocolloid skin barrier is designed to stay firmly in place, helping to minimize leakage.
- Comfort and discretion - with soft and quiet fabric, designed to allow for a secure and comfortable experience.


## Four Plateau Sizes

The Esteem® + Flex Convex system is available in 4 plateau sizes to address a range of stoma management issues.



## Esteem® + Flex Convex Drainable Pouch with Cut-to-Fit Hydrocolloid (HC) Skin Barrier, Safe Seal Clipless Closure and Filter

Opaque with 2-sided comfort panel with window, belt tabs, no tape collar.  
Pouch length: 28 cm (11"). Each box contains 10 filter covers and 1 tail closure.




CONVEX

Product Code	Plateau Size	Cutting Range	Box of
422159	V0 (Extra-Large)	20 mm - 57 mm (3/4" - 2 1/4")	10
421615	V1 (Large)	20 mm - 43 mm (3/4" - 1 2/3")	10
421619	V2 (Medium)	20 mm - 35 mm (3/4" - 1 3/8")	10
421622	V3 (Small)	20 mm - 25 mm (3/4" - 1")	10



## Esteem® + Flex Convex Urostomy Pouch with Cut-to-Fit Hydrocolloid (HC) Skin Barrier, Soft-Tap Closure

Opaque with 2-sided comfort panel with window, belt tabs, no tape collar.  
Pouch length: 20 cm (8"). Each box contains 2 night drainage adapters.



CONVEX

Product Code	Plateau Size	Cutting Range	Box of
421623	V1 (Large)	10 mm - 43 mm (3/8" - 1 2/3")	10
421626	V2 (Medium)	10 mm - 35 mm (3/8" - 1 3/8")	10
421628	V3 (Small)	10 mm - 25 mm (3/8" - 1")	10

# ESTEEM® Soft Convex


The flexibility of Esteem® + Soft Convex combined with a shallow convexity is designed to be gentle on the skin and deliver comfort and flexibility.

## Key benefits include:

- Peace of mind - Split fabric, easy-view window for accurate positioning and observation.
- Trust and reassurance - The filter helps minimize odours and ballooning.
- Simple, secure fastening - The Safe Seal clipless closure helps make cleaning and emptying easy.
- Support and security - Integral soft flexible belt loops provide added security.


## Three Plateau Sizes

The Esteem® + Soft Convex system is available in 3 plateau sizes to address a range of stoma management issues.



## Esteem® + Soft Convex Drainable Pouch with Cut-to-Fit Hydrocolloid (HC) Skin Barrier, Safe Seal Clipless Closure and Filter


Opaque with 2-sided comfort panel with window, belt tabs, no tape collar.  
Pouch length: 28 cm (11"). Each box contains 10 filter covers and 1 tail closure.



Product Code	Plateau Size	Cutting Range	Box of
422363	V1 (Large)	20 mm - 47 mm (3/4" - 1 7/8")	10
422364	V2 (Medium)	15 mm - 40 mm (3/5" - 1 3/5")	10
422365	V3 (Small)	10 mm - 28 mm (3/8" - 1 1/8")	10

## Esteem® + Soft Convex Urostomy Pouch with Cut-to-Fit Hydrocolloid (HC) Skin Barrier, Soft-Tap Closure

Opaque with 2-sided comfort panel with window, belt tabs, no tape collar.  
Pouch length: 20 cm (8"). Each box contains 2 night drainage adapters.



Product Code	Plateau Size	Cutting Range	Box of
422551	V1 (Large)	20 mm - 47 mm (3/4" - 1 7/8")	10
422557	V2 (Medium)	15 mm - 40 mm (3/5" - 1 3/5")	10
422562	V3 (Small)	10 mm - 28 mm (3/8" - 1 1/8")	10



# One-Piece Ostomy Systems

ActiveLife® – Combines the skin barrier  
and pouch in one unit

# ActiveLife® One-Piece Ostomy Systems




Lightweight, flexible, one-piece ActiveLife® ostomy system combines the skin barrier and pouch in one simple, easy-to-manage system. Pouches are available in transparent and opaque films.



**ActiveLife® Convex One-Piece Pre-Cut Drainable Pouch with Durahesive® Skin Barrier, no Filter, Clipped Closure**

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar. Includes belt tabs  
Pouch length: 30.5 cm (12"). 1 tail closure per box.



Product Code	Stoma Size	Box of
175777	19 mm (3/4")	5
175778	22 mm (7/8")	5
175779	25 mm (1")	5
175780	28 mm (1 1/8")	5
175781	32 mm (1 1/4")	5
175782	35 mm (1 3/8")	5
175783	38 mm (1 1/2")	5
175784	45 mm (1 3/4")	5



**ActiveLife® One-Piece Drainable Cut-to-Fit Pouch with Stomahesive® Skin Barrier, no Filter, Clipped Closure**

Transparent with 1-sided comfort panel  
Pouch length: 30.5 cm (12"). 1 tail closure per box.

Product Code	Cutting Range	Box of
125330	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 64 mm (2 1/2")	20
22771	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 64 mm (2 1/2")	10



**ActiveLife® One-Piece Drainable Cut-to-Fit Pouch with Durahesive Plus® Skin Barrier, no Filter, Clipped Closure**

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar  
Pouch length: 30.5 cm (12"). 1 tail closure per box.

Product Code	Cutting Range	Box of
400598	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 64 mm (2 1/2")	10

Opaque with 1-sided comfort panel and Acrylic (AC) tape collar  
Pouch length: 30.5 cm (12"). 1 tail closure per box.



Product Code	Cutting Range	Box of
400599	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 64 mm (2 1/2")	10



**ActiveLife® One-Piece Pre-Cut Drainable Pouch with Stomahesive® Skin Barrier, no Filter, Clipped Closure**

Opaque with 1-sided comfort panel and Acrylic (AC) tape collar  
 Pouch length: 25.4 cm (10"). 1 tail closure per box.

Product Code	Stoma Size	Box of
125331	19 mm (3/4")	20
125333	32 mm (1 1/4")	20
125334	38 mm (1 1/2")	20
125335	45 mm (1 3/4")	20
125336	50 mm (2")	20
22751	25 mm (1")	10



**ActiveLife® One-Piece Pre-Cut Drainable Pouch with Stomahesive® Skin Barrier, no Filter, Clipped Closure**

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar  
 Pouch length: 30.5 cm (12"). 1 tail closure per box.

Product Code	Stoma Size	Box of
125338	19 mm (3/4")	20
125339	25 mm (1")	20
125340	32 mm (1 1/4")	20
125341	38 mm (1 1/2")	20
125342	45 mm (1 3/4")	20
125343	50 mm (2")	20

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar  
 Pouch length: 30.5 cm (12"). 1 tail closure per box.

Product Code	Stoma Size	Box of
22764	19 mm (3/4")	10
22765	25 mm (1")	10
22766	32 mm (1 1/4")	10
22767	38 mm (1 1/2")	10
22768	45 mm (1 3/4")	10
22769	50 mm (2")	10





**ActiveLife® One-Piece Cut-to-Fit Drainable Pouch with Stomahesive® Skin Barrier, no Filter, Clipped Closure**

Transparent with 1-sided comfort panel and no tape collar  
Pouch length: 35.5 cm (14") 10 Duolock® tail closures per box

Product Code	Stoma Size Opening	Box of
650803	up to 100 mm (4")	10



**ActiveLife® One-Piece Stoma Cap with Stomahesive® Skin Barrier and Filter**

Opaque with no tape collar. 30 filter covers per box.

Product Code	Cutting Range	Box of
175611	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 50 mm (2")	30



**ActiveLife® Convex One-Piece Pre-Cut Urostomy Pouch with Durahesive® Skin Barrier**

Transparent with bendable tap and Acrylic (AC) tape collar. Includes belt tabs  
Pouch length: 23 cm (9"). Each box contains one universal night drainage adapter.

Product Code	Stoma Size	Box of
125361	13 mm (1/2")	10
125362	16 mm (5/8")	10
125363	19 mm (3/4")	10
125364	22 mm (7/8")	10
125365	25 mm (1")	10
125366	28 mm (1 1/8")	10
125367	32 mm (1 1/4")	10
125368	35 mm (1 3/8")	10
125369	38 mm (1 1/2")	10





**ActiveLife® One-Piece Pre-Cut Urostomy Pouch with Durahesive® Skin Barrier and Accuseal® Tap with Valve**

Transparent with 1-sided comfort panel and Acrylic (AC) tape collar. Includes belt tabs  
 Pouch length: 20.3 cm (8"). 2 Accuseal® night drainage adapters per box.

Product Code	Stoma Size	Box of
650828	19 mm (3/4")	10
650829	22 mm (7/8")	10
650830	25 mm (1")	10
650831	28 mm (1 1/8")	10
650832	32 mm (1 1/4")	10
650833	38 mm (1 1/2")	10



**ActiveLife® One-Piece Cut-to-Fit Urostomy Pouch with Stomahesive® Skin Barrier with Bendable Tap**

Transparent with 1-sided comfort panel and no tape collar. Includes belt tabs  
 Pouch length: 23 cm (9"). 1 night drainage adapter per box.

Product Code	Cutting Range	Box of
64927	Can be enlarged to accommodate stoma sizes 19 mm (3/4") to 45 mm (1 3/4")	10



# Little Ones® Pediatric Ostomy Systems

A product line designed to  
care for the youngest patients



# Little Ones® Pediatric Ostomy Systems



With all the challenges that go with raising a child, finding the right ostomy product needn't be one of them. That's because Convatec has a line of products designed just for kids. The line features products for children who have had colostomy, ileostomy, or urostomy surgery. As with our adult line, Convatec Little Ones® Ostomy products offer a full range of products to meet a child's need or parent's preference. Each has been scaled down in size to comfortably fit a child's unique body contour. Little Ones® products are designed with an active child in mind.

## Color Key and Sizing Guides

Match the color of the dot on the skin barrier package with the color of the dot on the pouch package.

Little Ones® Pediatric Ostomy System	
●	Stoma Size: 0-23 mm (0 - 9/10")
●	Stoma Size: 4.8-31 mm (1/5 - 1 1/4")
●	Flange Size: 32 mm (1 1/4")
●	Flange Size: 45 mm (1 3/4")



**Little Ones® One-Piece Extra Small Drainable Pouch with Cut-to-Fit Stomahesive® Skin Barrier and InvisiClose® Tail Closure System, no Filter**

Transparent with 1-sided comfort panel  
Pouch length: 11 cm (4")

Product Code	Cutting Range	Color Match	Box of
411633	0-23 mm (0-9/10")	●	10



**Little Ones® One-Piece Drainable Pouch with Cut-to-Fit Stomahesive® Skin Barrier, no Filter**

Transparent with 1-sided comfort panel and no tape collar  
Pouch length: 15 cm (6"). 15 soft, moldable tail closures per box.

Product Code	Cutting Range	Box of
20922	8-50 mm (5/16" - 2")	15



**Little Ones® One-Piece Urostomy Pouch with Cut-to-Fit Stomahesive® Skin Barrier**

Transparent and no tape collar  
Pouch length: 15 cm (6"). Each box contains two night drainage adapters.

Product Code	Cutting Range	Box of
20917	8-25 mm (5/16" - 1")	15



**Little Ones® Two-Piece Adhesive Coupling Technology™ Standard Flexible with Cut-to-Fit Stomahesive® Skin Barrier**

Hydrocolloid (HC) tape collar, white  
Overall dimensions: 8.5 cm x 8.5 cm (3 3/8" x 3 3/8")

Product Code	Cutting Range	Color Match	Box of
411642	4.8-31 mm (3/16" - 1 1/4")	●	5



**Little Ones® Two-Piece Adhesive Coupling Technology™ Standard Closed-End Pouch with Filter**

Pouch length: 15.2 cm (6"). 10 filter covers per box.

Product Code	Stoma Size	Color Match	Box of
Opaque 411639	4.8-31 mm (3/16" - 1 1/4")	●	10



**Little Ones® Two-Piece Adhesive Coupling Technology™ Standard Drainable Pouch with InvisiClose® Tail Closure System and Filter**  
 Pouch length: 20.3 cm (8"). 10 filter covers per box.

Product Code	Stoma Size	Color Match	Box of
Transparent 411637	4.8-31 mm (3/16" - 1 1/4")		10



**Little Ones® Two-Piece Standard Flexible Cut-to-Fit Stomahesive® Skin Barrier**  
 Acrylic (AC) tape collar, white  
 Overall dimensions: 8 cm x 8 cm (3 1/4" x 3 1/4")

Product Code	Cutting Range	Color Match	Flange Size	Box of
401925	12.7-19 mm (1/2" - 3/4")		32 mm (1 1/4")	5



**Little Ones® Two-Piece Standard Closed-End Mini Pouch, no Filter**  
 Opaque with 1-sided comfort panel  
 Pouch length: 13 cm (5")

Product Code	Flange Size	Color Match	Box of
401930	32 mm (1 1/4")		20



**Little Ones® Two-Piece Standard Drainable Pouch, no Filter**  
 Transparent with 1-sided comfort panel  
 Pouch length: 15 cm (6"). 10 soft, moldable tail closures per box—designed for pediatric use.

Product Code	Flange Size	Color Match	Box of
401927	32 mm (1 1/4")		10
401928*	45 mm (1 3/4")		10

\*Little Ones® pouch 401928 can be used with the Natura® Flexible Skin Barrier 125259, 45 mm (1 3/4").



**Little Ones® Two-Piece Standard Urostomy Pouch**  
 Transparent with 1-sided comfort panel  
 Pouch length: 13 cm (5"). Each box contains two night drainage adapters.

Product Code	Flange Size	Color Match	Box of
401929	32 mm (1 1/4")		10



# Eakin Wound Pouches™

A variety of products that provide complete and consistent wound and skin management



# Eakin® Fistula and Wound Pouches



Eakin Wound Pouches™ can help contain a fistula, handle an exuding wound and help protect vulnerable skin. They come in a variety of shapes and sizes to provide complete and consistent wound and skin management.





### Small Eakin® Wound Pouch with Fold and Tuck Closure

Suitable for cutting area of 45 mm x 30 mm (1.8" x 1.2")

Product Code	Size	Box of
839250	45 mm x 30 mm (1.8" x 1.2")	10



### Medium Eakin® Wound Pouch with Fold and Tuck Closure

Suitable for cutting area of 110 mm x 75 mm (4.3" x 3")

Product Code	Size	Box of
839251	110 mm x 75 mm (4.3" x 3")	10



### Large Eakin® Wound Pouch with Fold and Tuck Closure

Suitable for cutting area of 175 mm x 110 mm (6.9" x 4.3")

Product Code	Size	Box of
839252	175 mm x 110 mm (6.9" x 4.3")	10



### Small Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure

Suitable for cutting area of 45 mm x 30 mm (1.8" x 1.2")

Product Code	Size	Box of
839260	45 mm x 30 mm (1.8" x 1.2")	10



### Medium Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure

Suitable for cutting area of 110 mm x 75 mm (4.3" x 3")

Product Code	Size	Box of
839261	110 mm x 75 mm (4.3" x 3")	10



**Large Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure**

Suitable for cutting area of 175 mm x 110 mm (6.9" x 4.3")

Product Code	Size	Box of
839262	175 mm x 110 mm (6.9" x 4.3")	10



**Extra Large Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure**

Suitable for horizontal cutting area of 245 mm x 160 mm (9.7" x 6.3")

Product Code	Size	Box of
839263	245 mm x 160 mm (9.7" x 6.3")	5



**Extra Large Eakin® Wound Pouch with Tap Closure and Access Window**

Suitable for easy access to horizontal cutting area of 245 mm x 160 mm (9.7" x 6.3")

Product Code	Size	Box of
839264	245 mm x 160 mm (9.7" x 6.3")	5



**Extra Large Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure**

Suitable for vertical cutting area of 160 mm x 245 mm (6.3" x 9.7")

Product Code	Size	Box of
839265	160 mm x 245 mm (6.3" x 9.7")	5



**Extra Large Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure**

Suitable for vertical cutting area of 130 mm x 290 mm (5.1" x 11.4")

Product Code	Size	Box of
839266	130 mm x 290 mm (5.1" x 11.4")	5



### Small Plus Eakin® Wound Pouch with Remote Drainage Attachment and Tap Closure

Suitable for cutting area of 60 mm x 86 mm (2.4" x 3.4")

Product Code	Size	Box of
839267	60 mm x 86 mm (2.4" x 3.4")	10



### Extra Large Eakin® Wound Pouch with Tap Closure and Access Window

Suitable for vertical cutting area of 160 mm x 245 mm ( 6.3" x 9.7")

Product Code	Size	Box of
839268	160 mm x 245 mm (6.3" x 9.7")	5



### Eakin® Access Window

Add to Eakin® wound and fistula pouches to allow instant access to the wound (may be used with 839252, 839262, 839263, 839265, 839266)

Product Code	Size	Box of
839280	n/a	5

eakin  
**Cohesive**<sup>®</sup>  
seal

**No.1**

for absorbency  
to stop leaks  
occurring<sup>1</sup>

# Routinely Reassuring.

Make every day a good day with eakin **Cohesive**<sup>®</sup> seals and help your patients prevent leaks from messing with their routine.

For more information, call our Customer Relations Center at **1-800-465-6302** or visit [convatec.ca](http://convatec.ca)



<sup>1</sup> Mc Grogan G, Haughey S and McDowell K (2018) An absorbent, enzyme-inhibiting seal reduces peristomal skin complications. *Gastrointestinal Nursing* 16:1 42-4



# Critical Care

## Fecal Management System

### Flexi-Seal® PROTECT PLUS Fecal Management System and Flexi-Seal® PROTECT PLUS with ENFit™ Connector Fecal Management System



Fecal management systems that sense when the retention balloon is correctly filled and visually indicates over-inflation

- Medication/irrigation port with LUER or ENFit™ connection, allowing to deliver medication effectively via a secondary route
- Self-closing catheter designed to reduce the risk of cross-contamination during bag change
- Blue finger pocket allows easy and accurate placement of the soft silicone retention balloon
- Proven *in vitro* to contain *C. difficile* and its spread into the environment<sup>25</sup>
- Odor barrier runs the full length of the catheter, providing an end-to-end odor barrier<sup>28</sup>
- The blue sampling port allows safe and easy stool collection helping minimize the risk of fecal exposure
- Contains Protect Plus Privacy™ Collection Bag with APS filter with 6x more odor-capturing capabilities compared to the charcoal filtered collection bags<sup>29</sup>. This filter also eliminates the need of manually releasing the excess gas from the collection bag<sup>30</sup>. Includes 4 Convatec Diamonds™ gelling and odor control sachets

Product Code	Box of
422303	1 kit with LUER connection with one filtered bag
421703	1 kit with ENFit™ connection with one filtered bag
422291	Protect Plus Privacy Collection Bag with APS filter (5/box)



### Flexi-Seal® PROTECT Fecal Management System

A fecal management system that senses when the retention balloon is correctly filled and visually indicates over-inflation

- Blue finger pocket allows easy and accurate placement of the soft silicone retention balloon
- Proven *in vitro* to contain *C. difficile* and its spread into the environment<sup>25</sup>
- Odor barrier runs the full length of the catheter, providing an end-to-end odor barrier<sup>28</sup>
- The blue sampling port allows safe and easy stool collection helping minimize the risk of fecal exposure
- Contains Privacy™ Collection Bag with APS filter with 6x more odor-capturing capabilities compared to the charcoal filtered collection bags<sup>29</sup>. This filter also eliminates the need of manually releasing the excess gas from the collection bag<sup>30</sup>

Product Code	Box of
421630	1 kit with one filtered bag

ENFit™ is a trademark of Global Enteral Device Supplier Association (GEDSA).



### Flexi-Seal SIGNAL® Fecal Management System

- Proven *in vitro* to contain the spread of *C. difficile* and prevent *C. difficile* from spreading into the environment<sup>25</sup>
- The indicator is designed to prevent overfilling of the retention balloon
- Excess gas released through the collection bag filter minimizes the need for frequent bag changes, potentially reducing the aerosolized spread of *C. difficile* spores<sup>30,31</sup>
- The filter also reduces or eliminates the need of manually releasing the excess gas from the collection bag<sup>30</sup>
- Odor barrier runs the length of the catheter, providing an end-to-end odor barrier<sup>28</sup>

Product Code	Box of
418000	1 kit with 3 filtered bags

### Flexi-Seal® FMS Disposable Replacement Bags

Product Code	Description	Box of
411102	Disposable Replacement Collection Bags with Filters, 1000 mL capacity	10 bags
411108	Disposable Replacement Privacy™ Collection Bags with APS Filter, 1000 mL capacity	10 bags



# Continence Care



# GentleCath™ Uncoated Intermittent Urinary Catheters



GentleCath™ PVC intermittent urinary catheters are designed to make catheter insertion and removal as smooth and easy as possible.

- Uncoated intermittent catheter features polished, rounded eyelets designed to minimize friction and a soft, rounded tip to aid smooth insertion
- Available in male and female straight tip and male coude tip, for insertion through complicated passages
- Color coded connector for easy FR size identification
- DEHP-Free<sup>15,\*</sup>

Female, 6.5", Straight Tip

Product Code	French Size	Color Code	Box of
501019	08 FR	●	100
501020	10 FR	●	100
501021	12 FR	○	100
501022	14 FR	●	100
501023	16 FR	●	100

Male, 16", Straight Tip

Product Code	French Size	Color Code	Box of
419800	08 FR	●	100
501002	10 FR	●	100
501003	12 FR	○	100
501004	14 FR	●	100
501005	16 FR	●	100
501006	18 FR	●	100

Male, 16", Coude Tip

Product Code	French Size	Color Code	Box of
501011	08 FR	●	100
501012	10 FR	●	100
501013	12 FR	○	100
501014	14 FR	●	100
501015	16 FR	●	100
501016	18 FR	●	100

\*DEHP (di(2-ethylhexyl)phthalate) is used to soften PVC plastic and is suspected to reduce fertility.

# GentleCath™ Glide Hydrophilic Intermittent Urinary Catheters



GentleCath™ Glide with FeelClean™ Technology integrates hydrophilic properties<sup>21</sup> into the catheter and is designed to protect the urethra from damage.<sup>22</sup>

- Hydrophilic intermittent urinary catheter with water sachet
- FeelClean™ Technology doesn't contain PVP<sup>1</sup>, therefore it helps protect the urethra from damage that can result from sticky residues<sup>21</sup>
- A no-touch handling sleeve helps prevent contamination of the catheter with bacteria from the hands<sup>23</sup>
- Hydrophilic properties are instantly activated with purified water<sup>24</sup>
- Color coded connector for easy FR size identification
- DEHP-Free<sup>15,\*</sup>

Female, 7", Straight Tip

Product Code	French Size	Color Code	Box of
421570	08 FR	●	30
421571	10 FR	●	30
421572	12 FR	○	30
421573	14 FR	●	30
421574	16 FR	●	30

Male, 16", Straight Tip

Product Code	French Size	Color Code	Box of
421564	08 FR	●	30
421565	10 FR	●	30
421566	12 FR	○	30
421567	14 FR	●	30
421568	16 FR	●	30
421569	18 FR	●	30

Male, 16", Coudé Tip

Product Code	French Size	Color Code	Box of
421907	08 FR	●	30
421908	10 FR	●	30
421909	12 FR	○	30
421910	14 FR	●	30
421911	16 FR	●	30
421912	18 FR	●	30

\*DEHP (di(2-ethylhexyl)phthalate) is used to soften PVC plastic and is suspected to reduce fertility.

# Life is full of sticky situations

With advanced FeelClean™ Technology, the GentleCath™ Glide intermittent hydrophilic catheter won't be one of them.

No sticking, minimal mess and residue, and game-changing relief.



Feel the difference,  
For more information, visit [gentlecath.com](http://gentlecath.com)

In vitro data on file. AP-033147-MRL-CA






# Wound Management

**AQUACEL<sup>®</sup>** *Pro*  
*Foam*

# Protecting the delicate skin of your patient is paramount

AQUACEL<sup>®</sup> Foam Pro is the foam dressing that helps protect the skin from pressure injuries by creating a healthy skin microclimate with Hydrofiber<sup>®</sup> Technology, while protecting even the most fragile skin with a gentle silicone adhesive layer.



**convatec**  
— forever caring —

## Wound Management



### Wound Hydration



#### DuoDerm® Hydroactive™ Gel

- Designed for the hydration and management of partial and full-thickness wounds such as leg ulcers and pressure ulcers
- Creates a moist healing environment, which helps to promote the natural autolytic processes of debridement

Product Code	Size	Box of
187987	30 g	3
187990	15 g	10



### Specialty Primary Dressings



#### CarboFlex® Odor Control Dressing

- Multi-layered sterile dressing for malodorous wounds<sup>5</sup>
- Creates a gel in the presence of wound exudate
- Provides odor control and absorbency<sup>5</sup> in *in vitro* testing in one easy-to-use, non-adhesive dressing

Application tip: place the cream-colored fibrous side on the wound bed. The white surface (which features a water-resistant film) should be facing away from the wound

Product Code	Size	Shape	Box of
403202	10 cm x 10 cm (4" x 4")	Square	10
403203	7.5 cm x 15 cm (3" x 6")	Rectangle	5
403204	15 cm x 20 cm (6" x 8")	Rectangle	5

\*See package insert for complete Instructions for Use.



## Moisture-Retentive Dressings

### Hydrofiber® Technology Dressings



#### AQUACEL® Foam

Aquacel® Foam Adhesive Silicone Hydrofiber® Dressing and Aquacel® Foam Non Adhesive Hydrofiber® Dressing offer comfort and simplicity, plus the healing benefits of an AQUACEL® interface.

- The outer film layer provides a waterproof viral and bacterial barrier which protects the wound from external contaminants, reducing the risk of infection<sup>1</sup>
- The foam and Hydrofiber® materials within the pad absorb high amounts of wound fluid and bacteria<sup>1</sup>
- The Hydrofiber® wound contact layer creates a soft, cohesive gel that intimately conforms to the wound surface, maintains a moist environment in the wound and aids in the removal of non-viable tissue from the wound without damaging tissue<sup>1</sup>
- Gentle silicone adhesive designed to adhere to surrounding skin, not the wound bed<sup>42</sup>

Product Code	Size	Shape	Box of
<b>Adhesive</b>			
420619	12.5 cm x 12.5 cm (5" x 5")	Square	10
420621	17.5 cm x 17.5 cm (7" x 7")	Square	10
420623	21 cm x 21 cm (8.5" x 8.5")	Square	5
420624	25 cm x 30 cm (10" x 12")	Rectangle	5
420625	19.8 cm x 14 cm (8" x 5.5")	Heel	5
420626	20 cm x 16.9 cm (8" x 7")	Sacral	5
420680	10 cm x 10 cm (4" x 4")	Square	10
420804	8 cm x 8 cm (3.2" x 3.2")	Square	10
421151	10 cm X 20 cm (4" x 8")	Rectangle	10
421153	10 cm x 25 cm (4" x 10")	Rectangle	10
421155	10 cm x 30 cm (4" x 12")	Rectangle	10
<b>Non-Adhesive</b>			
420631	5 cm x 5 cm ( 2" x 2")	Square	10
420633	10 cm x 10 cm (4" x 4")	Square	10
420635	15 cm x 15 cm (6" x 6")	Square	5
420636	20 cm x 20 cm (8" x 8")	Square	5
420637	15 cm x 20 cm (6" x 8")	Rectangle	5
421156	10 cm x 20 cm (4" x 8")	Rectangle	5

\*See package insert for complete Instructions for Use.



## AQUACEL® Ag Foam

Aquacel® Ag Foam Adhesive Silicone Hydrofiber® Dressing with Silver and Aquacel® Ag Foam Non Adhesive Hydrofiber® with Silver offer the benefits of an AQUACEL® Ag interface, the comfort of foam and the bacteria-killing power of ionic silver.<sup>7,†</sup>

- Absorbs wound fluid and bacteria, providing a moist wound healing environment, aiding autolytic debridement and removing dead-space between the wound and dressing interface<sup>2</sup>
- Is designed to manage excess exudate levels which may further damage the wound bed and surrounding skin<sup>2</sup>
- Provides a physical microbial, viral, and waterproof barrier to protect the wound<sup>2</sup>
- Effectively kills bacteria, yeasts and moulds<sup>2</sup>

Product Code	Size	Shape	Box of
<b>Adhesive</b>			
420627	12.5 cm x 12.5 cm (5" x 5")	Square	10
420628	17.5 cm x 17.5 cm (7" x 7")	Square	10
420629	21 cm x 21 cm (8.5" x 8.5")	Square	5
420647	19.8 cm x 14 cm (8" x 5.5")	Heel	5
420648	20 cm x 16.9 cm (8" x 7")	Sacral	5
420681	10 cm x 10 cm (4" x 4")	Square	10
420805	8 cm x 8 cm (3.2" x 3.2")	Square	10
420807	25 cm x 30 cm (10" x 12")	Rectangle	5
<b>Non-Adhesive</b>			
420639	5 cm x 5 cm (2" x 2")	Square	10
420642	10 cm x 10 cm (4" x 4")	Square	10
420645	15 cm x 15 cm (6" x 6")	Square	5
420646	20 cm x 20 cm (8" x 8")	Square	5
420806	15 cm x 20 cm (6" x 8")	Rectangle	5



## AQUACEL® Foam Pro

AQUACEL® Foam Pro Adhesive Silicone Hydrofiber® Dressing is designed to protect the skin from breakdown caused by moisture, friction and shear forces when used as part of a protocol of care.<sup>†,\*</sup>

- Smooth, breathable, waterproof film helps minimize shear and friction<sup>\*,†</sup>, e.g., between bed linens and dressing
- Multilayered silicone foam powered by Hydrofiber® Technology helps provide a healthy skin microclimate by absorbing and locking away excess moisture from the skin<sup>\*,†, 20</sup>
- Perforated gentle silicone adhesive layer allows easy application and removal<sup>\*</sup>
- Adhesive does not stick to itself<sup>\*</sup>

Product Code	Size	Shape	Box of
421579	20 cm x 17.5 cm (8" x 7")	Sacral	5
421580	24 cm x 21 cm (9.5" x 8.5")	Sacral	5
422356	19.8 cm x 14 cm (8" x 5.5")	Rectangle	10
422357	10 cm x 10 cm (4" x 4")	Square	10
422358	15 cm x 15 cm (6" x 6")	Square	10
422359	8 cm x 8 cm (3.2" x 3.2")	Square	10

\*As demonstrated *in vitro*.

†When used as part of a protocol of care.

\*\*See package insert for complete Instructions for Use.



## Moisture-Retentive Dressings



### FoamLite™ ConvaTec\*

Designed specifically for the management of low to non-exuding chronic and acute wounds. FoamLite™ ConvaTec Adhesive Silicone Foam Dressing is a thin dressing with a breathable film and perforated, gentle, skin friendly silicone adhesive.<sup>†,40</sup> FoamLite™ ConvaTec Adhesive Silicone Foam Dressing is flexible and conforms<sup>†,40</sup> to the body.

- The outer film provides a waterproof, viral and bacterial barrier which helps protect the wound from external contamination<sup>†,3</sup>
- Perforated, gentle, skin-friendly silicone adhesive, easy to apply, re-position and atraumatic to the wound and periwound skin upon removal<sup>†,41</sup>
- Thin, flexible, soft foam layer designed to absorb low exudate levels and helping maintain a moist wound healing environment<sup>†,40</sup>

Product Code	Size	Shape	Box of
421557	8 cm x 8 cm (3.1" x 3.1")	Square	10
421559	10 cm x 10 cm (4" x 4")	Square	10
421561	15 cm x 15 cm (5.9" x 5.9")	Square	10
421563	5.5 cm x 12 cm (2.1" x 4.7")	Rectangle	10
421927	5 cm x 5 cm (approx. 2" x 2")	Square	10
421928	10 cm x 20 cm (approx. 4" x 8")	Rectangle	10

† As demonstrated *in vitro*

\*Please refer to package insert for complete instructions for use.



## Moisture-Retentive Dressings

### Other Moisture-Retentive Dressings

#### DuoDerm® Extra Thin\*

DuoDerm® Extra Thin CGF™ Dressing is easy to use<sup>43</sup>, moulds, and can be cut to shape\* to allow the convenient management of wounds in awkward areas.

- Covers and protects susceptible areas—can be used anywhere on the body to gently protect newly formed tissue or skin at risk of breakdown
- Supports moist wound healing<sup>37</sup>
- Creates a moist wound-healing environment that supports reepithelialisation<sup>38</sup>
- Thin design is flexible and conformable

Product Code	Size	Shape	Box of
187900	5 cm x 10 cm (2" x 4")	Rectangle	20
187901	7 cm x 7 cm (3" x 3")	Square	20
187902	10 cm x 15 cm (4" x 6")	Oval	10
187903	15 cm x 18 cm (6" x 7")	Triangle	10
187955	10 cm x 10 cm (4" x 4")	Square	10
187957	14 cm x 14 (6" x 6")	Square	10
187961	5 cm x 20 cm (2" x 8")	Rectangle	10
187932	4.4 cm x 3.8 cm (1.75" x 1.5")	Spots	20

#### DuoDerm Signal®\*

- DuoDerm Signal® A tapered Edge CGF™ Dressing with easy-to-use visual indicator signals when it's time to change the dressing - may help simplify dressing changes
- Tapered border edges—designed to reduce the potential for dressing edges to catch on bed linen and clothing and to enhance patient comfort
- Smooth, low-friction backing—designed to reduce shear that can cause dressing to prematurely dislodge
- Thin profile—makes the dressing flexible and conformable
- Translucent—enhances placement and initial monitoring of the wound
- For wounds with low-to-moderate exudate

Product Code	Size	Shape	Box of
403326	10 cm x 10 cm (4" x 4")	Square	5
403327	14 cm x 14 cm (6" x 6")	Square	5
403332	15 cm x 18 cm (6" x 7")	Triangle	5
403333	20 cm x 23 cm (8" x 9")	Triangle	5
410501	20 cm x 22 cm (8" x 9")	Butterfly	5
410510	11 cm x 19 cm (4.5" x 7.5")	Oval	5

\*See package insert for complete Instructions for Use.



### DuoDerm® CGF™

- DuoDerm® CGF™ Control Gel Formula Dressing is for the management of exuding wounds
- May be used for partial and full-thickness wounds with exudate including: dermal ulcers including full-thickness wounds, pressure injuries (Stage 2 - 4), leg ulcers, superficial wounds, partial thickness under the supervision of a Health Care Professional

Product Code	Size	Shape	Box of
187643	15 cm x 20 cm (6" x 8")	Rectangle	5
187644	20 cm x 29 cm (8" x 12")	Rectangle	5
187658	10 cm x 10 cm (4" x 4")	Square	20
187659	14 cm x 14 cm (6" x 6")	Square	20
187660	10 cm x 10 cm (4" x 4")	Square	5
187661	14 cm x 14 cm (6" x 6")	Square	5
187662	20 cm x 20 cm (8" x 8")	Square	5

\*See package insert for complete Instructions for Use.

The original



**AQUACELE<sup>®</sup>Extra<sup>™</sup>**

#realimpact

# REAL PROGRESS

Our innovative Hydrofiber<sup>®</sup> gives you the confidence that wounds are healing and your patients are making real progress.<sup>1</sup>

**88.1%**  
of wounds healed  
or improved in  
just 14 days<sup>1</sup>

Visit [www.convatec.ca](http://www.convatec.ca) for more information




1. Carrere C et al. Community Setting Survey Evaluating AQUACELE<sup>®</sup> Dressings. Journal of Wound Care VOL 30, NO 9. 2021 Sep 2;30(9):763-774.  
© 2022 Convatec @/™ indicate trademarks of the Convatec group of companies. AP-034369 - MRL- GBL

# Wound Management



## Exudate Management Dressings

### Hydrofiber® Technology




Locks in wound exudate<sup>11</sup> and traps bacteria<sup>8-9,†</sup>

- Helps protect periwound skin and reduce maceration<sup>10-11</sup>
- May help minimize cross-infection and risk of wound infection during removal<sup>9,12</sup>



Micro-contours to the wound bed<sup>13,†</sup>

- Minimizes “dead space” where bacteria can thrive<sup>13</sup>
- Maintains moisture balance in the wound bed<sup>14</sup>



Responds to wound conditions by forming a cohesive gel

- The addition of ionic silver provides rapid and sustained antimicrobial activity<sup>7, 9, 16,†</sup>
- Minimizes pain associated with dressing changes<sup>17-18</sup>

†As demonstrated *in vitro*  
All images are artist's impressions

### Hydrofiber® Technology Dressings



#### AQUACEL® EXTRA™

AQUACEL® EXTRA™ Hydrofiber® Dressing with Strengthening Fibre is a soft, sterile, non-woven pad dressing composed of Hydrofiber® (sodium carboxymethylcellulose). This conformable and highly absorbent dressing absorbs wound fluid and creates a soft cohesive gel, maintains a moist environment which supports the body's healing process and aids in the removal of necrotic material from the wound (autolytic debridement) without damaging newly formed tissue. The AQUACEL® EXTRA™ Hydrofiber® Dressing with Strengthening Fiber dressing is composed of two layers of Hydrofiber® Technology stitched together, and designed to make it 9X stronger<sup>19</sup> and provides an increased absorbency of 39%<sup>19</sup> higher than the original AQUACEL® wound dressing.

- Supports healing by providing a moist environment
- Soft, non-woven pad dressing can be used on acute and chronic wounds including: leg ulcers, diabetic foot ulcers, pressure injuries, post-op wounds or traumatic wounds left to heal by secondary intent, abrasions, lacerations

Product Code	Size	Shape	Box of
420671	5 cm x 5 cm (approx. 2" x 2")	Square	10
420672	10 cm x 10 cm (approx. 4" x 4")	Square	10
420673	15 cm x 15 cm (approx. 6" x 6")	Square	5

\*See package insert for complete Instructions for Use.



## AQUACEL®

Aquacel® Hydrofiber® Dressing with Strengthening Fiber is:

- Soft, non-woven ribbon dressing can be used on acute and chronic wounds including: leg ulcers, diabetic foot ulcers, pressure injuries, post-op wounds or traumatic wounds left to heal by secondary intent, abrasions, lacerations
- Gelling action is designed to prevent dressing from adhering to wounds

Product Code	Size	Box of
403770	2 cm x 45 cm (approx. 0.75" x 18")	5
420127	1 cm x 45 cm (approx. 0.39" x 18")	5

## Surgical Dressings

Aquacel® Surgical Hydrofiber® Cover Dressing and Aquacel® Ag Surgical Hydrofiber® Cover Dressing:

- Have been designed to reduce the incidence of superficial surgical site infections, a postoperative complication<sup>33,34</sup>
- Can be worn during rehabilitation after surgery, including during early mobilisation and enhanced recovery<sup>33</sup>
- Waterproof, bacterial and viral barrier, dressings are designed to conform to the incision during body movement<sup>33</sup>

### AQUACEL® Surgical



Product Code	Size	Shape	Box of
412017	9 cm x 10 cm (approx. 3.5" x 4")	Rectangle	10
412018	9 cm x 15 cm (approx. 3.5" x 6")	Rectangle	10
412019	9 cm x 25 cm (approx. 3.5" x 10")	Rectangle	10
412020	9 cm x 35 cm (approx. 3.5" x 14")	Rectangle	10
420669	9 cm x 30 cm (approx. 3.5" x 12")	Rectangle	10

### AQUACEL® Ag Surgical

Product Code	Size	Shape	Box of
412009	9 cm x 10 cm (approx. 3.5" x 4")	Rectangle	10
412010	9 cm x 15 cm (approx. 3.5" x 6")	Rectangle	10
412011	9 cm x 25 cm (approx. 3.5" x 10")	Rectangle	10
412012	9 cm x 35 cm (approx. 3.5" x 14")	Rectangle	10
420670	9 cm x 30 cm (approx. 3.5" x 12")	Rectangle	10

\*See package insert for complete Instructions for Use.

† Patients with infected wounds should always be managed with a comprehensive protocol of care to address wound infection.

## Other Exudate Management Dressings



### KALTOSTAT®

KALTOSTAT® Calcium-Sodium Alginate is a soft, sterile dressing made of calcium-sodium alginate fiber that converts to a firm gel/fiber mat on contact with fluid

- Designed to absorb exudate and to help control minor bleeding
- Supports a moist wound healing environment
- Non adherent to wound surface allowing for easy, non-traumatic removal

Product Code	Size	Shape	Box of
168117	2 g	Rope	5
168210	5 cm x 5 cm (2" x 2")	Square	10
168212	7 cm x 12 cm (3" x 5")	Rectangle	10
168214	10 cm x 19 cm (4" x 8")	Rectangle	10

## Exudate & Biofilm Management Dressings



### AQUACEL® Ag+ Extra™

Aquacel® Ag+ Extra™ Enhanced Hydrofiber® Dressing with Silver and Strengthening Fiber disrupts and destroys biofilm. Specifically developed to win the battle against biofilm, MORE THAN SILVER™ technology contains three components; ionic silver together with a surfactant and metal chelating agent, which work together to deliver superior anti-biofilm performance<sup>44</sup>

- Rapid, superior and sustained anti-biofilm activity against antibiotic-resistant biofilm and prevents reformation of biofilm<sup>32</sup>
- Absorbs wound fluid, bacteria and biofilm, providing a moist wound healing environment, aiding autolytic debridement and removing dead-space between the wound and dressing interface<sup>4</sup>
- Provides a barrier to protect the wound bed from contamination<sup>4</sup>
- Effectively kills bacteria, yeasts and moulds in both planktonic and biofilm form<sup>4</sup>

\*For instruction on how to use the dressing format on Partial Thickness (Second Degree) burns, please refer to the most recent Aquacel® Ag+ Extra™ dressing Instructions for Use

Product Code	Size	Shape	Box of
413566	5 cm x 5 cm (2" x 2")	Square	10
413567	10 cm x 10 cm (4" x 4")	Square	10
413568	15 cm x 15 cm (6" x 6")	Square	5
413569	20 cm x 30 cm (8" x 12")	Rectangle	5



### AQUACEL® Ag+ Extra™ Enhanced Hydrofiber® Ribbon Dressing with Silver and Strengthening Fibre

Product Code	Size	Box of
413570	1 cm x 45 cm (approx. 0.39" x 18")	5
413571	2 cm x 45 cm (approx. 0.75" x 18")	5

\*See package insert for complete Instructions for Use.

† As demonstrated *in vitro*

# YOU HAVE THE POWER TO DISRUPT AND DESTROY BIOFILM<sup>1</sup> TO ADVANCE HEALING<sup>2,3</sup>



## NON-HEALING WOUNDS NEED MORE THAN SILVER™

Specifically developed to win the battle against biofilm, **MORE THAN SILVER™** technology harnesses the synergistic power of three key components; ionic silver, a surfactant and a metal chelating agent working together to deliver superior<sup>1</sup> performance.\*



+



=

**AQUACEL<sup>®</sup> Ag+**  
Dressings

1. Bowler PG, et al. Parsons, Wound Medicine 14 (2016) 6-11. 2. Metcalf DG et al. J. Wound Care 2016; Vol25, No3. 3. Metcalf DG, et al. Int Wound J 2017; 14: 203-213. 4. Malone M et al. 2017. JWC; 20-25.

AQUACEL, Hydrofiber and MORE THAN SILVER are trademarks of Convatec Inc. ©2019 Convatec Inc. AP-019908-MM

\* When compared to AQUACEL<sup>®</sup> Ag Extra™ dressing and other silver-only competitor dressings: ACTICOAT™ 7 and SILVERCEL™ Non-Adherent dressings.





# Negative Pressure Wound Therapy

# Negative Pressure Wound Therapy



## Avelle™ Negative Pressure Wound Therapy System

The first system to combine negative pressure wound therapy with Hydrofiber® Technology. Negative Pressure Wound Therapy delivered in a simple, discreet, canister-less system: Easy to use - Small and portable - Quiet, no alarm - Disposable pump with 30 day lifetime† - Designed for flexibility - Suitable for use within and across different care settings - Separate pump and dressing packs - Supports patients resuming normal life.

### Feature description of Hydrofiber® Technology dressing:

- Gentle silicone adhesive border<sup>39\*</sup> secures the dressing in place.
- Stitch-bonded Hydrofiber® wound contact layer gels on contact with wound fluid and is designed to maintain dressing integrity on removal.<sup>39\*</sup>
- Hydrofiber® core gels on contact with wound fluid locking in exudate and the bacteria it contains.<sup>26,27,39\*</sup> Channels in the Hydrofiber® core are designed to allow fluid to move throughout the dressing ensuring exudate is locked away.<sup>39\*</sup>
- Foam layer is designed to aid distribution of negative pressure across the dressing.<sup>39\*</sup>
- Film layer provides a bacterial and viral barrier, is showerproof, and permits evaporation of exudate aiding overall fluid handling.<sup>39\*</sup>

Product Code	Description	Box of
422285	Pump unit	1
446650	Carry Case	1
422308	ConvaTec™ Hydrocolloid Sealing Strips 5 X 20 cm	10
421921	Avelle™ NPWT Pump Spare Tube Set	5

Product Code	Size	Inner Pad dimensions	Box of
421552	16 cm x 16 cm (approx. 6.3" x 6.3")	8 cm x 8 cm (approx. 3.2" x 3.2")	5
421553	16 cm x 21 cm (approx. 6.3" x 8.3")	8 cm x 13 cm (approx. 3.2" x 5.1")	5
421554	12 cm x 21 cm (approx. 4.7" x 8.3")	4 cm x 13 cm (approx. 1.6" x 5.1")	5
421555	12 cm x 31 cm (approx. 4.7" x 12.2")	4 cm x 23 cm (approx. 1.6" x 9.1")	5
422155	12 cm x 41 cm (approx. 4.7" x 16.1")	4 cm x 33 cm (approx. 1.6" x 13")	5
422156	21 cm x 26 cm (approx. 8.3" x 10.2")	13 cm x 18 cm (approx. 5.1" x 7.1")	5
422157	26 cm x 26 cm (approx. 10.2" x 10.2")	18 cm x 18 cm (approx. 7.1" x 7.1")	5

\* As demonstrated *in vitro*.

† Battery change may be required during pump lifetime

# Avelle™

Negative Pressure Wound Therapy System

NPWT  
POWERED BY

AQUACEL®



## A FLEXIBLE NPWT SYSTEM

Maximize healing in hard-to-heal wounds

Avelle™ NPWT System - combining  
AQUACEL® dressing technology with the  
healing power of NPWT for up to 30 days.



convatec  
— forever caring —

# References

1. Instruction for Use Aquacel® Foam
2. Instruction for Use Aquacel® Foam Ag
3. Instructions for Use Foam Lite™ ConvaTec
4. Instructions for Use Aquacel® Ag+ Extra™ dressing.
5. Williams C. Role of CarboFlex® in the nursing management of wound odour. *Br J Nurs.* 2001;10(2):122-125.
6. McNichol L, Lund C, Rosen T, Gray M. Medical Adhesives and Patient Safety: State of the Science; Consensus Statements for the Assessment, Prevention, and Treatment of Adhesive-Related Skin Injuries. *JWOCN.* 2013;40(4):365-380.
7. Jones SA, Bowler PG, Walker M, Parsons D. Controlling wound bioburden with a novel silver-containing Hydrofiber® dressing. *Wound Rep Reg.* 2004;12:288-294.
8. Waring MJ, Parsons D. Physico-chemical characterisation of carboxymethylated spun cellulose fibres. *Biomaterials.* 2001;22:903-912.
9. Walker M, Hobot JA, Newman GR, Bowler PG. Scanning electron microscopic examination of bacterial immobilisation in a carboxymethylcellulose (AQUACEL®) and alginate dressings. *Biomaterials.* 2003;24(5):883-890.
10. Coutts P, Sibbald RG. The effect of a silver-containing Hydrofiber® dressing on superficial wound bed and bacterial balance of chronic wounds. *Int. Wound J.* 2005; 2(4):348-356.
11. Robinson BJ. The use of a Hydrofiber® dressing in wound management. *J Wound Care.* 2000;9(1):32-34.
12. Bowler PG, Jones SA, Davies BJ, Coyle E. Infection control properties of some wound dressings. *J. Wound Care.* 1999;8(10):499-502.
13. Jones S, Bowler PG, Walker M. Antimicrobial activity of silver-containing dressings is influenced by dressing conformability with a wound surface. *Wounds.* 2005;17(9):263-270.
14. Bishop SM, Walker M, Rogers AA, Chen WYJ. Moisture balance: optimising the wound-dressing interface. *J Wound Care.* 2003;12:125-128.
15. GentleCath™ Product Data Sheet 2018. AP-019876-US and GentleCath™ Glide Product Data Sheet 2022. AP-58423-GBL-ENG-v1. Data on file.\*DEHP (di(2-ethylhexyl)phthalate) is used to soften PVC plastic and is suspected to reduce fertility.
16. Bowler PG, Jones SA, Walker M, Parsons D. Microbicidal properties of a silver-containing Hydrofiber® dressing against a variety of burn wound pathogens. *J Burn Care Rehabil.* 2004;25:192-196.
17. Vanscheidt W, Münter KC, Klövekorn W, Vin F, Gauthier JP, Ukat A. A prospective study on the use of a non-adhesive gelling foam dressing on exuding leg ulcers. *J Wound Care.* 2007;16(6):261-265.
18. Parish LC, Dryski M, Cadden S. Prospective clinical study of a new adhesive gelling foam dressing in pressure ulcers. *Int Wound J.* 2008;5:60-67.
19. Assessment of the Physical Properties of AQUACEL® EXTRA™ & AQUACEL® Dressings. Scientific Background Report WHRI3461 TA214. 2011 Data on File, ConvaTec.
20. Walker M, Hobot JA, Newman GR, Bowler PG. Scanning electron microscopic examination of bacterial immobilization in a carboxymethylcellulose (AQUACEL®) and alginate dressing. *Biomaterials.* 2003;24(5):883-890.
21. Pollard D, Allen D, Irwin N J, Moore J V, McClelland N, McCoy C P. Evaluation of an Integrated Amphiphilic Surfactant as an Alternative to Traditional Polyvinylpyrrolidone Coatings for Hydrophilic Intermittent Urinary Catheters. *Biotribology.* 2022;32(Dec 2022):100223
22. *In vitro* data on file
23. An-Sofie Goessaert, Stephanie Antoons, Melissa Van Den Driessche, Ali Tourchi, Ronny Pieters & Karel Everaert. No-touch intermittent catheterization: caregiver point of view on sterility errors, duration, comfort, and costs. Accepted for publication 10 October 2012
24. RPT-006447 Biological Evaluation Process. Data on file. ConvaTec Inc
25. Jones S, Towers V, Welsby S, et al. Clostridium difficile Containment Properties of a Fecal Management System.
26. Waring MJ, Parsons D. Physico-chemical characterisation of carboxymethylated spun cellulose fibres. *Biomaterials.* 2001;22:903-912.
27. Walker M, Hobot JA, Newman GR, Bowler PG. Scanning electron microscopic examination of bacterial immobilisation in a carboxymethylcellulose (AQUACEL®) and alginate dressings. *Biomaterials.* 2003;24(5):883-890.
28. Odor Barrier Testing. 130124-001. Data on file, ConvaTec Inc.
29. Flexi-Seal® Privacy Bag Filter Evaluation. 121412-001. Data on file, ConvaTec Inc.
30. Minimizing the spread of C. difficile spores from the release of gas. February 19, 2013. Data on file, ConvaTec Inc.
31. Roberts K, Smith CF, Snelling AM, et al. Aerial dissemination of Clostridium difficile spores. *BMC Infect Dis.* 2008;8:7.

## References

32. WHRI5860 MA322 - In vitro Antimicrobial Activity of AQUACEL® Ag + Extra™ and AQUACEL® Ag Extra™ against Acticoat™ 7, Silvercel™ NA and Urgotul™ Ag -V1. Convatec Data on File. 27 November 2018
33. WHRI6551 MS171 In vitro performance characteristics of AQUACEL® Ag+ SCD and SCD SP
34. JV Clarke, AH Deakin, JM Dillon, S Emmerson, AWG Kinninmonth. A prospective clinical audit of a new dressing design for lower limb arthroplasty wound. J Wound Care. Volume 18, No.1 January 2009
35. Roberts K, Smith CF, Snelling AM, et al. Aerial dissemination of Clostridium difficile spores. BMC Infect Dis. 2008;8:7.
36. Queen, D. Understanding Hydrofiber® Technology, Wounds International. 2010 1(5).
37. Greguric S, Dragomir B and Soldo-Belic A. Hydrocolloid dressing versus a conventional dressing using magnesium sulphate paste in the management of venous leg ulcers. Acta Dermatovenerol Croat. 1994;2(1): 65-71.
38. Greguric S, Budimcic D, Reepithelialisation, Soldo BA et al. Hydrocolloid Moist wound-healing dressing versus a environment, conventional dressing using Autolytic magnesium sulphate paste debridement in the management of venous leg ulcers. Acta Dermatovenerol Croat 1994; 2: 65-71.
39. HFM-2015/017. Data on file. 2015. ConvaTec Inc.
40. In vitro Performance of Foam Lite™ ConvaTec WHRI4680 MS132. 25th April 2016.
41. Meuleneire F, Rücknagel H. Soft silicones Made Easy. Wounds International 2013 (May). Available from: [www.woundsinternational.com](http://www.woundsinternational.com)
42. Healthy Volunteer Study Report, Protocol CW 0207 11 A736. Data on File, ConvaTec Inc. CW 0207 11 A734. Prospective, multi centre, performance evaluation of exudate management and adhesion of AQUACEL® foam adhesive dressing in pressure ulcers. Tickle, J., 2016. Positive clinical and patient outcomes with a next generation foam dressing. Wounds UK.
43. Forshaw A. Hydrocolloid Convenient, dressings in paediatric Mouldable, Cut to wound care. J Wound Care shape 1993; 2(4): 209-12.
44. Bowler PG, Parsons, D. Combatting wound biofilm and recalcitrance with a novel anti-biofilm Hydrofiber® wound dressing. Wound Medicine 14 (2016) 6-11

© 2023 Convatec Inc. ®/™ indicates trademarks of Convatec Inc. unless otherwise noted.

## Useful Resources

### Ostomy Canada Society

5800 Ambler Dr., Suite 210,  
Mississauga, ON L4W 4J4  
Tel.: (905) 212-7111  
Toll Free: 1 888 969-9698  
Fax: (905) 212-9002  
E-mail: info1@ostomycanada.ca  
www.ostomycanada.ca

### Canadian Cancer Society

55 St. Clair Avenue West, Suite 500  
Toronto, ON M4V 2Y7  
Tel.: (416) 961-7223  
Toll Free: 1 888 939-3333  
Fax: (416) 961-4189  
Email: info@cis.cancer.ca  
www.cancer.ca

### Crohn's & Colitis Canada

2110 - 439 University Ave.  
Toronto, Ontario M5G 1Y8  
Tel.: (416) 920-5035  
Toll Free: 1 800 387-1479  
Fax: (416) 929-0364  
E-mail: support@crohnsandcolitis.ca  
www.crohnsandcolitis.ca

### Nurses Specialized in Wound, Ostomy, & Continence Canada (NSWOCC) - Formerly CAET

66 Leopolds Drive  
Ottawa, Ontario  
K1V 7E3  
Tel.: (888) 739-5072  
Fax: (613) 523-1878  
E-mail: office@nswoc.ca  
www.nswoc.ca

### Colorectal Cancer Canada

1350 Sherbrooke St West, Suite 300  
Montreal, QC H3G 1J1  
E-mail: info@colorectalcancerCanada.com  
Tel.: (514) 875-7745  
Tel.: (877) 502 6566  
Fax: (514) 875-7746  
www.colorectalcancerCanada.com

### Friends of Ostomates Worldwide Canada (FOWC)

Tel.: (647) 951-3940  
Email: info@fowc.ca

### Wounds Canada

(The Canadian Association of Wound Care)  
P.O. Box 35569  
York Mills Plaza,  
North York, ON M2L 2Y4  
Tel.: (416) 485-2292  
Fax: (888) 710-2289  
www.woundscanada.ca

### Canadian Home Care Association

2000 Argentia Road  
Plaza 3, Suite 302  
Mississauga, ON L5N 1W1  
Tel.: (905) 567-7373  
Email: chca@cdnhomecare.ca  
www.cdnhomecare.ca

### Infection Prevention and Control Canada (IPAC Canada)

P.O. Box 46125 RPO Westdale  
Winnipeg, MB R3R 3S3  
Tel.: (866) 999-7111 or (204) 897-5990  
Fax: (204) 895-9595  
Email: info@ipac-canada.org  
www.ipac-canada.org

### Canadian Dermatology Association

320 March Road, Suite 400  
Ottawa, Ontario K2K 2E3  
Tel.: (800) 267-3376 or (613) 738-1748  
Fax: (866) 267-2178 or (613) 738-4695  
www.dermatology.ca

### Diabetes Canada

National Office  
1300-522 University Ave  
Toronto, ON M5G 2R5  
Tel.: (800) 226-8464 or (416) 363-3373  
Email: info@diabetes.ca  
www.diabetes.ca

### Health Canada

Address Locator 0900C2  
Ottawa, ON K1A 0K9  
Tel.: (866) 225-0709 or (613) 957-2991  
Fax: (613) 941-5366  
Email: hcinfo.infosc@canada.ca  
www.hc-sc.gc.ca  
www.canada.ca

## Ostomy Glossary of Terms

### **Adhesive Coupling Technology™ (ACT) –**

A two-piece ostomy pouching system that attaches with a tape-like adhesive rather than plastic rings snapping together.

**Closed-end pouch –** A pouch that has no opening at the bottom and is removed and thrown away after each use. This makes them ideal for special occasions, intimate moments, or for everyday use. Some are offered with a filter designed to reduce odor and gas buildup.

**Comfort panels –** A fine covering on one or both sides of the pouches, feels comfortable on the skin.

**Convexity –** The outward curving of a faceplate which begins at the aperture of the faceplate and extends outward.<sup>1</sup>

**Cut-to-Fit –** The opening of the skin barrier or appliance needs to be enlarged (using scissors) to match the stoma size.

**Drainable pouch –** Has an opening at the bottom which allows stool or urine to be drained and re-closed with a clip or self gripping fastener hook tape, thus can be re-used. Urostomy pouches are drained by the valve or tap.

**Durahesive® –** Is an extended wear barrier that provides the ultimate in security and adhesion for people with liquid to semi-liquid outputs.

**Durahesive Plus® –** Is an extended wear barrier available in a one-piece cut-to-fit pouch.

**Film –** All Convatec pouches have built in layers that help reduce the rustling noise under clothing, and a built in odor-barrier.

**Filter –** Some Convatec pouches have an integral filter that helps reduce gas buildup within the pouch.

**Flange –** The plastic connector ring on the ostomy skin barrier or pouch.

**InvisiClose® tail closure –** A clipless pouch that is closed via an integrated interlocking system.

**Moldable –** No scissors are needed to size a moldable skin barrier. The adhesive can be molded or shaped to the size of the stoma using one's fingers.

**One-piece ostomy pouching system –** The skin barrier and pouch are combined into one unit for a one-step application. It is flexible and offers a low profile. The pouch cannot be removed without removing the skin barrier/wafer.

**Peristomal skin –** The skin surrounding a stoma.

**Pre-cut –** The opening of the skin barrier or appliance is already cut to match the stoma size (No need of scissors).

**Skin barrier/wafer –** Round or square-shaped with a sticky side that is applied to the skin around the stoma (called peristomal skin). Helps protect the skin from stool or urine, which can be very irritating.

**Skin protective barrier –** A wipe or spray that forms a protective film on skin to help prevent skin breakdown by shielding it from irritating or caustic substances.

**Stomahesive® –** Is a skin barrier that provides excellent security for people with semi-formed to formed output because of its special formulation that enables it to adhere to both dry and wet skin.

**Tape collar –** Convatec skin barriers are available with or without tape collars. Most skin barriers with tape collars are available in white, dye-free tape or tan tape.

**Two-piece ostomy pouching system –** Skin barrier/wafer and pouch are separate and attached together with a round, plastic ring (called a flange). Pouch can be easily removed without also having to remove the skin barrier/wafer.

<sup>1</sup> Rolstad BS, Boarini J. Principles and techniques in the use of convexity. *Ostomy/Wound Management*; January/February 1996; (42)1: 26.

## Wound Care Glossary of Terms

**Bactericidal** – Bacteria will be killed and therefore growth will not restart when contact with the antimicrobial agent is stopped.

**Bacteriostatic** – Bacterial growth will be inhibited, but will start again when contact with the antimicrobial agent is stopped.

**Biofilm** – Biofilm describes surface-adherent bacteria, encased in and protected by an extracellular polymeric substance (EPS) composed of polysaccharides, DNA and proteins.

**Contamination** – Soiling by contact or introduction of organisms into a wound.<sup>1</sup>

**Hydrocolloid** – A substance that yields a gel with water.<sup>2</sup>

**Hydrofiber® Technology** – The Hydrofiber® Technology is exclusive to Convatec. It is composed of 100% pure sodium carboxymethylcellulose. It is the underlying technology behind AQUACEL®.

**Infection** – An overgrowth of microorganisms capable of tissue invasion and destruction. Local and systemic symptoms<sup>1</sup> may be observed.

**Microbicidal** – Bacteria, yeasts and fungi will be killed and therefore growth will not restart when contact with the antimicrobial agent is stopped.

**Partial-thickness burn** – (Second degree): A burn that damages the epidermal and some dermal tissues but does not damage the lower-lying hair follicles, sweat or sebaceous glands.<sup>3</sup>

<sup>1</sup> Barton P, Parslow N. Wound Care-A Comprehensive Guide for Community Nurses. Saint Elizabeth Health Care; 1996: 86-87.

<sup>2</sup> Webster's Ninth New Collegiate Dictionary. Philippines Copyright by Merriam-Webster Inc.; 1984: 589.

<sup>3</sup> Taber's Cyclopedic Medical Dictionary. Edition 19: F.A. Davis Company; 2001: 295.



# Product Code Index

Product Code	UPC Codes	Page	Product Code	UPC Codes	Page
20917	7 68455 10869 5	69	125339	7 68455 10200 6	64
20922	7 68455 10139 9	69	125340	7 68455 10201 3	64
21712	7 68455 10694 3	10	125341	7 68455 10202 0	64
21715	7 68455 10696 7	10	125342	7 68455 10203 7	64
22355	7 68455 11550 1	12	125343	7 68455 10204 4	64
22356	7 68455 11551 8	12	125361	7 68455 10214 3	65
22736	7 68455 10708 7	13	125362	7 68455 10215 0	65
22751	7 68455 10141 2	64	125363	7 68455 10216 7	65
22764	7 68455 10154 2	64	125364	7 68455 10217 4	65
22765	7 68455 10155 9	64	125365	7 68455 10218 1	65
22766	7 68455 10156 6	64	125366	7 68455 10219 8	65
22767	7 68455 10157 3	64	125367	7 68455 10220 4	65
22768	7 68455 10158 0	64	125368	7 68455 10221 1	65
22769	7 68455 10159 7	64	125369	7 68455 10222 8	65
22771	7 68455 10161 0	63	168117	7 68455 08116 5	95
25510	7 68455 10826 8	10	168210	7 68455 08114 1	95
25542	7 68455 10702 5	10	168212	7 68455 08115 8	95
27060	7 68455 10913 5	12	168214	7 68455 08111 0	95
27062	7 68455 10709 4	12	175507	7 68455 10718 6	12
37436	7 68455 10767 4	8	175611	7 68455 10236 5	65
37439	7 68455 10766 7	8	175635	7 68455 10719 3	11
37443	7 68455 10760 5	8	175650	7 68455 10707 0	13
37444	7 68455 10761 2	8	175652	7 68455 10867 1	11
62436	7 68455 10876 3	11	175777	7 68455 10251 8	63
64927	7 68455 10168 9	66	175778	7 68455 10252 5	63
125257	7 68455 10170 2	21	175779	7 68455 10253 2	63
125258	7 68455 10171 9	21	175780	7 68455 10254 9	63
125259	7 68455 10172 6	21	175781	7 68455 10255 6	63
125260	7 68455 10173 3	21	175782	7 68455 10256 3	63
125261	7 68455 10174 0	21	175783	7 68455 10257 0	63
125262	7 68455 10175 7	20	175784	7 68455 10258 7	63
125263	7 68455 10176 4	20	183910	7 68455 10690 5	10
125264	7 68455 10177 1	20	187643	7 68455 10269 3	91
125265	7 68455 10178 8	20	187644	7 68455 10670 7	91
125266	7 68455 10179 5	20	187658	7 68455 10270 9	91
125267	7 68455 10180 1	23	187659	7 68455 10271 6	91
125268	7 68455 10181 8	23	187660	7 68455 10697 4	91
125269	7 68455 10182 5	23	187661	7 68455 10272 3	91
125270	7 68455 10183 2	23	187662	7 68455 10698 1	91
125271	7 68455 10184 9	23	187900	7 68455 10277 8	90
125272	7 68455 10185 6	23	187901	7 68455 10278 5	90
125273	7 68455 10186 3	23	187902	7 68455 10279 2	90
125274	7 68455 10187 0	23	187903	7 68455 10280 8	90
125275	7 68455 10188 7	23	187932	7 68455 10281 5	90
125276	7 68455 10189 4	23	187955	7 68455 10691 2	90
125277	7 68455 10190 0	23	187957	7 68455 10692 9	90
125330	7 68455 10191 7	63	187961	7 68455 10282 2	90
125331	7 68455 10192 4	64	187987	7 68455 10706 3	86
125333	7 68455 10194 8	64	187990	7 68455 10864 0	86
125334	7 68455 10195 5	64	400598	7 68455 10289 1	63
125335	7 68455 10196 2	64	400599	7 68455 10290 7	63
125336	7 68455 10197 9	64	401500	7 68455 10319 5	33
125338	7 68455 10199 3	64	401501	7 68455 10321 8	33

# Product Code Index

Product Code	UPC Codes	Page	Product Code	UPC Codes	Page
401502	7 68455 10323 2	33	401912	7 68455 10400 0	13
401503	7 68455 10325 6	33	401913	7 68455 10401 7	13
401504	7 68455 10326 3	33	401914	7 68455 10402 4	13
401505	7 68455 10327 0	34	401918	7 68455 10405 5	13
401506	7 68455 10328 7	34	401919	7 68455 10406 2	13
401507	7 68455 10329 4	34	401925	7 68455 10412 3	70
401508	7 68455 10330 0	34	401927	7 68455 10414 7	70
401509	7 68455 10331 7	34	401928	7 68455 10415 4	70
401510	7 68455 10332 4	33	401929	7 68455 10416 1	70
401511	7 68455 10333 1	33	401930	7 68455 10417 8	70
401512	7 68455 10334 8	33	401989	7 68455 10425 3	13
401513	7 68455 10335 5	33	403202	7 68455 10428 4	86
401514	7 68455 10336 2	33	403203	7 68455 10429 1	86
401515	7 68455 10337 9	34	403204	7 68455 10430 7	86
401529	7 68455 10346 1	35	403326	7 68455 10725 4	90
401530	7 68455 10347 8	35	403327	7 68455 10723 0	90
401531	7 68455 10348 5	35	403332	7 68455 10724 7	90
401532	7 68455 10349 2	35	403333	7 68455 10727 8	90
401533	7 68455 10350 8	36	403770	7 68455 10677 6	94
401534	7 68455 10351 5	36	403948	7 68455 10434 5	45
401535	7 68455 10352 2	36	404006	7 68455 10710 0	24
401536	7 68455 10353 9	36	404007	7 68455 10711 7	24
401537	7 68455 10354 6	36	404008	7 68455 10712 4	24
401538	7 68455 10671 4	36	404009	7 68455 10713 1	24
401539	7 68455 10355 3	36	404010	7 68455 10714 8	24
401540	7 68455 10356 0	36	404011	7 68455 10715 5	24
401541	7 68455 10357 7	36	404012	7 68455 10716 2	24
401542	7 68455 10358 4	36	404013	7 68455 10717 9	24
401543	7 68455 10359 1	36	404014	7 68455 10437 6	33
401544	7 68455 10360 7	36	404015	7 68455 10438 3	33
401545	7 68455 10361 4	36	404016	7 68455 10439 0	33
401546	7 68455 10362 1	36	404017	7 68455 10440 6	33
401547	7 68455 10363 8	37	404018	7 68455 10441 3	33
401548	7 68455 10364 5	37	404025	7 68455 10442 0	33
401549	7 68455 10365 2	37	404026	7 68455 10443 7	33
401550	7 68455 10366 9	37	404027	7 68455 10444 4	33
401551	7 68455 10367 6	37	404028	7 68455 10445 1	33
401552	7 68455 10368 3	37	404029	7 68455 10446 8	33
401553	7 68455 10369 0	37	404030	7 68455 10447 5	34
401554	7 68455 10370 6	37	404031	7 68455 10448 2	34
401555	7 68455 10371 3	37	404032	7 68455 10449 9	34
401557	7 68455 10372 0	34	404033	7 68455 10450 5	34
401558	7 68455 10373 7	34	404034	7 68455 10451 2	34
401559	7 68455 10374 4	34	404592	7 68455 10452 9	18
401573	7 68455 10375 1	22	404593	7 68455 10453 6	18
401574	7 68455 10376 8	22	404594	7 68455 10454 3	18
401575	7 68455 10377 5	22	405451	7 68455 10491 8	47
401576	7 68455 10378 2	22	405456	7 68455 10496 3	44
401577	7 68455 10379 9	22	405457	7 68455 10497 0	44
401905	7 68455 10393 5	22	405458	7 68455 10498 7	44
401906	7 68455 10394 2	22	405468	7 68455 10500 7	45
401909	7 68455 10397 3	35	405474	7 68455 10505 2	45
401911	7 68455 10399 7	13	405475	7 68455 10507 6	45

# Product Code Index

Product Code	UPC Codes	Page	Product Code	UPC Codes	Page
405476	7 68455 10508 3	45	413160	7 68455 10620 2	21
405477	7 68455 10509 0	45	413161	7 68455 10621 9	21
405478	7 68455 10510 6	45	413162	7 68455 10622 6	21
405479	7 68455 10511 3	45	413163	7 68455 10623 3	21
405480	7 68455 10512 0	45	413164	7 68455 10624 0	20
405482	7 68455 10514 4	45	413165	7 68455 10625 7	20
405483	7 68455 10515 1	45	413166	7 68455 10626 4	20
409268	7 68455 10523 6	44	413167	7 68455 10627 1	20
409269	7 68455 10524 3	44	413168	7 68455 10628 8	20
409270	7 68455 10525 0	44	413177	7 68455 10637 0	19
410501	7 68455 10748 3	90	413178	7 68455 10638 7	19
410510	7 68455 10750 6	90	413179	7 68455 10639 4	19
411102	7 68455 10862 6	79	413180	7 68455 10640 0	19
411108	7 68455 12572 2	79	413181	7 68455 10641 7	19
411265	7 68455 10562 5	32	413182	7 68455 10642 4	19
411266	7 68455 10563 2	32	413183	7 68455 10643 1	19
411267	7 68455 10564 9	32	413184	7 68455 10644 8	19
411289	7 68455 10566 3	32	413185	7 68455 10645 5	19
411290	7 68455 10567 0	32	413186	7 68455 10646 2	19
411291	7 68455 10568 7	32	413187	7 68455 10647 9	19
411490	7 68455 10581 6	32	413188	7 68455 10648 6	19
411491	7 68455 10582 3	32	413417	7 68455 12391 9	18
411492	7 68455 10583 0	32	413418	7 68455 12392 6	18
411493	7 68455 10584 7	32	413419	7 68455 12393 3	18
411633	7 68455 10587 8	69	413420	7 68455 12394 0	18
411637	7 68455 10588 5	70	413421	7 68455 12395 7	18
411639	7 68455 10590 8	69	413422	7 68455 12396 4	18
411642	7 68455 10593 9	69	413423	7 68455 12397 1	18
411657	7 68455 11733 8	44	413424	7 68455 12398 8	18
411658	7 68455 11734 5	44	413437	7 68455 12885 3	30
411659	7 68455 11735 2	44	413438	7 68455 12887 7	30
411800	7 68455 10594 6	18	413509	7 68455 13005 4	55
411801	7 68455 10595 3	18	413510	7 68455 13006 1	55
411802	7 68455 10596 0	18	413511	7 68455 13007 8	55
411803	7 68455 10597 7	18	413515	7 68455 13017 7	51
411804	7 68455 10598 4	18	413516	7 68455 13018 4	51
411805	7 68455 10599 1	18	413521	7 68455 13023 8	51
411806	7 68455 10600 4	18	413522	7 68455 13024 5	51
411807	7 68455 10601 1	18	413523	7 68455 13025 2	51
412009	7 68455 11116 9	94	413566	7 68455 13206 5	95
412010	7 68455 11119 0	94	413567	7 68455 13207 2	95
412011	7 68455 11122 0	94	413568	7 68455 13208 9	95
412012	7 68455 11125 1	94	413569	7 68455 13209 6	95
412017	7 68455 11128 2	94	413570	7 68455 13210 2	95
412018	7 68455 11131 2	94	413571	7 68455 13261 4	95
412019	7 68455 11134 3	94	416400	7 68455 11318 7	29
412020	7 68455 11137 4	94	416403	7 68455 11321 7	29
413153	7 68455 10614 1	22	416406	7 68455 11324 8	29
413154	7 68455 10615 8	22	416409	7 68455 11327 9	29
413155	7 68455 10616 5	22	416412	7 68455 11330 9	29
413156	7 68455 10617 2	22	416415	7 68455 11333 0	27
413157	7 68455 10618 9	22	416416	7 68455 11334 7	27
413159	7 68455 10619 6	21	416417	7 68455 11335 4	27

# Product Code Index

Product Code	UPC Codes	Page	Product Code	UPC Codes	Page
416418	7 68455 11336 1	27	416798	7 68455 11484 9	46
416419	7 68455 11337 8	27	416799	7 68455 11485 6	46
416420	7 68455 11338 5	27	416900	7 68455 11486 3	46
416421	7 68455 11339 2	27	416901	7 68455 11487 0	46
416422	7 68455 11340 8	27	416902	7 68455 11488 7	46
416423	7 68455 11341 5	27	416903	7 68455 11489 4	46
416424	7 68455 11342 2	27	416908	7 68455 11812 0	53
416472	7 68455 11572 3	27	416915	7 68455 11649 2	58
416700	7 68455 11170 1	55	416975	7 68455 16962 7	51
416701	7 68455 11171 8	55	416976	7 68455 16964 1	52
416703	7 68455 11173 2	56	418000	7 68455 11717 8	79
416704	7 68455 11174 9	56	419800	7 68455 11139 8	81
416705	7 68455 11175 6	57	420127	7 68455 11821 2	94
416706	7 68455 11176 3	56	420619	7 68455 12287 5	87
416707	7 68455 11177 0	56	420621	7 68455 12291 2	87
416709	7 68455 11179 4	56	420623	7 68455 12294 3	87
416710	7 68455 11180 0	56	420624	7 68455 12298 1	87
416712	7 68455 11182 4	56	420625	7 68455 12302 5	87
416713	7 68455 11183 1	56	420626	7 68455 12306 3	87
416715	7 68455 11185 5	56	420627	7 68455 12311 7	88
416718	7 68455 11188 6	53	420628	7 68455 12315 5	88
416719	7 68455 11189 3	53	420629	7 68455 12318 6	88
416721	7 68455 11191 6	53	420631	7 68455 12331 5	87
416724	7 68455 11194 7	54	420633	7 68455 12336 0	87
416725	7 68455 11195 4	54	420635	7 68455 12341 4	87
416728	7 68455 11198 5	54	420636	7 68455 12345 2	87
416729	7 68455 11199 2	54	420637	7 68455 12349 0	87
416732	7 68455 11202 9	54	420639	7 68455 12354 4	88
416733	7 68455 11203 6	54	420642	7 68455 12359 9	88
416737	7 68455 11207 4	54	420645	7 68455 12364 3	88
416738	7 68455 11208 1	54	420646	7 68455 12368 1	88
416739	7 68455 11209 8	54	420647	7 68455 12322 3	88
416740	7 68455 11210 4	54	420648	7 68455 12326 1	88
416741	7 68455 11211 1	54	420669	7 68455 12507 4	94
416742	7 68455 11212 8	54	420670	7 68455 12511 1	94
416743	7 68455 11213 5	54	420671	7 68455 12559 3	93
416744	7 68455 11214 2	54	420672	7 68455 12665 1	93
416745	7 68455 11215 9	54	420673	7 68455 12564 7	93
416746	7 68455 11216 6	54	420680	7 68455 12712 2	87
416747	7 68455 11217 3	54	420681	7 68455 12716 0	88
416748	7 68455 11218 0	54	420695	7 68455 12761 0	34
416749	7 68455 11219 7	54	420696	7 68455 12762 7	34
416750	7 68455 11220 3	54	420697	7 68455 12763 4	34
416751	7 68455 11221 0	54	420787	7 68455 12983 6	8
416774	7 68455 11460 3	47	420788	7 68455 12984 3	8
416775	7 68455 11461 0	47	420789	7 68455 12985 0	9
416778	7 68455 11464 1	47	420790	7 68455 12986 7	9
416779	7 68455 11465 8	47	420791	5 06014 8230085	8
416780	7 68455 11466 5	47	420804	7 68455 12909 6	87
416781	7 68455 11467 2	47	420805	7 68455 12913 3	88
416783	7 68455 11469 6	47	420806	7 68455 12917 1	88
416795	7 68455 11481 8	46	420807	7 68455 12921 8	88
416796	7 68455 11482 5	46	421033	7 68455 13307 9	17

# Product Code Index

Product Code	UPC Codes	Page	Product Code	UPC Codes	Page
421034	7 68455 13308 6	17	421680	7 68455 14977 3	29
421035	7 68455 13309 3	17	421681	7 68455 14989 6	29
421039	7 68455 13319 2	17	421682	7 68455 14978 0	29
421040	7 68455 13320 8	17	421683	7 68455 14990 2	29
421041	7 68455 13321 5	17	421687	7 68455 14928 5	56
421151	7 68455 13642 1	87	421688	7 68455 14930 8	56
421153	7 68455 13644 5	87	421689	7 68455 14932 2	56
421155	7 68455 13646 9	87	421703	7 68455 14583 6	78
421156	7 68455 13647 6	87	421738	7 68455 15666 5	28
421454	7 68455 13826 5	20	421739	7 68455 15667 2	28
421456	7 68455 13830 2	19	421740	7 68455 15668 9	28
421458	7 68455 13834 0	20	421746	7 68455 15670 2	28
421460	7 68455 13838 8	19	421747	7 68455 15671 9	28
421462	7 68455 13842 5	20	421748	7 68455 15672 6	28
421464	7 68455 13846 3	19	421749	7 68455 15673 3	28
421552	7 68455 14887 5	98	421796	7 68455 15235 3	30
421553	7 68455 14888 2	98	421797	7 68455 15236 0	30
421554	7 68455 14889 9	98	421798	7 68455 14979 7	29
421555	7 68455 14890 5	98	421816	7 68455 15026 7	57
421557	7 68455 14227 9	89	421821	7 68455 14926 1	55
421559	7 68455 14230 9	89	421822	7 68455 14934 6	56
421561	7 68455 14233 0	89	421823	7 68455 14936 0	56
421563	7 68455 14236 1	89	421848	7 68455 15316 9	52
421564	7 68455 20531 8	82	421853	7 68455 15569 9	52
421565	7 68455 20532 5	82	421858	7 68455 15552 1	52
421566	7 68455 20533 2	82	421893	7 68455 14991 9	29
421567	7 68455 20534 9	82	421907	7 68455 20559 2	82
421568	7 68455 20535 6	82	421908	7 68455 20560 8	82
421569	7 68455 20536 3	82	421909	7 68455 20561 5	82
421570	7 68455 20537 0	82	421910	7 68455 20562 2	82
421571	7 68455 20538 7	82	421911	7 68455 20563 9	82
421572	7 68455 20539 4	82	421912	7 68455 20564 6	82
421573	7 68455 20540 0	82	421921	7 6845 515953 6	98
421574	7 68455 20541 7	82	421927	7 68455 15113 4	89
421579	7 68455 14053 4	88	421928	7 68455 15116 5	89
421580	7 68455 14052 7	88	422155	7 68455 15433 3	98
421615	7 68455 14321 4	59	422156	7 68455 15434 0	98
421619	7 68455 14329 0	59	422157	7 68455 15435 7	98
421622	7 68455 14335 1	59	422159	7 68455 15254 4	59
421623	7 68455 14337 5	59	422285	7 68455 15504 0	98
421626	7 68455 14343 6	59	422291	7 68455 16976 4	78
421628	7 68455 14347 4	59	422303	7 68455 15715 0	78
421630	7 68455 14301 6	78	422308	7 68455 15719 8	98
421631	7 68455 14609 3	57	422356	7 68455 15727 3	88
421638	7 68455 14374 0	19	422357	7 68455 15728 0	88
421639	7 68455 14376 4	19	422358	7 68455 15729 7	88
421640	7 68455 14378 8	19	422359	7 68455 15730 3	88
421641	7 68455 14380 1	19	422363	7 68455 15712 9	60
421642	7 68455 14382 5	19	422364	7 68455 15713 6	60
421643	7 68455 14384 9	19	422365	7 68455 15714 3	60
421676	7 68455 14975 9	29	422551	7 68455 16547 6	60
421678	7 68455 14976 6	29	422557	7 68455 16553 7	60
421679	7 68455 14988 9	29	422562	7 68455 16558 2	60

# Product Code Index

Product Code	UPC Codes	Page	Product Code	UPC Codes	Page
423281	7 68455 17001 2	8	839263	7 68455 12680 4	74
423282	7 68455 17002 9	9	839264	7 68455 12682 8	74
423288	7 68455 17012 8	9	839265	7 68455 12681 1	74
423289	7 68455 17013 5	8	839266	7 68455 12686 6	74
423548	7 68455 19812 2	38	839267	7 68455 12683 5	75
423549	7 68455 19813 9	39	839268	7 68455 12684 2	75
423550	7 68455 19814 6	39	839280	7 68455 12685 9	75
423551	7 68455 19815 3	39			
423552	7 68455 19816 0	40			
423553	7 68455 19817 7	40			
423554	7 68455 19818 4	38			
423555	7 68455 19819 1	38			
423556	7 68455 19820 7	38			
423557	7 68455 19821 4	39			
423558	7 68455 19822 1	39			
423559	7 68455 19823 8	39			
446650	7 68455 14907 0	98			
501002	7 68455 11140 4	81			
501003	7 68455 11141 1	81			
501004	7 68455 11142 8	81			
501005	7 68455 11143 5	81			
501006	7 68455 11144 2	81			
501011	7 68455 11149 7	81			
501012	7 68455 11150 3	81			
501013	7 68455 11151 0	81			
501014	7 68455 11152 7	81			
501015	7 68455 11153 4	81			
501016	7 68455 11154 1	81			
501019	7 68455 11157 2	81			
501020	7 68455 11158 9	81			
501021	7 68455 11159 6	81			
501022	7 68455 11160 2	81			
501023	7 68455 11161 9	81			
650803	7 68455 10875 6	65			
650828	7 68455 10661 5	66			
650829	7 68455 10662 2	66			
650830	7 68455 10663 9	66			
650831	7 68455 10664 6	66			
650832	7 68455 10665 3	66			
650833	7 68455 10666 0	66			
839001	7 68455 10872 5	11			
839002	7 68455 10873 2	11			
839003	7 68455 10933 3	11			
839004	7 68455 10934 0	11			
839005	7 68455 11825 0	11			
839006	7 68455 13917 0	11			
839010	(01)05060267416704	10			
839250	7 68455 12687 3	73			
839251	7 68455 12688 0	73			
839252	7 68455 12689 7	73			
839260	7 68455 12690 3	73			
839261	7 68455 12691 0	73			
839262	7 68455 12692 7	74			



# Autorisation de retour de marchandise Authorization to return merchandise

Numéro d'autorisation  
OBLIGATOIRE

Les retours seront acceptés seulement si ce numéro d'autorisation aura été obtenu préalablement du bureau des commandes de Convatec Canada au 1-866-987-1471. Returns will only be accepted if this Authorization Number has been obtained from Convatec Canada Order Desk at 1-866-987-1471.

N°

Authorization number  
MANDATORY

NOM DU CLIENT  
CUSTOMER NAME

N° DE COMPTE  
ACCOUNT N°.

ADRESSE  
ADDRESS

CODE POSTAL  
POSTAL CODE

N° DE CODE CODE N°.	QUANT.	DESCRIPTION	DATE DE PÉREMPTION EXP. DATE	N° DE LOT LOT N°.

**REMARQUE:** TOUTES LES MARCHANDISES RETOURNÉES DOIVENT ÊTRE PORT PAYÉES, SINON LES FRAIS DE TRANSPORT SERONT DÉDUITS DE LA NOTE DE CRÉDIT. **NOTE:** ALL RETURNED SHIPMENTS MUST BE PREPAID OTHERWISE TRANSPORTATION CHARGES WILL BE DEDUCTED FROM THE CREDIT NOTE. 

COMMENTAIRES:  
COMMENTS:

APPROUVÉ PAR REP:  
APPROVED BY REP:

(LETTRES MOULÉES) / (PLEASE PRINT)

No DE TERRITOIRE:  
TERRITORY No.

DATE:

SIGNATURE DU CLIENT / CUSTOMER SIGNATURE

ADRESSE DE RETOUR:  
RETURN ADDRESS:  
**Convatec Canada  
Bureau des commandes /  
Order Desk  
1375, Newton  
Boucherville, QC  
J4B 5H2**

À REMPLIR PAR LE SERVICE DE LA RÉCEPTION / TO BE COMPLETED BY RECEIVING DEPT.	
DATE	NOMBRE DE BOÎTES: NO. OF CTNS:
PAR: VIA:	
<input type="checkbox"/> PORT PAYÉ PREPAID:	<input type="checkbox"/> PORT DÛ COLLECT FRAIS DE TRANSPORT: TRANSPORT CHARGES:
ÉTAT DE LA MARCHANDISE: CONDITION OF GOODS:	
SIGNATURE: SIGNED:	N° DE LA NOTE: MEMO N°:



# Healthy bonds keep us strong

**The strength of our adhesives  
creates a healthy bond between  
your patient and their pouch**

**For our strongest bond, Esteem<sup>®</sup>+ pouches  
with Durahesive<sup>®</sup> Plus are designed to:**

- Bond to the skin with enhanced tack and cohesion
- Provide a long-lasting seal for extended wear time, even with high liquid output




**ESTEEM<sup>®</sup>** 

**One-piece Drainable Pouch**  
with Durahesive<sup>®</sup> Plus  
Adhesive Technology



# FlexiSeal<sup>®</sup> PROTECT



Designed to protect patients and clinicians against over-inflation risk



**convatec**  
— forever caring —

# Explanation of Convatec Package Code Numbering System

## LOT NUMBERS

Every Convatec product carries a lot number that serves to indicate when the item was produced, and from which a complete history as well as expiration date of the product can be derived.

**Most products use the following format:**

Example - Lot number **9J00978**

**9** = Year (2019)

**J** = Month (September)

**00978** = Lot number assignment for positive identification of each lot produced

### Expiration Designation:

- If an expiration date is noted on the box, it will appear below the lot number.
- If an expiration date is not noted on the box, it can be derived by adding five years to date of manufacture of product. See above for details on manufacture dates.

**Example:** If date of manufacture is determined as September 2009 (ex. 9J), the product's shelf life will be until September 2014; however, limit to return the product is six months prior to this date which would be: March 1, 2014 .

A - January	B - February	C - March	D - April
E - May	F - June	G - July	H - August
J - September	K - October	L - November	M - December

Please contact your Convatec Territory Manager or the **Convatec Customer Relations Center** at **1-800-465-6302** for clarification of control numbers or expiration dates.



# Product Catalog 2023 2024

For more information, please call our Customer Relations Center  
(Registered Nurses on staff) at **1-800-465-6302**, Monday through Friday,  
8:00 AM to 6:00 PM (EST), or visit **[www.convatec.ca](http://www.convatec.ca)**