



ConvaTec |

Ostomy Care

PRODUCTS, CARE *and* SUPPORT
SHAPED AROUND YOUR LIFE

MOLDABLE Skin Barriers

ConvaTec Moldable Technology™ Skin Barriers are designed to hug the stoma for a more custom and secure fit.

Features and Benefits

- ✓ Reduced application steps—fewer steps than traditional cut-to-fit skin barriers
- ✓ No cutting needed
- ✓ Memory – flexible materials gently expand and contract with the stoma
- ✓ Turtlenecking – the adhesive gently swells when in contact with output to help protect the skin and secure the seal.




95% of people who started on ConvaTec Moldable Technology kept their skin healthy¹


86% who switched to ConvaTec Moldable Technology improved their skin condition¹



Rolls, molds, and holds to the stoma's shape and size



Turtlenecks around the stoma to protect against leaks



Rebounding Memory Technology™ rolls back to achieve a custom fit

NATURA® DURAHEESIVE® SKIN BARRIERS

Moldable | Qty = 10

Flange Size	Stoma Size	Acrylic Tape Collar WHITE	Hydrocolloid Tape Collar WHITE
1 3/4" (45mm)	13-22mm (1/2 - 7/8")	413417	411800
1 3/4" (45mm)	22-33mm (7/8 - 1 1/4")	413418	411802
2 1/4" (57mm)	33-45mm (1 1/4 - 1 3/4")	413419	411804
2 3/4" (70mm)	33x45mm by 45x56mm (1 1/4 - 1 3/4" by 1 3/4 - 2 1/8") Oval shape	413420	411806

Moldable | Accordion Flange | Qty = 10

Flange Size	Stoma Size	Hydrocolloid Tape Collar TAN
2 1/4" (57mm)	13-22mm (1/2 - 7/8")	421039
2 1/4" (57mm)	22-33mm (7/8 - 1 1/4")	421040
2 3/4" (70mm)	33-45mm (1 1/4 - 1 3/4")	421041

Moldable | Convex | Qty = 10

Flange Size	Stoma Size	Acrylic Tape Collar WHITE
1 3/4" (45mm)	13-22mm (1/2 - 7/8")	404592
1 3/4" (45mm)	22-33mm (7/8 - 1 1/4")	404593
2 1/4" (57mm)	33-45mm (1 1/4 - 1 3/4")	404594

NATURA® STOMAHESIVE® SKIN BARRIERS

Moldable | Qty = 10

Flange Size	Stoma Size	Acrylic Tape Collar WHITE	Hydrocolloid Tape Collar WHITE
1 3/4" (45mm)	13-22mm (1/2 - 7/8")	413421	411801
1 3/4" (45mm)	22-33mm (7/8 - 1 1/4")	413422	411803
2 1/4" (57mm)	33-45mm (1 1/4 - 1 3/4")	413423	411805
2 3/4" (70mm)	33x45mm by 45x56mm (1 1/4 - 1 3/4" by 1 3/4 - 2 1/8") Oval shape	413424	411807

Moldable | Accordion Flange | Qty = 10

Flange Size	Stoma Size	Hydrocolloid Tape Collar TAN
2 1/4" (57mm)	13-22mm (1/2 - 7/8")	421033
2 1/4" (57mm)	22-33mm (7/8 - 1 1/4")	421034
2 3/4" (70mm)	33-45mm (1 1/4 - 1 3/4")	421035

ESTEEM SYNERGY® SKIN BARRIERS

Moldable | Stomahesive® | Qty = 10

Coupling Size	Stoma Size	Hydrocolloid Tape Collar WHITE
Small Stoma Opening	13-22mm (1/2 - 7/8")	411657
Medium Stoma Opening	22-33mm (7/8 - 1 1/4")	411658
Large Stoma Opening	33-45mm (1 1/4 - 1 3/4")	411659

Moldable | Convex | Durahesive® | Qty = 10

Coupling Size	Stoma Size	Acrylic Tape Collar WHITE
Small Stoma Opening	13-22mm (1/2 - 7/8")	409268
Medium Stoma Opening	22-33mm (7/8 - 1 1/4")	409269
Large Stoma Opening	33-45mm (1 1/4 - 1 3/4")	409270

ESTEEM®+ POUCHES (One-Piece Systems)

Moldable | Stomahesive® | Drainable | InvisiClose® Tail Closure | Qty = 10 | Pouch Length = 30.5cm (12")

No Filter

Stoma Size	Transparent
20-30mm (13/16 - 1 3/16")	413518
30-40 mm (1 3/16 - 1 9/16")	413519

Filter (10 filter covers per box)

Stoma Size	Transparent	Opaque w/ Easy-View Window
20-30mm (13/16 - 1 3/16")	413515	413521
30-40mm (1 3/16 - 1 9/16")	413516	413522

Moldable | Stomahesive® | Closed-End

Qty = 30 | Pouch Length = 20.3cm (8")

Filter (30 filter covers per box)

Stoma Size	Opaque w/ Easy-View Window
20-30mm (13/16 - 1 3/16")	413509
30-40mm (1 3/16 - 1 9/16")	413510
40-50mm (1 9/16 - 1 15/16")	413511



support

For more information, please call our Customer Relations Center (Registered Nurses on staff) at **1-800-465-6302**, Monday through Friday, 8:00 AM to 6:00 PM (EST), or visit our Web Site at **www.convatec.ca**

Reference 1. Maria Teresa Szweczyk, MD, PhD; Grazyna Majewska, RN, ETN; Mary V.Cabral, MS, FNP-BC, CWOCN-AP; and Karin HJzel-Piontek, RN; The Effects of Using a Moldable Skin Barrier on Peristomal Skin Condition in Persons with an Ostomy: Results of a Prospective, Observational, Multinational Study, Ostomy Wound Management 2014;60(12):16-26.

All trademarks are property of their respective owners.
©2020 ConvaTec Inc. AP-030196-CA V1603



ConvaTec

Ostomy Care