

# Take the pressure off of pouch application




Natura® Cut-to-Fit Accordion Flange lifts easily to provide generous finger room so coupling is easy and comfortable. Designed for tender abdomens or for anyone who wants everyday comfort.

Natura® Cut-to-Fit Accordion Flange

# NOW Natura® Cut-to-Fit Accordion Flange

Designed to make pouch application easy and comfortable



## Comfortable and discreet:


- In the “UP” position the Accordion Flange provides full finger width space so pouch can be attached with minimal pressure on the abdomen
- In the “DOWN” position it gives the wearer a flat low profile that is extremely flexible and conforms to body contours



“UP” position



“DOWN” position




## With the control of a cut-to-fit barrier:

- It's simple to measure and cut to create an opening tailored to each individual patient's stoma size
- Available in Stomahesive® and Durahesive® adhesive options



## Compatible with a range of Natura® pouch options:

- Pair the Accordion Flange with our state-of-the-art Natura®  comfort pouches in drainable and closed-end options

## Ordering Information

### Natura® Accordion Flange: Cut-to-Fit

| Product Code | Color Code | Description   | Stoma Size              | Flange Size   | Qty Per Box |
|--------------|------------|---|-------------------------|---------------|-------------|
| 421454       | ●          | Durahesive® Skin Barrier with Acrylic Collar 45mm (1 3/4")  | 13-21mm (1/2" - 3/4")   | 45mm (1 3/4") | 10          |
| 421456       | ●          | Stomahesive® Skin Barrier with Acrylic Collar 45mm (1 3/4") | 13-21mm (1/2" - 3/4")   | 45mm (1 3/4") | 10          |
| 421458       | ●          | Durahesive® Skin Barrier with Acrylic Collar 57mm (2 1/4")  | 13-33mm (1/2" - 1 1/4") | 57mm (2 1/4") | 10          |
| 421460       | ●          | Stomahesive® Skin Barrier with Acrylic Collar 57mm (2 1/4") | 13-33mm (1/2" - 1 1/4") | 57mm (2 1/4") | 10          |
| 421462       | ●          | Durahesive® Skin Barrier with Acrylic Collar 70mm (2 3/4")  | 13-45mm (1/2" - 1 3/4") | 70mm (2 3/4") | 10          |
| 421464       | ●          | Stomahesive® Skin Barrier with Acrylic Collar 70mm (2 3/4") | 13-45mm (1/2" - 1 3/4") | 70mm (2 3/4") | 10          |



For more information, please call our Customer Relations Center (Registered Nurses on staff) at **1-800-465-6302**, Monday through Friday, 8:00 AM to 6:00 PM (EST), or visit our Web Site at [www.convatec.ca](http://www.convatec.ca)