2020

OSTOMY Product Guide

Shaping Lives with Complete
Ostomy Products, Care and Support







Innovation begins by listening and working together.

ConvaTec is a global medical products and technologies company focused on therapies for the management of chronic conditions. Our portfolio of ostomy products, me+™ patient support program and collaboration with health care professionals help to improve the lives of the people we touch.

TABLE OF CONTENTS

me+™ SUPPORT	1
OSTOMY OVERVIEW	5
Introduction to Ostomy Care	5
Choosing the Right Ostomy System	7
GUIDE TO OSTOMY SYSTEMS	9
A Guide to Adhesives	9
A Guide to Skin Barriers	11
Moldable Skin Barriers	13
Accordion Skin Barriers	15
The Purpose of Convexity	17
Ostomy Pouch Options	21
Common Pouch Features	23
TWO-PIECE SYSTEMS	25
• Natura [™] , Natura [™] +, and Sur-FIT® Natura [™] Systems	27
• Esteem Synergy® Systems	41
ONE-PIECE SYSTEMS	45
• Esteem™+ Pouches	46
ActiveLife® Pouches	55

ACUTE KITS	60
PEDIATRIC PRODUCTS	71
ACCESSORIES	77
WOUND & FISTULA POUCHES	89
OSTOMY RESOURCES & SUPPORT	92
MEDICARE UTILIZATION CHART	93
PRODUCT CODE INDEX	94
REFERENCES	98

THERE IS A NEED FOR ONGOING PATIENT SUPPORT

The **2nd most common** reason for hospital readmissions is ileostomy and other enterostomy surgery¹



of readmissions following colorectal surgery may be preventable²



of ostomy patients exhibited peristomal skin complications within 90 days of surgery³



of respondents reported a need for adequate information and communication post-op⁴

Most common ostomy self-care challenges are related to leakage, skin issues, frequency of changes and difficulty caring for an ostomy⁵



THE **me+**™PROGRAM SUPPORTS PATIENTS

The me+[™] program is a comprehensive support program designed to complement clinical care by providing patients with support, inspiration, insights and ideas throughout their entire ostomy journey.

me+[™] program members can:

- Request free ostomy samples
- Speak to a WOC nurse or product specialist
- ♥ Connect to our me+ community
- Receive ongoing lifestyle education, support and much more

TO ENROLL YOUR PATIENTS in the me+ program

call us at 1-800-422-8811 or visit mepluscare.com



Our goal is to partner with you and your patients to deliver best-in-class care throughout the patient journey:

- **⊘** Pre-Surgery / Surgery
 - Access online resources at mepluscare.com
 - Request me+[™] patient Starter Kit and education materials
 - Enroll in me+ program
 - Learn about the me+ recovery series, developed to aid rehabilitation after surgery
- **⊘** Weeks Immediately After Surgery
 - Consult with me+ WOC nurses and product specialists
 - Personalized Starter Kit
 - me+ welcome email and support
 - me+ recovery phase 1 exercise booklet*
- First Year & Beyond
 - Explore online resources and tools
 - Get inspired by the me+ community
 - Sample new products
 - Continue me+ recovery series exercises*

Throughout their journey, your patients can call 1-800-422-8811 for live support from me+™ WOC nurses and product specialists

A me+™ Concierge team member can assist Healthcare professionals with:

- Personalized support
- Enrolling your patients into the me+[™] program
- Providing a me+ Starter Kit and free product samples for your patients

TO LEARN MORE contact your me+ concierge directly at 877-585-0470



"You guys have been a rock for me. Your expert insight and help picking the right product for my needs has been crucial at setting my mind at ease, and feeling that even with this significant life change that it's all going to be alright." - Sofia[†]

 $^{^\}dagger Model$ portrayal for illustrative purposes only.

An Introduction to OSTOMY CARE

There are 3 primary types of ostomy surgeries. Understanding the fundamental difference is important to select the right product and provide the right support.



COLOSTOMY (large intestine)

A colostomy is a surgical procedure that creates an opening on the abdomen through which the large intestine is redirected to bypass diseased or damaged parts of the large intestine. The large intestine is brought to the surface of the skin, and stitched to the abdomen to form a new exit for the stool from the bowel.

Things to know about a Colostomy

- Primarily formed due to colon or rectal cancer or trauma
- Can be temporary or permanent
- Stool tends to be semi-formed to formed and predictable 1-2 times per day

ILEOSTOMY (small intestine)

An ileostomy is a surgical procedure that creates an opening on the abdomen through which the end of your small intestine is redirected to bypass diseased or damaged parts of the digestive system. The small intestine is brought to the surface of the skin, and stitched to your abdomen to form a new exit for waste from your bowels.

Things to know about an Ileostomy

- Primarily formed due to inflammatory bowel disease (Crohn's, Ulcerative Colitis)
- Can be temporary or permanent
- Stool is liquid or pasty and constant throughout the day
- The small intestine contains digestive enzymes and acids, and the output may cause skin irritation if it comes in contact with the skin

UROSTOMY (urinary tract)

A urostomy is a surgical procedure that creates an opening on the abdomen through which urine is redirected away from a bladder that is damaged or diseased. After the bladder is removed, a piece of small intestine (ileum) is separated from the rest of the bowel. This piece is used as a conduit for urine to pass out of the body. The ureters are attached to this new passageway and one end is brought up to an opening made on your abdomen.

Things to know about a Urostomy

- Primarily formed due to cancer or inflammation of the bladder
- (Is always permanent
- Urine flows through the stoma as it is produced which means output is constant
- Urine may cause skin irritation if it comes in contact with the skin

Choosing the Right OSTOMY SYSTEM

The ostomy pouch collects stoma output. The pouch comes as part of a two-piece or one-piece system.



2 PIECE SYSTEM

The two-piece system includes a skin barrier and a separate pouch that securely join together. This option allows the pouch to be removed to empty or change, without removing the skin barrier as frequently.



Designed for patients

- who want to change their pouch more frequently than their skin barrier
- with an ileostomy who need to empty/ change the pouch frequently
- who exercise, travel and play sports

1 PIECE SYSTEM

The skin barrier and pouch are provided as an integrated solution that cannot be separated. It offers one simple and convenient step to attach or remove, providing a fresh start every time. This option is thin and flexible, and designed for more frequent system changes.



Designed for patients

- on the abdomen caused when connecting a skin barrier and pouch
- who prefer to change their ostomy system frequently
- who prefer a flexible and low profile ostomy system

A Guide to ADHESIVES

The adhesive is an important component of the skin barrier. Our product portfolio features a range of adhesive technologies to enable you to select an appropriate solution that meets the needs of your patients.



STOMAHESIVE®

Designed for people with mostly solid output, our unique Stomahesive® technology is a hydrocolloid adhesive that gently adheres to the skin, protecting it from the damaging effects of stoma effluent.

We also offer a modified Stomahesive® option, which features a touch of mineral oil that makes the barrier more flexible during wear, and removal easier on the skin.

Designed for patients

with semi-formed/formed output

who have sensitive skin

with frequent changes

who prefer a standard weartime

DURAHESIVE®

Designed for people with mostly liquid output, our unique Durahesive® technology is a hydrocolloid adhesive that swells up upon hydration to protect the stoma. This special effect called turtlenecking, helps create a more secure seal without harming the skin.

We also offer Durahesive® Plus, which features an enhanced tack, strength and cohesion designed for extended weartime or aggressive output.

Designed for patients

with semi-liquid/liquid output, including urine

who want or need longer weartime

ESTEEM™+ Soft & Flex CONVEX

Adhesive designed with a high initial tack, to aid in the initial attachment to the skin. The hydrocolloid barrier features a thin tapered edge to help prevent the edge catching on clothing, and a thicker center to help absorption around the stoma.⁶

Designed for patients

with both solid and liquid output

who prefer a standard weartime

TAPE COLLARS are available in two forms:

Acrylic: Tape collar designed with good tack **Hydrocolloid:** Designed for those with sensitivity to other tape collars

A Guide to SKIN BARRIERS

The skin barrier is the piece of the pouching system that sticks to the body. It should fit snugly around the stoma to help prevent leakage and have an opening customized to stoma size and shape.



CUT-TO-FIT

Cut-to-Fit barriers can be cut to match the shape and size of individual stomas. Using scissors, one can cut a circle, oval or any shape needed to fit the stoma.



Designed for patients

with a stoma size in-between pre-cut sizes

with an irregular stoma shape

PRE-CUT

Pre-cut barriers have a uniform, circular stoma hole and are ready to apply right most common stoma sizes.



out of the box. They are available for the



MOLDABLE

ConvaTec Moldable Technology™ skin barrier with Rebounding Memory Technology[™] creates an elastic-like seal that rebounds to fit most shapes and sizes, and then "turtlenecks" upward, to enhance the seal and protect the skin.7



Designed for patients

with an evenly round stoma that fits the available pre-cut size

with dexterity issues who find cutting challenging

Designed for patients

who have large, round, oval or irregular shaped stomas

with dexterity issues who find cutting challenging

MOLDABLE Skin Barriers

ConvaTec Moldable Technology[™] Skin Barriers are designed to hug the stoma opening for a more custom and secure fit.

CLINICALLY PROVEN


to prevent and improve peristomal skin complications⁷

Up to 55% of ostomates have skin issues8

- 95% of people who started on ConvaTec Moldable Technology kept their skin healthy⁷
- 86% of people who switched to ConvaTec Moldable Technology improved their skin condition⁷

Features and Benefits


- Reduced application steps—fewer steps than traditional cut-to-fit skin barriers
- No cutting needed
- Memory flexible materials gently expand and contract with the stoma
- Turtlenecking the adhesive gently swells when in contact with output & stool to help protect the skin and secure the seal.



Roll, mold, and hold to the stoma's shape and size



Turtlenecks around the stoma to protect against leaks



Rebounding Memory Technology™- rolls back to achieve a custom fit

ACCORDION Skin Barriers

The Accordion Flange is an option on a range of Natura® skin barriers.
The accordion flange lifts to provide generous finger room so coupling is easy and comfortable.



PRODUCT FEATURES


- 2X more finger space than any other lifting flange⁹
- · Simple to measure and cut or mold
- Range of flat, convex, Cut-to-Fit and Moldable options
- Compatible with SUR-FIT® Natura™ and Natura™+ pouch option

Designed for patients

- with dexterity issues as it provides finger room so coupling is easy and comfortable
- designed to prevent pressure or trauma if abdominal tenderness is present due to post-operative stoma surgery/incisions

HOW IT WORKS

- Designed to help minimize abdominal pressure and pain when attaching the pouch.
- Connects with an audible 'click' and physical 'snap' to confirm the pouch is firmly attached.



The Purpose of CONVEXITY



WHAT IS CONVEXITY?

All bodies are different. Convex wafers are designed to help manage irregular stomas and abdomens to help reduce leakage and simplify stoma management.

- Convexity is the outward curving of the baseplate or skin barrier.^{10,11}
- Convexity applies a uniform pressure to the immediate peristomal skin, improving the projection of the stoma for drainage.¹²
- Convexity helps maintain a secure seal between the pouch and the skin.¹³

GROWING CONVEXITY NEED

Data shows the need for convex skin barriers growing year on year.¹⁴ This is due to a number of factors including:

- Changing surgical techniques/outcomes of stoma shape & profile
- Increased obesity rates

WHEN TO USE CONVEXITY



Flush Stomas
Located at the same level as
the abdomen skin surface,
and does not protrude¹⁵



Skin Folds
A stoma which is deep in a skin fold/crease¹³



Retracted Stomas Stoma pulled into the abdomen below skin level or having little or no spout¹⁶



Soft or Flaccid Abdomens An abdomen requiring firm level of support to stabilize the pouching system¹⁵

The Characteristics of CONVEXITY


RIGIDITY

Firm/Rigid Convexity

Typically used on soft or flaccid abdomens.

Soft Convexity

Typically used on firm or toned abdomens.



DEPTH

Height of the convex plateau.








PLATEAU WIDTH

Distance across the widest flat part of the outward convex curve.



Why Plateau Width Matters

Successful convexity must press down on the peristomal skin enough to make a good seal around the base of the stoma and help with protrusion.¹³

- Convexity Pressure: A narrow plateau width may be used to create pressure directly around the stoma with lesser pressure as the plateau width gets wider.
- Plateau Width: A wider plateau width may be used to flatten out skin creases away from the stoma.

Choosing the right convexity solution for different body profiles can be challenging.

Please feel free to consult with one of our me+™ Specialists or WOC Nurses for further guidance.

Ostomy Pouch OPTIONS

ConvaTec offers a wide range of pouching options across both 2 Piece and 1 Piece Systems.



CLOSED

A closed pouch doesn't require draining. Once they are full, you simply remove and replace with a new pouch.

Additionally, ideal during physical activity, intimate moments, swimming, and travel.



DRAINABLE

Drainable pouches allow you to drain the output periodically throughout the day so that you do not have to change pouches as frequently.



Designed for patients

Designed for patients

(1-2 times per day)

who have loose, liquid stool or semi-formed stool that can be unpredictable or constant throughout the day

who have thick, formed stool that is predictable

UROSTOMY

A Urostomy pouch has a faucet tap, soft tap, or fold over tap at the bottom that allows for draining urine as frequently as needed.



Designed for patients

with a Urostomy


POUCH FEATURES



POUCH COLOR

DRAINABLE POUCH CLOSURE

Transparent (Clear front of pouch)






Soft and quiet material cover both the front and the body side of the pouch for more discretion and comfort.



Easy-View Window

Some options come with an easy-view window for accurate positioning and observation.







InvisiClose™ [™] Clipless Tail Closure

Self closing drainable pouch closure which provides an easy closure option that does not require a separate clip for added security.



Lock-it Pocket™

Folds inside out to hide away the secured drainage tail—providing the flexibility of a drainable pouch, with the feel and comfort of a closed-end pouch.

Traditional drainable pouches are closed with a separate clip.

UROSTOMY POUCH DRAINAGE



The Accuseal™ Tap with valve is a faucet type tap designed for easy manipulation and emptying. Other drainage options include a Soft Tap or Fold Over Tap.

All Urostomy pouches include an anti-reflux valve which is designed to help prevent backflow.

TWO-PIECE Systems

Two-piece systems include a **skin barrier and** a separate pouch that securely join together. There are two brands that make up ConvaTec's Two-Piece Ostomy Systems: Natura™ + and Esteem Synergy® +.





The Natura[™] brand (inclusive of Natura, Sur-Fit Natura, and Natura+) is an extensive range of skin barriers and pouches that connect together with a plastic ring called a flange in a simple snap-off, snap-on design.



ESTEEM SYNERGY® (





The Esteem Synergy®+ brand are skin barriers and pouches that connect together with Adhesive Coupling Technology™. This is a flat, flexible adhesive foam ring on the pouch that sticks to a flat 'landing zone' on the skin barrier.

Color Key & Sizing Guide

COLOR MATCH						
FLANGE	1 1/4"	1 1/2"	1 3/4"	2 1/4"	2 1/4"	4"
SIZE	(32mm)	(38mm)	(45mm)	(57mm)	(57mm)	(100mm)

Color Key & Sizing Guide

COLOR MATCH		•			
COUPLING	Small	Medium	Large	Extra-Large	
SIZE	Stoma Opening	Stoma Opening	Stoma Opening	Stoma Opening	



When ordering a Two-Piece System, make sure you select the same flange/coupling size for the skin barrier and pouch. Use the colored dots in the table to help.

NATURA™ DURAHESIVE® SKIN BARRIERS

Compatible with Natura™ + and SUR-FIT® Natura™ pouches







Cut-to-Fit

- Durahesive® (Extended Wear)
- Quantity = 10 (*Box of 5)

HCPCS A4409 (413157 and 401905 are HCPCS A4410)

	Flange Size	Cutting Range	Acrylic Tape Collar WHITE	Acrylic Tape Collar TAN	Without Tape Collar
•	1 1/4" (32mm)	1/2 - 3/4"	413159	413164	413153
•	1 1/2" (38mm)	1/2 - 1"	413160	413165	413154
•	1 3/4" (45mm)	1/2 - 1 1/4"	413161	413166	413155
•	2 1/4" (57mm)	1/2 - 1 3/4"	413162	413167	413156
•	2 3/4" (70mm)	1/2 - 2 1/4"	413163	413168	413157
•	4" (100mm)	1/2 - 3 1/2"	N/A	N/A	401905*





Cut-to-Fit | Accordion Flange

- Durahesive® (Extended Wear)
- Quantity = 10

	Flange Size	Cutting Range	Acrylic Tape Collar WHITE
•	1 3/4" (45mm)	1/2 - 13/16" (13mm-21mm)	421454
•	2 1/4" (57mm)	1/2 - 1 5/16" (13mm-33mm)	421458
•	2 3/4" (70mm)	1/2 - 1 3/4" (13mm-45mm)	421462



NATURA™ SKIN BARRIER

Moldable | Accordion Flange

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4409

	Flange Size	Stoma Size	Hydrocolloid Tape Collar TAN
•	2 1/4" (57mm)	1/2 - 7/8" (13-22mm)	421039
•	2 1/4" (57mm)	7/8 - 1 1/4" (22-33mm)	421040
•	2 3/4" (70mm)	1 1/4 - 1 3/4" (33-45mm)	421041

Recommended Product



NATURA™ SKIN BARRIER

Cut-to-Fit | Convex | Accordion Flange

- Durahesive® (Extended Wear)
- Quantity = 10

		HCPCS A4407		
	Flange Size	Cutting Range	Acrylic Tape Collar TAN	Hydrocolloid Tape Collar TAN
•	1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	421638	421639
•	2 1/4" (57mm)	1/2 - 1 3/8" (13-35mm)	421640	421641
•	2 3/4" (70mm)	1/2 - 1 7/8" (13-48mm)	421642	421643



NATURA™ SKIN BARRIER

Moldable

- Durahesive® (Extended Wear)
- Quantity = 10

Η	C	PC	S	Δ	44	10	9

	Flange Size	Stoma Size	Tape Collar WHITE	Tape Collar WHITE
•	1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	413417	411800
•	1 3/4" (45mm)	7/8 - 1 1/4" (22-33mm)	413418	411802
•	2 1/4" (57mm)	1 1/4 - 1 3/4" (33-45mm)	413419	411804
•	2 3/4" (70mm)	1 1/4 - 1 3/4" by 1 3/4 - 2 1/8" (33x45mm by 45x56mm) Oval shape	413420	411806



NATURA™ SKIN BARRIER

Moldable | Convex

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4407

	Flange Size	Stoma Size	Acrylic Tape Collar WHITE
•	1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	404592
•	1 3/4" (45mm)	7/8 - 1 1/4" (22-33mm)	404593
•	2 1/4" (57mm)	1 1/4 - 1 3/4" (33-45mm)	404594

SUR-FIT® NATURA™ SKIN BARRIER

Pre-cut | Convex

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4407

Ultrafina and Hadal

	Flange Size	Stoma Size	Acrylic Tape Collar WHITE
•	1 3/4" (45mm)	1/2" (13 mm)	413177
•	1 3/4" (45mm)	5/8" (16 mm)	413178
•	1 3/4" (45mm)	3/4" (19 mm)	413179
•	1 3/4" (45mm)	7/8" (22 mm)	413180
•	1 3/4" (45mm)	1" (25 mm)	413181
•	1 3/4" (45mm)	1 1/8" (28 mm)	413182
•	1 3/4" (45mm)	1 1/4" (32 mm)	413183
•	1 3/4" (45mm)	1 3/8" (35 mm)	413184
•	2 1/4" (57mm)	1 1/2" (38 mm)	413185
•	2 1/4" (57mm)	1 5/8" (41 mm)	413186
•	2 1/4" (57mm)	1 3/4" (45 mm)	413187
•	2 1/4" (57mm)	2" (50 mm)	413188

NATURA™ STOMAHESIVE® SKIN BARRIERS

Compatible with Natura™ + and SUR-FIT® Natura™ pouches





SUR-FIT® NATURA™ SKIN BARRIER

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Quantity = 10 (*Box of 5)

HCPCS A4414

(401577 and 401906 are HCPCS A4415)

	Flange Size	Cutting Range	Acrylic Tape Collar WHITE	Acrylic Tape Collar TAN	Without Tape Collar
•	1 1/4" (32mm)	1/2 - 3/4"	125257	125262	401573
•	1 1/2" (38mm)	1/2 - 1"	125258	125263	401574
•	1 3/4" (45mm)	1/2 - 1 1/4"	125259	125264	401575
•	2 1/4" (57mm)	1/2 - 1 3/4"	125260	125265	401576
•	2 3/4" (70mm)	1/2 - 2 1/4"	125261	125266	401577
•	4" (100mm)	1/2 - 3 1/2"	N/A	N/A	401906*



SUR-FIT® NATURA™ SKIN BARRIER

Pre-cut

- Stomahesive® (Standard Wear)
- Quantity = 10

	A44	

	Flange Size	Stoma Size	Acrylic Tape Collar TAN
•	1 3/4" (45mm)	1/2" (13 mm)	125267
•	1 3/4" (45mm)	5/8" (16 mm)	125268
•	1 3/4" (45mm)	3/4" (19 mm)	125269
•	1 3/4" (45mm)	7/8" (22 mm)	125270
•	1 3/4" (45mm)	1" (25 mm)	125271
•	1 3/4" (45mm)	1 1/8" (28 mm)	125272
•	1 3/4" (45mm)	1 1/4" (32 mm)	125273
•	1 3/4" (45mm)	1 3/8" (35 mm)	125274
•	2 1/4" (57mm)	1 1/2" (38 mm)	125275
•	2 1/4" (57mm)	1 5/8" (41 mm)	125276
•	2 1/4" (57mm)	1 3/4" (45 mm)	125277



MATURA™ SKIN BARRIER

Cut-to-Fit | Accordion Flange

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

	Flange Size	Cutting Range	Acrylic Tape Collar WHITE
•	1 3/4" (45mm)	1/2 - 13/16" (13-21mm)	421456
•	2 1/4" (57mm)	1/2 - 1 5/16" (13-33mm)	421460
•	2 3/4" (70mm)	1/2 - 1 3/4" (13-45mm)	421464





NATURA™ SKIN BARRIER

Moldable | Accordion Flange

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

	Flange Size	Stoma Size	Hydrocolloid Tape Collar TAN
•	2 1/4" (57mm)	1/2 - 7/8" (13-22mm)	421033
•	2 1/4" (57mm)	7/8 - 1 1/4" (22-33mm)	421034
•	2 3/4" (70mm)	1 1/4 - 1 3/4" (33-45mm)	421035



NATURA™ SKIN BARRIER

Moldable

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

	Flange Size	Stoma Size	Acrylic Tape Collar WHITE	Hydrocolloid Tape Collar WHITE
•	1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	413421	411801
•	1 3/4" (45mm)	7/8 - 1 1/4" (22-33mm)	413422	411803
•	2 1/4" (57mm)	1 1/4 - 1 3/4" (33-45mm)	413423	411805
•	2 3/4" (70mm)	1 1/4 - 1 3/4" by 1 3/4 - 2 1/8" (33x45mm by 45x56mm) Oval shape	413424	411807

NATURA™+ POUCHES

Compatible with SUR-FIT® Natura[™] and Natura[™] skin barriers





NATURA™+ POUCH

Drainable | InvisiClose™ Tail Closure

- Pouch Length = 12" (*pouch length = 14")
- Quantity = 10

HCPCS A5063 | No Filter

	Flange Size	Transparent	Opaque
•	1 3/4" (45mm)	416418	421738
•	2 1/4" (57mm)	416421	421739
•	2 3/4" (70mm)	416424	421740
•	4" (100mm)	416472*	N/A

HCPCS A4425 | Filter (Filter covers included)

	Flange Size	Transparent	Opaque w/ Easy-View Window	Opaque
•	1 1/2" (38mm)	416414	421746	416415
•	1 3/4" (45mm)	416416	421747	416417
•	2 1/4" (57mm)	416419	421748	416420
•	2 3/4" (70mm)	416422	421749	416423

SUR-FIT® NATURA™ POUCH

Drainable | InvisiClose™ Tail Closure

- Pouch Length = 12"
- Quantity = See table below

HCPCS A5063 | No Filter

	Flange Size	Transparent Qty = 10	Opaque Qty = 20
•	1 3/4" (45mm)	413312	413334
•	2 1/4" (57mm)	413313	413335
•	2 3/4" (70mm)	413314	413336

HCPCS A4425 | Filter (Filter covers included)

	Flange Size	Transparent Qty = 10	Opaque Qty = 20
•	1 1/2" (38mm)	N/A	411309
•	1 3/4" (45mm)	411360	411310
•	2 1/4" (57mm)	411361	411311
•	2 3/4" (70mm)	411362	411312

SUR-FIT® NATURA™ POUCH

Drainable | Clipped Closure

- Pouch Length = 10"
- 2 Tail Closures per box
- Quantity = See table below

HCPCS A5063 | No Filter

	Flange Size	Opaque Qty = 10	Opaque Qty = 20
•	1 1/4" (32mm)	401505	404030
•	1 1/2" (38mm)	401506	404031
•	1 3/4" (45mm)	401507	404032
•	2 1/4" (57mm)	401508	404033
•	2 3/4" (70mm)	401509	404034



Drainable | Clipped Closure

- Pouch Length = 12" (*pouch length = 14"; **pouch length = 10")
- 2 Tail Closures per box
- Quantity = See table below

HCPCS A5063 | No Filter

	Flange Size	Transparent Qty = 10	Opaque Qty = 10	Opaque Qty = 20	Opaque w/ 1 comfort panel Qty = 10	Opaque w/ 1 comfort panel Qty = 20
•	1 1/4" (32mm)	401510	401500	404025	401932	404014
•	1 1/2" (38mm)	401511	401501	404026	401933	404015
•	1 3/4" (45mm)	401512	401502	404027	401934	404016
•	2 1/4" (57mm)	401513	401503	404028	401935	404017
•	2 3/4" (70mm)	401514	401504	404029	401936	404018
•	4" (100mm)	401515*	N/A	N/A	N/A	N/A

HCPCS A4425 | Filter (Filter covers included)

	Flange Size	Transparent Qty = 10	Opaque Qty = 20	Opaque Qty = 20**
•	1 1/4" (32mm)	N/A	411489	411288**
•	1 1/2" (38mm)	411264	411490	411289**
•	1 3/4" (45mm)	411265	411491	411290**
•	2 1/4" (57mm)	411266	411492	411291**
•	2 3/4" (70mm)	411267	411493	411292**







Recommended Product

SUR-FIT® NATURA™ POUCH

Drainable | High-Output

- Anti-Reflux Valve
- Outlet with spout and cap for high output drainage
- Pouch Length = 14"
- 1 Tail Closure, 1 Night Drainage Adapter, and 5 Replacement Filters per box
- Quantity = 5

HCPCS A4413 | Filter

	Flange Size	Transparent	Opaque
•	1 3/4" (45mm)	420695	401557
•	2 1/4" (57mm)	420696	401558
•	2 3/4" (70mm)	420697	401559



NATURA™+ POUCH

Closed-End

- Pouch Length = 8"
- Quantity = 30

HCPCS A5054 | No Filter

	Flange Size	Opaque w/ Easy-View Window	Opaque	
•	1 1/4" (32mm)	421677	416401	
•	1 1/2" (38mm)	421679	416404	
•	1 3/4" (45mm)	421681	416407	
•	2 1/4" (57mm)	421683	416410	
•	2 3/4" (70mm)	421893	416413	

HCPCS A4419 | Filter (Filter covers included)

	Flange Size	Opaque w/ Easy-View Window	Opaque
•	1 1/4" (32mm)	421676	416400
•	1 1/2" (38mm)	421678	416403
•	1 3/4" (45mm)	421680	416406
•	2 1/4" (57mm)	421682	416409
•	2 3/4" (70mm)	421798	416412

NATURA™+ POUCH

Closed-End

- Pouch Length = 6"
- Quantity = 30

HCPCS A5054 | No Filter

	Flange Size	Opaque
•	1 1/4" (32mm)	416402
•	1 1/2" (38mm)	416405
•	1 3/4" (45mm)	416408
•	2 1/4" (57mm)	416411

HCPCS A4419 | Filter (Filter covers included)

	Flange Size	Opaque w/ Easy-View Window
•	1 3/4" (45mm)	421796
•	2 1/4" (57mm)	421797

SUR-FIT® NATURA™ POUCH

Closed-End

- Pouch Length = 8"
- Quantity = See table below

HCPCS A5054 | No Filter

	Flange Size	Opaque Qty = 30	Opaque Qty = 60
•	1 1/2" (38mm)	401521	413169
•	1 3/4" (45mm)	401522	413170
•	2 1/4" (57mm)	401523	413171
•	2 3/4" (70mm)	401524	413172

HCPCS A4419 | Filter (Filter covers included)

	Flange Size	Opaque Qty = 30	Opaque Qty = 60
•	1 1/2" (38mm)	401525	413173
•	1 3/4" (45mm)	401526	413174
•	2 1/4" (57mm)	401527	413175
•	2 3/4" (70mm)	401528	413176





SUR-FIT® NATURA™+ POUCH

Closed-End | Mini

- Pouch Length = 5"
- Quantity = 20

HCPCS A5054 | No Filter

	Flange Size	Opaque w/ 1 comfort panel
•	1 1/4" (32mm)	401529
•	1 1/2" (38mm)	401530
•	1 3/4" (45mm)	401531
•	2 1/4" (57mm)	401532

SUR-FIT® NATURA™+ POUCH

Closed-End | Flange Cap

• Quantity = 25

HCPCS A5055 | Filter (Filter covers included)

	Flange Size	Opaque w/ 1 comfort panel	
•	1 3/4" (45mm)	401909	



SUR-FIT® NATURA™+ POUCH

Urostomy | Accuseal™ Tap with Valve

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10

HCPCS A4432

	Flange Size	Transparent	Opaque	Opaque*
•	1 1/4" (32mm)	401542	401551	401547*
•	1 1/2" (38mm)	401543	401552	401548*
•	1 3/4" (45mm)	401544	401553	401549*
•	2 1/4" (57mm)	401545	401554	401550*
•	2 3/4" (70mm)	401546	401555	N/A

NATURA™ POUCH

Urostomy | Soft Tap

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Soft Tap Drainage Adapters per box
- Quantity = 10

HCPCS A5073

	Flange Size	Transparent	Transparent*
•	1 1/4" (32mm)	413435	421074*
•	1 1/2" (38mm)	413436	421075*
•	1 3/4" (45mm)	413437	421076*
•	2 1/4" (57mm)	413438	N/A
•	2 3/4" (70mm)	413439	N/A

SUR-FIT® NATURA™ POUCH

Urostomy | Fold Over Tab

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Universal Night Drainage Adapters per box
- Quantity = 10

HCPCS A5073

	Flange Size	Transparent	Transparent*
•	1 1/4" (32mm)	401533	401538*
•	1 1/2" (38mm)	401534	401539*
•	1 3/4" (45mm)	401535	401540*
•	2 1/4" (57mm)	401536	401541*
•	2 3/4" (70mm)	401537	N/A





ESTEEM SYNERGY® SKIN BARRIERS

Compatible with Esteem synergy®+ and Esteem synergy® Pouches

ESTEEM SYNERGY® SKIN BARRIER

0

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

	Coupling Size	Cutting Range	Acrylic Tape Collar WHITE
•	Small Stoma Opening	up to 1 3/8" (35mm)	405456
•	Medium Stoma Opening	up to 1 7/8" (48mm)	405457
•	Large Stoma Opening	up to 2 3/8" (61mm)	405458

HCPCS A4415 (*HCPCS is A4414)

	Coupling Size	Cutting Range	Without Tape Collar
•	Small Stoma Opening	up to 1 3/8" (35mm)	405467*
•	Medium Stoma Opening	up to 1 7/8" (48mm)	405468
•	Large Stoma Opening	up to 2 3/8" (61mm)	403948
•	Extra-Large Stoma Opening	up to 3 1/2" (89mm)	405483

ESTEEM SYNERGY® SKIN BARRIER



Moldable

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

	Coupling Size	Stoma Size	Hydrocolloid Tape Collar WHITE
•	Small Stoma Opening	1/2 - 7/8" (13-22mm)	411657
•	Medium Stoma Opening	7/8 - 1 1/4" (22-33mm)	411658
•	Large Stoma Opening	1 1/4 - 1 3/4" (33-45mm)	411659

ESTEEM SYNERGY® SKIN BARRIER



Moldable | Convex

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4407

	Coupling Size	Stoma Size	Acrylic Tape Collar WHITE
•	Small Stoma Opening	1/2 - 7/8" (13-22mm)	409268
•	Medium Stoma Opening	7/8 - 1 1/4" (22-33mm)	409269
•	Large Stoma Opening	1 1/4 - 1 3/4" (33-45mm)	409270

ESTEEM SYNERGY® SKIN BARRIER



Pre-Cut

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

	Coupling Size	Stoma Size	Acrylic Tape Collar WHITE
•	Small Stoma Opening	5/8" (16 mm)	405473
•	Small Stoma Opening	3/4" (19 mm)	405474
•	Small Stoma Opening	7/8" (22 mm)	405475
•	Small Stoma Opening	1" (25 mm)	405476
•	Small Stoma Opening	1 1/8" (28 mm)	405477
•	Small Stoma Opening	1 1/4" (32 mm)	405478
•	Small Stoma Opening	1 3/8" (35 mm)	405479
•	Medium Stoma Opening	1 1/2" (38 mm)	405480
•	Medium Stoma Opening	1 5/8" (41 mm)	405481
•	Medium Stoma Opening	1 3/4" (45 mm)	405482

ESTEEM SYNERGY® POUCHES

Compatible with Esteem Synergy® Skin Barriers

ESTEEM SYNERGY® + POUCH

Drainable | InvisiClose™ Tail Closure

- Pouch Length = 12"
- Quantity = 10

HCPCS A5063 | No Filter

	Coupling Size	Transparent
•	Small Stoma Opening	416797
•	Medium Stoma Opening	416900
•	Large Stoma Opening	416903

HCPCS A4425 | Filter (Filter covers included)

	Coupling Size	Transparent	Opaque	
•	Small Stoma Opening	416795	416796	
•	Medium Stoma Opening	416798	416799	
•	Large Stoma Opening	416901	416902	

ESTEEM SYNERGY® + POUCH

Closed-End

- Pouch Length = 8"
- Quantity = 30

HCPCS A4419 | Filter (Filter covers included)

	Coupling Size	Transparent	Opaque
•	Small Stoma Opening	416774	416775
•	Medium Stoma Opening	416778	416779
•	Large Stoma Opening	416782	416783

ESTEEM SYNERGY® + POUCH

Closed-End

- Pouch Length = 6"
- Quantity = 30

HCPCS A4419 | Filter (Filter covers included)

	Coupling Size	Transparent	Opaque
•	Small Stoma Opening	416776	416777
•	Medium Stoma Opening	416780	416781

ESTEEM SYNERGY® POUCH

Drainable | Clipped Closure

- Pouch Length = 14"
- 2 Tail Closures per box
- Quantity = 10

HCPCS A5063 | No Filter

	Coupling Size	Transparent
•	Extra-Large Stoma Opening	405406

ESTEEM SYNERGY® POUCH

Urostomy | Accuseal™ Tap with Valve

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10

HCPCS A4432

	Coupling Size	Transparent	Transparent*	Opaque	
•	Small Stoma Opening	405451	405446*	405448	
•	Medium Stoma Opening	405452	405447*	405449	
•	Large Stoma Opening	405453	N/A	405450	







ONE-PIECE Systems

TWO BRANDS make up ConvaTec's ONE-Piece Ostomy Systems



ESTEEM[™]

The Esteem™ + brand is an all-in-one solution combining the skin barrier and pouch in a single unit. Designed to be easy and convenient, it's one simple step to attach or remove.

Available in multiple adhesive options, growing convexity options, and with ConvaTec's latest pouch technologies.



ActiveLife

Active Life® is a legacy one-piece system brand that also combines the skin barrier and pouch in a single unit.

ESTEEM™ + POUCHES

Esteem $^{\text{TM}}$ + Pouches come without a tape collar unless otherwise noted.





SESTEEM™+ POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Durahesive® Plus (Extended Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- Quantity = 10

HCPCS A5056 | Filter (Filter covers included)

Cutting Range	Transparent	Opaque
3/4 - 2 1/2" (19-64mm)	416976	416975



SESTEEM™+ POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Modified Stomahesive® (Standard Wear)
- Pouch Length = 12"
- Quantity = 10

HCPCS A5061 | No Filter

Cutting Range	Transparent	Opaque w/ Easy-View Window
13/16 - 2 3/4" (20-70mm)	416721	421853

HCPCS A4424 | Filter (Filter covers included)

Cutting Range	Transparent	Opaque w/ Easy-View Window	Opaque
13/16 - 2 3/4" (20-70mm)	416718	421858	416719





SESTEEM™+ POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)*
- Pouch Length = 14"
- Quantity = 10

HCPCS A5061 | No Filter

Cutting Range	Transparent
Up to 4" (100mm)	416908
*Double-layer adhesive (Stomahesive body-side, E the skin friendly aspects of Stomahesive with the	



SESTEEM™+ POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Modified Stomahesive® (Standard Wear)
- Pouch Length = 12"
- 2 Tail Closures per box
- Quantity = 10

HCPCS A4424 | Filter (Filter covers included)*

Cutting Range	Transparent
13/16 - 2 3/4" (20-70mm)	421848



SESTEEM™+ POUCH

Pre-Cut | Drainable | InvisiClose™ Tail Closure

- Modified Stomahesive® (Standard Wear)
- Pouch Length = 12"
- Quantity = 10

HCPCS A5061 | No Filter

Stoma Size	Opaque w/ Easy-View Window
1" (25mm)	421854
1 3/16" (30mm)	421855
1 3/8" (35mm)	421856
1 9/16" (40mm)	421857

HCPCS A4424 | Filter (Filter covers included)

Stoma Size	Transparent	Opaque w/ Easy-View Window	Opaque
1" (25mm)	416724	421859	416725
1 3/16" (30mm)	416728	421860	416729
1 3/8" (35mm)	416732	421861	416733
1 9/16" (40mm)	416736	421862	416737



SESTEEM™+ POUCH

Moldable | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Pouch Length = 12"
- Quantity = 10

HCPCS A5061 | No Filter

	Stoma Size	Transparent	
13/16	- 1 3/16" (20-30mm)	413518	_
1 3/16	- 1 9/16" (30-40mm)	413519	
1 9/16	- 1 15/16" (40-50mm)	413520	

HCPCS A4424 | Filter (Filter covers included)

Stoma Size	Transparent	Opaque w/ Easy-View Window
13/16 - 1 3/16" (20-30mm)	413515	413521
1 3/16 - 1 9/16" (30-40mm)	413516	413522
1 9/16 - 1 15/16" (40-50mm)	413517	413523





ESTEEM™+ POUCH

Pre-Cut | Drainable | Convex | InvisiClose™ Tail Closure

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- Quantity = 10

HCPCS A5057 | Filter (Filter covers included)

Stoma Size	Transparent	Opaque
3/4" (19mm)	416738	416739
7/8" (22mm)	416740	416741
1" (25mm)	416742	416743
1 1/8" (28mm)	416744	416745
1 1/4" (32mm)	416746	416747
1 3/8" (35mm)	416748	416749
1 1/2" (38mm)	416750	416751



STEEM™+ FLEX CONVEX POUCH

Cut-to-Fit | Drainable | Convex

- Esteem[™] + *Flex* Convex Adhesive (*Standard Wear*)
- Convex wafter depth = 7mm
- Safe Seal Clipless Closure (1 Tail Closure per box)
- Pouch Length = 8 2/3"
- Quantity = 10

HCPCS A4389 | Filter (Filter covers included)


Cutting Range	Plateau Size	Opaque w/ Easy-View Window
3/4 - 1" (20-25mm)	V3	421622
3/4 - 1 3/8" (20-35mm)	V2	421619
3/4 - 1 2/3" (20-43mm)	V1	421615
3/4 - 2 6/25" (20-57mm)	VO	422159

Pre-Cut | Drainable | Convex

HCPCS A4389 | Filter (Filter covers included)

Stoma Size	Plateau Size	Opaque w/ Easy-View Window
1" (25mm)	V2	421620
1 1/6" (30mm)	V2	421621
1" (25mm)	V1	421616
1 1/6" (30mm)	V1	421617
1 3/8" (35mm)	V1	421618

Plateau Width Guidance for Esteem™+ Flex and Soft Convex





SESTEEM™+ SOFT CONVEX POUCH

Cut-to-Fit | Drainable | Convex

- Esteem[™] + *Soft* Convex Adhesive (*Standard Wear*)
- Convex wafter depth = 4mm
- Safe Seal Clipless Closure (1 Tail Closure per box)
- Pouch Length = 8 2/3"
- Quantity = 10

HCPCS A4389 | Filter (Filter covers included)

Cutting Range	Plateau Size	Opaque w/ Easy-View Window
3/8 - 1 1/8" (10-28mm)	V3	422365
3/5 - 1 3/5" (15-40mm)	V2	422364
3/4 - 1 7/8" (20-47mm)	V1	422363

Pre-Cut | Drainable | Convex

HCPCS A4389 | Filter (Filter covers included)

Stoma Size	Plateau Size	Opaque w/ Easy-View Window
1" (25mm)	V2	422545
1 1/6" (30mm)	V2	422546
1" (25mm)	V1	422539
1 1/6" (30mm)	V1	422540
1 3/8" (35mm)	V1	422541

Recommended Product



SESTEEM™+ POUCH

Cut-to-Fit | Closed-End

- Modified Stomahesive® (Standard Wear)
- Pouch Length = 8"
- Quantity = 30

HCPCS A4416 | Filter (Filter covers included)

Cutting Range	Transparent	Opaque w/ Easy-View Window	Opaque
13/16 - 2 3/4" (20-70mm)	416700	421821	416701

ESTEEM™+ POUCH

Cut-to-Fit | Closed-End

- Modified Stomahesive® (Standard Wear)
- Pouch Length = 6"
- Quantity = 30

HCPCS A4416 | Filter (Filter covers included)

Cutting Range	Opaque w/ Easy-View Window	Opaque
13/16 - 2 3/4" (20-70mm)	421816	416702

ESTEEM™+ POUCH

Pre-Cut | Closed-End

- Modified Stomahesive® (Standard Wear)
- Pouch Length = 8"
- Quantity = 30

HCPCS A4416 | Filter (Filter covers included)

Stoma Size	Transparent	Opaque w/ Easy-View Window	Opaque
1" (25mm)	416703	421687	416704
1 3/16" (30mm)	416709	421688	416707
1 3/8" (35mm)	416706	421689	416710
1 9/16" (40mm)	416712	421822	416713
2" (50mm)	N/A	421823	416715

ESTEEM™+ POUCH

Moldable | Closed-End

- Stomahesive® (Standard Wear)
- Pouch Length = 8"
- Quantity = 30

HCPCS A4416 | Filter (Filter covers included)

Stoma Size	Opaque w/ Easy-View Window
13/16 - 1 3/16" (20-30mm)	413509
1 3/16 - 1 9/16" (30-40mm)	413510
1 9/16 - 1 15/16" (40-50mm)	413511

ESTEEM™+ FLEX CONVEX POUCH

Cut-to-Fit | Closed-End | Convex

- Esteem[™] + *Flex* Convex Adhesive (*Standard Wear*)
- Convex wafter depth = 7mm
- Pouch Length = 8"
- Quantity = 30

HCPCS A4417 | Filter (Filter covers included)

Cutting Range	Plateau Size	Opaque w/ Easy-View Window
3/4 - 1" (20-25mm)	V3	421614
3/4 - 1 3/8" (20-35mm)	V2	421613
3/4 - 1 2/3" (20-43mm	V1	421612











SESTEEM™+ POUCH

Urostomy | Cut-to-Fit | Accuseal™ Tap with Valve

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Anti-Reflux Valve
- Pouch Length = 10"
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10

HCPCS A4428

Cutting Range	Transparent
1/2 - 1 3/4" (13-45mm)	421631



SESTEEM™+ FLEX CONVEX POUCH

Urostomy | Cut-to-Fit | Convex

- Esteem[™] + *Flex* Convex Adhesive (*Standard Wear*)
- Soft Tap Closure
- Convex wafter depth = 7mm
- Pouch Length = 7 2/3"
- 2 Esteem[™] + Flex Convex Night Drainage Adapters per box
- Quantity = 10

HCPCS A4392

	Cutting Range	Plateau Size	Opaque w/ Easy-View Window
Ī	3/8 - 1" (10-25mm)	V3	421628
	3/8 - 1 3/8" (10-35mm)	V2	421626
	3/8 - 1 2/3" (10-43mm)	V1	421623

Urostomy | Pre-Cut | Convex

HCPCS A4392

Stoma Size	Plateau Size	Opaque w/ Easy-View Window
9/16" (15mm)	V3	421629
9/16" (15mm)	V2	421627
9/16" (15mm)	V1	421624
4/5" (20mm)	V1	421625



SESTEEM™+ SOFT CONVEX POUCH

Urostomy | Cut-to-Fit | Convex

- Esteem[™] + *Soft* Convex Adhesive (*Standard Wear*)
- Soft Tap Closure
- Convex wafter depth = 4mm
- Pouch Length = 7 2/3"
- 2 Esteem[™] + Soft Convex Night Drainage Adapters per box
- Quantity = 10

HCPCS A4392

	Cutting Range	Plateau Size	Opaque w/ Easy-View Window
Ī	3/8 - 1 1/8" (10-28mm)	V3	422562
	3/5 - 1 3/5" (15-40mm)	V2	422557
	3/4 - 1 7/8" (20-47mm)	V1	422551

Urostomy | Pre-Cut | Convex

HCPCS A4392

Stoma Size	Plateau Size	Opaque w/ Easy-View Window
4/5" (20mm)	V3	422564
1" (25mm)	V3	422565
4/5" (20mm)	V2	422558
1" (25mm)	V2	422559
1 1/6" (30mm)	V2	422560
1" (25mm)	V1	422552
1 1/6" (30mm)	V1	422553

ACTIVELIFE® POUCHES

ACTIVELIFE® POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- 1 Tail Closure per box
- Quantity = See table below

HCPCS A5061 | No Filter

Cutting Range	Transparent	Transparent	Transparent
	Qty = 1	Qty = 10	Qty = 20
3/4 - 2 1/2" (19-64mm)	025433	022771	125330

ACTIVELIFE® POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Durahesive® Plus (Extended Wear)
- · Acrylic Tape Collar
- Pouch Length = 12"
- 1 Tail Closure per box
- Quantity = 10

HCPCS A4388 | No Filter

Cutting Range	Transparent	Opaque
3/4 - 2 1/2" (19-64mm)	400598	400599



Pre-Cut | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar
- Pouch Length = 12" (*pouch length = 10")
- 1 Tail Closure per box
- Quantity = See table below

HCPCS A5061 | No Filter

Stoma Size	Transparent Qty = 10	Opaque Qty = 10	Opaque Qty = 20	Opaque* Qty = 10	Opaque* Qty = 20
3/4" (19mm)	022764	022757	125338	022750*	125331*
1" (25mm)	022765	022758	125339	022751*	N/A
1 1/4" (32mm)	022766	022759	125340	022752*	125333*
1 1/2" (38mm)	022767	022760	125341	022753*	125334*
1 3/4" (45mm)	022768	022761	125342	022754*	125335*
2" (50mm)	022769	022762	125343	022755*	125336*
2 1/2" (64mm)	022770	022763	125344	022756*	125337*

ACTIVELIFE® POUCH

Pre-Cut | Drainable | Convex | **Clipped Closure**

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- 1 Tail Closure per box
- Quantity = See table below

HCPCS A4390 | No Filter

Stoma Size	Transparent Qty = 5	Transparent Qty = 10
3/4" (19mm)	175777	125352
7/8" (22mm)	175778	125353
1" (25mm)	175779	125354
1 1/8" (28mm)	175780	125355
1 1/4" (32mm)	175781	125356
1 3/8" (35mm)	175782	125357
1 1/2" (38mm)	175783	125358
1 3/4" (45mm)	175784	125359
2" (50mm)	175785	N/A





ACTIVELIFE® POUCH

Cut-to-Fit | Closed-End

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar
- Pouch Length = 8"
- Quantity = 30

HCPCS A4416 | Filter (Filter covers included)

Cutting Range	Transparent	Opaque
3/4 - 2" (19 - 50mm)	650422	650408

ACTIVELIFE® POUCH

Cut-to-Fit | Closed-End

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 8"
- Quantity = See table below

HCPCS A4416 | Filter (Filter covers included)

Cutting Range	Transparent	Transparent	Opaque	Opaque
	Qty = 15	Qty = 60	Qty = 15	Qty = 60
3/4 - 2 1/2" (19 - 64mm)	175761	413147	175768	413141

ACTIVELIFE® POUCH

Pre-Cut | Closed-End

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 8"
- Quantity = See table below

HCPCS A4416 | Filter (Filter covers included)

Stoma Size	Transparent Qty = 15	Transparent Qty = 60	Opaque Qty = 15	Opaque Qty = 60
1" (25mm)	175762	413148	175769	413142
1 1/4" (32mm)	175763	413149	175770	413143
1 1/2" (38mm)	175764	413150	175771	413144
1 3/4" (45mm)	175765	413151	175772	413145
2" (50mm)	175766	413152	175773	413146

Cut-to-Fit | Stoma Cap

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Quantity = 30

inci es Asoss i inter (i inter covers ii	iciaaca)
Cutting Range	Opaque
3/4 - 2" (19 - 50mm)	175611

ACTIVELIFE® POUCH

Urostomy | Pre-Cut | **Accuseal**[™] **Tap with Valve**

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Anti-Reflux Valve
- Pouch Length = 8"
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10

HCPCS A4428

Stoma Size Transparent
3/4" (19mm) 650828
7/8" (22mm) 650829
1" (25mm) 650830
1 1/8" (28mm) 650831
1 1/4" (32mm) 650832
1 1/2" (38mm) 650833









HCPCS A5055 | Filter (Filter covers included)



ACTIVELIFE® POUCH

Urostomy | Cut-to-Fit | Fold Over Tap

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 9"
- 1 Universal Night Drainage Adapter per box
- Quantity = 10

HCPCS A5071

Cutting Range	Transparent
3/4 - 1 3/4" (19-45mm)	064927

ACTIVELIFE® POUCH

Urostomy | Pre-Cut | Convex | Fold Over Tap

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Pouch Length = 8"
- 1 Universal Night Drainage Adapter per box
- Quantity = See table below

HCPCS A4393

Stoma Size	Qty = 5	Qty = 10
1/2" (13mm)	175790	125361
5/8" (16mm)	175791	125362
3/4" (19mm)	175792	125363
7/8" (22mm)	175793	125364
1" (25mm)	175794	125365
1 1/8" (28mm)	175795	125366
1 1/4" (32mm)	175796	125367
1 3/8" (35mm)	175797	125368
1 1/2" (38mm)	175798	125369



ConvaTec's range of Acute Kits includes Post-Operative/Surgical kits in Sterile and Non-Sterile combinations, and Unit Dose kits.

Information on the components included in the Acute Kits can be found below the product information tables.

Color Key & Sizing Guide

COLOR MATCH			•	•		•
FLANGE	1 1/4"	1 1/2"	1 3/4"	2 1/4"	2 1/4"	4"
SIZE	(32mm)	(38mm)	(45mm)	(57mm)	(57mm)	(100mm)

POST-OPERATIVE/ SURGICAL KIT TWO-PIECE SYSTEM

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Cut-to-Fit | Drainable | InvisiClose[™] Tail Closure

- Double Layer Skin Barrier[†] (Standard Wear)
- Acrylic Tape Collar White (*without tape collar)
- Pouch Color: Transparent
- Pouch Length = 12" (*pouch length is 14")
- 1 Low-Pressure Adaptor
- Quantity = 5 Kits

HCPCS N/A

	Flange Size	Non-Sterile Kit	Sterile Kit	Ostomy Loop Rod
•	1 3/4" (45mm)	416919	416920	None
•	2 1/4" (57mm)	416923	416924	None
•	2 3/4" (70mm)	416927	416928	65mm
•	4" (100mm)	416929*	416930*	90mm

†Double-layer adhesive (Stomahesive body-side, Durahesive top-side) which combines the skin friendly aspects of Stomahesive with the turtlenecking protection of Durahesive.

Above kits include the following components:

Kit	Comparable to Skin Barrier in Kit* Stomahesive®	Comparable to Skin Barrier in Kit* Durahesive ®	Pouch In Kit	Low Pressure Adaptor	Ostomy Loop Rod
416919/ 416920	401575	413155	416418	401993	None
416923/ 416924	401576	413156	416421	401994	None
416927/ 416928	401577	413157	416424	401995	022355
416929/ 416930	401906	401905	416472	401996	022356

^{*}The double-layer skin barrier in the kit is not available in the retail market. Comparable skin barriers are in the chart above.

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Moldable | Drainable | InvisiClose™ Tail Closure

- Durahesive® (Extended Wear)
- Acrylic Tape Collar (White)
- Pouch Color: Transparent
- Pouch Length = 12"
- 1 Low-Pressure Adaptor
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	1 3/4" (45mm)	416917
•	2 1/4" (57mm)	416921
•	2 3/4" (70mm)	416925

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit	Low Pressure Adaptor
416917	411802	416418	401993
416921	411804	416421	401994
416925	411806	416424	401995

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Moldable | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar (White)
- Pouch Color: Transparent
- Pouch Length = 12"
- 1 Low-Pressure Adaptor
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	1 3/4" (45mm)	416918
•	2 1/4" (57mm)	416922
•	2 3/4" (70mm)	416926

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit	Low Pressure Adaptor
416918	411803	416418	401993
416922	411805	416421	401994
416926	411807	416424	401995

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Moldable | Accordion Flange | Drainable | InvisiClose™ Tail Closure

- Durahesive® (Extended Wear)
- Acrylic Tape Collar (Tan)
- Pouch Color: Transparent
- Pouch Length = 12"
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	2 1/4" (57mm)	421045
•	2 1/4" (57mm)	421046
•	2 3/4" (70mm)	421047

Continued

Natura[™] Skin Barrier & Natura[™] + Pouch Kit Moldable | Accordion Flange | Drainable | InvisiClose[™] Tail Closure

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit
421045	421039	416421
421046	421040	416421
421047	421041	416424

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Moldable | Accordion Flange | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar (Tan)
- Pouch Color: Transparent
- Pouch Length = 12"
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

Elango Cizo

	Flange Size	Kit Product Code
•	2 1/4" (57mm)	421042
•	2 1/4" (57mm)	421043
•	2 3/4" (70mm)	421044

Vit Droduct Code

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit
421042	421033	416421
421043	421034	416421
421044	421035	416424

NATURA™ SKIN BARRIER & SUR-FIT® NATURA™ + POUCH KIT

Moldable | Urostomy | Accuseal™ Tap with Valve

- Durahesive® (Extended Wear)
- Acrylic Tape Collar (White)
- Anti-Reflux Valve
- Pouch Color: Transparent
- Pouch Length = 10"
- 1 Low Pressure Adaptor
- 1 Accuseal Night Drainage Adapter per box
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	1 3/4" (45mm)	416932
•	2 1/4" (57mm)	416940
•	2 3/4" (70mm)	416948

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit	Low Pressure Adaptor
416932	411802	401544	401993
416940	411804	401545	401994
416948	411806	401546	401995

NATURA™ SKIN BARRIER & SUR-FIT® NATURA™ + POUCH KIT

Moldable | Urostomy | Accuseal™ Tap with Valve

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar (White)
- Anti-Reflux Valve
- Pouch Color: Transparent
- Pouch Length = 10"
- 1 Low Pressure Adaptor
- 1 Accuseal Night Drainage Adapter per box
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	1 3/4" (45mm)	416934
•	2 1/4" (57mm)	416942
•	2 3/4" (70mm)	416950

Natura[™] Skin Barrier & SUR-FIT® Natura[™]+ Pouch Kit Moldable | Urostomy | Accuseal[™] Tap with Valve

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit	Low Pressure Adaptor
416934	411803	401544	401993
416942	411805	401545	401994
416950	411807	401546	401995

NATURA™ SKIN BARRIER & SUR-FIT® NATURA™ + POUCH KIT

Cut-to-Fit | Urostomy | Accuseal™ Tap with Valve

- Double Layer Skin Barrier[†] (Standard Wear)
- Acrylic Tape Collar (White)
- Anti-Reflux Valve
- Pouch Color: Transparent
- Pouch Length = 10"
- 1 Low Pressure Adaptor
- 1 Accuseal Night Drainage Adapter per box
- Quantity = 5 Kits

HCPCS N/A | Sterile

	Flange Size	Kit Product Code		
•	1 3/4" (45mm)	416938		
•	2 1/4" (57mm)	416946		
•	2 3/4" (70mm)	416954		

[†]Double-layer adhesive (Stomahesive body-side, Durahesive top-side) which combines the skin friendly aspects of Stomahesive with the turtlenecking protection of Durahesive.

Above kits include the following components:

Kit	Comparable to Skin Barrier in Kit* Stomahesive®	Comparable to Skin Barrier in Kit* Durahesive ®	Pouch In Kit	Low Pressure Adaptor
416938	401575	413155	401544	401993
416946	401576	413156	401545	401994
416954	401577	413157	401546	401995

^{*}The double-layer skin barrier in the kit is not available in the retail market. Comparable skin barriers are in the chart above.

ESTEEM SYNERGY® SKIN BARRIER & ESTEEM SYNERGY® POUCH KIT

Cut-to-Fit | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- No Tape Collar
- Pouch Color: Transparent
- Pouch Length = 5 14"
- 1 Tail Closure per box
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Coupling Size	Kit Product Code
•	6 x 6"	409275

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit	1 Tail Closure
409275	405483	405406	175635

CONVATEC LOOP OSTOMY ROD

Sterile

• Quantity = 10

HCPCS A4421

	Flange Size	Size	Product Code
•	2 3/4" (70mm)	65mm	022355
•	4" (100mm)	90mm	022356

POST-OPERATIVE/ SURGICAL KIT ONE-PIECE SYSTEM

ESTEEM™ + POUCH KIT

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Double Layer Skin Barrier† (Standard Wear)
- Acrylic Tape Collar
- · Color: Transparent
- Pouch Length = 14"
- Quantity = 5 Kits

HCPCS N/A

Cutting Range	Non-Sterile Kit	Sterile Kit
Can be enlarged to accommodate stoma sizes up to 4" (100mm)	416915	416916

†Double-layer adhesive (Stomahesive body-side, Durahesive top-side) which combines the skin friendly aspects of Stomahesive with the turtlenecking protection of Durahesive.

UNIT-DOSE KITS

SUR-FIT® NATURA™ SKIN BARRIER & SUR-FIT® NATURA™ + POUCH KIT

Cut-to-Fit | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar (White)
- Pouch Color: Transparent
- Pouch Length = 12"
- 1 Tail Closure per box
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	1 3/4" (45mm)	401920
•	2 1/4" (57mm)	402239
•	2 3/4" (70mm)	401921

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit	1 Tail Closure
401920	125259	401512	175635
402239	125260	401513	175635
401921	125261	401514	175635

SUR-FIT® NATURA™ SKIN BARRIER & SUR-FIT® NATURA™ + POUCH KIT

Cut-to-Fit | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- No Tape Collar
- Pouch Color: Transparent
- Pouch Length = 14"
- 1 Tail Closure per box
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	4" (100mm)	401922

Continued

SUR-FIT® Natura™ Skin Barrier & SUR-FIT® Natura™+ Pouch Kit Cut-to-Fit | Drainable | Clipped Closure

Above kits include the following components:

	Kit	Skin Barrier in Kit	Pouch In Kit	1 Tail Closure	
Ī	401922	401906	401515	175635	

SUR-FIT® NATURA™ SKIN BARRIER & SUR-FIT® NATURA™ + POUCH KIT

Cut-to-Fit | Urostomy | Accuseal™ Tap with Valve

- Durahesive® (Extended Wear)
- Acrylic Tape Collar (White)
- Pouch Color: Transparent
- Pouch Length = 10"
- 1 Accuseal Night Drainage Adapter per box
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

	Flange Size	Kit Product Code
•	1 3/4" (45mm)	401923
•	2 1/4" (57mm)	401924

Above kits include the following components:

Kit	Skin Barrier in Kit	Pouch In Kit
401923	413161	401544
401924	413162	401545

PEDIATRIC Products

LITTLE ONES is a pediatric Two-piece and One-piece system designed for babies smaller than 10 lbs, and children over 10 lbs.

The pouches have soft-welded edges and smaller pouch sizes designed to provide extra comfort for little ones.



Color Key & Sizing Guide



71

When ordering a Pediactric Two-Piece System, make sure you select the same flange/coupling size for the skin barrier and pouch. Use the colored dots in the table to help.

COLOR MATCH	•		•	
FLANGE/ COUPLING SIZE	0 - 9/10" (0-23mm)	1/5 - 1 1/4" (4.8-31mm)	1 1/4" (32mm)*	1 3/4" (45mm)*


^{*}Indicates Flange size measurement

LITTLE ONES® EXTRA SMALL TWO-PIECE OSTOMY SYSTEM

For babies smaller than 10lbs/4.5kg

LITTLE ONES® EXTRA SMALL SKIN BARRIER

Cut-to-Fit



- Adhesive Coupling Technology[™]
- No Starter Hole
- Quantity = 5

HCPCS A4414

	Cutting Range	Without Tape Collar	
•	0 - 9/10" (0-23mm)	411632	

LITTLE ONES® EXTRA SMALL POUCH

Drainable | InvisiClose™ Tail Closure

- Adhesive Coupling Technology[™]
- Pouch Length = 6"
- Quantity = 10

HCPCS A5063 | No Filter

Stoma Size		Transparent	
•	0 - 9/10" (0-23mm)	411631	



LITTLE ONES® EXTRA SMALL ONE-PIECE OSTOMY SYSTEM

For babies smaller than 10lbs/4.5kg

LITTLE ONES® EXTRA SMALL POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Without Tape Collar
- No Starter Hole
- Pouch Length = 6"
- Quantity = 10

HCPCS A5061 | No Filter

Stoma Size	Transparent
0 - 9/10" (0-23mm)	411633

LITTLE ONES® TWO-PIECE OSTOMY SYSTEM

Suggested for children over 10lbs/4.5kg

LITTLE ONES® SKIN BARRIER

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Adhesive Coupling Technology[™]
- Quantity = 5

HCPCS A4414

	Cutting Range	Hydrocolloid Tape Collar WHITE	No Tape Collar	
•	1/5 - 1 1/4" (4.8-31mm)	411642	411641	

LITTLE ONES® POUCH

Drainable | InvisiClose™ Tail Closure

- Adhesive Coupling Technology[™]
- Pouch Length = 8"
- Quantity = 10

HCPCS A4425 | Filter (Filter covers included)

	Stoma Size	Transparent	Opaque
•	1/5 - 1 1/4" (4.8-31 mm)	411637	411638

LITTLE ONES® POUCH

Closed-End

- Adhesive Coupling Technology[™]
- Pouch Length = 6"
- Quantity = 10

HCPCS A4419 | Filter (Filter covers included)

	Stoma Size	Transparent	Opaque	
•	1/5 - 1 1/4" (4.8-31 mm)	411640	411639	Ī

LITTLE ONES® SKIN BARRIER

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Flange System
- Quantity = 5

	Flange Size	Cutting Range	Acrylic Tape Collar WHITE
•	1 1/4" (32mm)	1/2 - 3/4" (13-19mm)	401925
•	1 3/4" (45mm)	1 1/8 - 1 1/4" (28-32mm)	401926





LITTLE ONES® POUCH

Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Flange System
- Pouch Length = 6"
- 10 Soft, Moldable Tail Closures per box (Designed for Pediatric Use)
- Quantity = 10

HCPCS A5063 | No Filter

	Flange Size	Transparent
•	1 1/4" (32mm)	401927
•	1 3/4" (45mm)	401928

LITTLE ONES® POUCH

Closed-End

- Flange System
- Pouch Length = 5"
- Quantity = 20

HCPCS A5054 | No Filter

	Flange Size	Opaque
•	1 1/4" (32mm)	401930
•	1 3/4" (45mm)	401931

LITTLE ONES® POUCH

Urostomy | Fold Over Tap

- Flange System
- Pouch Length = 6"
- 2 Universal Night Drainage Adapters per box
- Quantity = 10

HCPCS A5073 | No Filter

	Flange Size	Transparent
•	1 1/4" (32mm)	401929



Suggested for children over 10lbs/4.5kg

LITTLE ONES® POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 6"
- 15 Soft, Moldable Tail Closures per box (Designed for Pediatric Use)
- Quantity = 15

HCPCS A5061 I No Filter

Cutting Range	Transparent	
5/16 - 2" (8-50mm)	020922	

LITTLE ONES® POUCH

Urostomy | Cut-to-Fit | Fold Over Tap

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 5"
- 2 Universal Night Drainage Adapters per box
- Quantity = 15

Cutting Range	Transparent
5/16 - 1" (8-25mm)	020917





ACCESSORIES



Solutions that help promote healthy skin by reducing irritation and damage.

Choosing the right accessories

Whether you are looking for better skin protection or added security to your ostomy system, with the right combination of accessories, you can feel even more comfortable and at ease.



Skin **Preparation**



Added **Security**



Pouching Support



Pouch Removal

SENSI-CARE™

Sting-Free Barrier Wipe

- Designed to create a protective layer that shields skin from damage
- Helps reduce skin stripping and irritation from body waste
- May help increase the contact or adhesion to skin
- Silicone based
- Alcohol free
- CHG Compatible¹⁷

HCPCS A5120			
Liquid per Wipe	Product Code	Box of	Boxes per Case*
3 ml	413501	30	20

*Sold by the case

SENSI-CARE™

Sting-Free Barrier Spray

- Designed to create a protective layer that shields skin from damage
- Helps reduce skin stripping and irritation from body waste
- May help increase the contact or adhesion to skin
- Silicone based
- Alcohol free
- Non-Aerosol
- CHG Compatible¹⁷

LIC	DCC	A4369	١.
пС	rcs	A4305	,

Spray Bottle Size	Product Code	Box of	Boxes per Case*
28 mL	420797	1	12
50 mL	413502	1	12

*Sold by the case

SENSI-CARE™

Sting-Free Barrier Foam Applicator

- Designed to create a protective layer that shields skin from damage
- Helps reduce skin stripping and irritation from body waste
- May help increase the contact or adhesion to skin
- · Silicone based
- · Alcohol free
- CHG Compatible¹⁷

HCPCS A5120 | Sterile

Applicator Size	Product Code	Individually Packaged Applicators per Box	Boxes per Case*
1 mL	420793	5	20
1 mL	420794	25	4
3 mL	420795	5	20
3 mL	420796	25	4

^{*}Sold by the case

ALLKARE®

Protective Barrier Wipe

- Designed to create a protective layer that shields skin from damage
- Helps reduce skin stripping and irritation from body waste
- May help increase the contact or adhesion to skin
- · Contains Alcohol

HCPCS A5120

Product Code	Box of
037439	50
037444	100

STOMAHESIVE®

Protective Powder

- Helps protect irritated skin by absorbing excess moisture
- For use on weeping skin around a stoma
- Lightly dust affected area with powder, wiping off excess. May be followed up with a skin barrier wipe or spray

HCPCS A4371

Size	Product Code	Box of	
1 07	025510	1 hottle	



Products that help create the best fit around the stoma, prevent leakage and help level out the application surface area for better pouch performance.

EAKIN COHESIVE®

Seals

- Moldable barrier that helps seal the area around hard-to-fit stomas from any leakage under the skin barrier
- Acts as a filler on uneven skin surfaces to help create a flatter area for the skin barrier to adhere to
- Designed to absorb moisture and prevent it coming in contact with, and causing harm to the skin-forms a gel that protects the skin

HCPCS A4385

	Size	Product Code	Depth	Box of
Ī	Small 2"	839002	4.2 mm	20
	Slims 2"	839005	3 mm	10
	Large 4"	839001	2.1 mm	10
	StomaWrap 3 1/3"	839006	3 mm	10

STOMAHESIVE®

Paste

- Helps to protect exposed skin between the base of the stoma and the opening of the skin barrier. Can also be used to fill crevices, skin folds and scars to help create a flat pouching surface for optimal appliance attachment
- Helps to prevent leakage under the skin barrier that may result in skin breakdown
- · Pectin-based
- Contains alcohol

Size	Product Code	Box of
2 oz.	183910	1 tube





EAKIN COHESIVE®

Paste

- Use to fill crevices, skin folds and scars to help create an even pouching surface
- Leakage prevention by helping form a waterproof barrier around the stoma
- Alcohol free

HCPCS A4405

Size	Product Code	Box of
2.1 oz.	839010	1 tube

STOMAHESIVE®

Skin Barrier

- Helps fill uneven skin surfaces to create a flat pouching surface for optimal appliance attachment
- Designed to help protect the skin from contact with bodily fluids
- Large cutting surface to accommodate large and swollen stomas
- Sterile Stomahesive® skin barrier can be used post-operatively

HCPCS A4362

Size	Product Code	Box of
4 x 4"	021712	5
4 x 4 (Sterile)	021901	5

HCPCS A5122

Size	Product Code	Box of
8 x 8"	021715	3

EAKIN COHESIVE®

Skin Barrier

- Designed to help protect the skin from contact with body fluids
- Acts as a filler on uneven skin surfaces to help create a flat pouching surface for optimal appliance attachment
- · Moldable and moisture absorbing

HCPCS A4385

Size	Product Code	Box of
4 x 4"	839004	5
4 x 8"	839003	5

STOMAHESIVE®

Strips

- Helps to create a secure seal between the stoma and the appliance, to protect the skin against stoma output
- Can be cut, torn and shaped into a variety of shapes to help prevent leakage around the stoma by acting as a filler
- Designed to absorb moisture from the stoma, creating a snug fit around the stoma.

HCPCS A4362

	Size	Product Code	Box of	
Т	3 cards of 5 strips each	025542	15 strips	

SUR-FIT® NATURA™

Disposable Convex Inserts

- Convex inserts that can be added to a flat skin barrier to create a shallow convexity—insert into the inner diameter of the plastic flange on the pouch facing side of a SUR-FIT® Natura™ Durahesive® and Stomahesive® skin barrier
- Select an insert that is at least ¼" larger than the maximum diameter of the stoma
- For SUR-FIT® Natura™ skin barriers only. Not for use with ConvaTec Moldable Technology™ skin barriers, or any other product
- Quantity = 5

	For Use With Flange Size	Stoma Size no Wider Than	Convex Insert Opening Size	Product Code
•	1 1/2" (38mm)	1/2"	3/4" (19mm)	404006
•	1 1/2" (38mm)	5/8"	7/8" (22mm)	404007
•	1 1/2" (38mm)	3/4"	1" (25mm)	404008
•	1 3/4" (45mm)	7/8"	1 1/8" (29mm)	404009
•	1 3/4" (45mm)	1"	1 1/4" (32mm)	404010
•	1 3/4" (45mm)	1 1/8"	1 3/8" (35mm)	404011
•	2 1/4" (57mm)	1 1/4"	1 1/2" (38mm)	404012
•	2 1/4" (57mm)	1 3/8"	1 5/8" (41mm)	404013



SUR-FIT® NATURA™

Low-Pressure Adaptor

- For management of post-operative tenderness helps solve a temporary tenderness condition
- Optimal finger space
- Eliminate contact with abdomen
- Designed for low-pressure application when attaching pouch
- Tubular ring that is placed in between the skin barrier and pouch
- Compatible with SUR-FIT® Natura[™] skin barriers and pouches. Not for use with skin barriers with Accordion Flange
- Quantity = 10

HCPCS A4421

	Flange Size	Product Code
•	1 3/4" (45mm)	401993
•	2 1/4" (57mm)	401994
•	2 3/4" (70mm)	401995
•	4" (100mm)	401996

ADDED SECURITY

An extra layer of protection to keep your ostomy system in place.

EASE™

Strips

- Helps keep your skin barrier in place
- Can be applied to cover the edges of the skin barrier
- Designed to be thin and flexible
- Water-resistant

HCPCS $\triangle 4362$ (4 strips = 1 hillable unit of $\triangle 4362$)

1101 05 11 1302 (1301p3 = 15 maste and 51 11 1302)				
Size	Product Code	Quantity		
1 1/5 x 3 1/2" (3 x 9cm)	422163	80 strips		

CONVATEC®

Ostomy Appliance Belt

• Used to help keep your ostomy system in place

HCPCS A4367

Size	Product Code	Box of
Adjustable up to 42"	175507	1

DIAMONDS™

Gelling and Odor Control Sachets

- Gels liquid contents
- Helps reduce or eliminate excess gas that can lead to pouch ballooning
- Contains activated charcoal which helps to eliminate odors

HCPCS A4422

Product Code	Box of
TR105	100 Sachets

$\mathsf{DUOLOCK}^\mathsf{TM}$

Curved Tail Closure

• Helps secure the end of ConvaTec® Drainable pouch

HCDCS V1363

Product Code		Box of	
	175652	10	

CONVATEC®

Tail Closure

• Helps secure the end of ConvaTec® Drainable pouch

1101 03 74303		
Product Code	Box of	
175635	10	















CONVATEC®

Night Drainage Container Set*

- Plastic vessel to use for overnight urine collection
- Container capacity of 2000 cc/ml
- Reusable, can be washed out between uses

HCPCS A5102

Product Code	Quantity
027060	1 set

*Contains: 1 ConvaTec® Night Drainage Container with tubing, 1 ConvaTec® Night Drainage Container Cover and 1 Universal Adapter

CONVATEC®

Night Drainage Container Tubing

• Replacement tubing for the ConvaTec® Night Drainage Container

HCPCS A4331

110. 00 71 1331		
Product Code	Box of	
027062	1	



Solutions that help promote healthy skin by protecting it from irritation and damage caused by adhesives.

SENSI-CARE™

Sting-Free Adhesive Remover Wipe

- Designed to remove adhesive residue from the skin
- Usage with each appliance change helps prevent residue build-up
- Silicone based
- Alcohol free
- CHG Compatible¹⁷

HCPCS A4456

Liquid per Wipe	Product Code	Box of	Boxes per Case*
3 mL	413500	30	20

*Sold by the case

SENSI-CARE™

Sting-Free Adhesive Releaser Spray

- Tracks between adhesive & skin, gently releasing the adhesive, tape or appliance
- Helps minimize painful skin stripping on removal¹⁸
- Silicone based
- Alcohol free
- Non-Aerosol
- CHG Compatible¹⁷

HCPCS A4455

Spray Bottle Size	Product Code	Box of	Boxes per Case*
50 mL	413499	1	12
150 mL	420798 (Sterile)	1	12

*Sold by the case

ALLKARE®

Adhesive Remover Wipe



· Contains alcohol

HCPCS A4456

Product Code	Box of
037436	50
037443	100



IRRIGATION

VISI-FLOW®

Irrigator with Stoma Cone

HCPCS A4398 & A4399



VISI-FLOW®

Irrigation Sleeve

- Compatible with Natura[™] and SUR-FIT® Natura[™] skin barriers
- Color: Transparent
- 1 Tail Closure per box
- Quantity = 5

HCPCS A4397

	Flange Size	Product Code
•	1 1/2" (38mm)	401911
•	1 3/4" (45mm)	401912
•	2 1/4" (57mm)	401913
•	2 3/4" (70mm)	401914

VISI-FLOW®

Stoma Cone



VISI-FLOW®

Irrigation Sleeve Tail Closure

HCPCS A4363

11C1 C3 A4303	
Product Code	Box of
175650	10

VISI-FLOW®

Irrigation Adaptor Faceplate

- Compatible with Natura[™] and SUR-FIT® Natura[™] skin barriers
- Quantity = 1

HCPCS A4361

	Flange Size	Product Code
•	1 3/4" (45mm)	401918
•	2 3/4" (70mm)	401919


VISI-FLOW®

Irrigation Starter Set

- Contains:
 - 1 Visi-Flow® Irrigator with Stoma Cone
 - 2 SUR-FIT® Natura™ Stomahesive® skin barriers
 - 2 Visi-Flow® Irrigation Sleeves
 - 2 SUR-FIT® Natura™ Closed-End Pouches
 - 1 Stoma lubricant
 - 1 Ostomy belt
 - 1 Visi-Flow® Irrigation Sleeve Tail Closure
 - 1 Brush
 - 1 Carrying Case
- Quantity = 1

HCPCS A4398, A4399, A4397, A5054, A4367, A4414 (for 401916), A4415 (for 401917), A4402

	Flange Size	Product Code
•	1 3/4" (45mm)	401916
•	2 3/4" (70mm)	401917

















eakin® WOUND & FISTULA Pouches

Eakin wound pouches can help contain fistula and handle exuding wounds. Pouch options include the ability to connect tubing for remote drainage, or a fold & tuck option.





EAKIN®

Horizontal Cutting Wound Pouch

(9.7 x 6.3" Cutting Area)

• Quantity = 5

HCPCS A6154

Remote Drainage Cutting Range Attachment and Tap Closure

Remote Drainage Attachment and Tap Closure with Access Window

Fold and Tuck Closure

9.7 x 6.3" (24.6 x 16cm) 839263

839264 839253

EAKIN®

Vertical Cutting Wound Pouch

(6.3 x 9.7" Cutting Area)

• Quantity = 5

HCPCS A6154

Cutting Range

Remote Drainage Attachment and Tap Closure

Remote Drainage **Attachment and Tap Closure** with Access Window

839265 6.3 x 9.7" (16 x 24.6cm) 839268

EAKIN®

Vertical Cutting Wound Pouch

(5.1 x 11.4" Cutting Area)

• Quantity = 5

HCPCS A6154

Cutting Range

Remote Drainage Attachment and Tap Closure

5.1 x 11.4" (13 x 29cm)

839266

EAKIN®

Eakin Access Window

- Add to Eakin wound and fistula pouches to allow instant access to the wound (may be used with 839252, 839253, 839262, 839263, 839265, 839266)
- Quantity = 5

HCPCS N/A

Size **Remote Drainage Attachment and Tap Closure**

N/A 839280

OSTOMY RESOURCES AND SUPPORT

for Your Patient

ConvaTec

ConvaTec is a global medical products and technologies company developing ostomy products and solutions for almost 40 years.

Toll Free: 1-800-422-8811 www.convatec.com/ostomy

Ostomysecrets®

A ConvaTec company, offers a full line of stylish and functional undergarments, swim wear and accessories designed to support your ostomy solution.

Toll Free: 1-877-613-6246 www.ostomysecrets.com

me+™

Support, insights & product information for those living with an ostomy. Enrolling in me+ is fast, easy and free.

Toll Free: 1-800-422-8811 e-mail: CIC@convatec.com www.mepluscare.com

OTHER OSTOMY RESOURCES

The United Ostomy **Associations** of America (UOAA)

National organization providing support, information and advocacy.

Toll Free: 1-800-826-0826 www.ostomy.org

American Cancer Society

Nationwide communitybased, voluntary healthcare organization.

Toll Free: 1-800-ACS-2345 www.cancer.org

Wound Ostomy Continence Nurses Society™ (WOCN®)

Professional nursing society dedicated to providing expert health care to individuals with wound, ostomy and continence concerns.

Toll Free: 1-888-224-9626 www.wocn.org

Crohn's & Colitis Foundation of America

Non-profit, volunteer-driven organization dedicated to finding a cure for Crohn's Disease and ulcerative colitis.

Toll Free: 1-800-932-2423 www.ccfa.org

A4385 Skin Barrier, Solid 4x4 or Equivalent, Extended Wear As needed As needed Ostomy Pouch, Drainable, w/Extended Wear Barrier Attached (Ostomy Pouch, Drainable, w/Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4390 Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity Ostomy Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity (1-piece) A4391 Pouch, Drainable, w/Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) A4392 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) A4393 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) A4396 Ostomy Irrigation Bag 2 per 6 months A4397 Irrigation Sieeve 4 per month A4408 Ostomy Irrigation Sag 2 per 6 months A4409 Ostomy Skin Barrier, Non-Pectin Based, Paste 4 oz per month A4400 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4401 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4402 Skin Barrier, W/Flange, Extended Wear w/ Convexity As needed A4403 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4410 Skin Barrier, w/Flange, Extended Wear x= 4x4 As needed A4411 Skin Barrier, w/Flange, Extended Wear x= 4x4 As needed A4412 Skin Barrier, w/Flange, x= 4x4 20 per month A4415 Skin Barrier, w/Flange, x= 4x4 20 per month A4416 Pouch, Osed, w/Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4410 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4420 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4421 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4422 Pouch, Closed for use on Barrier w/Flange, w/Filter (1-piece) A4423 Pouch, Drainable, w/f Barrier Attached, w/ Filter (1-piece) A4424 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A5505 Pouch, Drainable, w/f Barrier Affange, w/Faucet-type Tap As needed A5505 Pouch, Drainable, w/f Barrier Affange, w/Faucet-type Tap As ne	HCPCS	HCPCS Descriptors	Medicare Usual Maximum Qty.
A4362 Skin Barrier, Solid, 4x4 or Equivalent A4363 Ostomy Clamp, Any Type, Replacement Only A567 Ostomy Belt 1 per month A4363 Ostomy Skin Barrier Liquid (spray or brush) 2 oz per month A4365 Ostomy Skin Barrier Liquid (spray or brush) 10 oz per 6 month A4368 Skin Barrier, Solid 4x4 or Equivalent, Extended Wear A4368 Pouch, Drainable, w/Extended Wear Barrier Attached Ostomy Pouch, Drainable, w/Extended Wear Barrier Attached Ostomy Pouch, Drainable, w/Barrier Attached, w/Built-In Cornexity, W/ Filter (1-piece) A4369 Pouch, Drainable, w/Barrier Attached, w/ Built-In Cornexity, W/ Filter (1-piece) A4390 Pouch, Drainable, w/Barrier Attached, w/ Built-In Cornexity, Ustomy Pouch, Urinary, w/ Barrier Attached, w/ Solomy Pouch, Urinary, w/ Barrier Attached, w/ Solomy Pouch, Urinary, w/ Barrier Attached, w/ Solomy Pouch, Urinary, w/ Extended Wear Barrier Attached, Ustomy Pouch, Urinary, w/ Extended Wear Barrier Attached, W/ Sullt-In Convexity (1-piece) A4393 Ustomy Prigation Cone/Catheter, and Brush A4399 Ostomy Irrigation Cone/Catheter, and Brush A4409 Ostomy Skin Barrier, Non-Pectin Based, Paste A4400 Ostomy Skin Barrier, Non-Pectin Based, Paste A4400 Ostomy Skin Barrier, Pectin-Based, Paste A4401 Ostomy Skin Barrier, Pectin-Based, Paste A4402 Skin Barrier, w/Flange, Extended Wear w/ Convexity A5 needed A4403 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4410 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4411 Skin Barrier, w/Flange, Set And 20 per month A4413 Pouch, Drainable, High Output, for Barrier w/Flange A5 needed A4414 Skin Barrier, w/Flange, > 4x4 A6 20 per month A4415 Skin Barrier, w/Flange, > 4x4 A6 20 per month A4416 Skin Barrier, w/Flange, > 4x4 A7 20 per month A4417 Ostomy Pouch, Closed, w/Barrier Attached, w/Filter A4410 Ostomy Pouch, Closed, w/Barrier Attached, w/Filter A4410 Ostomy Pouch, Drainable, w/Barrier Attached, w/Filter A4410 Ostomy Pouch, Drainable, w/Barrier Attached, w/Filter A4410 Ostomy Pouch, Drainable, w/Barrier Attached A4410 Ostomy Pouch, Drainable, w/Barrier Attached	A4331	Extension Drainage Tube	1 per 6 months
A4363 Ostomy Clamp, Any Type, Replacement Only A367 Ostomy Belt A4367 Ostomy Belt A4368 Ostomy Skin Barrier Liquid (spray or brush) 2 oz per month A4371 Ostomy Skin Barrier Liquid (spray or brush) 2 oz per 6 month A4367 Skin Barrier, Sloid 4x4 or Equivalent, Extended Wear A4368 Pouch, Drainable, wiExtended Wear Barrier Attached Ostomy Pouch, Drainable, wiExtended Wear Barrier Attached, wi Built-In Convexity,	A4361	Ostomy Faceplate	3 per 6 months
A4367 Ostomy Skin Barrier Liquid (spray or brush) 2 oz per month A4371 Ostomy Skin Barrier Liquid (spray or brush) 2 oz per month A4371 Ostomy Skin Barrier Powder 10 oz per 6 month A4385 Skin Barrier, Solid 4x4 or Equivalent, Extended Wear As needed A4388 Pouch, Drainable, w/Extended Wear Barrier Attached, w/ Built-in As needed A4389 Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity 20 per month A4390 Pouch, Urinary, w/ Barrier Attached, w/ Built-in Convexity 20 per month A4391 Crigation Sieeve 4 per month A4392 Irrigation Sieeve 4 per month A4393 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-in Convexity (1-piece) 20 per month A4402 Lubricant 4 per month A4403 Ostomy Irrigation Sieeve 4 per month A4402 Lubricant 4 oz per month A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste 4 oz per month A4406 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4407 Skin Barrier, W/Flange, Extended Wear > 4x4 <td>A4362</td> <td>Skin Barrier, Solid, 4x4 or Equivalent</td> <td>20 per month</td>	A4362	Skin Barrier, Solid, 4x4 or Equivalent	20 per month
A4369 Ostomy Skin Barrier Liquid (spray or brush) A4371 Ostomy Skin Barrier Powder A4385 Skin Barrier, Solid 4x4 or Equivalent, Extended Wear A588 Abin Barrier, Solid 4x4 or Equivalent, Extended Wear A588 Pouch, Drainable, WExtended Wear Barrier Attached A589 Pouch, Drainable, WExtended Wear Barrier, Attached, w/Convexity Wilter (1-piece) A4390 Pouch, Drainable, WExtended Wear Barrier, Attached, w/Convexity Corowexity, Wilter (1-piece) A4391 Pouch, Drainable, WExtended Wear Barrier, Attached, w/Convexity Costomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/Convexity Costomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/Convexity Costomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/Convexity A4393 Ostomy Irrigation Bag A4399 Ostomy Irrigation Cone/Catheter, and Brush A4409 Ostomy Skin Barrier, Non-Pectin Based, Paste A4400 Ostomy Skin Barrier, Poetin-Based, Paste A4400 Ostomy Skin Barrier, Pectin-Based, Paste A4401 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4410 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4411 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4412 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange A4414 Skin Barrier, w/Flange, > 4x4 A4416 Pouch, Closed, w/Barrier Attached, w/Filter A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, Wilter (1-piece) A4419 Pouch, Closed, w/ Barrier Attached, w/ Filter (1-piece) A4420 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4421 Ostomy Supply, Miscellaneous A5 needed A4422 Pouch, Drainable, with filter for use w/Flange A4423 Pouch, Drainable, w/Extended Wear Barrier w/Flange, w/Faucet-type Tap A5 needed A4428 Pouch, Drainable, w/Extended Wear Barrier w/Flange, w/Faucet-type Tap A5 needed A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Barrier Attached, w/ Filter (1-piece) A5067 Pouch, Drainable, w/Extended Wear Barrier w/Faucet-type Tap A5068 Pouch, Drainable, w/Barrier Attached	A4363	Ostomy Clamp, Any Type, Replacement Only	As needed
A4371 Ostomy Skin Barrier Powder A385 Skin Barrier, Solid 4x4 or Equivalent, Extended Wear As needed As nee	A4367	Ostomy Belt	1 per month
A4385 Skin Barrier, Solid 4x4 or Equivalent, Extended Wear A4388 Pouch, Drainable, w/Extended Wear Barrier Attached Ostorny Pouch, Drainable, w/Barrier Attached, w/ Built-In Convexity, W. Filter (1-piece) A4390 Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity Ostorny Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity (1-piece) A4391 Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity (1-piece) A4392 Ostorny Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) A4393 Ostorny Irrigation Bag A4399 Ostorny Irrigation Bag A4399 Ostorny Irrigation Cone/Catheter, and Brush A4402 Lubricant A4402 Lubricant A4405 Ostorny Skin Barrier, Non-Pectin Based, Paste A4406 Ostorny Skin Barrier, Pectin-Based, Paste A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity A4408 Skin Barrier, w/Flange, Extended Wear w/ Convexity A4409 Skin Barrier, w/Flange, Extended Wear w/ Ax44 Ax4410 Skin Barrier, w/Flange, Extended Wear w/ Ax44 Ax4411 Skin Barrier, w/Flange, Extended Wear w/ Ax44 Ax4412 Pouch, Drainable, High Output, for Barrier w/Flange Ax4413 Pouch, Drainable, High Output, for Barrier w/Flange Ax4414 Skin Barrier, w/Flange, x-4x4 Ax4415 Skin Barrier, w/Flange, x-4x4 Ax4416 Pouch, closed, w/Barrier Attached, w/Filter Ax4417 Ostorny Pouch, Closed, w/ Barrier Attached, w/ Filter Ax4418 Pouch, closed, w/Barrier Attached, w/Filter Ax4419 Pouch, closed for use on Barrier w/Flange, w/Filter Ax4420 Ostorny Supply, Miscellaneous Ax4420 Ostorny Absorbent Material (Sheet/Pad/Crystal Packet) Ax4421 Ostorny Supply, Miscellaneous Ax4422 Pouch, Drainable, w/f Barrier Attached, w/ Filter (1-piece) Ax4424 Ostorny Absorbent Material (Sheet/Pad/Crystal Packet) Ax4425 Pouch, Drainable, w/f Barrier Attached, w/ Filter (1-piece) Ax4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap Ax needed Ax4427 Pouch, Drainable, w/f Barrier Attached, w/ Filter (1-piece) Ax4428 Pouch, Drainable, w/f Barrier Attached Ax449 Pouch, Drainable, w/f Barrier Attached Ax440 Pouch, Drainable, w	A4369	Ostomy Skin Barrier Liquid (spray or brush)	2 oz per month
Pouch, Drainable, w/Extended Wear Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity w/ Filter (1-piece) Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity (1-piece) Pouch, Drainable, w/Extended Wear Barrier, Attached, w/ Built-In Convexity (1-piece) Pouch, Drainable, w/Extended Wear Barrier, Attached, w/ Built-In Convexity (1-piece) Pouch, Drainable, w/Extended Wear Barrier, Attached, w/ Built-In Convexity (1-piece) Pouch, Drainable, w/Extended Wear Barrier, Attached, w/ Built-In Convexity (1-piece) Pouch, Drainable, W/Extended Wear Barrier, Attached, w/ Built-In Convexity (1-piece) Pouch, Drainable, W/Extended Wear Barrier, Attached, w/ Built-In Convexity w/ Pouch, Drainable, W/ Pouch, Drainable, W/ Pouch, Drainable, W/ Pouch, Drainable, W/ Pouch, W/ Pouch, W/ Pouch, W/ Pouch, Drainable, W/ Pouch, W/ Pou	A4371	Ostomy Skin Barrier Powder	10 oz per 6 months
Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Built-In Convexity, Convexity, Convexity, Convexity, Constant Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity Ostomy Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity (1-piece) Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, W/ Built-In Convexity (1-piece) A4393 Ostomy Irrigation Bag A4399 Ostomy Irrigation Bag A4399 Ostomy Irrigation Cone/Catheter, and Brush A4401 Ostomy Skin Barrier, Non-Pectin Based, Paste A4406 Ostomy Skin Barrier, Non-Pectin Based, Paste A4407 Ostomy Skin Barrier, Pectin-Based, Paste A4408 A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4410 Skin Barrier, w/Flange, Extended Wear - 4x4 As needed A4411 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4413 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4414 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4415 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4416 Pouch, closed, w/Barrier Attached, w/Filter A4416 Pouch, closed, w/Barrier Attached, w/Filter A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4422 Ostomy Supply, Miscellaneous As needed A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Supply, Miscellaneous As needed A4424 Ostomy Supply, Miscellaneous As needed A4424 Ostomy Supply, Miscellaneous As needed A4425 Pouch, Urinary, w/ Extended Wear Barrier w/Flange Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4427 As needed A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4429 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4420 Ostomy Drainable, w/Extended Wear Barrier w/Faucet-type Tap As needed A5505 Pouch, Drainable, w/Extended Wear Barrier w/Faucet-type Tap As needed A6506 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier Attached Oper month A5063 Pouch,	A4385	Skin Barrier, Solid 4x4 or Equivalent, Extended Wear	As needed
Ad389 Convexity, W/ Filter (1-piece) Ad390 Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity Ostomy Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity (1-piece) Ad392 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) Ad393 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) Ad394 Irrigation Sleeve Ad395 Ostomy Irrigation Bag Ad399 Ostomy Irrigation Cone/Catheter, and Brush Ad402 Lubricant Ad405 Ostomy Skin Barrier, Non-Pectin Based, Paste Ad406 Ostomy Skin Barrier, Pectin-Based, Paste Ad407 Skin Barrier, w/Flange, Extended Wear w/ Convexity Ad408 Skin Barrier, w/Flange, Extended Wear w/ Convexity Ad409 Skin Barrier, w/Flange, Extended Wear w- Convexity Ad410 Skin Barrier, w/Flange, Extended Wear x = 4x4 As needed Ad411 Skin Barrier, w/Flange, Extended Wear x = 4x4 As needed Ad413 Pouch, Drainable, High Output, for Barrier w/Flange Ad414 Skin Barrier, w/Flange, > 4x4 Ad415 Skin Barrier, w/Flange, > 4x4 Ad416 Pouch, Closed, w/Barrier Attached, w/Filter Ad417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) Ad419 Pouch, closed for use on Barrier w/Flange, w/Filter Ad410 Ostomy Supphy, Miscellaneous Ad421 Ostomy Supphy, Miscellaneous Ad422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) Ad423 Pouch, Drainable, w/th filter for use w/Flange Ad424 Pouch, Drainable, w/th filter for use w/Flange Ad425 Pouch, Drainable, w/th filter for use w/Flange Ad426 Pouch, Drainable, w/th filter for use w/Flange Ad427 Pouch, Drainable, w/th filter for use w/Flange Ad428 Pouch, Drainable, w/th filter for use w/Flange Ad429 Pouch, Drainable, w/Extended Wear Barrier Attached, w/ Filter (1-piece) Ad420 Pouch, Drainable, w/Extended Wear Barrier Attached, w/ Filter Ad5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter Ad600 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter Ad600 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter	A4388	Pouch, Drainable, w/Extended Wear Barrier Attached	As needed
A4390 Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity Ostomy Pouch, Urinary, w/ Barrier Attached, w/ Built- In Convexity (1-piece) A4393 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built- In Convexity (1-piece) A4399 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built- In Convexity (1-piece) A4399 Ostomy Irrigation Bag A4399 Ostomy Irrigation Cone/Catheter, and Brush A4402 Lubricant A4402 Lubricant A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste A4406 Ostomy Skin Barrier, Pectin-Based, Paste A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity A5 needed A4409 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4410 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4411 Skin Barrier, w/Flange, Extended Wear - 4x4 A5 needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange A5 needed A4414 Skin Barrier, w/Flange, > 4x4 A4415 Skin Barrier, w/Flange, > 4x4 A4416 Pouch, closed, w/ Barrier Attached, w/ Built- In Convexity, w/ Filter (1-piece) A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built- In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4410 Ostomy Supply, Miscellaneous A4420 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4421 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4422 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4423 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4424 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4426 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4427 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4428 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4429 Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A5054 Pouch, Drainable, w/ Barrier Attached Pouch, Drainable, w/ Barrier Attached A5055 Pouch, Drainable, w/ Barrier Attached A5066 Pouch, Drainable,	A4389	Convexity,	20 per month
Ostomy Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity (1-piece) A4393 Virigation Sleeve 4 per month work Built-In Convexity (1-piece) A4397 Irrigation Sleeve 4 per month 4 per month 2 per 6 months 4 per month 5 per month 6 per month	A4390		As needed
(1-piece) Ad393 Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) Ad397 Irrigation Sleeve 4 per month Ad398 Ostomy Irrigation Bag 2 per 6 months Ad399 Ostomy Irrigation Bag 2 per 6 months Ad409 Lubricant 4 oz per month Ad400 Sotomy Skin Barrier, Non-Pectin Based, Paste 4 oz per month Ad400 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month Ad400 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed Ad400 Skin Barrier, w/Flange, Extended Wear v/ Convexity As needed Ad400 Skin Barrier, w/Flange, Extended Wear v/ Convexity As needed Ad410 Skin Barrier, w/Flange, Extended Wear v/ Convexity As needed Ad411 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed Ad413 Pouch, Drainable, High Output, for Barrier w/Flange As needed Ad414 Skin Barrier, w/Flange, <= 4x4 20 per month Ad415 Skin Barrier, w/Flange, > 4x4 20 per month Ad416 Pouch, closed, w/Barrier Attached, w/Filter Ad417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) Ad419 Pouch, closed for use on Barrier w/Flange, w/Filter Ad420 Ostomy Supply, Miscellaneous As needed Ad421 Ostomy Supply, Miscellaneous As needed Ad422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) Ad422 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) Ad423 Pouch, Urinary, w/ Extended Wear Barrier w/Flange 20 per month Ad424 Pouch, Urinary, w/ Extended Wear Barrier w/Flange 20 per month Ad425 Adhesive Remover Wipes, Any Type 50 per month Ad455 Adhesive Remover Wipes, Any Type 50 per month Ad456 Adhesive Remover Wipes, Any Type 50 per month Ad505 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5067 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5068 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5069 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5060 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrie	V4303		20 per month
A4397 Irrigation Sleeve 4 per month A4398 Ostomy Irrigation Bag 2 per 6 months A4399 Ostomy Irrigation Cone/Catheter, and Brush 2 per 6 months A4402 Lubricant 4 oz per month A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste 4 oz per month A4406 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4407 Skin Barrier, Wi-Range, Extended Wear w/ Convexity As needed A4409 Skin Barrier, wi-Range, Extended Wear v 4x4 As needed A4410 Skin Barrier, wi-Range, Extended Wear v 4x4 As needed A4411 Skin Barrier, wi-Range, Extended Wear v 4x4 As needed A4413 Pouch, Drainable, High Output, for Barrier wi-Range A4413 Pouch, Closed, wi-Barrier Attached, wi-Rilter A4415 Skin Barrier, wi-Range, v 4x4 20 per month A4416 Pouch, closed, wi-Barrier Attached, wi-Rilter A4417 Ostomy Pouch, Closed, wi-Barrier Attached, wi-Rilter A4419 Pouch, closed for use on Barrier wi-Range, wi-Filter A4410 Ostomy Supply, Miscellaneous A4420 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4421 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4423 Pouch, Drainable, with filter for use wi-Range A4424 Pouch, Drainable, with filter for use wi-Range A4428 Pouch, Urinary, wi-Extended Wear Barrier wi-Raucet-type Tap A4428 Pouch, Urinary, for use on Barrier wi-Range, wi-Raucet-type Tap A4428 Pouch, Urinary, for use on Barrier wi-Range, wi-Raucet-type Tap A4429 Pouch, Urinary, for use on Barrier wi-Range, wi-Raucet-type Tap A5054 Pouch, Urinary, for use on Barrier wi-Range A4432 Pouch, Urinary, for use on Barrier wi-Range A5054 Pouch, Drainable, wi-Extended Wear Barrier Attached, wi-filter A5057 Pouch, Drainable, wi-Extended Wear Barrier wi-Raucet-type Tap A5064 Pouch, Drainable, wi-Extended Wear Barrier wi-Raucet-type Tap A5065 Pouch, Drainable, wi-Extended Wear Barrier wi-Raucet-type Tap A5066 Pouch, Drainable, wi-Range A5077 Pouch, Drainable, wi-Range A5078 Pouch, Drainable, wi-Range A5079 Pouch, Drainable, wi-Range A5071 Pouch, Drainable, wi-Range A5071 Pouch, Drainable, wi-Range A5072 Pouch, Drainable, wi-		Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached,	·
A4398 Ostomy Irrigation Bag 2 per 6 months A4399 Ostomy Irrigation Cone/Catheter, and Brush 2 per 6 months A4402 Lubricant 4 oz per month A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste 4 oz per month A4406 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4407 Skin Barrier, WFlange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4411 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, < 4x4 20 per month A4415 Skin Barrier, w/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Flitter A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/Flitter A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4419 Pouch, Closed for use on Barrier w/Flange, w/Filter A4420 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4420 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4428 Pouch, Drainable, with filter for use w/Flange A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4428 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4428 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4429 Pouch, Urinary, for use on Barrier w/Flange A4430 Pouch, Urinary, for use on Barrier w/Flange A4431 Pouch, Urinary, for use on Barrier w/Flange A4432 Pouch, Urinary, for use on Barrier w/Flange A4433 Pouch, Urinary, for Barrier w/Flange A4445 Adhesive Remover Wipes, Any Type A5054 Pouch Closed, for Barrier w/Flange A5055 Oper month A5056 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5067 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5068 Pouch, Drainable, w/Extended Wear Barrier w/Flange A5071 Pouch, Urinary, for Barrier	A4207		A nor month
A4399 Ostomy Irrigation Cone/Catheter, and Brush A4402 Lubricant A4402 Lubricant A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste A4406 Ostomy Skin Barrier, Pectin-Based, Paste A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear - 4x4 As needed A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4411 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4412 Skin Barrier, w/Flange, < - 4x4 As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, < - 4x4 A4416 Pouch, Closed, w/Barrier Attached, w/Filter A4417 Ostomy Pouch, Closed, w/Barrier Attached, w/Filter A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4420 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap Adhesive Remover or Solvent A4427 Pouch Closed, for Barrier w/Flange, w/Faucet-type Tap Adhesive Remover Wipes, Any Type Douch Closed, for Barrier w/Flange Douch, Drainable, w/Extended Wear Barrier Attached, w/filter A5055 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5056 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5067 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5068 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5069 Ostomy Accessory, Convex Insert A5060 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5062 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, w/Extended Wear Barrier Attached A5064 Pouch, Drainable, w/Extended Wear Barrier Attached A5065 Pouch, Drainab		· ·	·
A4402 Lubricant A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste A4406 Ostomy Skin Barrier, Non-Pectin Based, Paste A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4411 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4412 Skin Barrier, w/Flange, < = 4x4 As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, < = 4x4 A20 per month A4415 Skin Barrier, w/Flange, > 4x4 A4416 Pouch, closed, w/Barrier Attached, w/Filter G50 per month A4417 Ostomy Pouch, Closed, w/Barrier Attached, w/ Built-In Convexity, w/Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4428 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4429 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5068 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5069 Ostomy Accessory, Convex Insert A5070 Pouch, Urinary, w/Barrier Attached A5090 Ostomy Accessory, Convex Insert A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 20 per month A5120 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month			
A4405 Ostomy Skin Barrier, Non-Pectin Based, Paste 4 oz per month A4406 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4411 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, <= 4x4 20 per month A4415 Skin Barrier, w/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter 60 per month A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Urinary, w/ Extended Wear Barrier w/Flange 20 per month A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4429 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4426 Adhesive Remover or Solvent 2 oz oz per month A4456 Adhesive Remover or Solvent 2 oz oz per month A456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5066 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5067 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5068 Pouch, Drainable, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, for Barrier w/Flange 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5073 Pouch, Urinary, for Barrie			·
A4406 Ostomy Skin Barrier, Pectin-Based, Paste 4 oz per month A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity As needed A4409 Skin Barrier, w/Flange, Extended Wear <= 4x4 As needed A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4411 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, <= 4x4 20 per month A4415 Skin Barrier, w/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter 60 per month A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Urinary, w/ Extended Wear Barrier w/Flange 20 per month A4428 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4425 Adhesive Remover or Solvent 2 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A4564 Pouch Orainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5063 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Urinary, w/Barrier Attached 20 per month A5063 Pouch, Urinary, w/Barrier Attached 20 per month A5064 Pouch Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5065 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5066 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5067 Pouch, Urinary, w/Barrier A			·
A4407 Skin Barrier, w/Flange, Extended Wear w/ Convexity A4409 Skin Barrier, w/Flange, Extended Wear <= 4x4 A4409 Skin Barrier, w/Flange, Extended Wear <= 4x4 A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 A5 needed A4411 Pouch, Drainable, High Output, for Barrier w/Flange A5 needed A4412 Skin Barrier, w/Flange, <= 4x4 A4415 Skin Barrier, w/Flange, <= 4x4 A4416 Pouch, closed, w/Barrier Attached, w/Filter A4416 Pouch, closed, w/Barrier Attached, w/Filter A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4421 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4427 As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4428 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4442 As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4445 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type 50 per month A4556 Pouch Closed, for Barrier w/Flange 60 per month A5057 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5056 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5064 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached			
A4409 Skin Barrier, w/Flange, Extended Wear <= 4x4 As needed A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, <= 4x4 20 per month A4415 Skin Barrier, w/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter 60 per month A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Drainable, with filter for use w/Flange 20 per month A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4455 Adhesive Remover or Solvent 2 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5063 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5063 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5063 Pouch, Drainable, w/extended Wear Barrier Attached 20 per month A5073 Pouch, Urinary, w/Barrier Attached 20 per month A5093 Ostomy Accessory, Convex Insert 10 pe		•	·
A4410 Skin Barrier, w/Flange, Extended Wear > 4x4 As needed A4413 Pouch, Drainable, High Output, for Barrier w/Flange As needed A4414 Skin Barrier, w/Flange, <= 4x4 20 per month A4415 Skin Barrier, w/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter 60 per month A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Drainable, with filter for use w/Flange 20 per month A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4455 Adhesive Remover Wipes, Any Type 50 per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5067 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5068 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, m/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Urinary, w/Barrier Attached 20 per month A5079 Pouch, Urinary, w/Barrier Attached 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5093 Ostomy Accessory, Convex Insert 10 per months A5120 Skin Barrier, Wipes, each 150 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months			
A4413 Pouch, Drainable, High Output, for Barrier w/Flange A4414 Skin Barrier, w/Flange, <= 4x4 20 per month A4415 Skin Barrier, w/Flange, >= 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4410 Ostomy Supply, Miscellaneous A4421 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4427 Adhesive Remover or Solvent A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4433 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4445 Adhesive Remover Wipes, Any Type A4456 Adhesive Remover Wipes, Any Type A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5066 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5067 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5068 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5069 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5093 Ostomy Accessory, Convex Insert A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable A5120 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent A5122 Skin Barrier, Solid, 8 X 8 or Equivalent A5122 Skin Barrier, Solid, 8 X 8 or Equivalent A5122 Skin Barrier, Solid, 8 X 8 or Equivalent A5122 Skin Barrier, Solid, 8 X 8 or Equivalent		, , ,	
A4414 Skin Barrier, w/Flange, <= 4x4 20 per month A4415 Skin Barrier, w/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter 60 per month A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Drainable, with filter for use w/Flange 20 per month A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4455 Adhesive Remover or Solvent 20 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5061 Pouch, Urinary, w/Barrier Attached 20 per month A5063 Pouch, Urinary, w/Barrier Attached 20 per month A5079 Pouch, Urinary, w/Barrier Attached 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5072 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, w/Barrier Attached 20 per month A5074 Pouch, Urinary, w/Barrier Attached 20 per month A5075 Pouch, Drainable, for use w/Flange 20 per month A5076 Pouch, Drainable, for use w/Flange 20 per month A5077 Pouch, Urinary, w/Barrier Attached 20 per month A5078 Pouch, Urinary, w/Barrier Attached 20 per month A5079 Pouch, Urinary, w/Barrier Attached 20 per month A5070 Pouch, Urinary, w/Barrier Attached 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5072 Pouch, Urinary, w/Barrier Attached 20			
A4415 Skin Barrier, W/Flange, > 4x4 20 per month A4416 Pouch, closed, w/Barrier Attached, w/Filter 60 per month A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Drainable, with filter for use w/Flange 20 per month A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4455 Adhesive Remover or Solvent 2 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5061 Pouch, Urinary, w/Barrier Attached 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, w/Barrier Attached 20 per month A5074 Pouch, Urinary, for Barrier w/Flange 20 per month A5075 Pouch, Urinary, w/Barrier Attached 20 per month A5076 Pouch, Urinary, w/Barrier Attached 20 per month A5077 Pouch, Urinary, w/Barrier Attached 20 per month A5078 Pouch, Urinary, for Barrier w/Flange 20 per month A5079 Pouch, Urinary, w/Barrier Attached 20 per month A5070 Pouch, Urinary, w/Barrier Attached 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5072 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, w/Barrier Attached 20 per month A5074 Pouch, Urinary, w/Barrier Attached 20 per month A5075 Pouch, Urinary, w/Barrier Attached 20 per month A5076 Pouch, Urinary, w/Barrier Attached 20 per month A5077 Pouch			
A4416 Pouch, closed, w/Barrier Attached, w/Filter A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4421 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4427 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4448 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4455 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Barrier Attached A5063 Pouch, Drainable, w/Barrier Attached A5064 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, for Barrier w/Flange A5074 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5074 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch,	A4414	Skin Barrier, w/Flange, <= 4x4	20 per month
A4417 Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter 60 per month A4421 Ostomy Supply, Miscellaneous As needed A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) No Fee Schedule A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) 20 per month A4425 Pouch, Drainable, with filter for use w/Flange 20 per month A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4455 Adhesive Remover or Solvent 2 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, w/Barrier Attached 20 per month A5074 Pouch, Urinary, w/Barrier Attached 20 per month A5075 Pouch, Urinary, m/Barrier Attached 20 per month A5076 Pouch, Urinary, m/Barrier Attached 20 per month A5077 Pouch, Urinary, m/Barrier Attached 20 per month A5078 Pouch, Urinary, m/Barrier Attached 20 per month A5079 Pouch, Urinary, m/Barrier Attached 20 per month A5070 Stomy Accessory, Convex Insert 10 per month A5071 Sedide Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months	A4415	Skin Barrier, w/Flange, > 4x4	20 per month
W/ Filter (1-piece) A4419 Pouch, closed for use on Barrier w/Flange, w/Filter A4421 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4429 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A44455 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, m/Barrier Attached A5064 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5074 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5074 Pouch, Urinary, w/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, w/Barrier Attached A5079 Pouch, Urinary, w/Barrier Attached A5070 Pouch, Urinary, w/Barrier At	A4416	Pouch, closed, w/Barrier Attached, w/Filter	60 per month
A4421 Ostomy Supply, Miscellaneous A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4426 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4427 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A44445 Adhesive Remover or Solvent A4455 Adhesive Remover Wipes, Any Type A4456 Adhesive Remover Wipes, Any Type A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Barrier Attached A5063 Pouch, Drainable, for use w/Flange A5063 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5072 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, for Barrier Attached A5074 Pouch, Urinary, for Barrier Attached A5075 Pouch, Drainable, for use w/Flange A5070 Pouch, Drainable, for use w/Flange A5071 Pouch, Urinary, m/Barrier Attached A5072 Pouch, Urinary, for Barrier Attached A5073 Pouch, Urinary, for Barrier Attached A5074 Pouch, Urinary, for Barrier w/Flange A5093 Ostomy Accessory, Convex Insert A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable A5120 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4417		60 per month
A4422 Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4455 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, w/Barrier Attached A5064 Pouch, Drainable, for use w/Flange A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, w/Barrier Attached A5078 Pouch, Urinary, for Barrier w/Flange A5079 Pouch, Urinary, for Barrier w/Flange A5070 Ostomy Accessory, Convex Insert A5070 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable A5070 Skin Barrier, Wipes, each A5120 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month 20 per month 20 per month	A4419	Pouch, closed for use on Barrier w/Flange, w/Filter	60 per month
A4424 Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) A4425 Pouch, Drainable, with filter for use w/Flange A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4455 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, m/Barrier Attached A5064 Pouch, Drainable, m/Barrier Attached A5075 Pouch, Urinary, w/Barrier Attached A5076 Pouch, Urinary, w/Barrier Attached A5077 Pouch, Urinary, m/Barrier Attached A5078 Pouch, Urinary, for Barrier w/Flange A5079 Pouch, Urinary, for Barrier w/Flange A5070 Ostomy Accessory, Convex Insert A5070 Skin Barrier, Wipes, each A5120 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month 20 per month 20 per month	A4421	Ostomy Supply, Miscellaneous	As needed
A4425 Pouch, Drainable, with filter for use w/Flange 20 per month A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap As needed A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap 20 per month A4455 Adhesive Remover or Solvent 2 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4422	Ostomy Absorbent Material (Sheet/Pad/Crystal Packet)	No Fee Schedule
A4428 Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4455 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, m/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5071 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5093 Ostomy Accessory, Convex Insert A5002 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable A5102 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month 20 per month 20 per month 20 per 6 months 20 per 6 months	A4424	Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece)	20 per month
A4432 Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap A4455 Adhesive Remover or Solvent A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange A5055 Stoma Cap, each A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5063 Pouch, Drainable, w/Barrier Attached A5063 Pouch, Drainable, for use w/Flange A5071 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, w/Barrier Attached A5093 Ostomy Accessory, Convex Insert A5093 Ostomy Accessory, Convex Insert A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable A5120 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4425	Pouch, Drainable, with filter for use w/Flange	20 per month
A4455 Adhesive Remover or Solvent 2 oz per month A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4428	Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap	As needed
A4456 Adhesive Remover Wipes, Any Type 50 per month A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4432	Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap	20 per month
A5054 Pouch Closed, for Barrier w/Flange 60 per month A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4455	Adhesive Remover or Solvent	2 oz per month
A5055 Stoma Cap, each 31 per month A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter 40 per month A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter 40 per month A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A4456	Adhesive Remover Wipes, Any Type	50 per month
A5056 Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Barrier Attached A5063 Pouch, Drainable, for use w/Flange A5071 Pouch, Urinary, w/Barrier Attached A5073 Pouch, Urinary, for Barrier w/Flange A5093 Ostomy Accessory, Convex Insert A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable A5120 Skin Barrier, Wipes, each A5122 Skin Barrier, Solid, 8 X 8 or Equivalent A5122 Skin Barrier, Solid, 8 X 8 or Equivalent A5075 Pouch, Urinary, for Barrier w/Flange A5120 Skin Barrier, Solid, 8 X 8 or Equivalent A5121 Schin Barrier, Solid, 8 X 8 or Equivalent A5122 Schin Barrier, Solid, 8 X 8 or Equivalent A5122 Schin Barrier, Solid, 8 X 8 or Equivalent A5125 Schin Barrier, Solid, 8 X 8 or Equivalent A5126 Schin Barrier, Solid, 8 X 8 or Equivalent A5127 Schin Barrier, Solid, 8 X 8 or Equivalent A5128 Schin Barrier, Solid, 8 X 8 or Equivalent A5129 Schin Barrier, Solid, 8 X 8 or Equivalent A5120 Schin Barrier, Solid, 8 X 8 or Equivalent A5120 Schin Barrier, Solid, 8 X 8 or Equivalent A5120 Schin Barrier, Solid, 8 X 8 or Equivalent A5120 Schin Barrier, Solid, 8 X 8 or Equivalent A5120 Schin Barrier, Solid, 8 X 8 or Equivalent A5120 Schin Barrier, Solid, 8 X 8 or Equivalent	A5054	Pouch Closed, for Barrier w/Flange	60 per month
A5057 Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A5055	Stoma Cap, each	31 per month
A5061 Pouch, Drainable, w/Barrier Attached 20 per month A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A5056	Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter	40 per month
A5063 Pouch, Drainable, for use w/Flange 20 per month A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A5057	Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter	40 per month
A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A5061	Pouch, Drainable, w/Barrier Attached	20 per month
A5071 Pouch, Urinary, w/Barrier Attached 20 per month A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month	A5063	Pouch, Drainable, for use w/Flange	20 per month
A5073 Pouch, Urinary, for Barrier w/Flange 20 per month A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month		Pouch, Urinary, w/Barrier Attached	20 per month
A5093 Ostomy Accessory, Convex Insert 10 per month A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month		•	20 per month
A5102 Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable 2 per 6 months A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month			·
A5120 Skin Barrier, Wipes, each 150 per 6 months A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month		31	•
A5122 Skin Barrier, Solid, 8 X 8 or Equivalent 20 per month			
	A6154	Wound Pouch	As needed

020922	77	76
	27	76
	28	
	58	
021901	8	31
	56	
	56	
	38	
022750)5	56
022751	5	56
	2	
	3	
022754	l5	56
022755	5	56
	55	
	7	
022758	35	56
022759	95	56
)5	
022762	<u> </u>	56
	35	
	15	
	5	
022766	5	56
	75	
	35	
	95	
022770)5	56
022771	5	5.5
	35	
)7	
025542	2	32
)8	
	2	
	3	
037439)7	79
037443	38	36
	17	
	75	
	73	
125258	33	31
)	
)3	
	13	
125262	23	31
	33	
	l3	
	53	31
125266	33	
125267		31
	7	
125268	73	32
	33	32
125269	3)	32 32 32
125269 125270	33)3)3	32 32 32
125269 125270	3)	32 32 32
125269 125270 125271	3	32 32 32 32
125269 125270 125271 125272	3	32 32 32 32 32
125269 125270 125271 125272 125273	3	32 32 32 32 32
125269 125270 125271 125272 125273 125274	3	32 32 32 32 32
125269 125270 125271 125272 125274 125274	3	32 32 32 32 32 32
125269 125270 125271 125272 125274 125274 125275	3	32 32 32 32 32 32 32
125269 125270 125271 125272 125274 125274 125275	3	32 32 32 32 32 32 32
125269 125270 125271 125273 125274 125274 125275 125276	3	32 32 32 32 32 32 32 32
125269 125270 125271 125273 125274 125275 125276 125277 125330	3	32 32 32 32 32 32 32 32 32 32 32
125269 125270 125271 125272 125273 125274 125275 125276 125277 125330 125331 125333	3	32 32 32 32 32 32 32 32 32 32 32 32 32
125269 125270 125271 125272 125273 125274 125275 125276 125277 125330 125331 125333	3	32 32 32 32 32 32 32 32 32 32 32 32 32
125269 125270 125271 125272 125273 125274 125276 125277 125331 125333 125333	3	32 32 32 32 32 32 32 32 32 32 32 32 32 3
125269 125270 125271 125272 125273 125274 125275 125276 125277 125331 125333 125334 125333	3	32 32 32 32 32 32 32 32 32 32 36 56
125269 125270 125271 125272 125273 125274 125276 125277 125330 125333 125334 125338 125336	S	32 32 32 32 32 32 32 32 32 32 32 32 32 3
125269 125270 125271 125272 125273 125274 125276 125277 125330 125331 125333 125334 125336 125336 125336	S	32 32 32 32 32 32 32 32 32 32 32 32 32 3
125269 125270 125271 125272 125273 125274 125276 125277 125377 125333 125333 125333 125335 125335 125335 125337 125337	3	32 32 32 32 32 32 32 32 32 32 32 32 32 3
125269 125270 125271 125272 125273 125274 125275 125277 125330 125331 125333 125334 125335 125336 125337 125338 125338	S	32332332332332332332333233323332333233
125269 125270 125271 125272 125273 125274 125277 125277 125330 125331 125336 125336 125336 125336 125336 125336 125336	3	32 33 32 33 33 33 33 33 33 33 33 33 33 3
125269 125270 125271 125272 125273 125274 125276 125277 125330 125333 125334 125335 125336 125337 125338 125338 125338 125338	S	32 32 33 32 33 32 33 33 33 33 33 33 33 3
125269 125270 125271 125273 125274 125275 125276 125277 125330 125331 125334 125336 125336 125336 125336 125336 125336 125336 125336 125336 125336	S	32 332 332 332 332 332 332 332 556 566 566 566 566 566
125269 125270 125271 125272 125273 125276 125276 125277 125330 125331 125336 125336 125336 125336 125336 125336 125336 125336 125336 125336 125336	3	32 332 332 332 332 332 332 332 332 332
125269 125270 125271 125272 125273 125274 125276 125276 125277 125330 125331 125333 125336 125336 125336 125336 125336 125336 125338 125338 125338 125338 125341 125341 125342 125343		32 332 332 332 332 332 332 332 332 332
125269 125270 125271 125272 125273 125276 125276 125277 125276 125277 125330 125331 125335 125336 125336 125336 125336 125336 125336 125340 125340 125342 125342 125342 125343	3	32 32 32 32 32 32 32 33 32 55 56 66 66 66 66 66 66 66 66 66 66 66
125269 125270 125271 125273 125274 125275 125276 125277 125330 125331 125333 125334 125335 125336 125337 125336 125340 125341 125342 125342 125343 125342 125343 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12534 12535 12534 12534 12534 12534 12535 12534 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535 12534 12535	3	32 32 32 33 32 33 32 33 32 33 32 33 32 33 32 33 32 33 33
125269 125270 125271 125272 125273 125274 125276 125277 125330 125331 125333 125336 125336 125340 125340 125342 125342 125342 125342 125342 125342 125342 125342 125342 125342 125344 125342 125344 125342 125344 125342 125344 125342 125344 125342 125344 125342 125344 125342 125344 125342 125344 12534 1		32 32 32 33 32 33 32 33 32 33 32 33 32 55 56 66 66 66 66 66 66 66 66 66 66 66
125268 125271 125271 125272 125272 125272 125277 125371 125331 125333 125333 125334 125341 125341 125341 125341 125342 12534 125342 125	3	32 32 32 32 32 32 32 32 32 32 32 55 56 66 66 66 66 66 66 66 66 66 66 66
125268 125271 125271 125272 125272 125272 125277 125371 125331 125333 125333 125334 125341 125341 125341 125341 125342 12534 125342 125		32 32 32 32 32 32 32 32 32 32 32 55 56 66 66 66 66 66 66 66 66 66 66 66

1253												
	57											56
1253												
1253												
1253												
1253												
1253												
1253												
1253	365.	•••••	•••••	• • • • •		 • • • • • •	• • • • • • • • • • • • • • • • • • • •	• • • • • •	••••	•••	• • •	59
1253												
1253												
1253												
1253												
1755	0/.					 • • • • • •	• • • • • • • • • • • • • • • • • • • •	• • • • • •		•••	• • •	84
1756												
1756												
1756												
1756	552.					 • • • • • •	• • • • • • • • • • • • • • • • • • • •	• • • • • •		•••	• • •	84
1757	61.					 				•••	• • •	5/
1757												
1757												
1757												
1757												
1757	66.	•••••	•••••	• • • • •		 • • • • • •	• • • • • •	• • • • • •	••••	•••	• • •	5/
1757 1757	'08.		••••			 	• • • • • •	• • • • • •	••••	•••	• • •	5/
1757	770	•••••				 ••••	••••	•••••		•••	• • •	57
1757 1757	/U.			••••	••••	 ••••	••••	•••••	••••	•••	• • •	5/
1757	/ I.			••••	••••	 ••••	••••	•••••	••••	•••	• • •	57
1757 1757	772.			••••	••••	 ••••	••••	•••••	••••	•••	• • •	57
1757	10. 777				••••	 ••••	••••	•••••	••••	•••	• • •	50
1757	11.				••••	 ••••	••••	•••••	••••	•••	• • •	50
1757												
1757												
1757	781			••••		 				•••	• • •	56
1757												
1757												
1757												
1757	85.					 						56
1757												
1757	791.					 						59
1757												
1757												
1757												
1757												
1757												
1757												
	98.											
						 						59
1839	10.					 						59 80
4005	910. 598.					 						59 80 55
4005	910. 598. 599.					 						59 80 55 55
4005 4005 4015	910. 598. 599. 500.					 					 	59 80 55 55 36
4005 4005 4015 4015	910. 598. 599. 500. 501.					 						59 80 55 55 36 36
4005 4005 4015 4015 4015	910. 598. 599. 500. 501.					 						59 80 55 55 36 36 36
4005 4005 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503.					 						59 80 55 55 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504.					 						59 80 55 55 36 36 36 36 36 35
4005 4005 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504.					 						59 80 55 55 36 36 36 36 35
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505.											59 80 55 55 36 36 36 36 35 35
4005 4005 4015 4015 4015 4015 4015 4015	910. 98. 999. 600. 601. 602. 603. 604. 605. 606. 607. 608.											59 80 55 55 36 36 36 36 35 35 35
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 506. 507. 508.											59 80 55 55 36 36 36 35 35 35 35 35
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 506. 507. 508. 509.											59 80 55 55 36 36 36 35 35 35 35 35
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 506. 507. 508. 510. 511.											59 80 55 55 36 36 36 35 35 35 35 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 506. 507. 508. 509. 511. 511.											59 80 55 36 36 36 35 35 35 35 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 506. 507. 508. 511. 512.											59 80 55 36 36 36 35 35 35 35 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 506. 506. 507. 511. 512. 513.											59 80 55 36 36 36 35 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 506. 511. 512. 513. 514.											59 80 55 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 506. 511. 512. 513. 514. 515.											59 80 55 36 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 511. 511. 511. 511. 511. 512. 521.											59 80 55 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 98. 998. 600. 601. 602. 603. 604. 605. 606. 611. 611. 611. 611. 611. 611. 612. 612. 613. 614. 615. 615. 615. 616. 617. 6											59 80 55 36 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4005 4015 4015 4015 4015 4015 4015	910. 98. 998. 600. 601. 602. 603. 604. 605. 606. 611. 611. 611. 611. 612. 611. 612. 612. 612. 613.											59 80 55 36 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4005 4015 4015 4015 4015 4015 4015	910. 998. 999. 600. 601. 602. 603. 603. 606. 607. 608. 611. 613. 614. 615. 621. 622. 632. 632. 632. 632. 633. 634. 644. 655. 665.											59 80 55 36 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 999. 600. 601. 602. 603. 604. 606. 607. 608. 611. 611. 612. 613. 621. 622. 622. 622. 622. 622. 622.											59 80 55 55 36 36 36 35 35 35 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 999. 600. 601. 602. 603. 604. 605. 606. 611. 612. 611. 612. 622. 622. 622. 622. 622. 622.											59 80 55 55 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 999. 600. 601. 602. 603. 604. 605. 606. 611. 612. 611. 612. 622. 622. 622. 622. 622. 622. 622.											59 80 55 55 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 999. 900. 901. 902. 903. 904. 905. 906. 907. 908. 907. 908. 909.											59 80 55 55 36 36 36 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 502. 503. 504. 505. 507. 511. 511. 512. 522. 522. 522. 522. 523. 524. 525. 527. 528. 530. 530.											59 59 55 55 56 36 36 36 36 36 36 36 36 36 3
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 503. 503. 504. 506. 511. 511. 511. 522. 522. 522. 522. 523. 524. 525. 527. 532. 533. 533.											59 80 55 55 56 36 36 36 36 36 36 36 36 36 3
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 501. 503. 503. 506. 507. 508. 511. 512. 513. 521. 522. 523. 524. 525. 526. 532. 533. 533. 533.											59 80 55 55 36 36 36 36 36 36 36 36 36 36
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 999. 900.											59 80 55 56 36 36 36 36 36 36 36 36 36 3
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 500. 500. 500. 500. 500. 500. 500. 501. 511. 511. 512. 522. 522. 522. 522. 532. 532. 532. 533. 533. 533. 533. 533.											5980 5555 3636 3635 3635 3636 3636 3636 363
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 501. 502. 503. 504. 506. 511. 512. 513. 522. 522. 522. 523. 523. 532. 532. 533. 533. 533. 533.											5980 5980 5555 5555 366 363 363 363 363 363 363 3
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 599. 5001. 5002. 5003. 5004. 5008. 5011. 5112. 5114. 5122. 5223. 5224. 5225. 523. 5331. 5327. 5331. 5331. 5331. 5331. 5331. 5331. 5331.											5980 5980 5555 5636 3636 3633 3633 3633 3633 363
4005 4015 4015 4015 4015 4015 4015 4015	010. 010. 0598. 0599. 0599. 0601. 0602. 0603. 0604. 0605. 0608. 0607. 0608. 0611. 0611. 0611. 0611. 0612. 0613. 06											598005555555555555555555555555555555555
4005 4015 4015 4015 4015 4015 4015 4015	010. 010. 0598. 0599. 0599. 0501. 0502. 0503. 0504. 0505. 0508. 0507. 0508. 0511. 0512. 0513. 0514. 0522. 0523. 0524. 0524. 0525. 0523. 0533. 0533. 0533. 0533.											598005555555555555555555555555555555555
4005 4015 4015 4015 4015 4015 4015 4015	910. 598. 598. 500. 500. 501. 502. 508. 508. 508. 511. 511. 512. 523. 524. 525. 526. 531. 532. 533. 533. 533. 533. 533. 533. 533. 533. 534. 533. 534. 535. 537. 538.											5980 5555 5555 3663 3636 3636 3636 3636 363
4005 4015 4015 4015 4015 4015 4015 4015	910. 988. 998. 900.											598055555366336336555555553663365555555555
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 998. 900.											598005555553663363636363363355555556336363636363636363636363636363636363636
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 998. 900. 901. 902. 903. 904. 905. 906. 907.											598005555553663365555555555555555555555555
4005 4015 4015 4015 4015 4015 4015 4015	910. 998. 999. 900. 901. 902. 903. 904. 905. 906. 907. 908.											5980555556363636363636335555555556363636363

401546	39	40545741
401547		40545841
401548		405467
401549		40546841
401550		40547342
401551		40547442
401552		40547542
401553	39	40547642
401554		40547742
401555		40547842
401557		40547942
401558		40548042
401559		40548142
401573		40548242
401574		40548341
401575	31	40926842
401576	31	40926942
401577	31	40927042
401905	27	40927567
401906	31	41126436
401909	39	41126536
401911		41126636
401912		41126736
401913		41128836
401914		41128936
401916		41129036
401917		41129136
401918		41129236
401919		41130935
401920		41131035
401921		41131135
401922 401923		411312
401924		411361
401925		411362
401926		411489
401927		411490
401928		411491
401929		411492
401930		411493
401931		411631
401932		41163272
401933	36	41163372
401934		41163774
401935		41163874
401936	36	41163974
401989	87	41164074
401993	83	41164173
401994	83	41164273
401995	83	41165741
401996		41165841
402239		41165941
403948		41180029
404006		41180133
404007		41180229
404008		41180333
404009		41180429
404011		41180533
404011		41180629
404012 404013		411807
404014		413142 57
404014		413143
404016		413144
404017		413145
404018		41314657
404025		41314757
404026	36	41314857
404027	36	41314957
404028		41315057
404029		41315157
404030		41315257
404031		41315327
404032		41315427
404033		41315527
404034		41315627
404592		41315727
404593		41315927
404594		41316027
405406 405446		41316127
405446		413162
405448		413164
405449		413165
405450		413166
405451		41316727
405452		41316827
405453		41316938
405456		41317038

41317138
41317238
41317338
41317438
41317538
41317638
41317730
41317830
41317930 41318030
413181
413182
413183
41318430
41318530
41318630
41318730
41318830
41331235
41331335
41331435
413334
413336
413417
41341829
41341929
41342029
41342133
41342233
413423
413424
413436
413437
413438
41343940
41349986
41350085
41350178
413502
413509
413510
413515
413516
41351748
41351848
41351948
41352048
41352148
413522
416400
416401
416402
41640337
41640437
41640538
41640637
416407
416408
416410
416411
41641237
41641337
41641434
41641534
41641634
416417
416418
416420
416421
416422
41642334
41642434
41647234
41670051
41670151
41670251 41670351
41670451
416706
41670751
41670951
41671051
41671251

41671351
41671551
416718
416721
416724
41672548
41672848
41672948
416732
416736
41673748
41673849
416739
416741
416742
41674349
416744
416745
416747
416748
41674949
416750
416751
416775
41677644
41677744
416778
416779
416781
41678243
41678343
416795
416797
41679843
41679943
416900
416901
416903
41690847
416915
416916
416918
41691961
41692061
41692162
416922
416924
41692562
416926
416927
416929
41693061
416932
416934
416940
41694265
416946
416948
416950
416975
41697646
420695
420696
420793
42079479
42079579
420796
420798
42103333
42103433
421035
421040
42104128

PRODUCT CODE INDEX REFERENCES

4210424
42104364
421044
421046
42104763
42107440
42107540
421076
421454
421458
42146032
42146228
421464
421612
421614
421615
42161649
42161749
421618
421619
421621
421622
42162353
42162453
42162553
421626
421628
42162953
42163153
42163829
421639
421641
421642
42164329
42167637
42167737
42167837 42167937
421680
42168137
42168237
42168337
421687
421689
421738
42173934
42174034
421746
421747
421749
42179638
42179738
421798
421816
421822
42182351
42184847
421853
421854
421856
42185748
42185846
42185948
421859
421859
421859
421859 48 421860 48 421861 48 421862 48 421893 37 422159 49
421859. 48 421860. 48 421861. 48 421862. 48 421893. 37 422159. 49 422163. 83
421859. 48 421860. 48 421861. 48 421862. 48 421862. 48 421893. 37 422159. 49 422163. 83
421859 48 421860 48 421861 48 421862 48 421893 37 422159 49 422163 83 422363 50 422364 50
421859. 48 421860. 48 421861. 48 421862. 48 421893. 37 422159. 49 422163. 83 422363. 50 422364. 50 422365. 50
421859 48 421860 48 421861 48 421862 48 421893 37 422159 49 422163 83 422364 50 422365 50 422599 50 42259 50 42259 50 42250 50
421859. 48 421860. 48 421861. 48 421862. 48 421862. 48 421893. 37 422159. 49 422163. 83 422364. 50 422364. 50 422539. 50 422540. 50
421859. 48 421860. 48 421861. 48 421862. 48 421862. 48 421863. 37 422159. 49 422163. 83 422364. 50 422364. 50 422364. 50 422540. 50 422540. 50 422541. 50
421859. 48 421860. 48 421861. 48 421862. 48 421862. 48 421893. 37 422159. 49 422163. 83 422364. 50 422364. 50 422539. 50 422540. 50

422552	
422553	.54
422557	.54
422558	
422559	.54
422560	.54
422562	
422564	
422565	
650408	.57
650422	.57
650828	
650829	
650830	.58
650831	
650832	
650833	
839001	
839002	
839003	
839004	
839005	
839006	
839010	
839250	
839251	
839252	
839253	
839260	
839261	
839262	
839263	
839264	
839265	
839266	
839267	
839268	
839280	
TR105	.84

- 1. Weiss, A., et. Al., Statistical Brief #154, Readmissions to U.S. Hospitals by Procedure, 2010, Healthcare Cost and Utilization Project, Agency for Healthcare Research and Quality, April 2013. Data derived from readmission rates of colostomy, Ileostomy, and other enterostomy procedures.
- 2. Hyde, L.Z., Al-Masrou, A.M., Kuritzkes, B.A., Suradkar, K., Valizadeh, N., Kiran, R.P. (2018). Readmissions after colorectal surgery: not all are equal. International Journal of Colorectal Disease, 33 (12), 1667-1674.
- 3. Taneja, C., Netsch, D., Rolstad, B.S., Inglese, G., Lamerato, L., Oster, G. (2017). Clinical and Economic Burden of Peristomal Skin Complications in Patients with Recent Ostomies. Journal of Wound, Ostomy Continence Nursing, 44 (4), 350-357.
- Ostomy, Advocacy page. United Ostomy Association of America website. https://www.ostomy.org/tag/survey/. Accessed August 8, 2018.
- Bulkeley, J.E., McMullen, C.K., Grant, M., Wendel, C., Hornbrook, M.C. Krouse, R.S. (2018). Ongoing Ostomy Self-Care Challenges of Long-Term Rectal Cancer Survivors. Supportive Care in Cancer, 26 (11), 3933-3939.
- 6. ConvaTec Internal in vitro lab testing. OCT 0112.
- 7. Maria Teresa Szewczyk, MD, PhD; Grazyna Majewska, RN, ETN; Mary V. Cabral, MS, FNP-BC, CWOCN-AP; and Karin Hölzel-Piontek, RN; The Effects of Using a Moldable Skin Barrier on Peristomal Skin Condition in Persons with an Ostomy: Results of a Prospective, Observational, Multinational Study, Ostomy Wound Management 2014;60(12):16–26.
- 8. Bosio G., Pisani F., Lucibello, L. A Proposal for Classifying Peristomal Skin Disorders: Results of a Multicenter Observational Study. OWM. 2007; 53(9): 38-43.
- 9. Clinical Evaluation of Sur Fit Natura Cut to Fit Skin Wafers Gillian Seaton. Technical data on file, reference document WHR13940 TA238 2014.
- 10. Rolstad, BS. Principles and techniques in the the use of convexity. Ostomy Wound Manage. 1996. 42(1); 24-32.
- 11. Pontieri, Lewis V. Basics of ostomy care. Medsurg Nurs. 2006; 15(4): 199-202.
- 12. Hanley, J. Convex stoma appliances: are we getting it right? British Journal of Nursing (Stoma Supplement). 2013; 22(16).
- 13. Turnbull, GB. The convexity controversy. Ostomy Wound Manage. 2003; 49(1): 16-www.o-wm.com/article/1203. Accessed 29 Feb, 2016.
- 14. DRG data 2015, in T5 markets: Japan, UK, USA, Germany, China .
- Boyd, K. Clinical Protocols for Stoma Care Use of convex appliances. Nursing Standard. 2004; 18(20): 37-38.
- 16. Burch, 2011: Management of Stoma Complications. NursingTimes.net. 10 Nov 2011.
- 17. Data on file. Trio Protocol. PMS 002. December 2007. ConvaTec Inc.
- 18. In-Vitro: Trio Products 110099. Test Report AC-0170. July, 2012. Data on file. ConvaTec Inc.



ENROLL YOUR PATIENT IN THE **me+**™ PROGRAM

The me+[™] program is designed to give your patient the support, insights and products they need to live life in a way that works for them.

It's free to sign up.

Call or e-mail us today:

1-800-422-8811 (*M-F*, 8:30am - 7:00pm EST) CIC@convatec.com

Visit us at mepluscare.com



Scan this code to enroll your patient in our me+ program



For general product information visit www.convatec.com/ostomy

Disclaimer:

The reimbursement information provided by ConvaTec is intended to provide general information relevant to coding and reimbursement of ConvaTec's products only. Coverage and payment policies for the same insurer can vary from one region to another and may change from time to time because of ongoing changes in government and insurance industry rules and regulations. Therefore, please confirm HCPCS codes with your local DME-MAC, private insurer, or Medicaid agency before processing claims. ConvaTec does not guarantee coverage or payment of its products listed herein.