



I CAN
live free

STOMA CARE
Product Guide



ConvaTec



ConvaTec Product Specialists

0800 28 22 54 (UK)
1800 721 721 (ROI)

ConvaTec Ltd Main Switchboard

01244 584 300

ConvaTec Website and Email

www.convatec.co.uk
stoma.webcare@convatec.com



Me+™

0800 467 866 (UK)
1800 721 721 (ROI)
www.convatec.co.uk/meplus
www.convatec.ie/meplus



Amcare™ Website
www.amcaregroup.co.uk



**WE EXIST TO IMPROVE
THE LIVES OF THE PEOPLE
WE TOUCH.**

**We drive for excellence in all we do – anticipating
and addressing our customers’ needs with
advanced technologies and best-in-class
products and services.**

CONTENTS

One Piece

Esteem[®]+, Esteem[®], Ileodress[®], Urodress[®], Colodress[®] and Stomadress[®]

Esteem [®] + Drainable and Closed	8 - 12
Esteem [®] Drainable, Closed and Urostomy	13 - 15
Ileodress [®] Drainable	16 - 17
Stomadress [®] Drainable	17
Urodress [®] Urostomy	18
Colodress [®] Closed	19

Two Piece

Esteem *synergy*[®]+

Flat Baseplates	22 - 24
Convex Baseplates	24
Drainable Pouches	25 - 26
Closed Pouches	27
Urostomy Pouches	28

Natura[®]

Flat Baseplates	29 - 32
Convex Baseplates	33
Natura [®] Accessories	34
Drainable Pouches	35 - 36
High Output	36
Closed Pouches & Stoma Cap	37
Urostomy	38

Two Piece continued

Consecura[®]

Flat Baseplates	39
Drainable Pouches	40
Closed Pouches	41
Urostomy Pouches	41

Accessories

Accessories	44 - 46
-------------	---------

Paediatric Products


Esteem [®]	48
Esteem <i>synergy</i> [®]	49
Natura [®]	50

CHOOSING WHAT'S RIGHT FOR YOUR PATIENT

The following pages detail the many types of products that are available.
Different products may be more appropriate at different stages or events in your patient's life.

Whether it is for a colostomy, ileostomy or urostomy, each ConvaTec stoma system is designed with your patient's /person's lifestyle in mind. The systems are made so that they are discreet even under leisure or sport clothing. They are secure and comfortable to the body and are designed to protect peristomal skin. They can be changed to fit in with the activity of the moment and are easy to use.

Look out for these symbols and abbreviations:



LOOK OUT
for this icon



PLEASE NOTE
Do not cut the
baseplates

Coupling size

This is the size of the attachment which connects our pouches to our baseplates. This makes it easy to identify pouches to the corresponding baseplate.

They are also colour co-ordinated as displayed on this table below.

Our baseplates are either cut-to-fit (CTF), pre-cut or ConvaTec Mouldable Technology™. This will be illustrated in this baseplate table.



Esteem synergy®+ Drainable Pouches with Lock-it Pocket™ Closure 10 units

Coupling	Opaque	Clear	Clear
	Standard	Standard	Standard
35mm	416796 ∇	416795 ∇	416797
48mm	416799 ∇	416798 ∇	416900
61mm	416902 ∇	416901 ∇	416903

The combination that you choose is based on your patient's needs and personal preference.

OUR ADHESIVES

Stomahesive®

Stomahesive® technology is a hydrocolloid adhesive that gently adheres to the skin, protecting it from the damaging effects of stoma effluent. Stomahesive® technology combines adhesive and hydrocolloid components for optimal performance under various conditions. Stomahesive® technology is a unique product, only from ConvaTec.

Durahesive®

Durahesive® skin baseplates are designed for people whose stoma output is mostly liquid. Unlike some other skin baseplates that can break down around liquid output, Durahesive® skin baseplates swell up to protect the stoma. This special effect, called turtlenecking, helps create a more secure seal without harming the skin. Durahesive® skin baseplates are also easy to fit, apply, and remove. They will adhere securely to the skin while showering, bathing and swimming.

ConvaTec Mouldable Technology™

Over 95% of new ostomates who started on ConvaTec Mouldable Technology™ maintained healthy skin¹. Our mouldable baseplates and pouches are designed to gently mould to any stoma shape or size for clinically proven skin protection¹.

This combination helps people feel in control and gives them the freedom to get on with the other important things in life.

- Secure - The elastic-like seal securely hugs any stoma. Designed to minimise leaks for healthier skin¹.
- Simple - Just roll, mould and hold to accommodate any stoma shape and size. No cutting, no stretching, no guesswork.
- Comfortable - Thin and flexible, designed to move with the body.

References 1. Clinical Summary: Osmose Study, Multinational Evaluation of the Peristomal Condition in New Ostomates Using Mouldable Skin Barriers. Data on file.

A photograph of two women exercising outdoors. The woman in the foreground has short, light-colored hair and is wearing a grey long-sleeved shirt. She is smiling and has her arms raised. The woman in the background is wearing a blue and white patterned shirt and also has her arms raised. The background is a bright, sunny outdoor setting, possibly a beach or park, with a blurred horizon.

I CAN
live free

ONE PIECE

ONE PIECE

Esteem®

Esteem®+ system is our one-piece solution, combining the baseplate and pouch in a single unit. Easy and convenient, it's one simple step to attach or remove and provides a fresh start every time. Our Esteem®+ drainable pouches feature an advanced InvisiClose™ secure closure and Lock-it Pocket™ for added comfort. The Esteem®+ range also benefits from our skin friendly modified Stomahesive® Technology with a touch of mineral oil in every baseplate.

We designed Esteem®+ one-piece with a range of benefits that provide security, confidence and freedom every day. It's available with closed and drainable pouches with a choice of pre-cut, cut-to-fit and mouldable. Our opaque Esteem®+ pouches also feature a split window cover.

Esteem® system with ConvaTec Mouldable Technology™ (CMT) Drainable Pouch



Esteem® System with ConvaTec Mouldable Technology™ 10 units

Hole Size	Opaque Split Cover	Clear
	Standard	Standard
20-30mm	413521 ✓	413515 ✓
30-40mm	413522 ✓	413516 ✓
40-50mm	413523 ✓	413517 ✓

MOULDABLE RANGE

ESTEEM®

The convenience of a one-piece with the security and clinically proven¹ skin protection of our Mouldable skin barriers.

ESTEEM SYNERGY®


Our Mouldable Esteem *synergy*® baseplates combined with our Esteem *synergy*® pouches are simple to use, helping people to feel in control of their stoma so they can get on with the important things in life.

NATURA® ACCORDION

The Natura® skin barrier with Accordion Baseplate takes the pressure off of pouch application for improved coupling and comfort. Designed for tender abdomens or for anyone who wants everyday comfort, the Accordion Baseplate lifts easily to provide full finger width space so the pouch can be attached with minimal pressure.

NATURA®

Mouldable barriers combined with our Natura® pouches are simple to use, helping people to feel in control of their stoma so they can get on with the important things in life.



References 1. Clinical Summary: Osmose Study, Multinational Evaluation of the Peristomal Condition in New Ostomates Using Mouldable Skin Barriers. Data on file.



Remove mouldable disc from pouch and warm in hands.



Remove liner marked **#1**. This is the moulding side.



Roll and mould the disc to the size of the stoma. **Do not stretch, pull or cut.** Pinch into place once desired size is reached.



Remove liner marked **#2**. This side will go on the skin.



Place the disk on the body. Push the roll of adhesive toward the stoma to ensure a close fit. Skin should be clean and dry.



Remove liner marked **#3**.



Remove pouch release liner.




Place pouch on the mouldable disc.



Esteem® Drainable Flat Pouch with Lock-It Pocket™ Closure 10 units

Hole size	Opaque	Clear	Clear
	Standard	Standard	Standard
20-70mm	416719 ✓	416718 ✓	416721
25mm	416725 ✓	416724 ✓	
30mm	416729 ✓	416728 ✓	
35mm	416733 ✓	416732 ✓	
40mm	416737 ✓		



Esteem® system with ConvaTec Mouldable Technology™ Closed Pouch 30 units

Hole size	Opaque Split Cover
	Standard
20-30mm	413509 ✓
30-40mm	413510 ✓
40-50mm	413511 ✓



Esteem® Closed Pouch 30 units

Hole size	Opaque		Clear
	Small	Standard	Standard
20-70mm	416702 ✓	416701 ✓	416700 ✓
25mm	416705 ✓	416704 ✓	416703 ✓
30mm	416708 ✓	416707 ✓	416706 ✓
35mm	416711 ✓	416710 ✓	416709 ✓
40mm	416714 ✓	416713 ✓	416712 ✓
50mm		416715 ✓	



Esteem® Convex (Depth 5.08mm) Lock-It Pocket™ Closure 10 units

Hole size	Opaque		Clear
		Standard	Standard
19mm	416739 ✓		
22mm	416741 ✓		416740 ✓
25mm	416743 ✓		416742 ✓
28mm	416745 ✓		416744 ✓
32mm	416747 ✓		416746 ✓
35mm	416749 ✓		416748 ✓
38mm	416751 ✓		

Suitable for use with Belt S210.

Esteem®

Esteem® system is our one-piece solution, combining the baseplate and pouch in a single unit. Our Esteem® drainable pouches are available with either an InvisiClose™ closure or clip closure. The Esteem® one-piece is available with closed and drainable pouches with a choice of pre-cut and cut-to-fit.



Esteem® Drainable Pouch with InvisiClose® Closure 30 units

Hole size	Opaque		Clear
	Small	Medium	Medium
20-70mm		S5079 Ț	S5078 Ț
40mm	S5076 Ț		



Esteem® Drainable Pouch with Straight Clip Closure 10 units

Hole size	Opaque			Clear	
	Small	Medium	Large	Medium	Large
20-70mm	S5041 Ț	S5048 Ț	S5062 Ț	S5047 Ț	S5055 Ț
25mm	S5049 Ț				
30mm	S5064 Ț				

Suitable for use with Curved Clip S202 and Soft Wire Tie S205.



Esteem® Closed Pouch 30 units

Hole size	Opaque			Clear	
	Small	Medium	Large	Medium	Large
20-70mm	S5001 ✓	S5013 ✓	S5027 ✓	S5006 ✓	S5020 ✓
25mm	S5002 ✓	S5014 ✓	S5028 ✓	S5007 ✓	
30mm	S5003 ✓	S5015 ✓	S5029 ✓	S5008 ✓	
35mm	S5011 ✓	S5018 ✓	S5034 ✓	S5012 ✓	
40mm	S5004 ✓	S5016 ✓	S5030 ✓	S5009 ✓	
50mm		S5017 ✓	S5031 ✓		
60mm			S5032 ✓		
70mm			S5033 ✓		



Esteem® Convex (Depth 5.08mm) Pouch Hydrocolloid Collar 10 units

Hole size	Opaque	Clear
	Standard	Standard
19mm	S5151 ✓	S5161 ✓
22mm	S5152 ✓	S5162 ✓
25mm	S5153 ✓	S5163 ✓
28mm	S5154 ✓	S5164 ✓
32mm	S5155 ✓	S5165 ✓
35mm	S5156 ✓	S5166 ✓
38mm	S5157 ✓	S5167 ✓
45mm	S5158 ✓	S5168 ✓
50mm	S5159 ✓	S5169 ✓

Suitable for use with Belt S210.



Esteem® Urostomy Pouch with Fold-up Tap
10 units

Hole size	Opaque		Clear	
	Small	Standard	Small	Standard
13-45mm	S5120	S5140	S5110	S5130
20mm		S5141		S5131
25mm		S5142	S5112	S5132
30mm		S5143	S5113	S5133
35mm		S5144	S5114	S5134

Suitable for use with Night Bag S320 and Connector 1019.

One Piece Heritage

One Piece Heritage

Colodress®, Ileodress® and Urodress® systems are our heritage one piece range. Drainable pouches within this range are used with a clip closure. All products within our heritage range feature our skin friendly Stomahesive® technology within every baseplate.



Ileodress® Plus with Straight Clip Closure 10 units			
Hole size	Opaque		Clear
	Small	Standard	Standard
19-70mm	S430	S411	S420
25mm		S412	
32mm		S413	
38mm		S414	S423
45mm		S415	
50mm		S416	

Suitable for use with Curved Clip S202 and Soft Wire Tie S205.

One Piece Heritage



Ileodress® with Stomahesive® Inner Section and Microporous Collar with Straight Clip Closure
10 units

Hole size	Opaque		Clear
	Small	Standard	Standard
19-64mm	S851	S831	S841
25mm	S852	S832	
32mm	S853	S833	
38mm	S855	S835	S845
45mm		S836	S846
50mm		S837	

Suitable for use with Curved Clip S202 and Soft Wire Tie S205.



Stomadress®
10 units

Hole size	Large Clear	
	Hospital	Community
8-100mm	S438	S439

Product S438 comes with 10 x S202 Curved Clips.
Product S439 comes with 1 x S202 Curved Clip.



Urodress® Urostomy Pouch with Fold-up Tap 10 units

Hole size	Clear
	Standard
19-45mm	S896
25mm	S897
32mm	S898
38mm	S899
45mm	S900

Suitable for use with Night Bag S320 and Connector 1019.




Urodress® Deep Convex Urostomy Pouch 5 units

Hole size	Clear
	Standard
16mm	S910
19mm	S911
22mm	S912
25mm	S913
28mm	S914
32mm	S915
35mm	S916

Suitable for use with Night Bag S320 and Connector 1019.

C



Colodress® Plus with Stomahesive® Technology
30 units

Hole size	Opaque		Clear
	Mini	Standard	Standard
19-45mm	S901 ✓		
19-64mm		S861 ✓	S871 ✓
25mm	S902 ✓	S862 ✓	S872 ✓
32mm	S903 ✓	S863 ✓	S873 ✓
38mm	S904 ✓	S864 ✓	S874 ✓
45mm		S865 ✓	S875 ✓
50mm		S866 ✓	S876 ✓
64mm		S867 ✓	

C



Stoma Cap
30 units

Hole size	Opaque
19-50mm	S821 ✓



I CAN
live free

TWO PIECE

TWO PIECE

Esteem synergy®

Esteem synergy® system is a two-piece stoma system with adhesive coupling technology for comfort, flexibility and discretion. Designed with security and confidence in mind, Esteem synergy® offers a choice of innovative features as well as our strong and skin-friendly adhesives. Our Esteem synergy® baseplates and Esteem synergy® pouches join securely together via a soft, foam adhesive ring on the pouch which attaches to a 'landing zone' on the baseplate. Flexible and convenient, it's perfect for frequent changing and ideal for people on the go. It has a low-profile that moves with the body so it's comfortable and discreet. Our Esteem synergy® drainable pouches feature our advanced InvisiClose™ Tail Closure and hide-away Lock-it Pocket™.

Esteem synergy® baseplates are available in cut-to-fit, pre-cut and mouldable. Esteem synergy® pouches are available in a range of closed, drainable and urostomy.



ConvaTec Mouldable Technology™ Baseplate with Flexible Hydrocolloid Collar
10 units

Hole Size		Coupling	
		48mm	61mm
13-22mm	411657		
22-33mm	411658		
33-45mm			411659



Esteem *synergy*® Stomahesive® (Inner Section) Baseplate with White Flexible Collar
10 units

Hole size	Coupling		
	35mm	48mm	61mm
13-35mm	S1100		
13-48mm		S1101	
13-61mm			S1102
25mm	S1144		
32mm	S1146		
35mm	S1147		



Esteem *synergy*® Stomahesive® (Inner Section) Baseplate with Flexible Hydrocolloid Collar
10 units

Hole size	Coupling		
	35mm	48mm	61mm
13-35mm	S1190		
13-48mm		S1191	
13-61mm			S1192



Esteem synergy® Fully Stomahesive® Baseplate 10 units

Hole size	Coupling			
	35mm	48mm	61mm	89mm
13-35mm	S1130			
13-48mm		S1131		
13-61mm			S1132	
13-89mm				S1133



Esteem synergy® ConvaTec Mouldable Technology™ Convex Baseplate (Depth 5.83mm) with White Flexible Collar 10 units

Hole size	Coupling		
	35mm	48mm	61mm
13-22mm	S1180		
22-33mm		S1181	
33-45mm			S1182



Esteem *synergy*® Drainable Pouches with Lock-it Pocket™ Closure
10 units

Coupling	Opaque	Clear	Clear
	Standard	Standard	Standard
35mm	416796 ✓	416795 ✓	416797
48mm	416799 ✓	416798 ✓	416900
61mm	416902 ✓	416901 ✓	416903



Esteem *synergy*® Drainable Pouch with InvisiClose® Closure
Right Sided Pouches
30 units

Coupling	Opaque	Opaque	Opaque
	Small	Standard	Small
35mm	S1230F ✓	S1210F ✓	S1235N*
48mm	S1231F ✓	S1211F ✓	S1236N*
61mm		S1212F ✓	

*10 units per box.



Esteem *synergy*® Drainable Pouches with InvisiClose® Closure Left Sided Pouches 30 units

Coupling	Clear
	Standard
35mm	S1220F ∇
48mm	S1221F ∇
61mm	S1222F ∇



Esteem *synergy*® Drainable Pouch 10 units

Coupling	Opaque	Opaque	Clear	Clear
	Standard	Standard	Standard	X Large
48mm	S1021F ∇			
61mm			S1007F ∇	
89mm				S1009N



Suitable for use with Curved Clip S202, Soft Wire Tie S205 and Straight Clip S207. Clip not included with this product.


C



Esteem synergy® Closed Pouch 30 units

Coupling	Opaque		Clear	
	Small	Standard	Small	Standard
35mm	416777 ✓	416775 ✓		416774 ✓
48mm	416781 ✓	416779 ✓	416780 ✓	416778 ✓
61mm		416783 ✓		416782 ✓

C



Esteem synergy® Closed Pouch 30 units

Coupling	Opaque		Clear		Opaque	
	Small	Standard	Small	Standard	Mini	Small
35mm	S1053F ✓	S1062F ✓		S1050F ✓	S1060N*	
48mm	S1054F ✓	S1063F ✓	S1057F ✓	S1051F ✓		S1054N
61mm	S1055F ✓	S1064F ✓		S1052F ✓		

*20 units per box.

U



Esteem *synergy*® Urostomy Pouch with Accuseal® Tap
10 units

Coupling	Opaque		Clear	
	Standard	Small	Standard	
35mm	S1080	S1075	S1085	
48mm	S1081	S1076	S1086	
61mm	S1082		S1087	

Suitable with – Night bag: S320, Connector: 1035 (2 x 1035 in box).

Natura®

Natura® system is a two-piece stoma solution, with a separate baseplate and pouch, great for busy, active lives. Designed with security and confidence in mind, the Natura® system is our latest two piece system, offering a choice of innovative features as well as our strong and skin friendly adhesives. With Natura® system, the pouch can be changed whilst the baseplate stays firmly in place. Its snap-lock fastening audibly clicks so you know that it's securely attached. Whether it's playing with the kids, gardening, dinner with friends or running, Natura® pouches offer convenience and security for people with busy lives.

Natura® baseplates are available in cut-to-fit, pre-cut and mouldable. Natura® pouches are available in a range of closed, drainable and urostomy.

The Natura® Accordion baseplate takes the pressure off the pouch application for improved coupling and comfort. Designed for tender abdomens or for anyone who wants everyday comfort, the Accordion Baseplate lifts easily to provide full finger width space so the pouch can be attached with minimal pressure. Once applied the baseplate gives a discreet, low profile that is flexible and conforms to body contours.



Natura® Accordion Baseplate with ConvaTec Mouldable Technology™ 5 units

Hole Size		Coupling
	57mm	70mm
13-22mm	421036	
22-33mm	421037	
33-45mm		421038

The system that allows coupling without pressure on the abdomen.



Natura® ConvaTec Mouldable Technology™ Baseplate with Flexible Hydrocolloid Collar
10 units

Hole size	Coupling		
	45mm	57mm	70mm
13-22mm	S7550		
22-33mm	S7551		
33-45mm		S7552	
45-56mm			S7553



Natura® Stomahesive® Baseplate with Flexible Hydrocolloid Collar
10 units

Hole size	Coupling				
	32mm	38mm	45mm	57mm	70mm
13-19mm	S7201				
13-25mm		S7202			
13-32mm			S7203		
13-45mm				S7204	
13-57mm					S7205



Natura® Stomahesive® (Inner Section) Baseplate with Microporous Collar 10 units

Hole size	Coupling				
	32mm	38mm	45mm	57mm	70mm
13-19mm	S7238				
13-25mm		S7239			
13-32mm			S7240		
13-45mm				S7241	
13-57mm					S7242



Natura® Stomahesive® Baseplate with Oval Microporous Collar 5 units

Hole size	Coupling				
	32mm	38mm	45mm	57mm	70mm
13-19mm	S7244				
13-25mm		S7245			
13-32mm			S7246		
13-45mm				S7247	
13-57mm					S7248



Natura® Flexible Fully Stomahesive® Baseplate 10 units

Hole size	Coupling				
	32mm	38mm	45mm	57mm	70mm
13-19mm	S7340				
13-25mm		S7341			
13-32mm			S7342		
13-45mm				S7343	
13-57mm					S7344



Natura® Fully Stomahesive® Baseplate 10 units

Hole size	Coupling					
	32mm	38mm	45mm	57mm	70mm	100mm (pack of 5)
13-19mm	S7294					
13-25mm		S7295				
13-32mm			S7296			
13-45mm				S7297		
13-57mm					S7298	
13-88mm						S7299B



Natura® Convex Baseplate with ConvaTec Mouldable Technology™
10 units

Hole Size		Coupling
45mm		57mm
13-22mm	S7304 (depth 5.46mm)	
22-33mm	S7305 (depth 5.83mm)	
33-45mm		S7306 (depth 6.17mm)



Natura® Convex Durahesive® Baseplate (Inner Section)
(Depth 5.08mm) with a Microporous Collar
5 units

Hole Size		Coupling
45mm		57mm
13mm	S7325	
16mm	S7326	
19mm	S7327	
22mm	S7328	
25mm	S7329	
28mm	S7330	
32mm	S7331	
35mm	S7332	
38mm		S7335
41mm		S7336
45mm		S7337
50mm		S7338



Natura® Convex Inserts 5 units

Coupling

Max Stoma Size	38mm	45mm	57mm
18mm	S7624		
25mm		S7626	
28mm		S7627	
35mm			S7629

Convex inserts are not suitable with ConvaTec Mouldable Technology™ baseplates.



Natura® Low Pressure Adapters 10 units

Coupling

45mm	401993
57mm	401994
70mm	401995



Natura® Drainable Pouch with Lock-it Pocket™ Closure 10 units

Coupling	Opaque	Clear	Clear
	Standard	Standard	Standard
38mm	416415 ✓	416414 ✓	
45mm	416417 ✓	416416 ✓	416418
57mm	416420 ✓	416419 ✓	416421
70mm	416423 ✓	416422 ✓	
100mm			416472



Natura® Drainable Pouch with InvisiClose® Closure 30 units

Coupling	Opaque	Opaque*	Clear
	Small	Standard	Standard
32mm	S7436F ✓	S7447F ✓	
70mm		S7435F ✓	S7445N

*Comfort backing on one side.



Natura® Drainable Pouch 10 units

Coupling	Opaque		Opaque		Clear
	Small	Standard	Small	Standard	Standard
32mm		S7400 ✓	S7279	S7269	
38mm		S7401 ✓	S7280	S7270	S7229
45mm	S7408 ✓	S7402 ✓	S7281	S7271	S7230
57mm	S7409 ✓	S7403 ✓	S7282	S7272	S7231
70mm		S7404 ✓	S7283	S7273	S7232
100mm					S7233

Suitable for use with Curved Clip S202,
Soft Wire Tie S205 and Straight Clip S207.
Clip not included with this product.



Natura® High Output Pouch 5 units

Coupling	Clear	Opaque
45mm	420695 ✓	S7467 ✓
57mm	420696 ✓	S7468 ✓
70mm	420697 ✓	S7469 ✓

All come with replacement filters,
clip and night drainage adaptor.

Suitable for use with 1044 Adaptor
Kit for high output.


C



Natura® Closed Pouch 30 units

Coupling	Opaque	Opaque	Opaque
	Standard	Small	Standard
32mm	416400 ∇	416402	416401
38mm	416403 ∇	416405	416404
45mm	416406 ∇	416408	416407
57mm	416409 ∇	416411	416410
70mm	416412 ∇		416413

C



Natura® Closed Pouch 30 units

Comfort Backing One Side				Comfort Backing Two Sides
Coupling	Opaque		Opaque	Opaque
	Small	Standard	Mini	Standard
32mm		S7254 ∇	S7290*	S7215
38mm	S7291F* ∇	S7255 ∇	S7291*	S7216
45mm	S7292F* ∇	S7256 ∇	S7292*	S7217
57mm	S7293F* ∇	S7257 ∇	S7293*	S7218
70mm		S7258 ∇		S7219

Natura® Stoma Cap 25 units

Coupling	
45mm	S7251
57mm	S7252

*20 units per box.



Natura® Urostomy Pouch with Fold-up Tap 10 units

Coupling	Clear	
	Small	Standard
32mm	S7360	S7350
38mm	S7361	S7351
45mm	S7362	S7352
57mm	S7363	S7353
70mm		S7354

Suitable for use with Night Bag S320 and Connector 1019.



Natura® Urostomy Pouch with Accuseal® Tap 10 units

Coupling	Opaque		Clear
	Small	Standard	Standard
32mm	S7390	S7370	S7380
38mm	S7391	S7371	S7381
45mm	S7392	S7372	S7382
57mm	S7393	S7373	S7383
70mm			S7384

Suitable for use with Night Bag S320 and Connector 1035.



ConSecura® Flexible Baseplate with Microporous Collar 5 units

Hole size	Coupling			
	35mm	45mm	57mm	70mm
13-23mm	S601LP			
13-32mm		S602LP		
13-45mm			S603LP	
13-57mm				S604LP



ConSecura® Flexible Fully Stomaheasive® Baseplate 5 units

Hole size	Coupling			
	35mm	45mm	57mm	70mm
13-23mm	S590LP			
13-32mm		S591LP		
13-45mm			S592LP	
13-57mm				S593LP



ConSecura® Drainable Pouch 10 units

Coupling	Opaque	Clear
	Standard	Standard
35mm	S606LP	
45mm	S607LP	S612LP
57mm	S608LP	S613LP
70mm	S609LP	S614LP

Suitable for use with Curved Clip S202.
Soft Wire Tie S205 and Straight Clip S207.
Clip not included with this product.

C



ConSecura® Closed Pouch 30 units

Coupling	Opaque	Opaque	Clear
	Standard	Standard	Standard
35mm	S616LP ✓	S606LP	S630LP ✓
45mm	S617LP ✓	S607LP	S631LP ✓
57mm	S618LP ✓	S608LP	S632LP ✓
70mm	S619LP ✓	S609LP	S633LP ✓

U



ConSecura® Urostomy Pouches 10 units

	Accuseal® Tap		Fold-up Tap
Coupling	Opaque	Clear	Clear
	Standard	Standard	Standard
35mm	S655LP	S640LP	S635LP
45mm	S656LP	S641LP	S636LP
57mm	S657LP	S642LP	S637LP
70mm	S658LP		

Accuseal® Tap Pouches are suitable for use with Night Drainage System S320.

Accuseal® Pouches are suitable with 1035 Connector.

Fold-up Tap Pouches are suitable with 1019 Connector.



ICAN
live free

ACCESSORIES

ACCESSORIES

ACCESSORIES/SUPPLEMENTARY PRODUCT



Stomahesive® Seals

Size	Code	Quantity
Standard Seal (48mm width, 4mm Depth)	413503	10
Standard Thin Seal (48mm width, 3mm Depth)	413504	10
Large Seal (98mm Width, 4mm Depth)	413505	10



Adhesive Removers

	Code	Quantity
Niltac™ Sting Free Adhesive Remover Wipes	TR102	30
Niltac™ Sting Free Adhesive Remover Spray	TR101	50ml
ConvaCare® Lemon Scented Adhesive Remover Wipes	S208	100



Solidifiers and Odour Control

	Code	Quantity
Diamonds™ Gelling and Odour Control Sachets	TR105	100
Diamonds™ Gelling and Odour Control Sachets (SAMPLE- not on FP10)	STR105	5

ACCESSORIES/SUPPLEMENTARY PRODUCT



Skin Fillers, Pastes, Barriers

	Code	Quantity
Sensi-Care® Sting Free Protective Skin Barrier Foam Applicator (Sterile)	420793	1ml x 5
Sensi-Care® Sting Free Protective Barrier Foam Applicator (Sterile)	420795	3ml x 5
Sillesse™ Sting Free Barrier Wipes	TR103	30
Sillesse™ Sting Free Barrier Gentle Spray	TR104	50ml
Stomahesive® Paste	S105	60G
Orabase® Paste	S103	30G

Skin Fillers, Pastes, Barriers

	Code	Quantity
Orahesive® Powder	S106	25G
ConvaCare® Barrier Wipes	S209	100
DuoDERM® Extra Thin dressing		
100mm x 100mm	S161	10
Stomahesive® Wafer		
200mm x 200mm	S101	3
100mm x 100mm	S100	5
Urihesive™ Strips	S120	3 x 5

ACCESSORIES/SUPPLEMENTARY PRODUCT



Clips (for Drainable Pouches)

	Code	Quantity
Curved Clip	S202	10
Straight Clip	S207	10
Soft Wire Tie	S205	50



Belt, Loops Rods and Urine Night Drainage Systems

	Code	Quantity
White Belt	S210	1
Urine Night Drainage System, 2 litre	S320	5
Accuseal® Connector	1035	
Urostomy Fold-up Tap Connector	1019	
Ostomy Loop Rod, 65mm	22355	10
Ostomy Loop Rod, 95mm	22356	10
Careline™ Night Drainage Stand	45-70-B	
Careline™ Night Bag Bed Hanger	45-70-H	

For further advice and product support call the ConvaTec Product Specialist Team on 0800 282 254.



I CAN
live free

PAEDIATRIC PRODUCTS

PAEDIATRIC

ESTEEM® PAEDIATRIC PRODUCTS

Little Ones™ Drainable Pouch



One Piece Drainable Pouch Clear without Filter

Code	Description	Hole size	Quantity
S880	Extra small. Clear pouch. No filter.	8-50mm	15



Pouches

Code	Description	Hole size	Qty
S5090	Cut-to-fit extra small, clear drainable pouch (comfort backing one sided) with InvisiClose® closure, no filter.	0-23mm	10

ESTEEM SYNERGY® PAEDIATRIC PRODUCTS



Baseplates

Code	Description	Hole size	Qty
S1171	Extra small. Full Stomahesive® Baseplate. No starter hole. Use with S1271.	0-23mm	5
S1172	Full Stomahesive® Baseplate. For use with S1272F and S1274F.	4.8-31mm	5
S1173	Stomahesive® Baseplate with flexible hydrocolloid collar. For use with S1272F and S1274F.	4.8-31mm	5



Pouches

Code	Description	Hole size	Qty
S1271	Extra small clear. Drainable pouch with InvisiClose® closure. One sided features comfort backing. No filter. For use with S1171.	4.8-31mm	10
S1272F	Clear closed pouch. One side features comfort backing. Filter. For use with S1172 or S1173.	31mm	10
S1273F	Opaque drainable pouch. Both sides feature comfort backing. Filter. For use with S1172 or S1173.	31mm	10
S1274F	Opaque closed pouch. Filter. For use with S1172 or S1173.	31mm	10
S1275F	Clear drainable pouch. One side features comfort backing. Filter. For use with S1172 or S1173.	31mm	10

NATURA® PAEDIATRIC PRODUCTS



Flexible Baseplates

Code	Coupling size	Hole size	Quantity
S7811	32mm	12-26mm	5
S7812	45mm	12-38mm	5



Drainable Pouch Clear without Filter

Code	Coupling size	Quantity
S7880	32mm	10
S7881	45mm	10



Closed Pouch Opaque without Filter

Code	Coupling size	Quantity
S7891	32mm	20
S7892	45mm	20



Urostomy Pouch Clear

Code	Coupling size	Quantity
S7850	32mm	5



Amcare™ Group

We deliver products from all manufacturers and have been awarded British standards ISO 9001 certification, so you have total confidence your patient's prescription will be delivered as promptly and discreetly as possible. When your patient joins our award winning Amcare™ Group home delivery service they will receive a range of complimentary items and commitment to support life with a stoma.

www.amcaregroup.co.uk



Amcare™

(England & Wales)

0800 88 50 50

enquiry@amcaregroup.co.uk

Amcare™

(Scotland)

0800 45 65 922

enquiry@amcaregroup.co.uk

AlphaMed™

0800 515 317

BCA Direct™

0800 854 753

Homestyle Positive™

0800 393 889

Trent Direct™

0800 028 4131



"My favorite part
of the day is taking
Ramsey to the dog park."

Rick



**me+™ supports
the 'whole'
person for life**

Healthcare professionals spend their days helping patients live their best lives. me+ complements this care by surrounding patients with the support, insights and products they need to live life with a stoma on their own terms.



**me+™
care**



**me+™
support**



**me+™
answers**



**me+™
style**



**me+™
community**

© 2016 ConvaTec Inc. ®/™ indicates a trade mark
of the ConvaTec group. AP-016769-GB OST056

Speak to our friendly me+™ advisor on 0800 467 866