

# OSTOMY

## Product Catalogue

Australia & New Zealand



Our Products

ONE PIECE

TWO PIECE  
Mechanical Coupling

TWO PIECE  
Adhesive Coupling

Accessories

Urine Drainage  
Systems



ConvaTec

Ostomy Care

Product  
Code Index

# Our PRODUCTS

Our ostomy care products are designed to help give people the confidence to achieve the things they want in life – whether it's going for a run, going out with friends or family, or going for a walk in the park.

Living with a stoma doesn't have to hold anyone back from where they want to go and we're here to help them get there.

Contact us.....	3
A guide to adhesives.....	4
A Guide to base plates.....	5
Common pouch features.....	6
Urostomy pouches .....	7

# *Support and advice EVERY STEP OF THE WAY*



For more information about ConvaTec products, Ostomate resources and ostomy management, go online or contact our ConvaTec Customer Care Centre.

## Australia

Freecall **1800 335 276**

or email [connection.au@convatec.com](mailto:connection.au@convatec.com)

**[www.convatec.com.au](http://www.convatec.com.au)**

## New Zealand

Freecall **0800 441 763**

or email [connection.nz@convatec.com](mailto:connection.nz@convatec.com)

**[www.convatec.co.nz](http://www.convatec.co.nz)**

# A guide to ADHESIVES



## STOMAHESIVE®

Stomahesive® technology is a hydrocolloid adhesive that gently adheres to the skin, protecting it from the damaging effects of stoma effluent. Stomahesive® technology combines adhesive and hydrocolloid components for optimal performance under various conditions. Stomahesive® technology is a unique product, only from ConvaTec.

Stomahesive® base plates are ideal for people whose stoma output is semi-formed to formed stool. This reliable base plate provides a special formula that enables the base plate to adhere to both dry and moist skin. It holds firmly to skin while helping to protect against irritating output. Another feature of the Stomahesive® base plates is an outer protective coating which helps keep it in place. Plus, Stomahesive® base plates are easy to fit, apply, and remove.

### Modified Stomahesive®

To further assist in the care of peristomal skin, ConvaTec have developed a modified Stomahesive® that allows a gentle adherence to the skin. Modified Stomahesive® base plates are more flexible and are easier to remove.

## DURAHESIVE®

Durahesive® skin base plates are designed for people whose stoma output is mostly liquid. Unlike some other skin base plates that can break down around liquid output, Durahesive® skin base plates swell up to protect the stoma.<sup>1</sup> This special effect, called turtle-necking, helps create a more secure seal without harming the skin. Durahesive® skin base plates are also easy to fit, apply, and remove. They will adhere securely to the skin while showering, bathing and swimming.

## TAPE COLLARS

*are available in two forms:*

**Acrylic:** Tape collar designed with good tack

**Hydrocolloid:** Designed for those with sensitivity to other tape collars

### Designed for patients

- with semi-formed/formed output
- who have sensitive intact skin
- with frequent changes
- who prefer a standard wear time

### Designed for patients

- with semi-liquid/liquid output, including urine
- who want or need longer wear time

# A guide to BASE PLATES



## CUT TO FIT

Cut-to-Fit base plates can be cut to match the shape and size of individual stomas. Using curved scissors, one can cut a circle, oval or any shape needed to fit the stoma.




### Designed for patients

- with a stoma size in-between pre-cut sizes
- with an irregular stoma shape

## PRE-CUT

Pre-cut base plates have a uniform, circular stoma hole and are ready to apply right out of the box. They are available for the most common stoma sizes.



### Designed for patients

- with an evenly round stoma that fits the available pre-cut size
- with dexterity issues who find cutting challenging

## MOULDABLE

ConvaTec Moldable Technology™ skin barrier with Rebounding Memory Technology™ creates an elastic-like seal that rebounds to fit most shapes and sizes, and then “turtlenecks” upward, to enhance the seal and protect the skin.<sup>2</sup>



### Designed for patients

- who have large, round, oval or irregular shaped stomas
- with dexterity issues who find cutting challenging



### Look for this icon

Please note: Do not cut the base plates

# Common pouch FEATURES



## DRAINABLE POUCH CLOSURE



### **InvisiClose® Clipless Tail Closure**

Self closing drainable pouch closure which provides an easy closure option that does not require a separate clip for added security.



### **Lock-it Pocket®**

Folds inside out to hide away the secured drainage tail—providing the flexibility of a drainable pouch, with the feel and comfort of a closed-end pouch.

*Traditional drainable pouches are closed with a separate clip.*

## POUCH COLOUR

### **Transparent (Clear front of pouch)**

Clear pouch front for easy output assessment and determining when a pouch needs to be changed or emptied. Soft and quiet material on the body side of the pouch.

### **Opaque (Tan-coloured pouch)**

Soft and quiet material cover both the front and the body side of the pouch for more discretion and comfort.



### **Easy-View Window**

Some options come with an easy-view window for accurate positioning and observation.

## FILTER OR NO FILTER



Filters help release and deodorize gas and keep the closed or drainable pouch flat.

Note that the filter must be covered with a filter cover while swimming, bathing or showering.



# Urostomy POUCHES

## UROSTOMY POUCH DRAINAGE

The ConvaTec Accuseal™ Tap quickly and easily drains liquid output. Its swivel mechanism glides easily into open or closed positions whilst offering extra security with the red marker that indicates the tap is in the open position. A stopper adds extra security. Accuseal™ is designed for adaptors on night drainage bags.

*Other drainage options include a Soft Tap or Fold Over Tap.*

*All Urostomy pouches include an anti-reflux valve which is designed to help prevent backflow.*



## URINE DRAINAGE BAG CONNECTORS

Each box of urostomy pouches contains two plastic adapters for attaching a urine/overnight drainage bag. Additional plastic adapters are available upon request for in-hospital use.

Telephone ConvaTec Customer Care for information or to order.

Freecall 1800 335 276 (Australia)

0800 441 763 (New Zealand)




# One Piece

<b>Esteem®+ Convexity Variant Guide</b> .....	9
<b>Esteem®+ Soft Convex</b> .....	10
Closed Pouches .....	11
Drainable Pouches.....	12
Urostomy Pouches.....	13
<b>Esteem®+ Flex Convex</b> .....	14
Closed Pouches.....	14
Drainable Pouches.....	15
Urostomy Pouches.....	15
<b>Esteem®+ Firm Convexity</b> .....	
Drainable Pouches.....	16
<b>ActiveLife® Firm Convexity</b> .....	17
Drainable Pouches.....	17
Urostomy Pouches.....	17
<b>Esteem®+ with Convatec Mouldable Technology™</b> .....	
Flat Closed Pouches.....	18
Flat Drainable Pouches.....	18
<b>Esteem®+</b> .....	
Flat Closed Pouches .....	19-20
Flat Drainable Pouches.....	21-22
<b>Esteem®</b> .....	
Flat Drainable Pouches.....	23
Flat Urostomy Pouches.....	23
<b>ActiveLife®</b> .....	
Flat Drainable Pouches.....	24
Flat Urostomy Pouches.....	24
<b>ActiveLife® Little Ones™</b> .....	25
Flat Drainable Pouches.....	25
Flat Urostomy Pouches.....	25



# CONVEXITY VARIANT GUIDE

The Esteem® portfolio offers a range of one-piece convex products with various depths to help provide comfort, freedom, security and reassurance to your patient, all designed for different stoma shapes and sizes.



Full assessment of the stoma and abdomen should be made when selecting convexity variant\*\*

## STOMA PROTRUSION

**Protruding:** use widest option available  
*Pressure distributed over a wider area around the stoma*

**Flush:** use narrowest option available  
*Assess the abdomen type and where pressure is required*

**Retracted:** use narrowest option available  
*Gives pressure closest to the stoma to help protrusion*

## ABDOMINAL TONE AND CONTOUR

**Skin folds/creases:** use widest option available  
*Help to move any skin folds out of the way*

**Rigid/rounded:** use widest option available  
*Pressure distributed over a wider area and can aid adhesion*

**Flaccid - Average tone:** use option most suitable to stoma protrusion

## PERISTOMAL SKIN

**Swollen/tender:** use widest option available to achieve desired protrusion  
*Pressure distributed over a wider area and can aid adhesion*

**Skin types:** use option most suitable to stoma protrusion and abdominal tone/contour

**Normal/healthy skin:** use option most suitable to stoma protrusion and abdominal tone/contour

\*Based on Convexity Plateau Variants Guidelines / Considerations selection guide AP-018355-MM.


\*\*Patients should be assessed by their Health Care Professional/Stoma Nurse for appropriate product choice and use of convexity.

# Introducing ESTEEM® Soft Convex

ONE PIECE



## ALL BODIES ARE DIFFERENT.

Esteem®  Soft Convex combines flexibility with a shallow convexity, and is designed to conform well to the body and be comfortable to wear.

- **Flexible Convexity:** conforms to the body and helps improve the pouch seal by minimising leaks
- **Easy-View Window:** for accurate positioning and observation
- **Filter:** designed to help minimise odours and ballooning\*
- **Safe Seal Clipless Closure:** designed to help make emptying and cleaning easy



**Designed  
to minimise  
pressure**



**Designed  
to be skin  
friendly**



**Flexible,  
conforms to  
body**



**Convex**

*\*Drainable options only.*

Closed Pouches with Soft Convexity



Hydrocolloid base plate with tapered edge, opaque, belt tabs, viewing window, with filter, (2 sided fabric)



Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
422360	20-47mm	V1 - 50mm	Cut to fit	30	80211P	60
422361	15-40mm	V2 - 40mm	Cut to fit	30	80211P	60
422362	10-28mm	V3 - 30mm	Cut to fit	30	80211P	60

						Pre-cut
422527	25mm	V1 - 50mm	Pre-cut	30	80211P	60
422528	30mm	V1 - 50mm	Pre-cut	30	80211P	60
422529	35mm	V1 - 50mm	Pre-cut	30	80211P	60
422530	40mm	V1 - 50mm	Pre-cut	30	80211P	60
422531	45mm	V1 - 50mm	Pre-cut	30	80211P	60
422532	20mm	V2 - 40mm	Pre-cut	30	80211P	60
422533	25mm	V2 - 40mm	Pre-cut	30	80211P	60
422534	30mm	V2 - 40mm	Pre-cut	30	80211P	60
422535	35mm	V2 - 40mm	Pre-cut	30	80211P	60
422536	15mm	V3 - 30mm	Pre-cut	30	80211P	60
422537	20mm	V3 - 30mm	Pre-cut	30	80211P	60
422538	25mm	V3 - 30mm	Pre-cut	30	80211P	60

Drainable Pouches with Soft Convexity



Hydrocolloid base plate with tapered edge, opaque, viewing window, with filter, integrated closure, belt tabs, (2 sided fabric)



Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
422363	20-47mm	V1 - 50mm	Cut to fit	10	80212Q	30
422364	15-40mm	V2 - 40mm	Cut to fit	10	80212Q	30
422365	10-28mm	V3 - 30mm	Cut to fit	10	80212Q	30

						Pre-cut
422539	25mm	V1 - 50mm	Pre-cut	10	80212Q	30
422540	30mm	V1 - 50mm	Pre-cut	10	80212Q	30
422541	35mm	V1 - 50mm	Pre-cut	10	80212Q	30
422542	40mm	V1 - 50mm	Pre-cut	10	80212Q	30
422543	45mm	V1 - 50mm	Pre-cut	10	80212Q	30
422544	20mm	V2 - 40mm	Pre-cut	10	80212Q	30
422545	25mm	V2 - 40mm	Pre-cut	10	80212Q	30
422546	30mm	V2 - 40mm	Pre-cut	10	80212Q	30
422547	35mm	V2 - 40mm	Pre-cut	10	80212Q	30
422548	15mm	V3 - 30mm	Pre-cut	10	80212Q	30
422549	20mm	V3 - 30mm	Pre-cut	10	80212Q	30
422550	25mm	V3 - 30mm	Pre-cut	10	80212Q	30

Urostomy Pouches with Soft Convexity




Hydrocolloid base plate with tapered edge, opaque, belt tabs, viewing window, anti-reflux secure valve, soft plug, (2 sided fabric)



Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
422551	20-47mm	V1 - 50mm	Cut to fit	10	80213R	30
422557	15-40mm	V2 - 40mm	Cut to fit	10	80213R	30
422562	10-28mm	V3 - 30mm	Cut to fit	10	80213R	30
						Pre-cut
422552	25mm	V1 - 50mm	Pre-cut	10	80213R	30
422553	30mm	V1 - 50mm	Pre-cut	10	80213R	30
422554	35mm	V1 - 50mm	Pre-cut	10	80213R	30
422555	40mm	V1 - 50mm	Pre-cut	10	80213R	30
422556	45mm	V1 - 50mm	Pre-cut	10	80213R	30
422558	20mm	V2 - 40mm	Pre-cut	10	80213R	30
422559	25mm	V2 - 40mm	Pre-cut	10	80213R	30
422560	30mm	V2 - 40mm	Pre-cut	10	80213R	30
422561	35mm	V2 - 40mm	Pre-cut	10	80213R	30
422563	15mm	V3 - 30mm	Pre-cut	10	80213R	30
422564	20mm	V3 - 30mm	Pre-cut	10	80213R	30
422565	25mm	V3 - 30mm	Pre-cut	10	80213R	30

# ESTEEM® Flex Convex



**Esteem®  Flex Convex system combines the comfort and freedom of flexibility with the firmness of convexity.**

Flex Convex adheres into creases, skin folds, and onto flaccid abdomens<sup>3</sup>. Its increased flexibility provides more comfortable wear compared to more rigid forms of convexity<sup>4</sup>. As Flex Convex creates less pressure on the peristomal skin<sup>6</sup> it is less likely to cause skin damage<sup>5</sup>.

## Closed Pouches with Flexible Convexity



Hydrocolloid base plate with tapered edge, opaque, belt tabs, viewing window, with filter, (2 sided fabric)



Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
422158	20-57mm	VO - 60mm	Cut to fit	30	80106D	60
421612	20-43mm	V1 - 50mm	Cut to fit	30	80106D	60
421613	20-35mm	V2 - 40mm	Cut to fit	30	80106D	60
421614	20-25mm	V3 - 30mm	Cut to fit	30	80106D	60



## Drainable Pouches with Flexible Convexity



Hydrocolloid base plate with tapered edge, opaque, viewing window, with filter, integrated closure, belt tabs, (2 sided fabric)



Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
422159	20-57mm	V0 - 60mm	Cut to fit	10	80113L	30
421615	20-43mm	V1 - 50mm	Cut to fit	10	80113L	30
421619	20-35mm	V2 - 40mm	Cut to fit	10	80113L	30
421622	20-25mm	V3 - 30mm	Cut to fit	10	80113L	30

Pre-cut

421616	25mm	V1 - 50mm	Pre-cut	10	80113L	30
421617	30mm	V1 - 50mm	Pre-cut	10	80113L	30
421618	35mm	V1 - 50mm	Pre-cut	10	80113L	30
421620	25mm	V2 - 40mm	Pre-cut	10	80113L	30
421621	30mm	V2 - 40mm	Pre-cut	10	80113L	30

## Urostomy Pouches with Flexible Convexity



Flex Convex hydrocolloid base plate with tapered edge, opaque, viewing window, anti-reflux secure valve, soft plug, (2 sided fabric)



Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
421623	10-43mm	V1 - 50mm	Cut to fit	10	80118R	30
421626	10-35mm	V2 - 40mm	Cut to fit	10	80118R	30
421628	10-25mm	V3 - 30mm	Cut to fit	10	80118R	30

Flex Convex hydrocolloid base plate with tapered edge, opaque, viewing window, with filter, anti-reflux secure valve, soft plug, (2 sided fabric)

Pre-cut

Product Code	Stoma Size	Plateau Size	Description	Pack Size	SAS Code	Monthly Allowance
421624	15mm	V1 - 50mm	Pre-cut	10	80118R	30
421625	20mm	V1 - 50mm	Pre-cut	10	80118R	30
421627	15mm	V2 - 40mm	Pre-cut	10	80118R	30
421629	15mm	V3 - 30mm	Pre-cut	10	80118R	30

Drainable Pouches with Firm Convexity and Lock-it Pocket™



Durahesive® base plate, standard, opaque with filter  
hydrocolloid collar, taped border and belt tabs (2 sided fabric) **Pre-cut**

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
416739	19mm	Pre-cut	10	5699Y	30
416741	22mm	Pre-cut	10	5699Y	30
416743	25mm	Pre-cut	10	5699Y	30
416745	28mm	Pre-cut	10	5699Y	30
416747	32mm	Pre-cut	10	5699Y	30
416749	35mm	Pre-cut	10	5699Y	30
416751	38mm	Pre-cut	10	5699Y	30

Durahesive® base plate, standard, transparent, with filter  
hydrocolloid collar, taped border and belt tabs (1 sided fabric) **Pre-cut**

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
416738	19mm	Pre-cut	10	5699Y	30
416740	22mm	Pre-cut	10	5699Y	30
416742	25mm	Pre-cut	10	5699Y	30
416744	28mm	Pre-cut	10	5699Y	30
416746	32mm	Pre-cut	10	5699Y	30
416748	35mm	Pre-cut	10	5699Y	30
416750	38mm	Pre-cut	10	5699Y	30

ONE PIECE

# ACTIVELife®



A lightweight, flexible, one-piece system, **ActiveLife®** combines the base plate and pouch in one simple, easy-to-manage system.

Pouches are available in both transparent and opaque.

## Drainable Pouches with Firm Convexity and Clip Closure



Durahesive® base plate, standard, transparent, *without filter belt tabs*, (1 sided fabric)

Pre-cut

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
175777	19mm	Pre-cut	5	3993E	30
175778	22mm	Pre-cut	5	3993E	30
175779	25mm	Pre-cut	5	3993E	30
175780	28mm	Pre-cut	5	3993E	30
175781	32mm	Pre-cut	5	3993E	30

## Urostomy Pouches with Firm Convexity



Durahesive® base plate, standard, transparent, *with bendable tap and belt tabs*

Pre-cut

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
175792	19mm	Pre-cut	5	3628Y	30
175793	22mm	Pre-cut	5	3628Y	30
175794	25mm	Pre-cut	5	3628Y	30
175795	28mm	Pre-cut	5	3628Y	30
175796	32mm	Pre-cut	5	3628Y	30
175797	35mm	Pre-cut	5	3628Y	30

Flat Closed Pouches with ConvaTec Mouldable Technology™



Mouldable base plate, standard, opaque, viewing window, with filter (2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
413509	20-30mm	Mouldable	30	80003Q	90
413510	30-40mm	Mouldable	30	80003Q	90
413511	40-50mm	Mouldable	30	80003Q	90

Flat Drainable Pouches with ConvaTec Mouldable Technology™



Mouldable base plate, standard, opaque, viewing window, with filter (2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
413521	20-30mm	Mouldable	10	80004R	30
413522	30-40mm	Mouldable	10	80004R	30
413523	40-50mm	Mouldable	10	80004R	30

Mouldable base plate, standard, transparent with filter (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
413515	20-30mm	Mouldable	10	80004R	30
413516	30-40mm	Mouldable	10	80004R	30
413517	40-50mm	Mouldable	10	80004R	30

Flat Closed Pouches



Full Stomahesive® base plate, small, opaque, viewing window, with filter (2 sided fabric)


Pre-cut

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421816	20-70mm	Cut to Fit	30	5677T	90
421817	25mm	Pre-cut	30	5677T	90
421818	30mm	Pre-cut	30	5677T	90
421819	35mm	Pre-cut	30	5677T	90
421820	40mm	Pre-cut	30	5677T	90



Full Stomahesive® base plate, standard, opaque, viewing window, with filter (2 sided fabric)


Pre-cut

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421686	20-70mm	Cut to Fit	30	5677T	90
421687	25mm	Pre-cut	30	5677T	90
421688	30mm	Pre-cut	30	5677T	90
421689	35mm	Pre-cut	30	5677T	90
421822	40mm	Pre-cut	30	5677T	90
421823	50mm	Pre-cut	30	5677T	90

Modified Stomahesive® base plate, standard, transparent, with filter (1 sided fabric)

Pre-cut
---------

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
416712	40mm	Pre-cut	30	5677T	90

Flat Closed Pouches

Full Stomahesive® base plate, large, opaque, viewing window, with filter (2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421825	20-70mm	Cut to Fit	30	5677T	90
					Pre-cut
421826	25mm	Pre-cut	30	5677T	90
421827	30mm	Pre-cut	30	5677T	90
421828	35mm	Pre-cut	30	5677T	90
421829	40mm	Pre-cut	30	5677T	90
421830	50mm	Pre-cut	30	5677T	90



Full Stomahesive® base plate, large, transparent, with filter (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421831	20-70mm	Cut to Fit	30	5677T	90
					Pre-cut
421832	25mm	Pre-cut	30	5677T	90
421833	30mm	Pre-cut	30	5677T	90
421834	35mm	Pre-cut	30	5677T	90
421835	40mm	Pre-cut	30	5677T	90
421836	50mm	Pre-cut	30	5677T	90



## Flat Drainable Pouches with InvisiClose® and Lock-it Pocket™

Modified Stomahesive® base plate, small, opaque with filter, (2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421849	20-70mm	Cut to Fit	10	80185G	30
					Pre-cut
421850	25mm	Pre-cut	10	80185G	30

Modified Stomahesive® base plate, standard, opaque with filter, viewing window, (2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421858	20-70mm	Cut to Fit	10	80185G	30
					Pre-cut
421859	25mm	Pre-cut	10	80185G	30
421860	30mm	Pre-cut	10	80185G	30
421861	35mm	Pre-cut	10	80185G	30
421862	40mm	Pre-cut	10	80185G	30

Modified Stomahesive® base plate, standard, transparent with filter, (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
416718	20-70mm	Cut to Fit	10	5678W	30
					Pre-cut
416724	25mm	Pre-cut	10	5678W	30
416728	30mm	Pre-cut	10	5678W	30
416732	35mm	Pre-cut	10	5678W	30
416736	40mm	Pre-cut	10	5678W	30



Flat Drainable Pouches with Clip Closure

Full Stomahesive® base plate, small, opaque, with filter  
(2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421837	20-70mm	Cut to Fit	10	80187J	30
					Pre-cut
421843	25mm	Pre-cut	10	80187J	30
421844	30mm	Pre-cut	10	80187J	30
421846	40mm	Pre-cut	10	80187J	30

Full Stomahesive® base plate, standard, opaque, with filter  
(2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421838	20-70mm	Cut to Fit	10	80187J	30
					Pre-cut
421839	25mm	Pre-cut	10	80187J	30
421840	30mm	Pre-cut	10	80187J	30

Full Stomahesive® base plate, standard, transparent  
with filter (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421848	20-70mm	Cut to Fit	10	80187J	30

Flat Drainable Pouches with InvisiClose®



Modified Stomahesive® base plate, standard, transparent *without filter*, Lock-it Pocket™ (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
416721	20-70mm	Cut to Fit	10	N/A	N/A

Modified Stomahesive® base plate, large, transparent *without filter* (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
416908	25-100mm	Cut to Fit	10	5678W	30

Flat Urostomy Pouches



Durahesive® base plate with Accuseal™ tap, transparent (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421631	13-45mm	Cut to fit	10	80128G	40

## Flat Drainable Pouches with Clip Closure



Stomahesive® base plate, standard, opaque  
without filter flexible white collar\* (2 sided fabric)

Pre-cut

Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
22757	19mm	Pre-cut	10	3699Q	60
22758	25mm	Pre-cut	10	3699Q	60
22759	32mm	Pre-cut	10	3699Q	60
22760	38mm	Pre-cut	10	3699Q	60
22761	45mm	Pre-cut	10	3699Q	60
22762	50mm	Pre-cut	10	3699Q	60

Stomahesive® base plate, standard, transparent  
without filter flexible white collar\* (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
22771	19-64mm	Cut to Fit	10	3688D	60

\*The outer white collar is sprayed with Stomahesive® technology to help keep the product skin friendly.

Stomadress® – Stomahesive® base plate, large  
transparent, without filter (2 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
650803	8-100mm	Cut to Fit	10	3636J	30

## Flat Urostomy Pouches



Stomahesive® base plate, standard, transparent with  
bendable tap and belt tabs



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
64927	19-45mm	Cut to Fit	10	3661Q	40
Pre-cut					
64928	25mm	Pre-cut	10	3661Q	40
64929	32mm	Pre-cut	10	3661Q	40
64930	38mm	Pre-cut	10	3661Q	40
64931	45mm	Pre-cut	10	3661Q	40

# ACTIVELife® LITTLE ONES™

Designed just for kids who have had colostomy, ileostomy, or urostomy surgery, **ActiveLife® Little Ones™** Ostomy System has been scaled down in size to comfortably fit a child's unique body contour. Little Ones® products have been created with an active child in mind.



## Paediatric Pouches



Stomahesive® base plate, drainable, transparent *without filter* (1 sided fabric)



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
20922	8-50mm	Drainable	15	3667B	60

## Paediatric Urostomy Pouches



Stomahesive® base plate, urostomy, transparent with bendable tap



Product Code	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
20917	8-25mm	Bendable Tap	15	3660P	60



# Two Piece


## Mechanical Coupling

<b>NATURA®</b> .....	27
Natura® Accordion Flange .....	28
Base Plates with Firm Convexity .....	29-30
Flat Base Plates .....	31-32
<b>NATURA® +</b> .....	
Closed Pouches .....	33
Drainable Pouches .....	34
Urostomy Pouches .....	35
<b>SUR-FIT® Plus</b> .....	36
Convex Base Plates .....	37
Flat Base Plates .....	37-38
Closed Pouches .....	39
Drainable Pouches .....	40
High Output Pouches .....	41
Urostomy Pouches .....	42
<b>SUR-FIT® AutoLock®</b> .....	43
Flat Base Plates .....	43
Closed Pouches .....	43
Drainable Pouches .....	44
Urostomy Pouches .....	44



# NATURA®

TWO PIECE  
Mechanical Coupling








ConvaTec developed these innovative pouches to be noticeably more comfortable and more discreet.

The Natura®<sup>+</sup> pouching system combined with ConvaTec advanced base plate adhesives offer security, comfort and discretion.

## Natura® Colour Key & Sizing Guide

Match the colour of the dot on the base plate box with the colour of the dot on the pouch box.

				
32 mm	38 mm	45 mm	57 mm	70 mm

# NATURA<sup>®</sup>

## Accordion Flange




### FEATURES

- The accordion flange lifts to provide generous finger room so coupling is easy and comfortable.
- Accordion Flange<sup>®</sup> has a low profile in the 'DOWN' position so it's not visible under clothing. The base plate moves as the body flexes
- Available in flat cut to fit, Convex cut to fit and flat ConvaTec Mouldable Technology<sup>™</sup>
- Connects with an audible 'click' and physical 'snap' to confirm the pouch is firmly attached

### Designed for patients

- with dexterity issues as it provides finger room so coupling is easy and comfortable
- to prevent pressure or trauma if abdominal tenderness is present due to post-operative stoma surgery/incisions

### HOW IT WORKS



## Base Plates with Firm Convexity and Accordion Flange®



Durahesive® base plate with flexible tape collar



Product Code		Flange Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421632	●	45mm	13-22mm	Cut to Fit	5	80129H	20
421634	●	57mm	13-35mm	Cut to Fit	5	80129H	20
421636	●	70mm	13-48mm	Cut to Fit	5	80129H	20



Durahesive® base plate with hydrocolloid collar



Product Code		Flange Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421633	●	45mm	13-22mm	Cut to Fit	5	80129H	20
421635	●	57mm	13-35mm	Cut to Fit	5	80129H	20
421637	●	70mm	13-48mm	Cut to Fit	5	80129H	20

## Base Plates with Firm Convexity and ConvaTec Mouldable Technology™



Durahesive® convex base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
404592	●	45mm	13-22mm	Mouldable	10	5682C	20
404593	●	45mm	22-33mm	Mouldable	10	5682C	20
404594	●	57mm	33-45mm	Mouldable	10	5682C	20

Base Plates with Firm Convexity



Durahesive® base plate with flexible white collar\*

Pre-cut

Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
413177	●	45mm	13mm	Pre-cut	10	5681B	20
413178	●	45mm	16mm	Pre-cut	10	5681B	20
413179	●	45mm	19mm	Pre-cut	10	5681B	20
413180	●	45mm	22mm	Pre-cut	10	5681B	20
413181	●	45mm	25mm	Pre-cut	10	5681B	20
413182	●	45mm	28mm	Pre-cut	10	5681B	20
413183	●	45mm	32mm	Pre-cut	10	5681B	20
413184	●	45mm	35mm	Pre-cut	10	5681B	20
413185	●	57mm	38mm	Pre-cut	10	5681B	20
413186	●	57mm	41mm	Pre-cut	10	5681B	20
413187	●	57mm	45mm	Pre-cut	10	5681B	20
413188	●	57mm	50mm	Pre-cut	10	5681B	20

\*The outer white collar is sprayed with Stomahesive® technology to help keep the product skin friendly.



## Flat Base Plates with Accordion Flange® and Convatec Mouldable Technology™

Stomahesive® Accordion Flange® base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421033	●	57mm	13-22mm	Mouldable	10	80025W	30
421034	●	57mm	22-33mm	Mouldable	10	80025W	30
421035	●	70mm	33-45mm	Mouldable	10	80025W	30

Duraheasive® Accordion Flange® base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421039	●	57mm	13-22mm	Mouldable	10	80025W	30
421040	●	57mm	22-33mm	Mouldable	10	80025W	30
421041	●	70mm	33-45mm	Mouldable	10	80025W	30

## Flat Base Plates with Accordion Flange®

Stomahesive® Accordion Flange® base plate with flexible white collar\*



Product Code		Flange Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421456	●	45mm	13-21mm	Cut to Fit	10	80057M	30
421460	●	57mm	13-33mm	Cut to Fit	10	80057M	30
421464	●	70mm	13-45mm	Cut to Fit	10	80057M	30

Duraheasive® Accordion Flange® base plate with flexible white collar\*



Product Code		Flange Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
421454	●	45mm	13-21mm	Cut to Fit	10	80057M	30
421458	●	57mm	13-33mm	Cut to Fit	10	80057M	30
421462	●	70mm	13-45mm	Cut to Fit	10	80057M	30

\*The outer white collar is sprayed with Stomahesive® technology to help keep the product skin friendly.



## Flat Base Plates with Convatec Mouldable Technology™



Stomahesive® base plate with hydrocolloid collar



Product Code		Base plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
411801	●	45mm	13-22mm	Mouldable	10	5680Y	20
411803	●	45mm	22-33mm	Mouldable	10	5680Y	20
411805	●	57mm	33-45mm	Mouldable	10	5680Y	20
411807	●	70mm	45-56mm	Mouldable	10	5680Y	20

## Flat Base Plates



Stomahesive® base plate with flexible hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
410293	●	32mm	13-19mm	Cut to Fit	10	5679X	30
410294	●	38mm	13-25mm	Cut to Fit	10	5679X	30
410295	●	45mm	13-32mm	Cut to Fit	10	5679X	30
410296	●	57mm	13-45mm	Cut to Fit	10	5679X	30
410297	●	70mm	13-57mm	Cut to Fit	10	5679X	30



Stomahesive® base plate with flexible tan collar



Product Code		Base plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
125262	●	32mm	13-19mm	Cut to fit	10	5679X	30
125264	●	45mm	13-32mm	Cut to fit	10	5679X	30
125265	●	57mm	13-45mm	Cut to fit	10	5679X	30
125266	●	70mm	13-57mm	Cut to fit	10	5679X	30



Full Stomahesive® base plate



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
401574	●	38mm	13-25mm	Cut to Fit	10	5679X	30
401575	●	45mm	13-32mm	Cut to Fit	10	5679X	30
401576	●	57mm	13-45mm	Cut to Fit	10	5679X	30
401577	●	70mm	13-57mm	Cut to Fit	10	5679X	30



## Closed Pouches



Small, opaque, viewing window, with filter (2 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
421794	●	32mm	Small with Filter	30	5683D	90
421795	●	38mm	Small with Filter	30	5683D	90
421796	●	45mm	Small with Filter	30	5683D	90
421797	●	57mm	Small with Filter	30	5683D	90



Standard, opaque, viewing window, with filter (1 sided fabric)

Product Code		Base plate Size	Description	Pack Size	SAS Code	Monthly Allowance
421676	●	32mm	Standard with Filter	30	5683D	90
421678	●	38mm	Standard with Filter	30	5683D	90
421680	●	45mm	Standard with Filter	30	5683D	90
421682	●	57mm	Standard with Filter	30	5683D	90
421798	●	70mm	Standard with Filter	30	5683D	90



Large, opaque, with filter, viewing window, with filter (2 sided fabric)

Product Code		Base plate Size	Description	Pack Size	SAS Code	Monthly Allowance
421799	●	38mm	Large with Filter	30	5683D	90
421800	●	45mm	Large with Filter	30	5683D	90
421801	●	57mm	Large with Filter	30	5683D	90
421802	●	70mm	Large with Filter	30	5683D	90

Standard, transparent, with filter (1 sided fabric)

Product Code		Base plate Size	Description	Pack Size	SAS Code	Monthly Allowance
421803	●	45mm	Standard with Filter	30	5683D	90
421804	●	57mm	Standard with Filter	30	5683D	90
421805	●	70mm	Standard with Filter	30	5683D	90

## Drainable Pouches with InvisiClose<sup>®</sup>

Standard, opaque, viewing window, with filter (2 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
421746	●	38mm	Standard with Filter	10	5685F	60
421747	●	45mm	Standard with Filter	10	5685F	60
421748	●	57mm	Standard with Filter	10	5685F	60
421749	●	70mm	Standard with Filter	10	5685F	60



Standard, transparent, with filter (1 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416414	●	38mm	Standard with Filter	10	5685F	60
416416	●	45mm	Standard with Filter	10	5685F	60
416419	●	57mm	Standard with Filter	10	5685F	60
416422	●	70mm	Standard with Filter	10	5685F	60

Urostomy Pouches



Small, transparent, with flexible soft tap (1 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
421074	●	32mm	Small with flexible soft tap	10	80024T	60
421075	●	38mm	Small with flexible soft tap	10	80024T	60
421076	●	45mm	Small with flexible soft tap	10	80024T	60

Standard, transparent, with flexible soft tap (1 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
413437	●	45mm	Standard with flexible soft tap	10	80024T	60
413438	●	57mm	Standard with flexible soft tap	10	80024T	60
413439	●	70mm	Standard with flexible soft tap	10	80024T	60

# SUR-FIT® Plus




TWO PIECE  
Mechanical Coupling



**SUR-FIT® Plus** Ostomy System is a comprehensive two piece product line with a plastic coupling ring allowing the pouch and base plate to snap together for audible security.

## SUR-FIT® Plus Colour Key & Sizing Guide

Match the colour of the dot on the base plate box with the colour of the dot on the pouch box.

				
32 mm	38 mm	45 mm	57 mm	70 mm

## Convex Base Plates with Firm Convexity and ConvaTec Mouldable Technology™



Durahesive® base plate with hydrocolloid collar



Product Code		Base plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
411449	●	45mm	13-22mm	Mouldable	10	9986N	20
411450	●	45mm	22-33mm	Mouldable	10	9986N	20
411451	●	57mm	33-45mm	Mouldable	10	9986N	20

Durahesive® base plate with flexible tape collar

Pre-cut

Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
125028	●	45mm	19mm	Pre-cut	5	3630C	20
125029	●	45mm	22mm	Pre-cut	5	3630C	20
125030	●	45mm	25mm	Pre-cut	5	3630C	20
125031	●	45mm	28mm	Pre-cut	5	3630C	20
125032	●	45mm	32mm	Pre-cut	5	3630C	20
125033	●	45mm	35mm	Pre-cut	5	3630C	20
125034	●	57mm	38mm	Pre-cut	5	3630C	20
125035	●	57mm	41mm	Pre-cut	5	3630C	20
125036	●	57mm	45mm	Pre-cut	5	3630C	20
125037	●	57mm	50mm	Pre-cut	5	3630C	20

## Flat Base Plates with ConvaTec Mouldable Technology™



Stomahesive® base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
411817	●	45mm	13-22mm	Mouldable	10	9877W	20
411819	●	45mm	22-33mm	Mouldable	10	9877W	20
411821	●	57mm	33-45mm	Mouldable	10	9877W	20
411823	●	70mm	45-56mm	Oval, Mouldable	10	9877W	20



## Flat Base Plates



Stomahesive® base plate with flexible white collar\*



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
125132	●	32mm	12-26mm	Cut to Fit	5	3752L	30
125133	●	38mm	12-32mm	Cut to Fit	5	3752L	30
125134	●	45mm	12-39mm	Cut to Fit	5	3752L	30
125135	●	57mm	12-50mm	Cut to Fit	5	3752L	30
125136	●	70mm	12-64mm	Cut to Fit	5	3752L	30

\*The outer white collar is sprayed with Stomahesive® technology to help keep the product skin friendly.



Full Stomahesive® base plate



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
125142	●	32mm	12-26mm	Cut to Fit	5	3763C	30
125143	●	38mm	12-32mm	Cut to Fit	5	3763C	30
125144	●	45mm	12-39mm	Cut to Fit	5	3763C	30
125145	●	57mm	12-50mm	Cut to Fit	5	3763C	30
400949	●	70mm	12-64mm	Cut to Fit	5	3763C	30



Flexible Stomahesive® base plate



Product Code		Base plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
125137	●	32mm	12-26mm	Cut to Fit	5	3763C	30
125138	●	38mm	12-32mm	Cut to Fit	5	3763C	30
125139	●	45mm	12-39mm	Cut to Fit	5	3763C	30
125140	●	57mm	12-50mm	Cut to Fit	5	3763C	30
400965	●	70mm	12-64mm	Cut to Fit	5	3763C	30



Full Durahesive® base plate



Product Code		Base plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
125022	●	32mm	12-26mm	Cut to Fit	5	3680Q	20
125023	●	38mm	12-32mm	Cut to Fit	5	3680Q	20
125024	●	45mm	12-39mm	Cut to Fit	5	3680Q	20
125025	●	57mm	12-50mm	Cut to Fit	5	3680Q	20



## Closed Pouches



Mini, opaque, *without filter*

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402510	●	38mm	Without Filter	20	3753M	90
402511	●	45mm	Without Filter	20	3753M	90



Standard, opaque, *without filter*

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402516	●	38mm	Without Filter	30	3758T	90
402517	●	45mm	Without Filter	30	3758T	90

Standard, opaque, *with filter*

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402521	●	32mm	With Filter	30	3649C	90
402522	●	38mm	With Filter	30	3649C	90
402523	●	45mm	With Filter	30	3649C	90
402524	●	57mm	With Filter	30	3649C	90
402525	●	70mm	With Filter	30	3649C	90

## Drainable Pouches with Integrated Closure

Small, opaque, with filter and InvisiClose®

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
411401	●	45mm	With Filter	10	9987P	60
411402	●	57mm	With Filter	10	9987P	60
411403	●	70mm	With Filter	10	9987P	60

Standard, opaque, with filter and InvisiClose®

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
411385	●	38mm	With Filter	10	9987P	60
411386	●	45mm	With Filter	10	9987P	60
411387	●	57mm	With Filter	10	9987P	60
411388	●	70mm	With Filter	10	9987P	60

## Drainable Pouches with Clip Closure

Small, opaque, *without filter*

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402526	●	32mm	Without Filter	10	3764D	60
402527	●	38mm	Without Filter	10	3764D	60
402528	●	45mm	Without Filter	10	3764D	60
402529	●	57mm	Without Filter	10	3764D	60

Standard, opaque, *without filter*

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402531	●	32mm	Without Filter	10	3759W	60
402532	●	38mm	Without Filter	10	3759W	60
402533	●	45mm	Without Filter	10	3759W	60
402534	●	57mm	Without Filter	10	3759W	60
402535	●	70mm	Without Filter	10	3759W	60





High Output Pouches

High output, opaque, with filter and cap

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
125303	●	45mm	With Filter	5	355OW	10
125304	●	57mm	With Filter	5	355OW	10
125305	●	70mm	With Filter	5	355OW	10



## Urostomy Pouches



Standard, opaque, with Accuseal® tap

Product Code	Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402553	● 32mm	Accuseal® Tap	10	3756Q	60
402554	● 38mm	Accuseal® Tap	10	3756Q	60
402555	● 45mm	Accuseal® Tap	10	3756Q	60
402556	● 57mm	Accuseal® Tap	10	3756Q	60



Small, transparent, with bendable tap

Product Code	Base plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402543	● 32mm	Bendable Tap	10	3728F	60
402544	● 38mm	Bendable Tap	10	3728F	60
402545	● 45mm	Bendable Tap	10	3728F	60
402546	● 57mm	Bendable Tap	10	3728F	60



Standard, transparent, with bendable tap

Product Code	Base plate Size	Description	Pack Size	SAS Code	Monthly Allowance
402548	● 32mm	Bendable Tap	10	3727E	60
402549	● 38mm	Bendable Tap	10	3727E	60
402550	● 45mm	Bendable Tap	10	3727E	60
402551	● 57mm	Bendable Tap	10	3727E	60

TWO PIECE  
Mechanical Coupling

# SUR-FIT® AutoLock®

SUR-FIT® AutoLock® is a locking two piece system which locks with a light touch and an audible click. The pouch can be rotated whilst locked in place.



## Flat Base Plates



Durahesive® flexi base plate with flexible tape collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
401306	●	45mm	12-39mm	Cut to Fit	5	3779X	20
401307	●	57mm	12-51mm	Cut to Fit	5	3779X	20

Stomahesive® flexi base plate with flexible tape collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
401429	●	35mm	12-29mm	Cut to Fit	5	3777T	30
401430	●	45mm	12-39mm	Cut to Fit	5	3777T	30
401431	●	57mm	12-51mm	Cut to Fit	5	3777T	30

## Closed Pouches



Standard, opaque, *without filter*



Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
401677	●	35mm	With Filter	15	3780Y	90
401678	●	45mm	With Filter	15	3780Y	90
401679	●	57mm	With Filter	15	3780Y	90

Colour  
Key

## Drainable Pouches with Clip



Standard, opaque, *without filter*

Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
401434	●	35mm	12-29mm	Without Filter	10	3781B	60
401435	●	45mm	12-39mm	Without Filter	10	3781B	60
401436	●	57mm	12-51mm	Without Filter	10	3781B	60



Standard, transparent, *without filter*

Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
401440	●	45mm	12-39mm	Without Filter	10	3781B	60
401441	●	57mm	12-51mm	Without Filter	10	3781B	60

## Urostomy Pouches



Standard, opaque, with Accuseal® tap

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
401462	●	35mm	Accuseal® Tap	10	3747F	60
401463	●	45mm	Accuseal® Tap	10	3747F	60
401464	●	57mm	Accuseal® Tap	10	3747F	60



Standard, transparent, with Accuseal® tap

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
401452	●	35mm	Accuseal® Tap	10	3747F	60
401453	●	45mm	Accuseal® Tap	10	3747F	60
401460	●	57mm	Accuseal® Tap	10	3747F	60



# Two Piece

Adhesive Coupling

**Esteem synergy®** ..... 46

**Esteem synergy® with Convatec Mouldable Technology™**

Base Plates with Firm Convexity ..... 47

Flat Base Plates ..... 47

**Esteem synergy®**

Flat Base Plates ..... 48

Closed Pouches ..... 49

Drainable Pouches ..... 50

Our Products

ONE PIECE

TWO PIECE  
Mechanical Coupling

TWO PIECE  
Adhesive Coupling

Accessories

Urine Drainage  
Systems

Product  
Code Index

# ESTEEM *synergy*®

TWO PIECE  
Adhesive Coupling



Esteem *synergy*® utilises an innovative Adhesive Coupling Technology™ that connects the pouch to the base plate offering an alternative to mechanical coupling systems.

Esteem *synergy*® is a ConvaTec innovation that is a skin friendly and versatile two-piece system, the provides the low profile, comfort, and flexibility of a conventional one-piece system.

Esteem *synergy*® has a flat, flexible adhesive foam ring on the pouch that securely adheres to a clear, flat, flexible plastic “Landing Zone” flange on the base plate.

## Esteem *synergy*® Colour Key & Sizing Guide

Match the colour of the dot on the base plate box with the colour of the dot on the pouch box.

		
35 mm	48 mm	61 mm

## Base Plates with Firm Convexity and ConvaTec Mouldable Technology™



Durahesive® base plate with flexible tape collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
409268	●	35mm	13-22mm	Mouldable	10	9922F	20
409269	●	48mm	22-33mm	Mouldable	10	9922F	20
409270	●	61mm	33-45mm	Mouldable	10	9922F	20

## Flat Base Plates with ConvaTec Mouldable Technology™



Stomahesive® base plate with hydrocolloid collar



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
411657	●	48mm	13-22mm	Mouldable	10	9786C	20
411658	●	48mm	22-33mm	Mouldable	10	9786C	20
411659	●	61mm	33-45mm	Mouldable	10	9786C	20

## Flat Base Plates



Stomahesive® base plate with flexible hydrocolloid collar



Product Code		Base plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
409288	●	35mm	13-35mm	Cut to Fit	10	9908L	30
409289	●	48mm	13-48mm	Cut to Fit	10	9908L	30
409290	●	61mm	13-61mm	Cut to Fit	10	9908L	30



Stomahesive® base plate with flexible hydrocolloid collar

Pre-cut

Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
403964	●	35mm	19mm	Pre-cut	10	9908L	30
403965	●	35mm	22mm	Pre-cut	10	9908L	30
403966	●	35mm	25mm	Pre-cut	10	9908L	30
403967	●	35mm	28mm	Pre-cut	10	9908L	30
403968	●	35mm	32mm	Pre-cut	10	9908L	30
403969	●	35mm	35mm	Pre-cut	10	9908L	30
403970	●	48mm	38mm	Pre-cut	10	9908L	30
403971	●	48mm	41mm	Pre-cut	10	9908L	30
403972	●	48mm	45mm	Pre-cut	10	9908L	30



Full Stomahesive® base plate



Product Code		Base Plate Size	Stoma Size	Description	Pack Size	SAS Code	Monthly Allowance
405467	●	35mm	13-35mm	Cut to Fit	10	3877C	30
405468	●	48mm	13-48mm	Cut to Fit	10	3877C	30
403948	●	61mm	13-61mm	Cut to Fit	10	3877C	30



## Closed Pouches



Small, opaque, with filter (2 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416777	●	35mm	Small Opaque	30	5684E	90
416781	●	48mm	Small Opaque	30	5684E	90



Standard, opaque, with filter (2 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416775	●	35mm	Standard Opaque	30	5684E	90
416779	●	48mm	Standard Opaque	30	5684E	90
416783	●	61mm	Standard Opaque	30	5684E	90

Small, transparent, with filter (1 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416776	●	35mm	Small Transparent	30	5684E	90
416780	●	48mm	Small Transparent	30	5684E	90

Standard, transparent, with filter (1 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416774	●	35mm	Standard Transparent	30	5684E	90
416778	●	48mm	Standard Transparent	30	5684E	90
416782	●	61mm	Standard Transparent	30	5684E	90

Drainable Pouches with InvisiClose® and Lock-it Pocket™



Standard, opaque, with filter (2 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416796	●	35mm	Standard Opaque	10	5686G	30
416799	●	48mm	Standard Opaque	10	5686G	30
416902	●	61mm	Standard Opaque	10	5686G	30



Standard, transparent, with filter (1 sided fabric)

Product Code		Base Plate Size	Description	Pack Size	SAS Code	Monthly Allowance
416795	●	35mm	Standard Transparent	10	5686G	30
416798	●	48mm	Standard Transparent	10	5686G	30
416901	●	61mm	Standard Transparent	10	5686G	30



## Accessories and Supplementary Products

Niltac™ Adhesive Remover .....	52
ConvaCare™ Adhesive Remover .....	52
Sillesse™ Protective Barrier Wipe and Spray.....	53
ConvaCare™ Protective Barrier Wipes .....	53
Diamonds™ Super Absorbent Gelling Sachets.....	54
ease™ Strips .....	54
Stomahesive® Seal .....	55
Full Stomahesive® Skin Barrier.....	56
Stomahesive® Protective Powder .....	56
Stomahesive® Paste .....	56
Orabase® Protective Paste.....	56
Belt .....	57
Clips .....	57

# SKIN Care

Looking after peristomal skin is an essential part of stoma management. The ConvaTec Skincare range includes second generation technology in Niltac™ Adhesive Remover and Silesse™ Liquid Skin Barrier.

## Niltac™ Adhesive Remover

Niltac™ is a Sting Free, fast drying adhesive remover that does not require any extra cleansing prior to the application of the next device, tape or adhesive dressing.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
420788	TR102 Adhesive Remover Wipes	30 wipes	80006W	60
420787	TR101 Adhesive Remover Spray	50ml spray can	80005T	2
420798	TR1015 Adhesive Remover Spray – Sterile	150ml spray can	N/A	N/A

## ConvaCare® Adhesive Remover

ConvaCare® Adhesive Remover Contains coconut oil and lemon essence which effectively removes adhesive residue and assists in the easy removal of adhesive devices, tapes and dressings. ConvaCare® moisturises the skin and does not contain alcohol. Skin must be thoroughly washed with water and dried before next device, tape or dressing application.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
37443	ConvaCare® Adhesive Remover Wipes	100	3522J	60



## Sillesse™ Sting Free Barrier Wipe and Spray

Sillesse™ provides a flexible, breathable skin barrier that protects the skin from chemical and mechanical damage commonly caused by adhesives and/or corrosive body fluids.

Product Code		Pack Size	SAS Code	Monthly Allowance
420789	TR103 Sting Free Barrier Wipe	30	N/A	N/A
420790	TR104 Sting Free Barrier Spray	50ml spray can	80009B	6 Annually
420797	TR109 Sting Free Barrier Spray	28ml pump spray	80031E	1
420793	TR107-5 Sting Free Barrier Foam Applicator 1ml	5	N/A	N/A
420795	TR108-5 Sting Free Barrier Foam Applicator 3ml	5	N/A	N/A

## ConvaCare® Protective Barrier Wipe

ConvaCare® Protective Barrier Wipe is a non-water soluble skin barrier that forms a thin protective film on the skin which allows better adhesion of adhesive devices, tapes and dressings. ConvaCare® protective barrier formula is designed not to peel or chip. Contains alcohol.

Product Code		Pack Size	SAS Code	Monthly Allowance
37444	ConvaCare® Barrier Wipes	100	3502H	700 Annually



Diamonds™ Super Absorbent Gelling Sachets



Diamonds™ contains ActiveOne™, a unique odour control formula intended for patients with loose stool to use every day. These super absorbent gelling sachets also provide discretion for Ileostomy patients looking to reduce noise and control the profile of their pouch by being able to manipulate the pouch contents into a flattened shape. Gelling output with Diamonds™ also reduces the possibility of back-flow around the edge of the base plate.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
420791	Diamonds™ Gelling and Odour Control Sachets	100 Sachets/jar	80103Y	7 per year with an authority

ease™ Strips



Thin flexible strips designed to adapt to your body and move with you. ease™ strips overlap the edge of the skin barrier for extra security. The corner or edge of the wafer could lose stickiness, when worn for long periods. The strips help maintain the contact with the skin. They can also be used to prevent the edge of the wafer sticking to clothing after swimming or bathing.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
422160	ease™ Thin strips 3cm x 11 cm	60	80188K	60
422164	ease™ Thin curve 3cm x 9 cm	60	80188K	60

## Stomahesive® Seal



ConvaTec Stomahesive® Seal is a mouldable, moisture absorbing skin barrier designed to protect peristomal skin from effluent. Designed to aid the adherence of a pouch/dressing to the skin around the stoma, Stomahesive® seals are suitable for patients with an ileostomy, urostomy or colostomy.

The seals are single use devices and should be replaced on application of each new ostomy device or dressing.

Product	Description for Use
48mm (Thin)	Ideal skin protection from leakage for all types, shapes and sizes of stoma.
48mm (Standard)	Thicker seal, ideal for extra skin protection from significant leakage problems.
98mm (Large)	Large seal designed to cover a greater area around the stoma for added skin protection from leakage.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
413503	48mm Standard	10	80008Y	30
413504	48mm Thin	10	80008Y	30
413505	98mm Large	10	80008Y	30



## Full Stomahesive® Skin Barrier

Used to protect peristomal skin, Stomahesive® Skin Barrier can be cut and moulded for use in difficult to manage stomas, provide a base for tape to adhere to and can also be used to 'window frame' large wounds.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
21712	10 x 10cm	5	3532X	30
21715	20 x 20cm	3	3580K	30



## Stomahesive® Protective Powder

Protective barrier against excoriating discharge on weeping areas. Increases wear time by absorbing moisture helping to prevent skin irritation.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
25510	28gm Bottle	1	3511T	2



## Stomahesive® Paste

Used as a 'crease or gap filler' Stomahesive Protective Paste® provides a flat pouching surface with the aim of preventing leakage under the base plate.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
183910	57gm Tube	1	3503J	3



## Orabase® Protective Paste

To protect and soothe the minor irritations of the skin.

Product Code	Description	Pack Size	SAS Code	Monthly Allowance
149710	Protective Paste 5gm	1	N/A	N/A
149715	Protective Paste 15gm	1	3787H	2



Belt

For use with convex appliances to assist in improving protrusion of flush/retracted stomas. ConvaTec Ostomy Belts are also designed to assist in providing Ostomates with extra security.



Product Code	Description	Pack Size	SAS Code	Monthly Allowance
175507	White Belt, up to 106cm	1	3792N	4 per year

Clips

For use with drainable appliances.



Product Code	Description	Pack Size	SAS Code	Monthly Allowance
175635	Straight Clip	10	3760X	1 Box 10 Clips per year



# Urine Drainage Systems

Unomedical® Urine Drainage Bags.....	59
Unomedical® Intermittent /Nelaton Catheters.....	60

Unomedical® Urine Drainage Bags

Product Code	Description	Pack Size (units)	SAS Code	Monthly Allowance
370050C	Urine Drainage Bag With Non-Return Valve, Drainable, 2 Litre Capacity, 100cm Inlet Tube, Push/Pull Tap, Non-Sterile	250	3652F	5
373050C	Urine Drainage Bag With Non-Return Valve, Drainable, 2 Litre Capacity, 90cm Inlet Tube, Cross Tap, Non-Sterile	250	3653G	5
373130C	Urine Drainage Bag With Non-Return Valve, Drainable, 2 Litre Capacity, 120cm Inlet Tube, Cross Tap, Non-Sterile	250	3653G	5
38 155 1S-20	A4 Closed System Urine Drainage Bag, 2 Litre Capacity With Non Return Valve and Kombikon™ Needle-Free Sample Port, 120cm Inlet Tube, Cross Outlet Tap, Sterile	50	3800B	5



## Unomedical® Intermittent / Nelaton Catheters – Male Length (40cm)

Product Code	Description	Pack Size (units)	SAS Code	Monthly Allowance
001.06.040	Rounded Tip, Firm, 6 FG	100	3671F	100
001.08.040	Rounded Tip, Firm, 8 FG	100	3671F	100
001.10.040	Rounded Tip, Firm, 10 FG	100	3671F	100
001.12.040	Rounded Tip, Firm, 12 FG	100	3671F	100
001.14.040	Rounded Tip, Firm, 14 FG	100	3671F	100
001.16.040	Rounded Tip, Firm, 16 FG	100	3671F	100
001.18.040	Rounded Tip, Firm, 18 FG	100	3671F	100
001.20.040	Rounded Tip, Firm, 20 FG	100	3671F	100
003.10.040	Rounded Tip, Soft, 10 FG	100	3671F	100
003.12.040	Rounded Tip, Soft, 12 FG	100	3671F	100
003.14.040	Rounded Tip, Soft, 14 FG	100	3671F	100
003.16.040	Rounded Tip, Soft, 16 FG	100	3671F	100
003.18.040	Rounded Tip, Soft, 18 FG	100	3671F	100
003.20.040	Rounded Tip, Soft, 20 FG	100	3671F	100
002.06.040	Rounded Tip, Firm, 6 FG, Double Wrapped	100	3671F	100
002.08.040	Rounded Tip, Firm, 8 FG, Double Wrapped	100	3671F	100
002.10.040	Rounded Tip, Firm, 10 FG, Double Wrapped	100	3671F	100
002.12.040	Rounded Tip, Firm, 12 FG, Double Wrapped	100	3671F	100
002.14.040	Rounded Tip, Firm, 14 FG, Double Wrapped	100	3671F	100
002.16.040	Rounded Tip, Firm, 16 FG, Double Wrapped	100	3671F	100
002.18.040	Rounded Tip, Firm, 18 FG, Double Wrapped	100	3671F	100
002.20.040	Rounded Tip, Firm, 20 FG, Double Wrapped	100	3671F	100
006.10.040	Rounded Tip, Soft, 10 FG, Double Wrapped	100	3671F	100
006.12.040	Rounded Tip, Soft, 12 FG, Double Wrapped	100	3671F	100
006.14.040	Rounded Tip, Soft, 14 FG, Double Wrapped	100	3671F	100
006.16.040	Rounded Tip, Soft, 16 FG, Double Wrapped	100	3671F	100
006.18.040	Rounded Tip, Soft, 18 FG, Double Wrapped	100	3671F	100
006.20.040	Rounded Tip, Soft, 20 FG, Double Wrapped	100	3671F	100
03007182*	Tiemann Tip (Angled, Tapered Tip) 10 FG	100	9962H	100
03008182*	Tiemann Tip (Angled, Tapered Tip) 12 FG	100	9962H	100
03009182*	Tiemann Tip (Angled, Tapered Tip) 14 FG	100	9962H	100
03011182*	Tiemann Tip (Angled, Tapered Tip) 16 FG	100	9962H	100

\*Requires authorisation form signed by health care professional or Stomal Therapy Nurse.

# Product Code Index

03007182	20917	125139	401435	402549	411449
03008182	20922	125140	401436	402550	411450
03009182	21712	125142	401440	402551	411451
03011182	21715	125143	401441	402553	411657
001.06.040	22757	125144	401452	402554	411658
001.08.040	22758	125145	401453	402555	411659
001.10.040	22759	125262	401460	402556	411801
001.12.040	22760	125264	401462	403948	411803
001.14.040	22761	125265	401463	403964	411805
001.16.040	22762	125266	401464	403965	411807
001.18.040	22771	125303	401574	403966	411817
001.20.040	25510	125304	401575	403967	411819
002.06.040	37443	125305	401576	403968	411821
002.08.040	37444	149710	401577	403969	411823
002.10.040	64927	149715	401677	403970	413177
002.12.040	64928	175507	401678	403971	413178
002.14.040	64929	175635	401679	403972	413179
002.16.040	64930	175777	402510	404592	413180
002.18.040	64931	175778	402511	404593	413181
002.20.040	125022	175779	402516	404594	413182
003.10.040	125023	175780	402517	405467	413183
003.12.040	125024	175781	402521	405468	413184
003.14.040	125025	175792	402522	409268	413185
003.16.040	125028	175793	402523	409269	413186
003.18.040	125029	175794	402524	409270	413187
003.20.040	125030	175795	402525	409288	413188
006.10.040	125031	175796	402526	409289	413437
006.12.040	125032	175797	402527	409290	413438
006.14.040	125033	183910	402528	410293	413439
006.16.040	125034	370050C	402529	410294	413503
006.18.040	125035	373050C	402531	410295	413504
006.20.040	125036	373130C	402532	410296	413505
38 155 1S-20	125037	400949	402533	410297	413509
	125132	400965	402534	411385	413510
	125133	401306	402535	411386	413511
	125134	401307	402543	411387	413515
	125135	401429	402544	411388	413516
	125136	401430	402545	411401	413517
	125137	401431	402546	411402	413521
	125138	401434	402548	411403	413522

# Product Code Index

413523	416902	421626	421822	422528
416414	416908	421627	421823	422529
416416	420787	421628	421825	422530
416419	420788	421629	421826	422531
416422	420789	421631	421827	422532
416712	420790	421632	421828	422533
416718	420791	421633	421829	422534
416721	420793	421634	421830	422535
416724	420795	421635	421831	422536
416728	420797	421636	421832	422537
416732	420798	421637	421833	422538
416736	421033	421676	421834	422539
416738	421034	421678	421835	422540
416739	421035	421680	421836	422541
416740	421039	421682	421837	422542
416741	421040	421686	421838	422543
416742	421041	421687	421839	422544
416743	421074	421688	421840	422545
416744	421075	421689	421843	422546
416745	421076	421746	421844	422547
416746	421454	421747	421846	422548
416747	421456	421748	421848	422549
416748	421458	421749	421849	422550
416749	421460	421794	421850	422551
416750	421462	421795	421858	422552
416751	421464	421796	421859	422553
416774	421612	421797	421860	422554
416775	421613	421798	421861	422555
416776	421614	421799	421862	422556
416777	421615	421800	422158	422557
416778	421616	421801	422159	422558
416779	421617	421802	422160	422559
416780	421618	421803	422164	422560
416781	421619	421804	422360	422561
416782	421620	421805	422361	422562
416783	421621	421816	422362	422563
416795	421622	421817	422363	422564
416796	421623	421818	422364	422565
416798	421624	421819	422365	650803
416799	421625	421820	422527	

**References:**

1. ConvaTec data on file. In vitro testing.
2. Maria Teresa Szewczyk, MD, PhD; Grazyna Majewska, RN, ETN; Mary V. Cabral, MS, FNP-BC, CWOCN-AP; and Karin Hölzel-Piontek, RN; The Effects of Using a Moldable Skin Barrier on Peristomal Skin Condition in Persons with an Ostomy: Results of a Prospective, Observational, Multinational Study, *Ostomy Wound Management* 2014;60(12):16–26.
3. Hoeflok J. Use of Convexity in Pouching - A Comprehensive Review. *J Wound Ostomy Continence Nurs.* 2012; 40(5):506-512
- 4 Kruse and Storling, Considering the benefits of a new stoma appliance: a clinical trial, *BJN*, Vol 24, issue sup 22, pS12-18
5. Boyd. K. Clinical Protocols for Stoma Care - Use of convex appliances. *Nursing Standard.* 2004; 18(20): 37-38