



Visit our web site:
b2b.convatec.com

Infusion Care Product Range

Conventional Infusion Sets


Steel needle and soft cannula
for manual insertion




comfort™



comfort™ short



contact™ detach



Quick-set™



inset™



inset™ 30



Mio™ Advance



Medtronic Extended

Single Use Sarter Technology

Soft cannula infusion sets with integrated
insertion device

Injection port

Soft catheter with automatic insertion











I-port Advance™

Product Overview

Infusion Sets

Injection Ports

	comfort™ / comfort™ short	Quick-set™	inset™	inset™ 30	Mio™ Advance	Medtronic Extended	contact™ detach	i-port Advance™
								
Cannula								
Soft	●	●	●	●	●	●		●
Stainless steel							●	
Insertion								
90 degrees		●	●		●	●	●	●
Angled insertion	●			●				
Disconnect								
At site	●	●	●	●	●	●		
Separate site							●	
Integrated inserter			●	●	●	●		●
Tubing length								
18 inches (46 cm)	●	●	●		●		●	
23 inches (60 cm)	●	●	●	●	●	●	●	
31 inches (80 cm)	●	●			●	●	●	
43 inches (110 cm)	●	●	●	●	●	●	●	
Cannula length (mm)	17 and 13	6 and 9	6 and 9	13	6 and 9	6 and 9	6, 8, 10	6 and 9

Infusion Care

ConvaTec Infusion Care, is a global leader within design and manufacturing of infusion set technologies and insertion devices for insulin pump treatment.

ConvaTec has researched and developed insulin infusion sets for the treatment of diabetes since the early 1980s, and has production facilities in Denmark and Mexico.

At ConvaTec we strive to create better tailor-made solutions for the treatment of diabetes. Both improvement of current and development of future infusion sets are essential parts of our business.

Our success is based on experience, knowledge and a close cooperation with a global network of customers, professionals and the people using our products.

Unomedical a/s is the legal manufacturer of the infusion care product range.

For further information about ConvaTec and our products, please visit our web site at b2b.convatec.com or contact your local infusion set supplier.

