



ConvaTec

Ostomy Care

PRODUCTS, CARE *and* SUPPORT  
SHAPED AROUND YOUR LIFE

# MOLDABLE Skin Barriers

ConvaTec Moldable Technology™ Skin Barriers are designed to hug the stoma for a more custom and secure fit.

## Features and Benefits

- ✓ Reduced application steps—fewer steps than traditional cut-to-fit skin barriers
- ✓ No cutting needed
- ✓ Memory – flexible materials gently expand and contract with the stoma
- ✓ Turtlenecking – the adhesive gently swells when in contact with output to help protect the skin and secure the seal.




**95%** of people who started on ConvaTec Moldable Technology kept their skin healthy¹


**86%** who switched to ConvaTec Moldable Technology improved their skin condition¹



Roll, mold, and hold to the stoma's shape and size



Turtlenecks around the stoma to protect against leaks



Rebounding Memory Technology™ - rolls back to achieve a custom fit

## NATURA® DURAHEESIVE® SKIN BARRIERS

### Moldable | Qty = 10 | HCPCS A4409

Flange Size	Stoma Size	Acrylic Tape Collar WHITE	Hydrocolloid Tape Collar WHITE
1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	413417	411800
1 3/4" (45mm)	7/8 - 1 1/4" (22-33mm)	413418	411802
2 1/4" (57mm)	1 1/4 - 1 3/4" (33-45mm)	413419	411804
2 3/4" (70mm)	1 1/4 - 1 3/4" by 1 3/4 - 2 1/8" (33x45mm by 45x56mm) Oval shape	413420	411806

### Moldable | Accordion Flange | Qty = 10 | HCPCS A4409

Flange Size	Stoma Size	Hydrocolloid Tape Collar TAN
2 1/4" (57mm)	1/2 - 7/8" (13-22mm)	421039
2 1/4" (57mm)	7/8 - 1 1/4" (22-33mm)	421040
2 3/4" (70mm)	1 1/4 - 1 3/4" (33-45mm)	421041

### Moldable | Convex | Qty = 10 | HCPCS A4407

Flange Size	Stoma Size	Acrylic Tape Collar WHITE
1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	404592
1 3/4" (45mm)	7/8 - 1 1/4" (22-33mm)	404593
2 1/4" (57mm)	1 1/4 - 1 3/4" (33-45mm)	404594

## NATURA® STOMAHESIVE® SKIN BARRIERS

### Moldable | Qty = 10 | HCPCS A4414

Flange Size	Stoma Size	Acrylic Tape Collar WHITE	Hydrocolloid Tape Collar WHITE
1 3/4" (45mm)	1/2 - 7/8" (13-22mm)	413421	411801
1 3/4" (45mm)	7/8 - 1 1/4" (22-33mm)	413422	411803
2 1/4" (57mm)	1 1/4 - 1 3/4" (33-45mm)	413423	411805
2 3/4" (70mm)	1 1/4 - 1 3/4" by 1 3/4 - 2 1/8" (33x45mm by 45x56mm) Oval shape	413424	411807

### Moldable | Accordion Flange | Qty = 10 | HCPCS A4414

Flange Size	Stoma Size	Hydrocolloid Tape Collar TAN
2 1/4" (57mm)	1/2 - 7/8" (13-22mm)	421033
2 1/4" (57mm)	7/8 - 1 1/4" (22-33mm)	421034
2 3/4" (70mm)	1 1/4 - 1 3/4" (33-45mm)	421035



Enroll your patients to get the support and resources they need.

Contact your me+™ concierge at 1-877-585-0470 or visit [mepluscare.com](https://mepluscare.com)

## ESTEEM SYNERGY® SKIN BARRIERS

### Moldable | Stomahesive® | Qty = 10 | HCPCS A4414

Coupling Size	Stoma Size	Hydrocolloid Tape Collar WHITE
Small Stoma Opening	1/2 - 7/8" (13-22mm)	411657
Medium Stoma Opening	7/8 - 1 1/4" (22-33mm)	411658
Large Stoma Opening	1 1/4 - 1 3/4" (33-45mm)	411659

### Moldable | Convex | Durahesive® | Qty = 10 | HCPCS A4407

Coupling Size	Stoma Size	Acrylic Tape Collar WHITE
Small Stoma Opening	1/2 - 7/8" (13-22mm)	409268
Medium Stoma Opening	7/8 - 1 1/4" (22-33mm)	409269
Large Stoma Opening	1 1/4 - 1 3/4" (33-45mm)	409270

## ESTEEM +® POUCHES (One-Piece Systems)

### Moldable | Stomahesive® | Drainable | InvisClose™ Tail Closure | Qty = 10 | Pouch Length = 12"

#### HCPCS A5061 | No Filter

Stoma Size	Transparent
13/16 - 1 3/16" (20-30mm)	413518
1 3/16 - 1 9/16" (30-40mm)	413519
1 9/16 - 1 15/16" (40-50mm)	413520

#### HCPCS A4424 | Filter (Filter covers included)

Stoma Size	Transparent	Opaque w/ Easy-View Window
13/16 - 1 3/16" (20-30mm)	413515	413521
1 3/16 - 1 9/16" (30-40mm)	413516	413522
1 9/16 - 1 15/16" (40-50mm)	413517	413523

### Moldable | Stomahesive® | Closed-End

Qty = 30 | Pouch Length = 8"

#### HCPCS A4416 | Filter (Filter covers included)

Stoma Size	Opaque w/ Easy-View Window
13/16 - 1 3/16" (20-30mm)	413509
1 3/16 - 1 9/16" (30-40mm)	413510
1 9/16 - 1 15/16" (40-50mm)	413511

