

Instruction For Use

Vacuum Control Fingertip

Product Description

A fingertip connects a suction catheter to the connecting tube. An open suction catheter is introduced into the trachea for removal of secretions from the upper and lower respiratory tract by suctioning, thus preventing stagnation of secretions. The fingertip allows for vacuum control. Fingertips are sold sterile.

Intended Use

The fingertip allows for vacuum control and used with Suction Catheter.

Suction catheters with fingertip are intended to be connected to a connecting tube for suction (which is applied between the suction catheter and the vacuum source/mobile unit). The intended use is oro-, nasopharyngeal and tracheobronchial suction of the upper and lower airways to remove excessive secretion. Suction Catheter with Fingertip is an invasive device with respect to body orifices and is indirectly connected to an active medical device.

Instruction for use

Suction Catheters with fingertip are not supplied with an IFU, since they are used by trained healthcare professionals, so can be safely used without instruction.

Precautions

Suction catheters with fingertip are intended for single patient/procedure use only.

Storage Condition

Keep Dry. Store at room temperature

Packing Quantity

100 pieces/box

Legal manufacturer

ConvaTec Limited, Flintshire, United Kingdom

Manufactured by

FE Unomedical Ltd., Minsk Region, Belarus

Imported by/Contact:

ConvaTec (Thailand) Co., Ltd.

87 M.Thai Tower, All Season Place, 9th Floor, Wireless Road, Lumpini, Phatumwan, Bangkok Tel. 02-6147000