



Which infusion set would you use for a SuperHero?



Making infusion easy to handle

neria™ guard infusion set is

all-in-one: soft cannula and insertion aid in one device

simple to use: few user steps and fully automatic insertion.

needle safety: no visible needle and automatic needle retraction



NO VISIBLE NEEDLE



Steps for insertion of neria™ guard*

Follow the below user steps *



1. Remove paper



2. Remove safeguard



3. Position on skin and press button



4. Connect tubing

For more information please visit
www.convatec.com/infusion-care

*Please note that this guide only illustrates some of the steps to follow when inserting the neria™ guard infusion sets and it is not a substitute for reading the instructions for use which comes with the product.

neria is a trademark of Unomedical a/s © ConvaTec Inc. 2018 AP-019665-MM